RAPTOR

SERVER ANALYSIS SOFTWARE

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RAPTOR

SERVER ANALYSIS SOFTWARE

Welcome to RAPTOR SERVER Analysis Software, the analysis software for the [ALEX multiplex test](#). Please log in below to continue.

Log in with Azure AD / Office 365

If your organization integrates the software into its Microsoft Azure/Office 365 user directory (single-sign-on) click the button below. You will be redirected to the Microsoft login screen.

 Log in

Log in with Username & password

Please enter your username and password.

 Log in

Forgot your password? No problem, you can [request a new one](#).

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Contract approval required

Please review the following documents. You are required to confirm that you have read and agree with each document in order to use RAPTOR SERVER:

- ☐ [Bitte lesen](#) (last updated 23.11.2023 18:00)
- ☐ [RAPTOR SERVER IFU](#) (last updated 04.12.2023 11:05)
- ☐ [Terms of Service](#) (last updated 22.11.2023 01:00)

Confirm

By clicking the button, I confirm that I have reviewed these documents and agree.

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
New Measurements					Last 7 days
QC	Product type / Source	QR-Code	Sample Code	Analysis Date	
	Research MAX 45k	99AFG0F7	NK FOX	30/03/2023 20:50	Details Report
	Research MAX 45k	99AFG0F8	N.G.	30/03/2023 20:50	Details Report
	Research MAX 45k	99AFG0F9	E18735	30/03/2023 20:50	Details Report
	Research MAX 45k	99AFG0FA	20220325w56	30/03/2023 20:50	Details Report
	Research MAX 45k	99AFG0EE	HM100038	30/03/2023 20:50	Details Report
	Research MAX 45k	99AFG04D	786118	30/03/2023 20:37	Details Report
	Research MAX 45k	99AFG016	HM100225	30/03/2023 20:37	Details Report
	Research MAX 45k	99AFG0CA	E17685_M	30/03/2023 20:37	Details Report
	Research MAX 45k	99AFG0CC	E17685_M	30/03/2023 20:37	Details Report
	Research MAX 45k	99AFG015	E17685_F	30/03/2023 20:37	Details Report
In total there are 270 new measurements in the selected time period.					
View all new measurements					
Approved Measurements					
QR-Code	Product type / Source	Sample Code	Analysis Date		
02AVJ227	ALEX2	Anna LMacey	28/03/2023 21:17	Details	Report
02ARA379	ALEX2 IX	SM_28032023	28/03/2023 14:29	Details	Report
02ARA393	ALEX2 IX	RG_28032023	28/03/2023 14:26	Details	Report
02AYA31F	ALEX2 MAX 45k	PR_Mar2023_11	21/03/2023 15:05	Details	Report
02AVI3B9	ALEX2 MAX 45k	KF	17/03/2023 14:46	Details	Report
02AY1337	ALEX2 MAX 45k	2846846+CCD	14/03/2023 19:06	Details	Report
02ATL349	ALEX2 IX		22/02/2023 14:21	Details	Report
02AVE383	ALEX2 MAX 45k	AO	03/02/2023 09:41	Details	Report
02AXE29A	ALEX2 MAX 9k	OW	19/01/2023 19:12	Details	Report
02AUY3AE	ALEX2 MAX 45k	BSZ_Ok2022_8	18/10/2022 18:12	Details	Report
View all approved measurements					

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9K-Test-Instrument

Not initialized

[Initialize](#)
[Shutdown](#)
[Monthly Maintenance](#)

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E-Mail address:

Language: English (en)

The language is used to determine the display language.

Culture: English (United States) (en)

The culture is used to format numbers and dates.

Timezone: (UTC+01:00) Amsterdam, Berlin, Bern, Rome, Stockholm, Vienna

Configure your timezone to see date and times in your local time instead of UTC.

Multi-Factor Authentication

[Manage MFA](#)

Use multi-factor authentication (MFA) to significantly reduce the risk of account compromise. It is highly recommended to enable MFA; this adds a second layer of security to your account against phishing attacks or weak passwords.


Configuration

Research spots: Show research spots

If enabled, "research use only" arrays. This does not affect the

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Configure "MADx"

General configuration	
* Name:	MADx
Default timezone:	(UTC+01:00) Amsterdam, Berlin, Bern, R...
Configure a timezone to see date and times in (instead of UTC). Please note that users can choose to use a different timezone.	
Default report template (ALEX2):	ALEX2
Default report template (ALEX3, ALEX Air, ALEX Food):	check 1.21
Default report template (FOX):	FOX (2025)
Default report template (PAX):	PAX
Default report language:	English (en)
Default ALEX MOD:	< No MOD >
Default FOX MOD:	< No MOD >
Default PAX MOD:	< No MOD >
Force MOD:	False
Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.	
Hide Global MOD:	False
If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.	
Display MOD Excluded Spots As Research Use Only:	False
If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".	
MFA:	☐ Force MFA
If this setting is enabled, users are forced to set up an MFA on sign-in, in case they haven't yet.	

MAX device default settings	
* Default Minimum RunId:	10
When a new MAX device is created, this value will be used for its Minimum Run Id	

Tenant-specific batch releases	
Restrict this tenant to only allow the following batches to be scanned or to restrict certain batches. Enter one per line.	
Use:	
☒ Whitelist	
☐ Blacklist	
Whitelisted Batches:	
Enter one batch code including product identifier per line, for example "02AAA"	
Blacklisted Batches:	
Enter one batch code including product identifier per line, for example "02AAA"	

Company Details

* Legal Entity Name:	MacroArray Diagnostics
* Street and Number:	Street 4
* Zip Code:	1230
* City:	New York
* Country:	United States of America
* Phone:	+4318652573
* E-Mail:	office@macroarraydx.com
* Legal Representative:	Dr. Christian Harwanegg
Website:	
Vat Number:	
Company Register Number:	

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Tenant-specific batch releases

Restrict this tenant to only allow the following batches to be scanned or to restrict certain batches. Enter one per line.

Use:

☒ **Whitelist**

☐ **Blacklist**

Whitelisted Batches:

02AAA
02ABC
02BAA

Enter one batch code including product identifier per line, for example "02AAA"

Blacklisted Batches:

Enter one batch code including product identifier per line, for example "02AAA"

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You requested a new measurement from ImageExplorer "IX Trainingsraum". Please wait until the measurement is completed ...

Invalid batch detected! The measurement can't be processed because its batch is blocked.

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
Login Credentials

Type: Username and Password

Please select if the user should be validated against Azure AD (Office 365), or if you want to specify username and password yourself. It is recommended to use Azure AD whenever possible. This allows for increased security by supporting two-factor authentication. Also, it is more convenient for the user, because he or she does not have to remember a new username and password (single-sign-on).

Username:

It is a good idea to use the users e-mail-address as his username, because the username needs to be unique across the entire system.

Password:

Automatically send invitation e-mail

Access Rights

Preconfigured groups

Use these pre-configured groups to automatically assign a set of access rights to this user. You can still add or remove specific access rights.

Reset Lab User Super Lab User Support/Distributor User System Administrator

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Tenants

This user is connected to the following tenants and has the following access rights:

☐ Tenants	Can share measurements	Can edit measurements	Can edit personal data	Can approve measurements	Can manage basic tenant configuration	Can manage API keys	Can manage ImageXplorers	Can view agent connections	Can manage users	Can manage users of Result Sharing Portal	Can manage report templates	Can manage measurement data completion	Can unapprove tenant specific batch releases	System Administration	Can use a custom calibration string	Can manage Run Configuration	Can skip
☐ Another MADx Tenant	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
☐ MADx Tenant	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒

Additionally this user can be added to the following tenants by checking the tenant name and hitting save afterwards. User rights w applied from current open tenant user.

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Additionally this user can be added to the following tenants by checking the tenant name and hitting save afterwards. User right applied from current open tenant user.

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Tenants

This user is connected to the following tenants and has the following access rights:

 Edit selected tenants rights

☐ Tenants		Can share measurements	Can edit measurements	Can edit personal data	Can approve measurements	Can manage basic tenant configuration	Can manage API keys	Can manage Imagey
☒	Another MADx Tenant	☒	☒	☒	☒	☒	☒	☒

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ImageExplorers

Use this to configure ImageExplorers, or add new ones. Every ImageExplorer needs to be known to Raptor Server, otherwise it won't be able to connect.

Manage ImageExplorers

Add new ImageExplorer

General information

Name:

Key:

This key is required when installing an Agent.

ImageExplorer Serialnumber:

Serialnumber of connected ImageExplorer.

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Details


Barcode: 02AE3488
Sample Code: 12
Requisition Code:
Sampling Date: 01.10.2021
Assay Date: 04.10.2021
Patient Code: TE
Patient Name: Test
Date of Birth: 17.05.2010
Info: Info
Doctor Info: Dr. A
Choose a doctor... Insert doctor info
Tags: Tag Info
MOD: < No MOD >
Report Template: Version 2 (ALEX2 only) - Deutsch (de)

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



Analysis Log

```
'PreProcessedPixels' took 128ms.
Using definition '02_ALEXIL_IVD_v1.6.0'.
'ApplyArrayDefinition' took 74ms.
Parsing calibration string ...
'ApplyDefinitions' took 75ms.
'FindGrid' took 1ms.
'AlignFeatures' took 862ms.
'CropImage' took 1ms.
'QuantifyFeatures' took 231ms.
'SubtractBackground' took 681ms.
Parsing calibration string ...
Using calibrators from QR-Code.
Total IgG (tIgG) normalized by calculating 8,330999755859375E-07x*+0,05697999954223633x+-38,69269921875
(x=8515), which results in 515.
'ReQuantify' took 57ms.
```

QC Log

These are the values used for QC (ALEX v2): Curve Correction=0,9997370867894242, tIgG=8515, GD1=28085 (508=185.38), GD2=29369 (508=185.28), GD3=28569 (508=147.53), GD4=28305 (508=128.1). Reference values are as follows: Guide Dots must have a value above 5000; each Guide Dot spot must exist and be unique; Guide Dot spot must not be set to bad. Average background must be between 4000 and 10000. QC (ALEX v2) succeeded.

Background Log

Calculated background statistics for 164 pixels (from a total of 182, skipped smallest 9, and highest 9 pixels). The average is 9.419, the minimum 7973 and the maximum 11792 (minimum and maximum is calculated including skipped pixels). This results in an average difference of 628, which is an SD of 49,17 and a CV of 0,52%. For comparison, there is an overall total of 282944 pixels. The overall average is 14,071, the minimum 6,012 and the maximum 64,965.

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
	All allergens	☑ All allergens
	Only positive allergens	
Tissue &	ire	kUA/L
Pollen Inhalation	E Aca m	≤ 0.10
Mites	E Aca s	7.11

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
New Measurements

Search by QR-Code, patient, sample, tag, analysis date, approval date ...

Selected measurements

☐	Source	QC	Tags	Product type
☐	IX	✓ OK		ALEX2
☒	IX	✓ OK		ALEX2
☒	IX	✓ OK		ALEX2
☒	IX	✓ OK		ALEX2

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Delete

Re-process measurements

You are about to re-process the following measurements. undone, the existing analysis result will be irreversibly over

QR-Code	Sample Code
02APK3C8	Dummy8
02APK258	Dummy1
02APK07F	Dummy13

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New Measurements

Search by QR-Code, patient, sample, tag, analysis date, approval date ...

	Source	QC	Tags	Product type
☐	Recaptured	✓ OK	recaptured	ALEX1
☐	⚠ Recaptured	✓ OK	recaptured	ALEX1
☐	QR-Code was scanned more than once.		recaptured	ALEX1

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RAPTOR New measurements Approved measur

New Measurements

Search by QR-Code, patient, sample, tag, analysis date, approval date ...

Selected measurements

	Source	QC	Tags	Product type	QR-Code
☒	IX	✓ OK		ALEX2	02AMC180
☒	IX	✓ OK		ALEX2	02AMC078
☒	IX	✓ OK		ALEX2	02AKV2D9

Selected measurements

- Download CSV format
- Download Excel format
- Approve**
- Update interpretation
- Re-Process
- Move
- Delete

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
Approve measurements

You are about to approve the following measurements.

There are measurements with one or more warnings. Please make sure to double-check these measurements before approval.

QR-Code	Sample Code	Analysis Date
02AKV035		28.04.2022 11:14

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New measurements
Approved measurements
Start m

Approved Measurements

Search by QR-Code, patient, sample, tag, analysis date ...

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New measurements
Approved measurements
Start measurement
Automats
MADx Testtenant

Exports

Certain types of data exports (for example, the export of many PDF reports at once) take a few minutes to prepare for download. Here download your requested exports. Please note that you can keep working, all exports are done in the background and you can download when they are completed.

Export requested	# of measurements included	Export completed
8/16/2021 11:07 AM	4	0 / 4
8/16/2021 11:06 AM	4	8/16/2021 11:06 AM
8/16/2021 11:04 AM	3	8/16/2021 11:04 AM

Download

Administration
Tenant Admin
Exports
Personal Settings
Log out
Support
Imprint
Download
Download

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MADx Tenant

Administration

- Tenant Admin
- System Admin
- Exports**
- Personal Settings
- System Notifications
- Log out

Support

- Imprint

Exports

Certain types of data exports (for example, the export of many PDF reports at a time) are possible. Please note that you can keep working, all exports are done in the background.

Export multiple measurements

Export requested	# of measurements included
06.07.2022 08:52	3

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Export multiple measurements

Please note that **Report** and **ASTM** export is only available for certain product types.

Product Type

ALEX2

Analysis Date From

TT . MM . JJJJ

Tag

☐ **Export CSV**

QR Code

Analysis Date

TT . MM . J

Approved or not

New measurements

☐ **Export Excel**

Search Measurements

Cancel

Export multiple measurements

Please note that **Report** and **ASTM** export is only available for certain product types.

Product Type

ALEX2

QR Code

Analysis Date From

TT . MM . JJJJ

Analysis Date To

TT . MM . JJJJ

Tag

Approved by

Approved by

use , to separate tags

☐ Export CSV

☐ Export Report

 Start export ...

← Back

122 Measurements found that match your criteria

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Edit report template

Your changes have been saved successfully.

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
Configure report template "New report template"

Field	Value	Description
Name:	New report template	This is the name of the template which shows up in the dropdown list while generating a report.
Base Template:	Version 2 (ALEX2 only)	A base template which specifies basic functionality must be chosen.
Preferred language:	English (United States) (en-us)	A default language can be specified which will be used if there isn't yet a language stored to a measurement's report.

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PDF naming

Setting Type: Manual

The filename of the report PDF can be customized. Either built automatically based on your selection below or manually adjusted in the order of your needs.

Name Pattern: Report_{PatientId}_{MeasurementQrCode}

Available Elements: {PatientId} ; {PatientName} ; {SampleId} ; {MeasurementQrCode} ; {DateAndTime}

--

PDF naming

Setting Type: Automatic

The filename of the report PDF can be customized. Either built automatically based on your selection below or manually adjusted in the order of your needs.

Include Patient ID: ☒

Include Patient Name: ☒

Include Sample ID: ☐

Include Measurement QR-Code: ☐

Include Date & Time: ☒

Footer

Footer text:

Will be printed in the footer if provided and rendered as HTML.

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Configure report template "New report template"

The screenshot shows a configuration window for a report template. It has four tabs: General, Design, Header, and Content. The 'Design' tab is active. Under 'PDF format', 'A4' is selected with a radio button, and 'Letter' is an option. A description states: 'Specifies the paper size of the output PDF.' Under 'Font', a dropdown menu shows 'Arial'. A description states: 'The font family of the whole report texts can be adjusted.' Under 'Background image', there is an empty text field and a 'Browse' button. A description states: 'If set, this image is used as background on the first page of the PDF.' Below the field, it says 'Maximum file size 1 MB'.

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Logo (left): Browse

Maximum file size 1 MB

This logo will appear on the **top left corner** of the first page of the report.

Text (center):

This text will appear either on the left, right, or center of the top of the first page (depending on your configured logos). Supports the placeholders `{{ PATIENT_PORTAL_QRCODE }}` and `{{ PATIENT_PORTAL_URL }}` to dynamically insert measurement information.

Logo (right): Browse

Maximum file size 1 MB

This logo will appear on the **top right corner** of the first page of the report.

Small logo: Browse

This logo is used on all detail pages. If not set, the default logo is used.

Signatures

Left: Browse

Maximum file size 1 MB

This signature or image will appear on the last page, left aligned.

Center: Browse

Maximum file size 1 MB

This signature or image will appear on the last page, center aligned.

Right: Browse

Maximum file size 1 MB

This signature or image will appear on the last page, right aligned.

Indicators

Custom intensity indicator:

A custom SVG for displaying allergen intensity can be pasted here. Will be styled according to specified colors.

The SVG element must contain the attributes **width** and **height** with a minimum dimension of 16 for correct display in PDF.

SVGs paths which should be colored must require the **colorize** class as well as a fill defined in the **defs** segment of the SVG.

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General	Design	Header	Content
Allergens to include:	All allergens		Specifies which kind of allergens should be included in the report.
Allergens sort order:	Alphabetically by allergen		Specifies the sort order of the allergens.
Include Overview:	☒		Whether the allergen overview should be printed or not.
Include product type in overview:	☒		Whether the product type should be printed in the header of allergen overview.
Include Total IgE:	On the front page		Whether the information of total IgE should be included in the report and where.
Cross-reactive information:	Include text		This section will be displayed underneath the allergen list itself.
Include Interpretation:	Include manually		It can be specified if the Raven report interpretation information should be included in the report or not. This section is printed underneath the allergens overview.

i Additional information (below the list of results)		Those informations are integrated into the report in an extra section underneath the allergen list.
Approval date:	☐	
Assay date:	☐	
Doctor info:	☐	The information about the doctor which accompanied the measurement/patient.
Print date:	☒	When the report was generated.
Report info:	☐	The „additional info“ field belonging to the measurement.
Sample date:	☒	

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Name:	New report template - Duplicate
Support products:	ALEX<sup>2</sup>
Preferred language:	Deutsch (de)

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
General configuration

* Name:	MADx
Default timezone:	(UTC+01:00) Amsterdam, Berlin, Bern, Rome, Stockholm, Vienna
Configure a timezone to see date and times in (instead of UTC). Please note that users can choose to use a different timezone.	
Default report template (ALEX 1):	Version 1 (ALEX1 & ALEX2)
Default report template (ALEX 2):	Version 2 (ALEX2 only)
Default report language:	English (United Kingdom) (en-gb)
Default MOD:	< No MOD >
Require ImageXplorer testarray scan:	Warn after 60 days, enforce after 90 days.

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
Raptor New measurements Approved measurements **New measurements** Add new

Report

Template: Version 2 (ALEX2 only)

Language: English (en)

Usually this tenant's default template and language will be used when sharing or exporting measurements. You may change the template and language associated with this measurement by saving your current configuration.

You can download this report in PDF format once the measurement has been approved.

You have access to manage the report templates and settings via the tenant administration. Please note that any changes will affect all users.

MADX
MACRO ARRAY DIAGNOSTICS

MADX
MACRO ARRAY DIAGNOSTICS

PATIENT ID
PATIENT NAME
DATE OF BIRTH
SAMPLE ID
BARCODE
ANALYZED ON
TESTED ALLERGENS
ZDS
TESTMETHODE
ALEX2

REFERRING PHYSICIAN

NOTE
The internal QC (Feasibility check for GD) was within acceptance range.

Lab report: Summary on detectable sensitisations

POLLEN
Grass Pollen

MICROORGANISMS
Fungal Spores & Yeast

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The screenshot shows the RAPTOR software interface. At the top, there are tabs for 'New measurements', 'Approved measurements', and 'Start m...'. The 'Approved Measurements' tab is active. Below the tabs, there is a search bar with the text 'Search by QR-Code, patient, sample, tag, analysis date ...'. Below the search bar, there is a dropdown menu labeled 'Download selected measurements'. The dropdown menu is open, showing options: 'Report (PDF format)', 'CSV format', 'Excel format', and 'ASTM format'. The 'Report (PDF format)' option is highlighted with a yellow circle. Below the dropdown menu, there is a table with columns: 'Source', 'Tags', 'Product type', 'QR-Code', and 'Sample C'. The table contains three rows of data, all with 'IX' in the 'Source' column, 'ALEX2' in the 'Product type' column, and '02AKV39D' in the 'QR-Code' column. The 'Sample C' column is empty. To the right of the table, there is a vertical menu with options: 'Update interpretation', 'Choose Reporttemplate', 'Share', 'Move', and 'Unapprove'.

	Source	Tags	Product type	QR-Code	Sample C
☒	IX		ALEX2	02AKV39D	
☒	IX		ALEX2	02AKV39D	
☒	IX		ALEX2	02AKV39D	

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Tags:

MOD:

Report Template: Version 2 (ALEX2 only) - Deutsch (de)

 Download ▾

 Report

 Save

» Move to other tenant

[New measurements](#)
[Approved measurements](#)
[Start measurement](#)
[Automats](#)

Report

Template: Default
Language: English

Save
Use defaults

Usually this tenants default template and language will be used when sharing or exporting measurements. You may change the template and language associated with this measurement by saving your current configuration.

Download PDF

Manage report template and settings

You have access to manage the report templates and settings via the tenant administration. Please note that any changes will affect all users.

Patient ID:	
Patient Name:	
Date of Birth:	
Sample ID:	
Barcode:	02AAE
Sampled on:	
Tested on:	12/18/
Printed on:	12/19/
Note:	The internal QC (Plausibility)

Lab report: Summary on de


Pollen	Grass Pollen
	Tree Pollen
	Weed Pollen
Mites	House Dust M
Microorganisms	Fungal Spores
	Legumes
	Grain
Plant-Based Food	Spices
	Fruits
	Vegetables
	Nuts & Seeds
	Milk

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Rule set

RAVEN Ruleset A for ALEX

A rule set defines the RAVEN or RAVEN² engine that is used to interpret results.

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AIT

Medication

MADx data declaration for smartRaven

☒ I accept to enter clinical data for using the full interpretation function of smartRAVEN

Patient has reported any symptoms

☐ Yes ☐ No

Additional information about the patient

☐ Atopic Dermatitis

☐ Allergen-specific immunotherapy (AIT) - ongoing, completed or discontinued

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[Request new interpretation](#)

This measurement already has an interpretation stored. Usually, it should not be necessary to request those again, but it could be necessary when the patients information has changed.

Interpretation

Raven Interpretation Summary

Sample Information

The sample was tested on ALEX² Barcode 02AAB019, interpretation date 12/19/2019.

Of the tested 295 allergens, 25 were above the detection limit of 0.3 kU_L/L. A sensitization can be an indicator of an allergy. For selected allergens, comments for interpretation are listed below.

Total IgE: 512 kU/L

The measured total IgE was 512 kU/L. Above a total IgE of 100 kU/L, allergy is likely.

Cross-Reactive Allergen Sensitization detected

Sensitization against allergen molecules which are markers of (broad) cross-reactivity between different allergen sources.

Detected cross-reactive allergen sensitizations:

- PR-10s: Aln g 1, Bet v 1, Cor a 1.0103, Fag s 1
- Polcalcins: Aln g 4, Phi p 7

PR-10 Proteins (PR10)

PR-10 inhalative: The major birch pollen allergen, Bet v 1, represents the prototype of all PR-10 allergens and is the primary sensitizer in regions with birch pollen exposure. The presence of PR-10 allergens in Fagales tree pollen explains IgE cross-reactivity between pollen from hazel, alder, beech, oak and hornbeam. PR-10 nutritive: PR-10 allergens in raw fruits, nuts, vegetable and legumes can induce oral allergy syndrome and sometimes severe allergic reactions in sensitized individuals. PR-10 allergens are not stable to heat and digestion.

Polcalcins (Ca++ binding proteins, PC)

Members of the PC family can cause inhalative symptoms. PC allergens can be found in pollen from grasses, trees and weeds. Inhalative symptoms manifest as allergic rhino-conjunctivitis and/or allergic asthma.

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
Tags:

MOD:

Report Template: Version 2 (ALEX2 only) - English (United States) (en-us)

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
Type:

Deactivated:

You cannot delete an MOD once it is in use. But you can deactivate it, so that it's not available for new measurements anymore.

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* Name:

Default timezone:

Configure a timezone to see date and times in (instead of UTC). Please note that users can choose to use a different timezone.

Default report template (ALEX 1):

Default report template (ALEX 2):

Default report template (FOX):

Default report template (PAX):

Default report language:

Default ALEX2 MOD:

Default FOX MOD:

Default PAX MOD:

Force MOD:

Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.

Hide Global MOD:

If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.

Display MOD Excluded Spots As Research Use Only:

If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".

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MOD: < No MOD >

Report Template: < No MOD >
Global_MOD_ALEX_v1.17
MOD_ALEX_v1.17

Force MOD: True
Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.

Hide Global MOD: True
If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.

Display MOD Excluded Spots As Research Use Only: False
If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".

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Force MOD: False
Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.


Hide Global MOD: False
If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.

Display MOD Excluded Spots As Research Use Only: True
If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".

Configuration

Calibration string:
Allows to override the calibration information from the cartridge's QR-Code. Leave empty to use the calibration information from the QR-Code on the cartridge (the usual behaviour).

Research spots: Show research spots
If enabled, "research use only" spots will be displayed on measurement details and arrays. This does not affect the export behaviour or the reports of measurements.

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

New Measurements

Search by barcode, patient, sample, tag, analysis date, a

Selected measurements ▾

- Download CSV format
- Download Excel format
- Approve
- Update interpretation
- Re-Process
- » Move
- Delete**

Tags	Product type
	FOX
	FOX
	ALEX2
	FOX
	ALEX2
	ALEX1

Delete measurements

You are about to delete the following measurements. Please note that this action is irreversible and cannot be undone.

Barcode	Analysis Date
@2AAB16@	16.12.2019 13:41

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Plant Food Food	M	Jug r 3	0.40
Plant Food Food	M	Jug r 1	0.38

Plant Food	M	Jug r 3	0.40
Food			
Plant Food	M	Jug r 1	0.38
Food			

Delete request scheduled now. [Cancel delete request](#)

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ign in. Can be

Background jobs

View details about running background jobs and manually schedule runs.

[View dashboard](#)

Measurement delete / unapproval log & requests

Show logging and pending requests for measurement delete and unapproval.

[Show log & requests](#)

Measurement unapproval / delete requests

Open requests

[Log & Finished requests](#)

Requested UTC	Requesting User (Email)	QR-Code	Request Type	Action
11.10.2023 13:05:28	Test Test (Test@test.com)	02AYQ145	Unapprove	Accept Reject
11.10.2023 12:59:27	Test Test (Test@test.com)	02AYQ14C	Delete	Accept Reject

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

Add new user on Result Sharing Portal

General configuration

First name:

Last name:

E-Mail address:

Password: Send an invitation e-mail to user automatically

Username:

It is a good idea to use the users e-mail-address as username as well.

Status: Locked (user cannot log in)

Enabled (user can log in)

Locked (user cannot log in)


Doctor Info Text:

This information will be printed on a patients report.

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≡ Share

You are about to share a measurement with a user of the User Portal. Please note that the doctor will have access to the results of the measurement. Once you share a measurement you'll be provided with a link to access the shared information via the encrypted User F

Barcode: 02AAB021

Share with: X Dr. (testuser2@macroarraydx.com)

If you can't find a specific user of the User Portal, go to the [tenant administration](#) to configure the list of available users of the User Portal.

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RAPTOR
New measurements
Approved measurements
Start measurement with ImageXplorer
Automats
MADx Testtenant

Share

Share this measurement


User	Link	Shared
RP Bak	https://resultshare.raptor-server.com/s/07PXEPC7ZNYthM4PdRfmc4NR5pUtWq5FJcqOY7Ry6oEx08U5j	7 months ago

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 Result Sharing Portal English (en) ▾

Welcome to the Result Sharing Portal

The Raptor-Server Result Sharing Portal allows users to view shared measurements. You will need a username and password to log in to view all measurements that have been shared with you. Also, a sharing link will allow to view specific shared measurements.

Please enter your username and your password:

Did you forget your password? No problem, you can [request a new one](#).

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Result Sharing Portal
English (en) ▾
Change password
Log out

Measurements

These are the measurements that are available for you. You may view the details of each measurement, or download the report PDF file.

Measurement	Product	Patient	Sample	Date
01AAK239	ALEX1		11852162	7/12/2019
01AAK281	ALEX1		04824821	12/18/2019
01AAK891	ALEX1		66560260	1/28/2020
02AAA015	ALEX2			2/21/2020

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Measurement 02AYQ143

Product: ALEX2
Patient:
Date of Birth:
Sample:
Info:
Doctor Info:
Sampling Date:
Assay Date: 10/12/2023
Processed Date: 10/12/2023
Approval Date: 10/12/2023

Report Language:

English (United Kingdom) ▼

You can request a new RAVEN interpretation that will be included in the report.

[Request new RAVEN interpretation](#)

You can download a PDF report that includes all information and results of this measurement.

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© Copyright by MacroArray Diagnostics