# RAPTOR

SERVER ANALYSIS SOFTWARE

## POKYNY NA POUŽITIE

### OBSAH

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Welcome to RAPTOR SERVER Analysis Software, the analysis software for the [ALEX multiplex test](#). Please log in below to continue.

<p>Log in with Azure AD / Office 365</p> <p>If your organization integrates the software into its Microsoft Azure/Office 365 user directory (single-sign-on) click the button below. You will be redirected to the Microsoft login screen.</p> <p> Log in</p>	<p>Log in with Username &amp; password</p> <p>Please enter your username and password.</p> <p>Your username</p> <p>Your password </p> <p> Log in</p> <p>Forgot your password? No problem, you can <a href="#">request a new one</a>.</p>
--	---

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Log in with Username & password

Please enter your username and password.

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### Contract approval required

Please review the following documents. You are required to confirm that you have read and agree with each document in order to use RAPTOR SERVER:

- [Bitte lesen](#) (last updated 23.11.2023 18:00)
- [RAPTOR SERVER IFU](#) (last updated 04.12.2023 11:05)
- [Terms of Service](#) (last updated 22.11.2023 01:00)

**Confirm**

By clicking the button, I confirm that I have reviewed these documents and agree.

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New Measurements					Last 7 days
QC	Product type / Source	QR-Code	Sample Code	Analysis Date	
❌ Failed	Research MAX 45k	99AFG0F7	NK_FOX	30/03/2023 20:50	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG0F8	N.G.	30/03/2023 20:50	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG0F9	E18735	30/03/2023 20:50	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG0FA	20220325w56	30/03/2023 20:50	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG0EE	HM100038	30/03/2023 20:50	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG04D	786118	30/03/2023 20:37	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG016	HM100225	30/03/2023 20:37	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG0CA	E17685_M	30/03/2023 20:37	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG0CC	E17685_M	30/03/2023 20:37	<a href="#">Details</a> <a href="#">Report</a>
❌ Failed	Research MAX 45k	99AFG015	E17685_F	30/03/2023 20:37	<a href="#">Details</a> <a href="#">Report</a>

In total there are 270 new measurements in the selected time period.

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
#### 9K-Test-Instrument

Not initialized

[Initialize](#) [Shutdown](#) [Monthly Maintenance](#)

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
The screenshot shows the RAPTOR user interface. At the top, there is a navigation bar with the RAPTOR logo and buttons for "New measurements", "Approved measurements", "Start measurement", and "Automats". Below this, the "Personal settings" section includes fields for Name, E-Mail address, Language (set to English (en)), Culture (set to English (United States) (en)), and Timezone (set to UTC+01:00). To the right, the "Multi-Factor Authentication" section has a "Manage MFA" button and explanatory text. Below that, the "Configuration" section has a "Research spots" dropdown menu set to "Show research spots".

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**General configuration**

**\* Name:**

**Default timezone:**    
Configure a timezone to see date and times in (instead of UTC). Please note that users can choose to use a different timezone.

**Default report template (ALEX2):**

**Default report template (ALEX3, ALEX Air, ALEX Food):**

**Default report template (FOX):**

**Default report template (PAX):**

**Default report language:**

**Default ALEX MOD:**

**Default FOX MOD:**

**Default PAX MOD:**

**Force MOD:**    
Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.

**Hide Global MOD:**    
If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.

**Display MOD Excluded Spots As Research Use Only:**    
If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".

**MFA:**  **Force MFA**   
If this setting is enabled, users are forced to set up an MFA on sign-in, in case they haven't yet.

**MAX device default settings**

**\* Default Minimum RunId:**    
When a new MAX device is created, this value will be used for its **Minimum Run Id**

**Tenant-specific batch releases**

Restrict this tenant to only allow the following batches to be scanned or to restrict certain batches. Enter one per line.

Use:

**Whitelist**

**Blacklist**

**Whitelisted Batches:**    
Enter one batch code including product identifier per line, for example "02AAA"

**Blacklisted Batches:**    
Enter one batch code including product identifier per line, for example "02AAA"

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#### Tenant-specific batch releases

Restrict this tenant to only allow the following batches to be scanned or to restrict certain batches. Enter one per line.

Use:

**Whitelist**

**Blacklist**

**Whitelisted Batches:**

02AAA  
02ABC  
02BAA

Enter one batch code including product identifier per line, for example "02AAA"

**Blacklisted Batches:**

Enter one batch code including product identifier per line, for example "02AAA"

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You requested a new measurement from ImageExplorer "IX Trainingsraum". Please wait until the measurement is completed ...

Invalid batch detected! The measurement can't be processed because its batch is blocked.

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
#### Login Credentials

Type:

Please select if the user should be validated against Azure AD (Office 365), or if you want to specify username and password yourself. It is recommended to use Azure AD whenever possible. This allows for increased security by supporting two-factor authentication. Also, it is more convenient for the user, because he or she does not have to remember a new username and password (single-sign-on).

Username:

Password:

It is a good idea to use the users e-mail-address as his username, because the username needs to be unique across the entire system.

#### Access Rights

##### Preconfigured groups

Use these pre-configured groups to automatically assign a set of access rights to this user. You can still add or remove specific access rights.

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Tenants

This user is connected to the following tenants and has the following access rights:

Tenants	Can share measurements	Can edit measurements	Can edit personal data	Can approve measurements	Can manage basic tenant configuration	Can manage API keys	Can manage ImageExplorers	Can manage MAX devices	Can view agent connections	Can manage users	Can manage users of Result Sharing Portal	Can manage report templates	Can manage measurement data completion	System Administration	Can unapprove measurements	Can use a custom calibration string	Can manage Run Configuration	Can skip r
<input type="checkbox"/> Another MADx Tenant	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				
<input type="checkbox"/> MADx Tenant	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				

Additionally this user can be added to the following tenants by checking the tenant name and hitting save afterwards. User rights w applied from current open tenant user.

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Additionally this user can be added to the following tenants by checking the tenant name and hitting save afterwards. User right applied from current open tenant user.

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### Tenants

This user is connected to the following tenants and has the following access rights:

[Edit selected tenants rights](#)

<input type="checkbox"/> Tenants	Can share measurements	Can edit measurements	Can edit personal data	Can approve measurements	Can manage basic tenant configuration	Can manage API keys	Can manage Imagey
<input checked="" type="checkbox"/> Another MADx Tenant	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				

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General information

Name:

Key:   
This key is required when installing an Agent.

ImageXplorer Serialnumber:   
Serialnumber of connected ImageXplorer.

Notes:

Internal Notes:   
For internal use only.

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
In order to exclude the influence of artifacts in the automated image analysis (satellite spots, sample contaminations, dust, smeared spots, ...), the images must be checked by a trained operator before the results are approved in order to exclude false results. In case of discrepancies between the processed array and the image acquired by the RAPTOR SERVER please consult your local distributor or MADX Support.

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

**Show Analysis & QC Logs**

**Reset everything and reprocess**

This discards all modifications to the array, and re-runs the complete analysis process from the original image data. Please note that any manual changes will be irreversible lost.

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#### Analysis Log

```
'PreProcessedPixels' took 128ms.
Using definition '02_ALEXII_IVD_v1.6.0'.
'ApplyArrayDefinition' took 74ms.
Parsing calibration string ...
'ApplyDefinitions' took 75ms.
'FindGrid' took 1ms.
'AlignFeatures' took 862ms.
'CropImage' took 1ms.
'QuantifyFeatures' took 231ms.
'SubtractBackground' took 681ms.
Parsing calibration string ...
Using calibrators from QR-Code.
Total Igt (tIgt) normalized by calculating 0,330999755859375E-07x*0,05697999954223633x+-30,69269921875
(x=8515), which results in 515.
'ReQuantify' took 57ms.
```

#### QC Log

These are the values used for QC (ALEX v2): Curve Correction=0,9997370867894242, tIgt=8515, GD1=28085 (SNR=105.38), GD2=29369 (SNR=185.28), GD3=28569 (SNR=147.53), GD4=28385 (SNR=128.1)  
Reference values are as follows: Guide Dots must have a value above 5000; each Guide Dot spot must exist and be unique; Guide Dot spot must not be set to bad. Average background must be between 4000 and 16000. QC (ALEX v2) succeeded.

#### Background Log

Calculated background statistics for 164 pixels (from a total of 182, skipped smallest 9, and highest 9 pixels). The average is 9.419, the minimum 7973 and the maximum 11792 (minimum and maximum is calculated including skipped pixels). This results in an average difference of 628, which is an SD of 49,17 and a CV of 0,52%.  
For comparison, there is an overall total of 282944 pixels. The overall average is 14.871, the minimum 6.012 and the maximum 64.965.

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**New Measurements**

Search by QR-Code, patient, sample, tag, analysis date, approval date ...

Selected measurements

	Source	QC	Tags	Product type
<input type="checkbox"/>	IX	✓ OK		ALEX2
<input checked="" type="checkbox"/>	IX	✓ OK		ALEX2
<input checked="" type="checkbox"/>	IX	✓ OK		ALEX2
<input checked="" type="checkbox"/>	IX	✓ OK		ALEX2

- Download CSV format
- Download Excel format
- Approve
- Update interpretation
- Choose Reporttemplate
- Modify Tags**
- Move
- Delete

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**Re-process measurements**

You are about to re-process the following measurements. If you click OK and the re-processing is undone, the existing analysis result will be irreversibly overwritten.

QR-Code	Sample Code
02APK3C8	Dummy8
02APK258	Dummy1
02APK07F	Dummy13

Re-process these measurements Cancel

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New Measurements

Search by QR-Code, patient, sample, tag, analysis date, approval date ...

	Source	QC	Tags	Product type	QR Code	Sample Code
<input type="checkbox"/>	IX	✓ OK		ALEX2	02AGT957	
<input type="checkbox"/>	IV	✓ OK				

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## New Measurements

<input type="checkbox"/>	Source	QC	Tags	Product type
<input type="checkbox"/>	Recaptured	✓ OK	recaptured	ALEX1
<input type="checkbox"/>	Recaptured	✓ OK	recaptured	ALEX1
<input type="checkbox"/>	QR-Code was scanned more than once.		recaptured	ALEX1

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The screenshot shows the RAPTOR software interface. At the top, there are navigation tabs for 'New measurements' and 'Approved measurements'. Below this is a search bar and a 'Selected measurements' dropdown menu. The main area contains a table with columns: Source, QC, Tags, Product type, and QR-Code. Three rows are visible, each with a blue checkmark in the first column, indicating they are selected. A context menu is open on the right side of the table, listing several actions: 'Download CSV format', 'Download Excel format', 'Approve' (highlighted with a yellow circle), 'Update interpretation', 'Re-Process', 'Move', and 'Delete'.

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
### Approve measurements

You are about to approve the following measurements.

There are measurements with one or more warnings. Please make sure to double-check these measurements before approval.

QR-Code	Sample Code	Analysis Date
02AKV035		28.04.2022 11:14

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The screenshot shows the 'Approved Measurements' section of the software. A table lists three measurements with source 'IX' and product type 'ALEX2'. A dropdown menu is open, showing options: Report (PDF format), CSV format, Excel format, ASTM format, Update interpretation, Choose Reporttemplate, Share, Move, and Unapprove.

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The screenshot shows the 'Exports' section. A table lists export requests with columns: Export requested, # of measurements included, and Export completed. A user menu is open on the right with options: Administration, Tenant Admin, Exports, Personal Settings, Log out, Support, and Imprint.

Export requested	# of measurements included	Export completed
8/16/2021 11:07 AM	4	0 / 4
8/16/2021 11:06 AM	4	8/16/2021 11:06 AM
8/16/2021 11:04 AM	3	8/16/2021 11:04 AM

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MADx Tenant

- Administration
  - Tenant Admin
  - System Admin
  - Exports**
  - Personal Settings
  - System Notifications
  - Log out
- Support
- Imprint

### Exports

Certain types of data exports (for example, the export of many PDF reports at a time) are processed in the background. Please note that you can keep working, all exports are done in the background.

**Export multiple measurements**

Export requested	# of measurements included
06.07.2022 08:52	3

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### Export multiple measurements

Please note that **Report PDF** and **ASTM** export is only available for certain reports.

QR-Code/Product/Batch  Sa

Analysis Date (from)

Analysis Date (to)

Tag

Status

Export CSV  Export Excel

**Search Measurements**

## Export multiple measurements

Please note that **Report PDF** and **ASTM** export is only available for single measurements.

QR-Code/Product/Batch

02aaa

Analysis Date (from)

TT . MM . JJJ

Analysis Date (to)

TT . MM . JJJ

Tag

Status

All

Export CSV

Export Excel

 Start export ...

← Back

66 Measurements found that match your criteria

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## Edit report template

Your changes have been saved successfully.

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
### Configure report template "New report template"

Field	Value	Description
Name:	New report template	This is the name of the template which shows up in the dropdown list while generating a report.
Base Template:	Version 2 (ALEX2 only)	A base template which specifies basic functionality must be chosen.
Preferred language:	English (United States) (en-us)	A default language can be specified which will be used if there isn't yet a language stored to a measurement's report.

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PDF naming

Setting Type: Manual

The filename of the report PDF can be customized. Either built automatically based on your selection below or manually adjusted in the order of your needs.

Name Pattern: Report\_{PatientId}\_{MeasurementQrCode}

Available Elements: {PatientId} ; {PatientName} ; {SampleId} ; {MeasurementQrCode} ; {DateAndTime}

---

PDF naming

Setting Type: Automatic

The filename of the report PDF can be customized. Either built automatically based on your selection below or manually adjusted in the order of your needs.

Include Patient ID:

Include Patient Name:

Include Sample ID:

Include Measurement QR-Code:

Include Date & Time:

Footer

Footer text:

Will be printed in the footer if provided and rendered as HTML.

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### Configure report template "New report template"

General	Design	Header	Content
<b>PDF format:</b> <input checked="" type="radio"/> A4 <input type="radio"/> Letter		Specifies the paper size of the output PDF.	
<b>Font:</b> <input type="text" value="Arial"/>		The font family of the whole report texts can be adjusted.	
<b>Background image:</b> <input type="text"/> <input type="button" value="Browse"/>		If set, this image is used as background on the first page of the PDF.	
Maximum file size 1 MB			

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---

### Logos

**Logo (left):**   Maximum file size 1 MB  
This logo will appear on the **top left corner** of the first page of the report.

**Text (center):**   
This text will appear either on the left, right, or center of the top of the first page (depending on your configured logos). Supports the placeholders `{{ PATIENT_PORTAL_QRCODE }}` and `{{ PATIENT_PORTAL_URL }}` to dynamically insert measurement information.

**Logo (right):**   Maximum file size 1 MB  
This logo will appear on the **top right corner** of the first page of the report.

**Small logo:**   Maximum file size 1 MB  
This logo is used on all detail pages. If not set, the default logo is used.

---

### Signatures

**Left:**   Maximum file size 1 MB  
This signature or image will appear on the last page, left aligned.

**Center:**   Maximum file size 1 MB  
This signature or image will appear on the last page, center aligned.

**Right:**   Maximum file size 1 MB  
This signature or image will appear on the last page, right aligned.

---

### Indicators

**Custom intensity indicator:**   
A custom SVG for displaying allergen intensity can be pasted here. Will be styled according to specified colors.

- ⚠ The SVG element must contain the attributes **width** and **height** with a minimum dimension of 16 for correct display in PDF.
- 🌈 SVGs paths which should be colored must require the **colorize** class as well as a fill defined in the **defs** segment of the SVG.

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Configure report template "New report template"

Header content:		Left half	Right half
		Custom headline	Custom headline
<input type="checkbox"/>	Patient Name	<input type="checkbox"/>	Doctor Information
<input checked="" type="checkbox"/>	Patient Name	<input type="checkbox"/>	None
	Custom headline		
<input type="checkbox"/>	Date of Birth	<input type="checkbox"/>	None
<input type="checkbox"/>	Sample ID	<input type="checkbox"/>	None
<input type="checkbox"/>	QR-Code	<input type="checkbox"/>	Additional Info text
<input type="checkbox"/>	Analysis Date	<input type="checkbox"/>	None
<input type="checkbox"/>	Amount of tested allergens	<input type="checkbox"/>	None
<input type="checkbox"/>	Test method	<input type="checkbox"/>	Quality Control text

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General	Design	Header	Content
<b>Allergens to include:</b>	All allergens		Specifies which kind of allergens should be included in the report.
<b>Allergens sort order:</b>	Alphabetically by allergen		Specifies the sort order of the allergens.
<b>Include Overview:</b>	<input checked="" type="checkbox"/>		Whether the allergen overview should be printed or not.
<b>Include product type in overview:</b>	<input checked="" type="checkbox"/>		Whether the product type should be printed in the header of allergen overview.
<b>Include Total IgE:</b>	On the front page		Whether the information of total IgE should be included in the report and where.
<b>Cross-reactive information:</b>	Include text		This section will be displayed underneath the allergen list itself.
<b>Include Interpretation:</b>	Include manually		It can be specified if the Raven report interpretation information should be included in the report or not. This section is printed underneath the allergens overview.
<b>Additional information (below the list of results)</b>		Those informations are integrated into the report in an extra section underneath the allergen list.	
<b>Approval date:</b>	<input type="checkbox"/>		
<b>Assay date:</b>	<input type="checkbox"/>		
<b>Doctor info:</b>	<input type="checkbox"/>		The information about the doctor which accompanied the measurement/patient.
<b>Print date:</b>	<input checked="" type="checkbox"/>		When the report was generated.
<b>Report info:</b>	<input type="checkbox"/>		The „additional info“ field belonging to the measurement.
<b>Sample date:</b>	<input checked="" type="checkbox"/>		

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
**General configuration**

**\* Name:** MADx

**Default timezone:** (UTC+01:00) Amsterdam, Berlin, Bern, Rome, Stockholm, Vienna  
Configure a timezone to see date and times in (instead of UTC). Please note that users can choose to use a different timezone.

**Default report template (ALEX 1):** Version 1 (ALEX1 & ALEX2)

**Default report template (ALEX 2):** Version 2 (ALEX2 only)

**Default report language:** English (United Kingdom) (en-gb)

**Default MOD:** < No MOD >

**Require ImageExplorer testarray scan:** Warn after 60 days, enforce after 90 days.

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
The screenshot shows the RAPIR web interface. At the top, there are navigation tabs: "New measurements", "Approved measurements", and "Start m". Below this is a search bar with the text "Search by QR-Code, patient, sample, tag, analysis date ...". A dropdown menu "Download selected measurements" is open, showing options: "Report (PDF format)", "CSV format", "Excel format", "ASTM format", "Update interpretation", "Choose Reporttemplate", "Share", "Move", and "Unapprove". Below the menu is a table of "Approved Measurements".

	Source	Tags	Product type	QR-Code	Sample C
<input checked="" type="checkbox"/>	IX		ALEX2	02AKV39D	
<input checked="" type="checkbox"/>	IX		ALEX2	02AKV39D	
<input checked="" type="checkbox"/>	IX		ALEX2	02AKV39D	

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### Details

The screenshot shows the "Details" page for a measurement. It contains the following fields:

- QR-Code: 02AKV39D
- Sample Code:
- Requisition Code:
- Sampling Date:
- Assay Date: 27.09.2021
- Patient Code:
- Patient Name:
- Date of Birth:
- Info:
- Doctor Info:
- Analysis: 27.09.2021 15:39
- Approved: 27.09.2021 15:39
- ImageExplorer: 23167844
- Tags:
- MOD:
- Report Template: Version 2 (ALEX2 only) - Deutsch (de)

At the bottom, there is a row of buttons: "Download", "Report" (circled in yellow), "Save", and "Move to other tenant".

X Raptor
 New measurements
Approved measurements
Start measurement
Automats

### Report

Template:

Language:

Usually this tenants default template and language will be used when sharing or exporting measurements. You may change the template and language associated with this measurement by saving your current configuration.

Patient ID:	
Patient Name:	
Date of Birth:	
Sample ID:	
Barcode:	02AAE
Sampled on:	
Tested on:	12/18/
Printed on:	12/19/
Note: The internal QC (Plausibility)	

**Lab report: Summary on de**

<b>Pollen</b>	Grass Pollen
	Tree Pollen
	Weed Pollen
<b>Mites</b>	House Dust M
<b>Microorganisms</b>	Fungal Spores
	Legumes
	Grain
	Spices
<b>Plant-Based Food</b>	Fruits
	Vegetables
	Nuts & Seeds
	Milk


You have access to manage the report templates and settings via the tenant administration. Please note that any changes will affect all users.

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
AIT

Medication

### MADx data declaration for smartRaven

I accept to enter clinical data for using the full interpretation function of smartRAVEN

### Patient has reported any symptoms

Yes  No

### Additional information about the patient

Atopic Dermatitis

Allergen-specific immunotherapy (AIT) - ongoing, completed or discontinued

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[Request new interpretation](#)

This measurement already has an interpretation stored. Usually, it should not be necessary to request those again, but it could be necessary when the patients information has changed.

**Interpretation**

**Raven Interpretation Summary**

**Sample Information**

The sample was tested on ALEX<sup>2</sup> Barcode 02AAB019, interpretation date 12/19/2019.

Of the tested 295 allergens, 25 were above the detection limit of 0.3 kU<sub>L</sub>/L. A sensitization can be an indicator of an allergy. For selected allergens, comments for interpretation are listed below.

**Total IgE: 512 kU/L**

The measured total IgE was 512 kU/L. Above a total IgE of 100 kU/L, allergy is likely.

**Cross-Reactive Allergen Sensitization detected**

Sensitization against allergen molecules which are markers of (broad) cross-reactivity between different allergen sources.

Detected cross-reactive allergen sensitizations:

- PR-10s: Aln g 1, Bet v 1, Cor a 1.0103, Fag s 1
- Polcalcins: Aln g 4, Phi p 7

**PR-10 Proteins (PR10)**

PR-10 inhalative: The major birch pollen allergen, Bet v 1, represents the prototype of all PR-10 allergens and is the primary sensitizer in regions with birch pollen exposure. The presence of PR-10 allergens in Fagales tree pollen explains IgE cross-reactivity between pollen from hazel, alder, beech, oak and hornbeam. PR-10 nutritive: PR-10 allergens in raw fruits, nuts, vegetable and legumes can induce oral allergy syndrome and sometimes severe allergic reactions in sensitized individuals. PR-10 allergens are not stable to heat and digestion.

**Polcalcins (Ca++ binding proteins, PC)**

Members of the PC family can cause inhalative symptoms. PC allergens can be found in pollen from grasses, trees and weeds. Inhalative symptoms manifest as allergic rhino-conjunctivitis and/or allergic asthma.

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Raptor | New measurements | Approved measurements | Start measurement | Automate

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

Raptor New measurements Approved measurements Start measurement Automats

### Create new MOD

Name:

Type:

Deactivated:

You cannot delete an MOD once it is in use. But you can deactivate it, so that it's not available for new measurements anymore.

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\* Name:

Default timezone:

Configure a timezone to see date and times in (instead of UTC). Please note that users can choose to use a different timezone.

Default report template (ALEX 1):

Default report template (ALEX 2):

Default report template (FOX):

Default report template (PAX):

Default report language:

Default ALEX2 MOD:

Default FOX MOD:

Default PAX MOD:

Force MOD:

Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.

Hide Global MOD:

If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.

Display MOD Excluded Spots As Research Use Only:

If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".

Require ImageExplorer ConfigExplorer scan: Warn after 60 days, enforce after 90 days.

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
<b>Force MOD:</b>	False	◆
	Requires a measurement to always have an MOD assigned. If this setting is enabled, a measurement cannot be approved without an MOD assigned.	
<b>Hide Global MOD:</b>	False	◆
	If this setting is enabled, only the tenant-specific MOD will be available to select, and the global system-wide MOD will be hidden.	
<b>Display MOD Excluded Spots As Research Use Only:</b>	True	◆
	If this setting is enabled, spots excluded by an MOD will not be hidden from the measurement image, but displayed as a "RUO spot".	

### Configuration

<b>Calibration string:</b>	<input type="text"/>	
	Allows to override the calibration information from the cartridge's QR-Code. Leave empty to use the calibration information from the QR-Code on the cartridge (the usual behaviour).	
<b>Research spots:</b>	Show research spots	◆
	If enabled, "research use only" spots will be displayed on measurement details and arrays. This does not affect the export behaviour or the reports of measurements.	

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
**New Measurements**

Search by barcode, patient, sample, tag, analysis date, a

Selected measurements ▾

- Download CSV format
- Download Excel format
- Approve
- Update interpretation
- Re-Process
- » Move
- Delete**

Tags	Product type
	FOX
	FOX
	ALEX2
	FOX
	ALEX2
	ALEX1

### Delete measurements

You are about to delete the following measurements. Please note that this action is irreversible and cannot be undone.

Barcode	Analysis Date
02AAB160	16.12.2019 13:41

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Plant Food Food	M Jug r 3	0.40
Plant Food Food	M Jug r 1	0.38

Request delete

Plant Food Food	M Jug r 3	0.40
Plant Food Food	M Jug r 1	0.38

Delete request scheduled now. [Cancel delete request](#)

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Background jobs  
View details about running background jobs and manually schedule runs.

[View dashboard](#)

---

**Measurement delete / unapproval log & requests**  
Show logging and pending requests for measurement delete and unapproval.

[Show log & requests](#)

Measurement unapproval / delete requests

Requested UTC	Requesting User (Email)	QR-Code	Request Type	Action
11.10.2023 13:05:28	Test Test (Test@test.com)	02AYQ145	Unapprove	<a href="#">Accept</a> <a href="#">Reject</a>
11.10.2023 12:59:27	Test Test (Test@test.com)	02AYQ14C	Delete	<a href="#">Accept</a> <a href="#">Reject</a>

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


Add new user on Result Sharing Portal

General configuration

<p><b>First name:</b></p> <input type="text"/>	<p><b>Last name:</b></p> <input type="text"/>
<p><b>E-Mail address:</b></p> <input type="text"/>	<p><b>Status:</b></p> <p>Enabled (user can log in) <span>⇩</span></p>
<p><b>Username:</b></p> <input type="text"/>	<p><b>Password:</b></p> <p>Send an invitation e-mail to user automatically <span>⇩</span></p>
<p><input checked="" type="checkbox"/> <b>Automatically send an email summary in case new measurement(s) have been shared</b>  <small>This email is sent once per day as a summary, just in case a new measurement was shared.</small></p>	
<p><b>Doctor Info Text:</b></p> <div style="border: 1px solid #ccc; height: 40px; width: 100%;"></div> <p><small>This information will be printed on a patients report.</small></p>	

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### Share

You are about to share a measurement with a user of the User Portal. Please note that the doctor will have access to the results of the measurement. Once you share a measurement you'll be provided with a link to access the shared information via the encrypted User F

Barcode: 02AAB021

Share with: X Dr. (testuser2@macroarraydx.com)


If you can't find a specific user of the User Portal, go to the [tenant administration](#) to configure the list of available users of the User Portal.

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
### Welcome to the Result Sharing Portal

The Raptor-Server Result Sharing Portal allows users to view shared measurements. You will need a username and password to log in to view all measurements that have been shared with you. Also, a sharing link will allow to view specific shared measurements.

Please enter your username and your password:

Did you forget your password? No problem, you can [request a new one](#).

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✕ Result Sharing Portal English (en) ▼ [Change password](#) [Log out](#)

### Measurements

These are the measurements that are available for you. You may view the details of each measurement, or download the report PDF file.

🔍 Search by barcode, patient or sample.

Measurement	Product	Patient	Sample	Date
<a href="#">01AAK239</a>	ALEX1		11852162	7/12/2019
<a href="#">01AAK281</a>	ALEX1		04824821	12/18/2019
<a href="#">01AAK891</a>	ALEX1		66560260	1/28/2020
<a href="#">02AAA015</a>	ALEX2			2/21/2020

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## Measurement 02AYQ143

**Product:** ALEX2  
**Patient:**  
**Date of Birth:**  
**Sample:**  
**Info:**  
**Doctor Info:**  
**Sampling Date:**  
**Assay Date:** 10/12/2023  
**Processed Date:** 10/12/2023  
**Approval Date:** 10/12/2023

**Report Language:**

English (United Kingdom) ▼

You can request a new RAVEN interpretation that will be included in the report.


[Request new RAVEN interpretation](#)

You can download a PDF report that includes all information and results of this measurement.

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



## Spustenie merania nie je možné:

The ImageXplorer IX training cannot be used, because it requires a ConfigXplorer scan.

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
ImageXplorer: 40287249

Automatic QC: **× Failed**

Analysis Error: **Insufficient data transfer - check connection.**

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© Copyright by MacroArray Diagnostics

MacroArray Diagnostics (MADx)

Lemböckgasse 59, Top 4

1230 Vienna, Austria

+43 (0)1 865 2573

[www.madx.com](http://www.madx.com)

Version number: 22-IFU-01-SK-09 (EN-26)

Date of Issue: 2026-02