



# Digital panels enable targeted enrichment of haematological cancer-associated genomic regions during sequencing

Adaptive Sampling enables targeted enrichment of digitally defined genomic regions and supports interrogation of haematological cancer-associated genes from blood-derived and bone marrow-derived DNA




Fig. 1 a) Library preparation with barcoding, b) Adaptive Sampling schematic, c) Digital panel design.

## Enrichment of genomic regions associated with pan-haematological cancers using Adaptive Sampling

Adaptive Sampling enables efficient enrichment of regions of interest during sequencing using a digital panel, without additional sample preparation. High molecular-weight DNA is fragmented to ~8 kb and end-prepared before barcode and adapter ligation using the Native Barcoding Kit 24 V14 (Fig. 1a). MinkNOW uses a reference FASTA file and a BED file defining the target regions to identify reads for continued sequencing, while off-target reads are rejected after approximately 400 bp (Fig. 1b). To demonstrate this approach, we used a BED file targeting 144 genes and regions associated with acute myeloid leukaemia (AML), acute lymphoblastic leukaemia (ALL), multiple myeloma (MM), and lymphoma (Fig. 1c).




Fig. 3 a) Sequencing yield, b) on-target lengths, c) on-target coverage, and d) background coverage.

## Adaptive Sampling of blood-derived and bone marrow-derived DNA provides sufficient enrichment for analysis

Genomic DNA extracted with the Puregene kit (n=3) was sequenced on PromethION Flow Cells (R10.4.1) using Adaptive Sampling with the pan-haematological cancer panel BED file for 72 hours. Comparable sequencing yields were obtained for blood-derived and bone marrow-derived DNA, at  $65 \pm 4$  Gb and  $66 \pm 1$  Gb, respectively (Fig. 3a), corresponding to an average of 112 million reads with an on-target read length N50 of 7.7 kb (Fig. 3b). Downsampling showed that 24 hours of sequencing was sufficient to achieve  $>100\times$  coverage at 95% of target sites for both sample types. After 72 hours, mean on-target coverage reached  $\sim 198\times$  for blood and  $\sim 180\times$  for bone marrow (Fig. 3c), with background coverage of  $\sim 14\times$  for each (Fig. 3d).




Fig. 2 a) DNA yield from blood and bone marrow, b) Fragment lengths, and c) Sample input volumes.

## DNA extraction from whole blood and bone marrow research samples using common extraction kits

Using 200  $\mu$ l of blood (n=3), the highest gDNA yield was obtained with GenFind ( $5.8 \pm 0.5$   $\mu$ g), followed by Puregene ( $3.6 \pm 1.1$   $\mu$ g) and QIAamp DSP ( $3.0 \pm 0.2$   $\mu$ g). For bone marrow (n=1), GenFind and Puregene produced comparable DNA yield per input volume (25.3  $\mu$ g/ml and 26.9  $\mu$ g/ml) (Fig. 2a). Femtopulse analysis showed that Puregene gave the longest DNA fragments for both sample types, with mean fragment lengths of 115 kb  $\pm$  25 kb for blood and 108 kb  $\pm$  6 kb for bone marrow. GenFind produced shorter fragments, while QIAamp DSP generated the shortest DNA from blood (Fig. 2b). Reduced inputs of 60  $\mu$ l for blood and 30  $\mu$ l for bone marrow each yielded more than 1  $\mu$ g of DNA, supporting downstream sequencing workflows (Fig. 2c).




Fig. 4 a) *De novo* SNV calling, b) Horizon Mimix Myeloid Cancer Panel *FLT3* ITD insertion.

## Benchmarking of methods to detect single nucleotide variants (SNVs) with digital panel data

For *de novo* SNV analysis, on-target COLO829 reads were basecalled using the super accuracy (SUP) model, then ClairS-TO and DeepSomatic were compared across binned variant allele frequencies (VAFs) for increasing sequence depth. ClairS-TO showed more balanced performance, whereas DeepSomatic recovered more true variants at 5–10% VAF, but with a higher false-positive rate (Fig. 4a). To assess structural variant detection, combined on-target and off-target data from the Horizon Mimix Myeloid Cancer Panel gDNA Reference Standard were examined for the *FLT3* internal tandem duplication (*FLT3*-ITD) 300 bp insertion present at 5% VAF. The insertion was visible in the IGV alignment (Fig. 4b).