

## Nanopore-only microbial isolate sequencing solution (NO-MISS): flexible, rapid bacterial and fungal isolate whole-genome sequencing

Across public health, clinical microbiology research, food safety, and the biopharma sector, whole-genome microbial isolate sequencing is a valuable tool. It is especially important during outbreak investigations, where it delivers crucial information on antimicrobial resistance genes, plasmids and other mobile genetic elements, and further features that may affect strain phenotypes. For the biopharma sector, whole-genome sequencing of production strains ensures that they are well characterised, stable, and free from mutations and contamination.

Complete, contiguous, high-quality microbial genome assembly, however, is limited with legacy short-read sequencing technology. Oxford Nanopore sequencing overcomes these limitations, generating reads of unrestricted length to deliver comprehensive genome assemblies. Plus, high-accuracy nanopore sequencing shows no bias in GC-rich regions<sup>1</sup> and spans repeat-rich sequences and structural variants that are inaccessible to short-read methods. Oxford Nanopore sequencing offers scalability, flexible sample batching, and rapid turnaround time from sample to answer. You can now generate complete, reference-quality microbial genome sequences in-house, in one go, without needing to outsource or use multiple validation techniques.

Here we present a rapid end-to-end workflow for whole-genome sequencing of microbial isolates, compatible with any of our sequencing devices and including analysis with the EPI2ME™ platform.

### Extraction:

lysing cells and obtaining high-quality DNA

View extraction protocol recommendations for your sample type, plus guidance on DNA storage and contaminants:  
[nanoporetech.com/documentation/prepare](https://nanoporetech.com/documentation/prepare)

To ensure high outputs of long reads in Oxford Nanopore sequencing, it is important to select an extraction method that yields high-quality DNA. Depending on your research requirements and sample, we recommend different cell lysis and extraction methods. Our workflow is suitable for multiple microbial isolate types.

We provide four options for extracting DNA from your isolates, comprising a universal bead-beating method for high-throughput requirements and manual column-based techniques for bacteria, hard-to-lyse bacteria, and fungi. Our protocols also include guidance on how to perform quality checks.

### Library preparation: multiplexing samples

Find out more about library preparation:  
[nanoporetech.com/prepare](https://nanoporetech.com/prepare)

To prepare your library for sequencing and downstream analysis, you can choose from either the 24- or 96-plex **Rapid Barcoding Kits**. These PCR-free kits use a transposase to fragment and attach barcodes to your microbial isolate DNA before adding a sequencing adapter. We recommend multiplexing 4–24 microbial isolate genomes per MinION™ Flow Cell or 4–96 genomes per PromethION™ Flow Cell.

Through multiplexed sequencing on a single flow cell, the cost per isolate can be considerably reduced. You can wash and reuse flow cells that are not run at full sample capacity to ensure efficient sample batching while maintaining a low cost per sample. The **Flow Cell Wash Kit** provides a cost-effective method to wash and re-run your flow cell multiple times.



