Oxford Nanopore  
Technologies

# Safety and regulatory information

for Oxford Nanopore Q-Line products – EN

## Use (and misuse)

---

Oxford Nanopore Q-Line products must be exclusively used with other Q-Line products that are on the same product version.

The Oxford Nanopore Q-Line product must only be used with the official Oxford Nanopore Q-Line documentation at **[nanoporetech.com/documentation](https://nanoporetech.com/documentation)**. The instructions in this safety and regulatory document, and the accompanying quick start guide, are to be followed together to use the product as intended by Oxford Nanopore Technologies plc. Notably, the technical specification in the quick start guide contains environmental conditions, which must be adhered to. Safety, regulatory, and performance risks may occur with the product and to the user if the Oxford Nanopore Q-Line documentation is not followed correctly.

Oxford Nanopore Q-Line products have expiry dates and their use once expired is not supported or warranted under the Oxford Nanopore Technologies plc Terms and Conditions.

## Oxford Nanopore Q-Line online

---

Visit: **[nanoporetech.com/documentation](https://nanoporetech.com/documentation)** for Oxford Nanopore Q-Line product documentation.

This document allows safe operation of the Oxford Nanopore Q-Line products when access to the Oxford Nanopore Technologies plc webpages is restricted (e.g. offline access).

## Manufacturing information

---

Oxford Nanopore Q-Line products are made by Oxford Nanopore Technologies plc, Gosling Building, Edmund Halley Road, Oxford Science Park, OX4 4DQ, UK. The product is protected by patents and patents pending: [nanoporetech.com/patents](https://nanoporetech.com/patents)

## Product description

---

The Q-Line GridION is intended for sequencing libraries prepared from native or amplified nucleic acids. It is a long- and short-read nanopore-based sequence analyser platform intended to generate sequence data from DNA and RNA.

DNA and RNA samples can be amplified by PCR or maintained in their native state. A range of fragment lengths can be used, and they can be single or multiplexed. Signals are measured and interpreted by accompanying Q-Line GridION software into DNA and RNA sequence information.

The Q-Line GridION provides all components needed to produce an industry-standard data format (FASTQ files) that contains sequences for downstream bioinformatic analysis.

The Q-Line GridION is not intended for clinical use and should not be used in any diagnostic procedure.

The Q-Line GridION is intended for use by trained laboratory personnel. For Research Use Only.

## Protection and prevention from hazards – categories

---

The safety information provides the details needed to use Oxford Nanopore Q-Line products safely and effectively.

**WARNING:** Indicates a hazardous situation which, if not avoided, could result in death or serious injury. It is important not to proceed until all stated conditions are met and clearly understood.

**CAUTION:** Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury. It is important not to proceed until all stated conditions are met and clearly understood.

**ADVISORY:** Indicates instructions that must be followed to avoid damage to the product or other equipment. The safety notices below are intended to complement and not supersede the normal safety requirements stipulated in the country of use.

## Lifting

---

**CAUTION:** Please refer to the product weight stated in the quick start guide and follow appropriate manual handling procedures based on local Health and Safety requirements.

**CAUTION:** Do not operate any Oxford Nanopore Q-Line product that has been damaged during transit or use.

Contact: **support@nanoportech.com** if necessary.

## Installation

---

**ADVISORY:** If the Q-Line GridION is to be moved or stored, it should be done within its original packaging, within the environmental temperature and humidity conditions stated in the quick start guide. Follow installation instructions in the quick start guide when unpacking.

**ADVISORY:** After transport or storage, the Q-Line GridION should be inspected. If there is any visible damage, the Q-Line GridION must not be used. Prior to sequencing, perform a Position check, as detailed in the quick start guide.

**ADVISORY:** Place the Q-Line GridION on a flat and stable surface.

**CAUTION:** Ensure any cables that are to be plugged into the Q-Line GridION are not stretched. Taut cables may pose electrical and physical risks.

**CAUTION:** Only use the bundled power cable shipped with the Q-Line GridION. Do not use any other power cable that has not been supplied by Oxford Nanopore Technologies plc to power the Q-Line GridION.

**CAUTION:** Allow 30 cm clearance at the rear and sides of the Q-Line GridION.

**ADVISORY:** Avoid exposing Oxford Nanopore Q-Line products to direct sources of heat or sunlight, unless otherwise stated in the product documentation.

## Power supplies

---

**WARNING:** The power cable must connect to the mains/Universal Power Supply (UPS). Only use the mains cable supplied, or a suitable, approved, and rated power cable provided by Oxford Nanopore Technologies plc. Ensure the Q-Line GridION is connected to a protective earth mains supply; only then, turn the power source on.

## General precautions

---

**CAUTION:** Observe the safety regulations of your region when handling toxic, radioactive, or pathogenic samples, as defined by the WHO Laboratory Biosafety Manual (available here: [www.who.int/publications/i](http://www.who.int/publications/i) by searching "WHO Laboratory Biosafety Manual").

**ADVISORY:** The Q-Line GridION is designed to operate up to an altitude of 2000 m above sea level.

## Use and performance

---

**ADVISORY:** Oxford Nanopore Q-Line products must be shipped, stored, and used within the correct temperature and humidity ranges (if stated); these can be found using the information below. Products must be kept in these temperature and humidity ranges, as applicable, to ensure performance risks are minimised.

You will find the following Oxford Nanopore Q-Line products' environmental conditions:

- In the technical specification of the accompanying quick start guide (for Q-Line GridION).
- On the product packaging labels, which use the symbols at the end of this document (for Oxford Nanopore Q-Line reagents).
- The flow cell pouch, which uses the symbols at the end of this document (Q-Line GridION Flow Cells).

**CAUTION:** The Oxford Nanopore Q-Line products are designed for use in a professional laboratory facility. Oxford Nanopore Q-Line products are to be used by trained and authorised personnel only. All Oxford Nanopore Q-Line products are not intended for use in locations where children may be present.

**CAUTION:** If it is suspected that performance is affected by electromagnetic interference, correct operation may be restored by increasing the distance between the equipment and the source of the interference.

**ADVISORY:** The total storage capacity of the Q-Line GridION is 7 TB. Ensure the assay definition configuration files contain an accurate estimation of output for a given assay to receive appropriate warnings about storage limitations before starting.

**ADVISORY:** Oxford Nanopore Q-Line reagents must be stored under the indicated storage conditions before and after use.

**CAUTION:** Do not use Oxford Nanopore Q-Line reagents or flow cells if the package is damaged and/or the reagent tubes or flow cells are leaking (this may mean reagents could be compromised and may not meet performance claims). Contact: **support@nanoporetech.com**

**CAUTION:** Do not use Oxford Nanopore Q-Line reagents or flow cells after their expiration dates.

**CAUTION:** The Q-Line GridION must be checked before use and must not be operated if visibly damaged.

**CAUTION:** This Q-Line GridION is not to be opened or tampered with, unless by authorised and trained engineers. Contact Oxford Nanopore Technologies plc support team at: **support@nanoporetech.com**

**CAUTION:** Opening the Q-Line GridION can expose users to components that may cause burns and other physical harm leading to serious injury. Do not open the equipment. Contact your local Oxford Nanopore Technologies Field Service Engineer.

**CAUTION:** All Q-Line GridION lids must be fully closed before operation.

**ADVISORY:** Please contact Oxford Nanopore Technologies plc support team for further information on the Q-Line GridION preventative maintenance and health checks.

**ADVISORY:** Ensure no flow cells are sequencing, or reagents are nearby when cleaning the Q-Line GridION. Clean the Q-Line GridION with 70% ethanol in water on a lint-free cloth. Use deionised water and a tissue to wipe down the surface(s) afterwards. Dispose of any cleaning materials according to your local waste disposal regulations.

## Personal protection

---

**WARNING: IN EMERGENCY** – Immediately disconnect the Q-Line GridION from the mains power supply if the device is emitting more heat or hot air than usual and this is noticed without touching its surface. Contact [support@nanoporetech.com](mailto:support@nanoporetech.com).

**CAUTION:** The adjacent symbol indicates a hot surface. Do not touch. Surface will remain hot during and after operation.

**WARNING:** Wear gloves (compliant with EN374 or equivalent), safety glasses (compliant with EN166 or equivalent), and a laboratory coat when handling chemicals. Ensure you have read the SDS for the chemical(s) being handled and adhere to your facility's chemical management guidance before their use.

**CAUTION:** Specimens and reagents containing materials from humans should be treated as potentially infectious. Use laboratory safety procedures to handle human materials according to your national, provincial, and local regulations and laws for biosafety, such as: Biosafety in Microbiological and Biomedical Laboratories at [www.cdc.gov/labs/bmbl/index.html](http://www.cdc.gov/labs/bmbl/index.html)

**CAUTION:** Take all necessary actions to avoid spreading hazardous agents in the vicinity of the Q-Line GridION. The facility should follow national, provincial, and local code of practice for biosafety regulations and laws.

**CAUTION:** Samples being loaded into the flow cell should be used, stored, and disposed of according to required safety regulations. Consult the responsible authority in your laboratory to comply with your national, provincial, and local regulations and laws.

**CAUTION:** Samples containing infectious agents should be handled with the greatest of care and in accordance with the required national, provincial, and local regulations and laws.

**ADVISORY:** It is good laboratory practice to always wear safety glasses, gloves (2 pairs, compliant with EN374 or equivalent), if working with infectious agents, and a laboratory coat. There may be other locally advised items which can be added to this recommended list. Consult the responsible delegate for safety in your laboratory for processes adhering to national, provincial, and local regulations and laws.

## Safety Data Sheets (SDS)

---

**ADVISORY:** The SDSs are available from the Oxford Nanopore Store at: [store.nanoporetech.com](https://store.nanoporetech.com)

**CAUTION:** Read and understand the SDSs before handling, operating, or storing the substances and mixtures being used within the Q-Line GridION.

**CAUTION:** Follow appropriate H&S guidelines: protective clothing, safety glasses, and gloves (compliant with EN374 or equivalent). SDSs will carry specific requirements.

**CAUTION:** Use of appropriate ventilation is highly recommended.

**CAUTION:** Check for any spills or leakages during operation. Follow the clean-up guidelines provided in Section-6 of the SDS if required.

**CAUTION:** The components of the Oxford Nanopore Q-Line products should be handled, stored, and disposed of in accordance with national, provincial, and local regulations and laws.

## Security

---

**ADVISORY:** For security purposes, sample sheets created and used with the Q-Line GridION should not contain any personal identification information.

**ADVISORY:** The user is responsible for ensuring the accuracy of sample sheets before use to mitigate risk of assigning results to the incorrect sample.

**ADVISORY:** When accessing data through the file browser or terminal, it is the user's responsibility to ensure that data copied or moved is stored securely.

**ADVISORY:** Store and operate the Q-Line GridION in a physically secure environment. When leaving the Q-Line GridION unattended, it should be locked, and user accounts should be logged out.

## Maintenance

---

**CAUTION:** Repairs must only be performed by an Oxford Nanopore Technologies plc Field Service Engineer, and components should not be replaced by the user. Contact: **support@nanoporetech.com** if the Oxford Nanopore Q-Line GridION or the flow cells become damaged.

**ADVISORY:** If in doubt, contact Oxford Nanopore Technologies plc to ensure that your intended decontamination method will not damage the Oxford Nanopore Q-Line GridION and/or flow cell before using cleaning agents, except those stipulated by Oxford Nanopore Technologies plc.

## Disposal and recycling instructions

---

**ADVISORY:** Visit [nanoporetech.com/support/customer-support/flow-cell-and-devices-returns/flow-cell-returns](https://nanoporetech.com/support/customer-support/flow-cell-and-devices-returns/flow-cell-returns) for information on how to flush and return your used flow cells.

**ADVISORY:** Visit [nanoporetech.com/support/customer-support/flow-cell-and-devices-returns/device-returns](https://nanoporetech.com/support/customer-support/flow-cell-and-devices-returns/device-returns) and follow the instructions on how to return the Oxford Nanopore Q-Line GridION to Oxford Nanopore Technologies plc.

**CAUTION:** Used plasticware which includes reagents, tubes, and pipette tips, must be collected, and disposed of properly, in accordance with national, provincial, and local regulations and laws.

## Flow cell and reagent use

---

**CAUTION:** Non-prescribed use of the Oxford Nanopore Q-Line flow cell(s) may result in leakage of hazardous substances.

**CAUTION:** Observe the specific local regulations when handling toxic, radioactive, or pathogenic samples.

**CAUTION:** Loading excess buffers, samples, or deionised water to the flow cell will cause an overflow of the waste compartment. Use absorbent material to absorb samples and buffers which will come out through the waste port. Dispose of all material in line with local regulations for biological waste.

**CAUTION:** The flow cell buffer, the sample preparation kit buffers, and flow cell wash kit buffers must not be mixed in ways other than those stated in the Oxford Nanopore Q-Line protocols and must be kept away from strong acids and alkalis.

**CAUTION:** The flow cell buffer, flow cell wash kit buffers, and sample preparation kit buffers must be disposed of according to the local regulations. The used liquid product must not be disposed of through domestic water waste infrastructure (e.g. a sink).

**CAUTION:** The Terms and Conditions for use of an Oxford Nanopore Q-Line product stipulate any flow cells that have been used with or otherwise been in contact with materials of Biohazard Level 3 and above (correlating to equivalent global biohazard categories), referenced as: contaminated flow cells, must not be returned. Biohazard Level 2-exposed flow cells (or equivalent) should be labelled with a biohazard symbol when returned to Oxford Nanopore Technologies plc. Proof of legal and appropriate destruction of any contaminated flow cells will be required.

**CAUTION:** The Q-Line GridION shall be decontaminated before decommissioning, and all local regulations for electronic and electrical waste shall be followed regarding disposal of the components if they are not being returned to Oxford Nanopore Technologies plc.

**CAUTION:** Flow cells containing human DNA/RNA shall not be returned to Oxford Nanopore Technologies plc and shall be disposed of as hazardous biological waste following national, provincial, and local regulations and laws.

## **Regulatory – Electromagnetic Compatibility (EMC), radio, and electrical safety**

---

The Q-Line GridION is occasionally referred to as the 'device' for ease and consistency in the regulatory section of this document.

## **Declaration of Conformity**

---

See our Declarations of Conformity (DoC) in the Oxford Nanopore Q-Line GridION documentation at [nanoporetech.com/documentation](https://nanoporetech.com/documentation) to download the Q-Line GridION DoC.

## General – requirements

---

1. The Q-Line GridION emissions may exceed the limits stated in BS EN 61326 when connected to auxiliary equipment.
2. Ethernet cables connected to the Ethernet port must not be longer than 100 m. Cables from other signal ports, other wired network ports, control ports, and DC power ports of products containing radios and auxiliary equipment must not be longer than 3 m.
3. The Q-Line GridION is intended to be operated under a controlled electromagnetic (EM) environment, such as a laboratory. Do not use the Q-Line GridION in proximity to sources of strong electromagnetic radiation (e.g. unshielded intentional RF sources), as these can interfere with proper operation.
4. The Q-Line GridION has minimum performance criteria in Electromagnetic Interference (EMI). The Q-Line GridION shall:
  - Resume operation after powering the Q-Line GridION off and on again.
  - Continue to acquire data in the Sequencing Software.
  - Save the sequencing files (for example, FASTQ) to the hard drive on the Q-Line GridION.
5. Use of the Q-Line GridION in a dry environment, especially if synthetic materials are present (e.g. synthetic clothing, etc.) can cause electrostatic discharge damage and errors on your Q-Line GridION.

## **EMC – US**

---

The device is tested and certified to Class A limits, pursuant to Part 15 of the FCC rules. See 'Declaration of Conformity'.

The device is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

## **FCC Part 15 Statement – unintentional transmitter Class A**

---

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their own expense.

## **Electrical Safety – US**

---

The device is tested and certified to UL/CSA C22.2 No 61010-1 3rd ed. See 'Declaration of Conformity'.

## EMC – EU

---

The device complies with the CISPR 11 limits and has been tested to EN 61326-2-6. See 'Declaration of Conformity'.

## Electrical Safety – EU

---

See the 'Declaration of Conformity'.

## Considerations

---

**WARNING:** Any changes or modifications to the device not expressly approved by Oxford Nanopore Technologies plc could void your authority to operate the equipment.

### Manufactured by:

Oxford Nanopore Technologies plc  
Gosling Building, Edmund Halley Road  
Oxford Science Park OX4 4DQ  
United Kingdom

## Use symbols

---

Catalogue  
number

Caution

Caution,  
hot surface

CE mark

Consult product  
manual

Serial  
Number

Date of  
manufacture

Distributor

Do not  
use if packaging  
has been  
damaged

EU  
importer

Fragile  
contents

Should not be  
discarded as  
unsorted waste

## Use symbols

---

Store and handle  
this way up

Keep dry

Lot number

Manufacturer

Temperature  
limit

Transport/storage  
humidity limits

Use by date

# Oxford Nanopore Technologies

**phone** +44 (0)845 034 7900

**email** support@nanoporetech.com

**X** @nanopore

Gosling Building, Edmund Halley Road,  
Oxford Science Park, Oxfordshire, OX4 4DQ, UK

**www.nanoporetech.com**

Oxford Nanopore Technologies, the Wheel icon, and GridION, are registered trademarks or the subject of trademark applications of Oxford Nanopore Technologies plc in various countries. Information contained herein may be protected by copyright, patents or patents pending of Oxford Nanopore Technologies plc. All other brands and names contained are the property of their respective owners. © 2026 Oxford Nanopore Technologies plc. All rights reserved. Oxford Nanopore Technologies products are RUO. Products labelled/branded as Oxford Nanopore Diagnostics may be RUO or may be regulated as in-vitro diagnostic devices in some jurisdictions, please check individual product labelling.

Safety and regulatory information for Oxford Nanopore Technologies Q-Line products.

ONT-08-00898-00 Rev. 2 | SD\_1142(EN)\_V2\_30Apr2026

# Safety and regulatory information

for Oxford Nanopore Technologies Q-Line products – Multilingual

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**[www.nanoporetech.com](http://www.nanoporetech.com)**

Oxford Nanopore Technologies, the Wheel icon, and GridION, are registered trademarks or the subject of trademark applications of Oxford Nanopore Technologies plc in various countries. Information contained herein may be protected by copyright, patents or patents pending of Oxford Nanopore Technologies plc. All other brands and names contained are the property of their respective owners. © 2026 Oxford Nanopore Technologies plc. All rights reserved. Oxford Nanopore Technologies products are RUO. Products labelled/branded as Oxford Nanopore Diagnostics may be RUO or may be regulated as in-vitro diagnostic devices in some jurisdictions, please check individual product labelling. Safety and regulatory information for Oxford Nanopore Technologies Q-Line products.

ONT-08-01278-00 Rev. 2 | SD\_1283(ML)\_V2\_30Apr2026