

Hematoxylin and Eosin (H&E) Staining – Manual Protocol

(From Baylor College of Medicine)

Protocol for H&E Staining

- Place slides containing paraffin sections in a slide holder (glass or metal)

- Deparaffinize and rehydrate sections:

3 x 3' Xylene (*blot excess xylene before going into ethanol*)
3 x 3' 100% ethanol
1 x 3' 95% ethanol 1 x 3' 80%
ethanol
1 x 5' deionized H₂O

- While sections are in water, skim surface of hematoxalin with a Kimwipe to remove oxidized particles. Blot excess water from slide holder before going into hematoxalin.

- Hematoxalin staining:

1 x 3' Hematoxalin
Rinse deionized water
1 x 5' Tap water (*to allow stain to develop*) Dip 8-
12x (fast) Acid ethanol (*to destain*) Rinse 2 x 1' Tap
water
Rinse 1 x 2' Deionized water (*can leave overnight at this stage*)
(Blot excess water from slide holder before going into eosin)

- Eosin staining and dehydration:

1 x 30 sec Eosin (*up to 45 seconds for an older batch of eosin*)
3 x 5' 95% ethanol
3 x 5' 100% ethanol (*blot excess ethanol before going into xylene*) 3
x 15' Xylene

- Coverslip slides using Permunt (xylene based).
- Place a drop of Permunt on the slide using a glass rod, taking care to leave no bubbles.
- Angle the coverslip and let fall gently onto the slide. Allow the Permunt to spread beneath the coverslip, covering all the tissue.
- Dry overnight in the hood.

Reagents for H&E Staining

- **Acid Ethanol:** 1 ml concentrated HCl + 400 ml 70% ethanol
- **Hematoxylin:** Poly Scientific (Bayshore, NY) #s212A
Harris hematoxylin with glacial acetic acid
- **Eosin:** Poly Scientific (Bayshore, NY) #s176
Eosin Phloxine stain, working
- **Permunt:** Fisher Scientific #SP15-100
Histological mounting medium