Ultrasonic Bone Surgical System

VarioSurg 4

OPERATION
MANUAL

CE 0197 MADE IN JAPAN

English	1
Deutsch.....	53
Français.....	109
Español	165
Italiano.....	221
Português.....	277

TABLE OF CONTENTS

English

1	Preface	3
1.1	User and intended purpose.....	3
1.2	Related manuals	3
1.3	Precautions for handling and operation	3
1.4	Symbol.....	8
1.5	Disposing of products	11
1.6	Warranty.....	11
2	Product Description	12
2.1	Package contents.....	12
2.2	Applicable products	13
2.3	Part names.....	14
2.4	Control panel	17
3	Preparation for Use	20
3.1	Connecting AC power cord	20
3.2	Installing batteries (Only for wireless foot control).....	20
3.3	Installing foot control hanger (wireless foot control).....	21
3.4	Connecting foot control	21
3.4.1	Wireless foot control.....	21
3.4.2	Wired foot control (Option).....	21
3.5	Mounting coolant solution hanger post	22
3.6	Installing irrigation tube.....	22
3.7	Mounting irrigation tube onto bag / bottle	23
3.8	Connecting handpiece cord.....	24
3.9	Attaching tip.....	24
3.10	Mounting tube holder.....	24
3.11	Handpiece stand	25
3.12	Check before treatment.....	26
4	Operation Procedure	27
4.1	Setting the power level.....	27
4.2	Activation	27
4.3	Sleep mode (Only for wireless foot control).....	28
4.3.1	About sleep mode	28
4.4	Memory function	28
5	Post-use Maintenance	29
5.1	Preparation	30
5.2	Cleaning handpiece irrigation path.....	30
5.3	Cleaning and disinfection	31
5.4	Sterilization	31
5.5	Storage	32
6	Maintenance.....	33

6.1	Periodical maintenance checks.....	33
6.2	Foot control calibration.....	34
7	Setting.....	35
7.1	Setting key.....	35
7.2	Pairing with wireless foot control	38
7.3	Programming system operation	39
7.4	Link function.....	40
8	Troubleshooting	42
8.1	Error code.....	42
8.2	Problems and solutions.....	43
9	Spare Parts and Option List.....	45
10	Specification.....	46
10.1	System specifications.....	46
10.2	Bluetooth specifications.....	47
10.3	Classification of equipment	47
10.4	Operation principle	47
10.5	Factory settings	48
11	EMC Information (Electromagnetic Compatibility Information)	49

1 Preface

In this operation manual, VarioSurg 4 is called “this product”.

1.1 User and intended purpose

- User

Qualified professionals (dentist, oral surgeon, hygienist)

- Intended purpose

The dental surgical device is intended for the following application(s):

Removal of accretions and deposits from the surface of teeth and dental restorations (fillings and prostheses), root canal preparation for retrograde endodontics (including removal of root canal filling materials), cutting, removal, incision, and perforation of jawbones, detachment of periodontal membrane due to tooth extraction, and detachment of mucosal membrane due to implant surgery.

1.2 Related manuals

For proper use with the combined items, refer to the following operation manuals:

Operation manual title	Operation manual No.
Ultrasonic Bone Surgery System VarioSurg 4	OM-DE1165EN
VarioSurg 4 Handpiece REPROCESSING MANUAL	OM-DE1136MA
Varios Tips VarioSurg Tips OPERATION MANUAL	OM-Z0288E
VarioSurg Tips OPERATION MANUAL	OM-Z0315E
Link Stand3 OPERATION MANUAL	OM-DE1158MA
Surgic Pro2 OPERATION MANUAL	OM-E0858E

1.3 Precautions for handling and operation

- Read these safety cautions thoroughly before use and operate the product properly.
- These indicators are to allow you to use the product safely, prevent danger and harm to you and others. These are classified by degree of danger, damage and seriousness. All indicators concern safety, be sure to follow them.

Class	Degree of risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

WARNING

Patients and user this product is inapplicable to

- Do not use on the following patients.
 - Patients with medical complications.
 - Patients who have pre existing conditions (e.g., cardiac, pulmonary, renal disturbance, or high blood pressure).
 - Patients who are pregnant, potentially pregnant, or breastfeeding.
 - Patients (including infants) and users with external/implantable pacemakers and implantable cardioverter defibrillators.
- It may affect the operation of external/implantable pacemakers and implantable cardioverter defibrillators.

Environment

- Do not use in rooms with a risk of explosion or near flammable substances.
Also, do not use on or near patients under flammable anesthesia (nitrous oxide).
- Use this product in a cool and dry environment that meets the requirements for the use environment. Otherwise, it may cause malfunction. => “■ Use, transportation and storage environment” in “10.1 System specifications”.

Designated items

- Use NSK's genuine tip (designated item). The use of an item other than those designated may cause the following faults or unexpected accidents:
 - Inappropriate vibration caused by forcible attachment of another manufacturers tip.
 - Accidental ingestion of a broken tip by the patient.
 - Damage or premature wear of the handpiece threads.
- Problems including handpiece damage, failure, and accidents resulted from the use of non-NSK tips are not covered in the warranty.
- Do not use non-NSK genuine irrigation tube. The use of an item other than those designated may cause the following faults or unexpected accidents:
 - Failure to ensure the amount of water injected as stated in the output setting may result in overheating of the handpiece, which may cause burns. => “4.1 Setting the power level”.
 - Water leak from the irrigation connection part or irrigation tube.
- Be sure to plug the AC power cord into a wall outlet specified for medical use with a dedicated ground terminal. If a short circuit occurs inside this product, it may cause an electric shock and/or fire.
- Grounding reliability can only be achieved when the equipment is connected to an equipment receptacle marked “Hospital Only” or “Hospital Grade”. If there is no “Hospital Only” or “Hospital Grade” receptacle, make sure that the equipment is grounded.
- Do not handle the AC power cord and the electric device with wet hands. Doing so may cause an electric shock.

Before use

- About sterilizable products:
 - After purchasing or repairing this product, clean, disinfect, and sterilize it following the instructions in this manual. => “5 Post-use Maintenance”.
- About irrigation tube (sterile item)
 - The irrigation tube is a single-use item. Do not reuse nor re-sterilize it as product breakage or infection could occur.

⚠ WARNING

- Check before use that the sterile package is not open or damaged.
- If the package is open or damaged, do not use irrigation tube as product breakage or infection could occur.
- Observe the use-by date written on the irrigation tube package. Do not use expired irrigation tube as product breakage or infection could occur.
- This product is a medical device. Observe this manual's instructions concerning the user, intended purpose, and use method when using this product.
- For the sake of safety of the patient and medical professionals, wear eye protection, a mask, and gloves during the use of this product.
- LED handpiece with cord, tip and tip wrench require reprocessing prior to first use. For reprocessing, refer to the OM-DE1136MA / OM-Z0288E / OM-Z0315E OPERATION MANUAL listed in “1.2 Related manuals”.

During use

- If the product produces smoke or smells of burning resin, immediately turn off the power, unplug the AC power cord, and then contact your authorized NSK dealer.
- Do not turn on and off the main power switch repeatedly. Doing so may blow a fuse.
- If water droplets, chemicals, solvents, liquid disinfectant, or dirt, etc. come into contact with the control unit, handpiece, handpiece cord, etc., wipe them off immediately after the use. Failure to do so may result in discoloration or deformation.
- Be careful not to get water droplets, chemicals, solvents, or liquid disinfectant in the control unit and the AC power cord. This could cause a short circuit leading to an electrical shock or fire. If it is splashed on the control unit, turn off the power supply, disconnect the power cord, and wipe it off immediately. Failure to do so may result in discoloration or deformation.
- The handpiece may get hot when some types of tips are used at high power level or in an overloaded state. (The temperature of the handpiece surface might exceed 41°C.) If such abnormality is found, stop using for about 5 minutes to allow the handpiece to cool.
- The temperature of the tip surface may be more than 65°C when used without coolant solution. To avoid this event, be sure to irrigate with coolant solution.
- Be sure to irrigate with coolant solution during the usual use. If the coolant water is insufficient, the handpiece may be overheated or the bones and tissues may be damaged.
- Do not aim the light of the handpiece at the eye of the patient, operator, or assistant. Doing so may damage eyesight.
- This product is a medical device. Observe this manual's instructions concerning the user, intended purpose, and use method when using this product.

⚠ CAUTION

Environment

- This product is for indoor use only.

Installation

- Do not drop, hit, or subject to excessive shock. Especially for the handpiece, the glass rod and handpiece may be damage.

CAUTION

- This product may malfunction when used in the presence of an electromagnetic interference wave. Do not install this product in the vicinity of any device which emits electromagnetic waves. Turn off this product when using an ultrasonic vibration device or an electrode knife near this product.
- This product can be stacked when connected to Surgic Pro2 using Link Stand 3 (option).
Should not be used adjacent to, or stacked with, other equipment.
If adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.
- Portable and mobile RF (Radio Frequency) communications equipment can affect this product. Do not use RF equipment in the vicinity of the product.
- For safety reasons, install the control unit in a place providing space of approximately 10cm around the product and outlet for where the AC power cord can be easily removed. (It is possible to disconnect the control unit from the power source by removing the AC power cord.)
- Install the control unit in a level and stable location.

Connect

- Before connecting or disconnecting the handpiece cord or irrigation tube make sure to remove the tip first to avoid injury by the tip.
- Only connect the VarioSurg 4 handpiece to this product. Do not connect VarioSurg 4 handpiece to VarioSurg 3.

Before use

- When operating the product, always consider the safety of the patient.
- The end user shall be responsible for any judgment that relates to the application of this product to a patient.
- This product does not consider patient's age (except infants), gender, weight or nationality.
- This product does not consider operator's age (mature person), height, weight, gender, or nationality.
- The user is responsible for the operational control, maintenance, and continual inspection of this product.
- Read this operation manual before use to fully understand the product functions and file for future reference.
- Make sure each component is securely connected before turning on the power.
- If the product has not been used for a long period, check that it is functioning correctly and safely before using on a patient.

During use

- Should the product function abnormally, stop using the product immediately and contact your authorized NSK dealer.
- If the irrigation pump or the irrigation tube get wet, wipe it dry. If it remains wet, the roller in the pump may slip and the pump may not operate properly.
- Do not bend or fold the irrigation tube while the irrigation pump is operating as the irrigation tube may be damaged or may detach.
- If abnormality is detected in irrigation, it may be due to wear of the irrigation tube or leaking of saline from the tube. In such cases, replace the irrigation tube.
- Make sure that cords and other parts do not touch the control panel while the control unit power is on. It may result in malfunction.

⚠ CAUTION

- Never use the device continuously for over 1 minute. Between use, please allow a period of 5 minute rest time. Continuous use over a long time may lead to overheating of the handpiece, which in turn may result in burns. (The maximum surface temperature might reach 50°C depending on by the situation in which this product is used.) In such cases, stop using the handpiece until the cools down.

Operating Mode

Non-continuous Operation	
ON	OFF
1 min	5 min

- Do not point the handpiece LED light directly to the eyes of the patient, or the operator, as it may be harmful to eyes.
- If any abnormalities (dark, does not light, flashes, etc.) are found with the LED, stop using it immediately and contact an authorized NSK dealer.
- Do not press any of the keys on the control panel with a pointed object (ball-point pen, etc.). Doing so may result in malfunction.

After use

- After using this product, turn off the power switch. When you intend to leave this product without using for a long time, unplug the AC power cord from the control unit. This will reduce power consumption and prevent accidents due to unintended operation.

Post-use Maintenance

- Immediately after using on a patient, clean, disinfect, and sterilize this product with the correct method. => “5 Post-use Maintenance”.

Maintenance

- Do not attempt to disassemble the product nor tamper with the mechanism except as recommended by NSK in this operation manual.
- To avoid clinical downtime, it is recommended that a spare be kept on hand in case of a breakdown during surgery.

Check

- Be sure to conduct periodical maintenance check on the device and parts.
- This product is rated medical electrical equipment. EMC (Electromagnetic compatibility) is described in the documentation included.
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information. => “11 EMC Information (Electromagnetic Compatibility Information)”.
- The use of ACCESSORIES such as handpieces and cables, other than those specified by the manufacturer, with the exception of handpieces and handpiece cord sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- Do not use parts other than the components for this product listed in this operation manual.

Report

- Please report any serious incident that may occur in connection with this product and its accessories to the manufacturer as well as to your national authority.

*Emission: Unnecessary energy generated by the device.

1 Preface

*Immunity: Tolerance to unnecessary energy received by the device.

NOTICE

- No special training is required for this device.
 - The handpiece or handpiece cord during operation could interfere with computers, LAN cables in close vicinity, or could cause noise in radio receivers nearby.
 - For service requirements and spare parts, contact an authorized NSK dealer.
-

1.4 Symbol

Battery cover closed

Battery cover open

The instrument can be used with a washer-disinfector for thermal disinfection

This product can be sterilized in a steam sterilizer at 135°C

Consult operation instructions

Caution

Type BF applied part

Marking on the outside of equipment or equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment

Protected against the effects of continuous immersion in water

This product is designed not to become the ignition source in air and flammable anesthetic gas

Dispose of this device and its accessories via methods approved for electronic device and in compliance with the Directive 2012/19/EU

Conforms to European "Medical Device Regulation"

Mandatory conformance mark on many products placed on the market in the European Economic Area

TUV Rheinland of North America is a Nationally Recognized Testing Laboratory (NRTL) in the United States and is accredited by the Standards Council of Canada to certify electro-medical products with Canadian National Standards

Regulatory compliance mark (RCM) indicates compliance with electrical safety, EMC, EME and telecommunications requirements for Australia and New Zealand

Korean KC compliance marking for wireless product

Brazil ANATEL compliance marking for wireless product

National Communications Commission (NCC)
Taiwan NCC compliance marking for wireless equipment

Manufacturer

Importer

Authorized representative in the European Community

Authorised representative for Switzerland

Medical Device

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician

Catalog number (Order code)

Serial No.

Batch code

Date of manufacture

Unique Device Identifier

GS1 DataMatrix for Unique Device Identifier

1 Preface

Temperature limitation

Humidity limitation

Atmospheric pressure limitation

Keep dry

Fragile / handle with care

This is the correct upright position of the distribution packages for transport and/or storage

MODEL

SUPPLY VOLTAGE

FREQUENCY

POWER INPUT

MODE OF OPERATION

Contains FCC ID

Contains IC

MADE IN JAPAN

NON-CONTINUOUS OPERATION

ON: 1min

OFF: 5min

1.5 Disposing of products

When your medical device has reached the end of its service life, it must be disposed of via a specialist collection, removal, recycling or destruction channel. Prior to doing so, in order to avoid the health risks to operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, you need to confirm the equipment is decontaminated.

Therefore contact your nearest NSK dealer, or one of the NSK branches to find out how to proceed. Refer to the link below:


1.6 Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model. Contact an authorized NSK dealer if repairs are necessary.

2 Product Description

2.1 Package contents

No.	Part name	Quantity
1	AC power cord	1
2	Wireless foot control ^{*1}	1
3	Control unit	1
4	Coolant solution hanger post	1
5	LED handpiece (with 2m handpiece cord) ^{*2}	1
6	Wired foot control connector cap	1
7	Handpiece stand	1
8	AAA Alkaline battery 3 pack	1
9	Tube holder	7
10	Tip wrench	1

No.	Part name	Quantity
11	Irrigation tube	5
12	SG3	1
	SG5	1
	SG6D	1
	SG7D	1
	SG11	1
	H-SG1	1
13	Sterilization cassette	1
14	Tip holder	1

*1: Wireless foot control is not included in the “without foot control” set.

*2: LED handpieces are referred to as handpieces in the operation manual.


2.2 Applicable products

Products that can be used in combination with the device


Product name
VarioSurg Tip
Irrigation Tube
Link Stand3
Surgic Pro2

2.3 Part names


■ Control unit

■ LED handpiece with handpiece cord


■ Wireless foot control


2 Product Description

■ Wired foot control (Option)

2.4 Control panel

No.	Key	Name	Display	Function
1		POWER key	A	Set the power level. Power level range [PERIO] : 5 – 100 % [ENDO] : 5 – 100 % [SURG] : 10 – 150 %
2		PROGRAM key	B	Select the available programs. Recalls the setting values stored in the program. Program No. [PERIO] : 8 and 9 [ENDO] : 6 and 7 [SURG] : 1 to 5
3		MEMORY key	-	Memorize the program parameters. (E.g. mode, program No., power level, burst output level, coolant solution flow level, and light level)

2 Product Description

No.	Key	Name	Display	Function
4		BURST key	D	Select the burst output level (Only SURG mode). No display: BURST off  : BURST1 (Frequency 10 Hz)  : BURST2 (Frequency 30 Hz)  : BURST3 (Frequency 60 Hz)
5		MODE key	E	Selects the mode to be used.  mode: Maintenance  mode: Root canal treatment  mode: Osteotomy
6		CLEANING key	-	Pre-cleans the inside of the handpiece. => "5.2 Cleaning handpiece irrigation path".
7		Setting key	-	Select the settings of the control unit, wireless foot control or a Bluetooth connection. => "7.1 Setting key".
8		Enter key	-	
9		Cancel key	-	
10		COOLANT key	H	Select the coolant solution flow level. Max.  Min.
11		LIGHT key	I	Select the light level 
12	-	-	F	The relevant power level is indicated by the number of indicator bars illuminated. For each 10% 1 bar  mode: Max.10 bars:  For each 50% 1 bar, Max 2 bars:   mode: Max 15 bars:  For each 50% 1 bar, Max 3 bars: 
13	-	-	G	No display: No vibration  : Vibration

No.	Key	Name	Display	Function
14	-	-	J	Displays the battery level of the wireless foot control. Check the battery level before and after treatment. Replace the batteries when the battery display level is flashing.  : 50-100%  : 10-50%  (flashing): Less than 10%  : 0%
15	-	-	J,K	Status of connection with the wireless foot control No display: Not connected  : Connected  (flashing): Not connected, paired with the wireless foot control
16	-	-	L	Status of connection with the wired foot control (option) No display: Not connected  : Connected
17	-	-	C	Used during wireless foot control calibration. => "6.2 Foot control calibration".
18	-	-	M	 LINK: Connected with Surgic Pro2. Surgic Pro2 operates when depressing the foot control.  (flashing): Connected with Surgic Pro2. VarioSurg 4 operates when depressing the foot control. No display: Not connected


NOTICE

- When the transparent protective sheet on the control panel is peeled off or when an object charged with static electricity is placed near the LCD, fine lines may appear in the display area. This is normal and the lines will disappear after a short time.

3 Preparation for Use

3.1 Connecting AC power cord

Correctly align then insert the AC power cord into the AC power cord connector socket at the back of the control unit.


⚠ CAUTION

- Hold the plug when pulling out the AC power cord. Pulling the cord may cause a broken wire or damage.
- Be sure to stop the device and turn off the main power switch before plugging or unplugging the AC power cord and the handpiece cord.

3.2 Installing batteries (Only for wireless foot control)

1. Turn the cover to the OPEN  and remove the battery cover.
2. Insert three AAA batteries.
3. Put the battery cover back and turn in the CLOSE  direction to close the cover.

⚠ WARNING

- When replacing batteries during use, keep away from the patient. Touching batteries while contacting the patient may cause electric shock to the operator and patient.

⚠ CAUTION


- When inserting batteries, check the marking on the battery case and make sure batteries are inserted in the correct orientation.
- Use only high-quality disposable AAA (Micro/LR03) 1.5V alkaline batteries. Using the wrong type of battery may lead to product malfunction.
- Do not mix new and old batteries or different types of batteries. Do not use rechargeable batteries. Otherwise, the product may not operate properly or malfunction. If this occurs, contact your authorized NSK dealer.

⚠ CAUTION

- If the device will not be used for a long period of time, remove the batteries for safety reasons. Failure to do so poses a risk of product malfunction due to battery heating or leakage.
- Make sure that the O-ring attached to the cover is free of debris before closing the battery cover.

3.3 Installing foot control hanger (wireless foot control)

1. Insert the hanger into the hanger holes.
2. Push the hanger until it stops.

To remove, pull out the hanger from the hanger holes.

If it is difficult to remove the hanger as it is stiff, pull it out alternating between right and left.

⚠ CAUTION

- When holding the wireless foot control, do not hold with the hanger retracted. It may cause the hanger to come out and the wireless foot control to fall.
- If the hanger is pulled forcibly, the hanger insertion part may be damaged.

3.4 Connecting foot control


3.4.1 Wireless foot control

The wireless foot control is paired at the factory.

If disconnected perform pairing again: => "7.2 Pairing with wireless foot control".

3.4.2 Wired foot control (Option)

1. Remove the connector cap using tweezers.
2. Turn the wired foot control cord plug so that the screw faces down.
3. Insert the plug so that it correctly fits into the wired foot control connector socket of the control unit.
4. Secure the lock nut of the wired foot control cord plug by turning clockwise.


NOTICE

- Operation method of the wired foot control is the same as wireless foot control. When the wired foot control is connected, the wireless foot control will not operate.

3.5 Mounting coolant solution hanger post

Mount the coolant solution hanger post onto the holder of the control unit.


⚠ CAUTION


- Insert the coolant solution hanger post all the way down firmly. If the post guide and the groove guide are not aligned, the coolant solution hanger post cannot be inserted all the way down.
- Use a pack of solution with the capacity of less than 500mL. Maximum pack size is 800g for the coolant solution hanger post.

3.6 Installing irrigation tube


1.

3.

2.


⚠ CAUTION

- Make sure that the irrigation tube is securely set on the rollers when closing the pump cover. If the irrigation tube is not correctly positioned on the rollers and the cover is closed, the irrigation tube could be cut or sheared.
- Only use genuine NSK irrigation tube.
- When mounting the irrigation tube, be careful not to get your hand caught in the moving part of the irrigation pump.
- Do not mount the irrigation tube while the main power switch is ON.



⚠ CAUTION

- If the irrigation tube is not installed in an appropriate location, it may break and cause the following risks:
 - Damage
 - Unable to irrigate
 - Leakage
- Do not use the irrigation tube if the pouch is punctured or the seal is broken. It is not sterile.



3.7 Mounting irrigation tube onto bag / bottle

1. Close the tube clamp between the irrigation tube needle and the irrigation pump.
3. Insert the other irrigation tube tip into the handpiece irrigation nozzle.

2. Insert the irrigation tube needle into the bottle cap.
4. Open the cap to let air into the pack.

5. Open the tube clamp.

⚠ CAUTION

- Do not operate the irrigation pump if the irrigation tube is bent or the tube clamp is in the closed position. Doing so could cause the tube to break or slip out of the bottle.
- Always use a coolant solution pack with a sufficient quantity of saline.


3 Preparation for Use

! CAUTION

- Before opening the tube clamp, always close the pump cover. If the tube clamp is opened while the pump cover is open, saline will flow out of the irrigation tube tip.
- When the coolant solution runs short, stop using the system and replace the coolant solution pack. Be sure to check that the coolant solution is flowing properly before using the system again.

3.8 Connecting handpiece cord

1. Align the [●] mark on the handpiece cord plug with the [◐] mark on the control unit.
2. Push the base of handpiece cord plug into the socket until the lock joint 'clicks' into place.

! CAUTION


- When disconnecting the handpiece cord, hold and pull out the lock joint straight toward from the control unit. Pulling the cord may cause a broken wire or damage.
- Insert the handpiece cord straight and do not put excessive load on the handpiece cord socket as it may cause breakage or deformation.
- Make sure that the handpiece cord plug and the handpiece cord socket on the control unit are free of rust and foreign matter. This may cause the handpiece not to vibrate or the LED not to illuminate.
- Be sure to only connect the handpiece with handpiece cord, that comes with this product.

3.9 Attaching tip

For mounting and removing of the tip, refer to the OM-Z0288E / OM-Z0315E OPERATION MANUAL listed in “1.2 Related manuals”.

3.10 Mounting tube holder


1. Insert the tube holder into the handpiece cord.
2. Push the irrigation tube into the tube holder.

To remove the tube holder, pull out the irrigation tube from the tube holder. Then remove the tube holder from the handpiece cord.

⚠ CAUTION

- Connect the handpiece cord and the irrigation tube at a total of 7 locations.



3.11 Handpiece stand

When the handpiece is not in use, place the handpiece horizontally on the handpiece stand.

⚠ CAUTION

- Make sure that the tip does not make contact the handpiece stand.
- Make sure that the handpiece rests horizontally on the handpiece stand, and not at an angle.

3.12 Check before treatment

Perform the following checks before using the product on a patient to confirm that there are no abnormalities. If abnormal vibration, noise, or overheating is detected, stop using the product and contact your authorized NSK dealer.

1. Make sure each component is securely connected.
2. Confirm that the irrigation tube correctly connected to the coolant solution pack, the irrigation pump and the handpiece.
3. Make sure the control unit main power switch is in the OFF position (O side), and then connect to a medical outlet.
4. Put the control unit main power switch in the ON position (I side).
5. “*” flashes and “CONNECT” is illuminated on the LCD, and the wireless foot control is automatically connected to the paired device.
6. Open the tube clamp.
7. Run the handpiece for approximately more than 10 seconds by depressing the wireless foot control and check the following.
 - Is the tip vibrating correctly?
 - Are there any abnormalities at the tip, e.g. vibration, noise and heat?
 - Can irrigation be performed and is the coolant solution flow level appropriate?
 - Does irrigation solution correctly flow from the tip and is the coolant solution flow level appropriate?
*When a new irrigation tube is used, it may take several seconds before irrigation solution flow commences.
 - Does the handpiece light illuminate?
8. Stop the handpiece and make sure there is no abnormal heating of the tips and the handpiece surface.

*When the wireless foot control is connected, the normal screen is displayed. If it is not switched to the normal screen, => “8 Troubleshooting”.


CAUTION

- Touching the control panel may increase the risk of cross-infections. In order to minimize such a risk, it is possible to use a barrier film in combination with VarioSurg 4. Do not use a film with poor operability and/or visibility.
- Confirm a beep sounds when operating control keys and that a notification tone sounds when starting the operation.

4 Operation Procedure

4.1 Setting the power level

1. Press **MODE** to select the mode.
2. Press **▽PROGRAM△** to select the program No.
3. Press **▽POWER△** to adjust the power level.
4. Press **○** to adjust the coolant flow level.
5. Press **ⓘ** to adjust the handpiece illumination intensity.
6. Press **■** to select the burst output level. (Only SURG mode).

Power outputs at various modes are compared in the above figure.

⚠ CAUTION

- Use the tip **ONLY** in the correct mode and never exceed the maximum power defined for the tip by the tip operation manuals. Refer to OM-Z0288E / OM-Z0315E OPERATION MANUAL listed in “1.2 Related manuals”.
- If used outside the prescribed power range the tip may break, and the tooth surface or soft tissue may be damaged.

4.2 Activation

The tip will commence vibrating when the main pedal (Ultrasonic ON-OFF) of the foot control is depressed.

⚠ CAUTION

- Start with a low power setting within the power range, observing the patient's reaction and according to the treatment plan.
- Always use the irrigation. If the irrigation is not sufficient, the handpiece may excessively heat, or oral injury may occur.
- Do not place the active tip anywhere other than the specific treatment area, as the ultrasonic vibration may cause oral injury.
- Even if the power level is correctly set for the tip in use application of excessive pressure to the tip may cause tip damage. Should tip breakage occur, immediately evacuate the oral cavity of any tip pieces using strong suction.
- If tip vibration reduces during operation, release the main pedal (Ultrasonic ON-OFF) of the foot control, and then depress it again. If the vibration level still does not return to normal check that the tip is correctly and firmly positioned. Incorrectly inserted tips may loosen during operation and power loss may be experienced.
- Do not connect or disconnect the AC power or the handpiece cord while foot control is depressed.
- While the handpiece is in operation the settings cannot be changed.

4.3 Sleep mode (Only for wireless foot control)

4.3.1 About sleep mode

This mode saves the battery of the wireless foot control.

You can select one of 8 times of inactivity to enter sleep mode.

OFF, 5 min, 10 min, 15 min, 20 min, 30 min, 45 min, and 60 min

*If "OFF" selected, sleep mode is disable. Factory default setting is disabled.

*The battery of the foot control continues to be consumed even when there are no operations.

We recommend that you turn on the sleep mode when you want to save battery of the foot control.

■ Sleep mode setting

1. Press and hold  for 2 seconds, "BLE SET" is displayed.
2. Press  to select "FC SET".
3. Press .
4. Press  to select "FC MODE".
5. Press  to select "FC".
6. Press .
7. Press . "FC SET" is displayed.
8. Press .
9. Press  to select "SLEEP".
10. Press  to select from 8 types.
11. Press . Make sure that beep sounds and selected enter sleep type is displayed on the display.
12. Press  twice to end setting mode.

■ Release the sleep mode

When the foot control enters sleep mode, "SLEEP" is illuminate and  flashes on the displayed.

To release the sleep mode, keep depressing the main pedal of the wireless foot control until it returns to the normal screen.

*When the control unit is powered on, "PUSH FC" is displayed on the screen and press the main pedal of the foot control.

4.4 Memory function

The set values displayed in the control panel (E.g. mode, program No., power level, burst output level, coolant solution flow level, and light level) can be memorized in each program. Once memorized, the programs can be called up whenever required.

1. Press  and select the program No. for memory.
2. Check the values of the parameters to be stored.
3. Hold down  for approximately one second. When the notification tone beeps, the setting is memorized.

5 Post-use Maintenance

After each patient maintain the product as follows.

WARNING

- For cleaning and disinfection, use disinfectant ethanol or disinfectant isopropanol that does not contain any additives (any quaternary ammonium salts). Use of other disinfectants may cause discoloration or cracking.
*For details, contact the manufacturer of the disinfectant.
*In this operation manual, disinfectant ethanol or disinfectant isopropanol that does not contain any additives (any quaternary ammonium salts) is called "disinfectant alcohol".
- Use appropriate disinfectant alcohol according to the regulations and guidelines of each country.
- Always obey the local directives, standards, guidelines, etc., as to the cleaning, disinfection, and sterilization of the product.
- Wear eye protection, protective mask, and gloves for safety and infection prevention before starting the maintenance.
- Use only the cleaning and disinfection agent products specified in this operation manual.
- For maintenance of the handpiece with handpiece cord, tip, tip wrench, E tip wrench, VS tip wrench, tip holder, tube holder, handpiece stand, refer to the OM-DE1136MA / OM-Z0288E / OM-Z0315E OPERATION MANUAL listed in "1.2 Related manuals".
- Do not use solvents such as electrolyzed-oxidizing water (strong acidic water or super acidic water), strong acidic or alkaline chemical, chlorine-containing solvent, benzene and thinner when cleaning, immersing, or wiping.
- To remove residues, conduct the following maintenance immediately (within 30 minutes) after the treatment. Failure to carry out this maintenance may cause infection, premature failure, and thermal injury resulting from overheating.
- After thermo disinfection, take the product out of the washer-disinfector immediately (within an hour). Failure to do so causes corrosion.
- Use clean and dry compressed air for cleaning.
- Before disinfection, be sure to clean the product to remove any residues. If blood remains in the product, it coagulates causing a malfunction.
- After the thermo disinfection, dry the product completely. Remaining moisture causes internal corrosion.
- The control unit and foot control are cannot be sterilized.
- Do not sterilize using ultraviolet light. Doing so may result in discoloration.
- The product lifetime and product's ability to operate correctly are mainly determined by mechanical stresses during use and the chemical influences due to processing. If you detect any abnormality such as heat generation, vibration or noise, or if the exterior markings (serial number or model number) have become unreadable before or during use, stop using the product immediately and request your nearest NSK dealer, or one of the NSK branches for a periodic inspection*.
- Disappearance of exterior markings is not covered under warranty for repair.

5 Post-use Maintenance

⚠ WARNING

- When the product has been reprocessed 250 times or has been in use for one year, it is recommended that a periodic inspection be performed by an authorized NSK service*.

*Refer to the link below:

⚠ CAUTION

- Failure to properly maintain the products may cause infection, product failure, or overheating and injuries. Be sure to clean and disinfect the products immediately (within 1 hour) after use to remove residue.
- The following products require reprocessing after each patient treatment. (Handpiece with handpiece cord, tip, tip wrench, E tip wrench, VS tip wrench, tip holder, tube holder, handpiece stand)
- The handpiece cannot be removed from the handpiece cord.
- Removing or twisting the handpiece could cause damage to the internal parts and failure may occur.

5.1 Preparation

- Wear eye protection, a mask, and gloves to prevent infection.

5.2 Cleaning handpiece irrigation path

■ Auto-cleaning

- Remove the irrigation needle from the solution pack.
- Place the irrigation needle in distilled or deionized water in an open container.
- Place the tip of the handpiece into the water.
- Start cleaning by holding down  for approximately 2 seconds.
While cleaning is active “CLEAN” is displayed.
- The bar graph indicates the time remaining (Cleaning lasts for approximately 30 seconds).
*Cleaning can be stopped anytime by pressing .
- Once cleaning has completed, the display returns to normal screen.
- Turn off the power switch of the control unit.
- Remove the tip from the handpiece.
- Remove the handpiece with handpiece cord from the control unit.

*For maintenance of the handpiece with handpiece cord and tips, refer to the OM-DE1136MA / OM-Z0288E / OM-Z0315E OPERATION MANUAL listed in “1.2 Related manuals”.

⚠ WARNING

- Perform cleaning after every use. If cleaning is omitted, the Handpiece may clog or malfunction.

5.3 Cleaning and disinfection

- Control unit and foot control

Wipe the surface of the products with a cloth moistened with disinfectant alcohol.

- Handpiece stand and tube holders

Cleaning

1. Clean the surface of the products under running water with a soft bristled brush for 15 seconds.
Water temperature: $\leq 38^{\circ}\text{C}$
Water quality: Drinking water
Water flow rate: $\geq 3.5\text{ L/min}$
2. Remove any remaining moisture from the product using a dry cloth or filtered clean and dry compressed air ($\leq 0.35\text{ MPa}$) ensuring it is dry internally and externally.
3. Under appropriate lighting ($\geq 500\text{ lx}$), inspect the products for contamination. If any visible contamination remains, repeat the process until it is visually clean following step 1.

Disinfection

Handpiece stand:

Wipe the surface of the products with disinfectant wipes (MinutenWipes from ALPRO) for 2 minutes each.

Tube holders:

Spray the WL-cid (ALPRO) on to the products.

Wait for the WL-cid (ALPRO) to take effect following the instructions provided by the manufacturer.

Proceed to the section "5.4 Sterilization".

5.4 Sterilization

After each patient treatment, sterilize the products as follows.

Sterilizable equipment:

Handpiece stand and tube holders

Insert the products into a sterilization pouch that conforms to EN ISO 11607-1 and seal the pouch.

Sterilize using the following conditions.

	Gravity displacement	Pre-vacuum (dynamic air removal)
Temperature	132 °C	134 °C
Holding time	15 min	3-18 min

⚠ CAUTION

- Do not use sterilization pouches containing water-soluble adhesive components such as PVA (polyvinyl alcohol). Adhesive components washed out during sterilization may enter the product, causing problems such as poor rotation, sticking and failure to operate correctly.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Be sure to use sterilizers that can perform sterilization up to 135 °C. In some sterilizers, the chamber temperature may exceed 135 °C. Do not use these sterilizers as failure of the product could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Do not heat or cool the product too quickly. Rapid change in temperature will cause damage to the product.
- Do not touch the product immediately after autoclaving as it will be very hot.

5 Post-use Maintenance

CAUTION

- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.
 - Immediately after sterilization is complete (within 1 hour), remove the products from the sterilizer. Failure to do so may cause corrosion.
-

NOTICE

- NSK recommends use of a Class B sterilizer as stated in EN 13060.
 - For details, confirm the manufacturer's sterilizer operation manuals.
-

5.5 Storage

Store the product in a place where it is kept clean until it is used next.

CAUTION

- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in "10 Specification". The air should be free from dust, salt and sulphur.
-

6 Maintenance

6.1 Periodical maintenance checks

Points to check	Details
Control unit operation	There should be no abnormal mechanical sound. LCD is displayed appropriately. The control panel does not operate abnormally.
Wireless foot control operation	• Check the battery level displayed on the control unit display. — If battery level is low, replace the batteries. => “3.2 Installing batteries (Only for wireless foot control)”. • Check for the tip vibrates, when depressing the main pedal. • Depress each button and make sure that function corresponding to each button can be confirmed. If the wireless foot control does not function properly, perform calibration. => “6.2 Foot control calibration”. • Remove the wireless foot control battery cover and make sure its O-ring is not damaged or worn. If damage or wear is seen, contact an authorized NSK dealer for replacement.
Wireless foot control hanger	Make sure the wireless foot control hanger does not easily come off the wireless foot control.
Coolant solution flow	Coolant solution flows appropriately with no leakage.
Handpiece	• Operate for 1 minute and make sure there is no abnormal heating, sound, or vibration. • Make sure the handpiece LED illuminates normally.
LCD	Make sure all the displays are lit once immediately after the power is turned on.

Additional safety information


Periodical safety testing shall be performed every two years following IEC 62353 (EN 62353), in accordance with test procedure OM-DE1029EN. As periodical safety testing is limited exclusively to professional dental installers who have been trained and qualified under the control of NSK, contact your authorized NSK dealer.

6.2 Foot control calibration

Main pedal of the wireless foot control and button response may become poor due to wear of parts over time. In that case, perform calibration as described below.


1. Press and hold for 2 seconds or more to switch to setting mode.
2. Press to select "FC SET".
3. Press .
4. Press to select "FC CAL".
5. Press to select "ON".
6. Press , and confirm that "----" is displayed.
7. Fully depress the main pedal and all of the buttons of the wireless foot control.

8. When "0000" is displayed, calibration is successful. Press .
9. Make sure that beep sounds and "OK" is displayed.
10. Press twice to end setting mode.

CAUTION

- Press the enter key only when "0000" is displayed.
If calibration is not performed correctly, "NG" is displayed. In this case, start again from step 1.

NOTICE

- Wired foot control does not require calibration.

7 Setting

7.1 Setting key

Control unit and wireless foot control operational settings can be changed.

For setting items and factory settings. => "<Setting value>".

1. Press and hold for 2 seconds or more to switch to setting mode, and "BLE SET" is displayed.
2. Press to select main menu.
BLE SET : Settings for Bluetooth connection with control unit
UNT SET : Settings for control unit
FC SET : Settings for wireless foot control
3. Press .
4. Press to select sub menu.
5. Press to change the setting value.
6. Press to finalize the setting value. To change other sub menu setting items, repeat procedure from step 4 to 6.
7. Press to return to main menu selection screen.
8. To change other main menu setting items, repeat procedure from step 2 to 6.
9. Press , twice to end setting mode.

<Setting value>

Main menu	Sub menu	Setting details	Setting range	Factory setting
BLE SET	PAIRING	Pair the wireless foot control with the control unit.	FC	-
		Surgic Pro2 with the control unit	LK	-
	BLE OFF	Unpair the wireless foot control from the control unit	FC	-
		Surgic Pro2 with the control unit	LK	-
UNT SET	BZ VOL	Volume for beep that sounds at key operations and when activating.	H: high L: low	H
	LCD BL	Control panel brightness	1 - 10 Larger value: Brighter	10

7 Setting

Main menu	Sub menu	Setting details	Setting range	Factory setting
UNT SET	TS LV	Touch sensor sensitivity level (When sensitivity is poor such as when using barrier film or surgical gloves)	1 - 3 Larger value: More sensitive	1
	F-RESET	Initializing settings OFF: Return to the previous screen without initializing settings. ON: Restore to factory settings.	OFF ON	OFF
FC SET	FC CAL	Calibration of main pedal on the wireless foot control. Button response may become poor due to wear to parts over time. In that case perform calibration. OFF: Return to the previous screen without performing calibration. ON: Perform calibration. => "6.2 Foot control calibration".	OFF ON	-
	FC MODE	AUT: Sets the wireless foot control to automatic connection when the control unit power is turned on. FC: Enable sleep mode. (It may be necessary to keep pressing the main pedal of the wireless foot control to connect the control unit and the foot control automatically.) => "4.3 Sleep mode (Only for wireless foot control)".	AUT FC	AUT
	SLEEP	Set the time until the wireless foot control shifts to sleep mode. *This can be set only when FC MODE is changed from AUT to FC. => "4.3 Sleep mode (Only for wireless foot control)".	OFF 5 m 10 m 15 m 20 m 30 m 45 m 60 m	5 m

⚠ CAUTION

- When setting values are changed, press the  to finalize.
If you press the  POWER  or  without finalizing with the  or if you turn the power off, setting details will not be applied.

NOTICE

- Initializing settings returns “UNT SET” in setting mode and program settings to factory default. “BLE SET” and “FC SET” in setting mode are not initialized.
- Finalized setting details will be retained even if the power is turned off.

7.2 Pairing with wireless foot control


1. Press and hold for 2 seconds, "BLE SET" is displayed.
2. Press .
3. The following screen will be displayed.
"PAIRING" is illuminated.
"FC" and flash.

4. Press .

5. and flash. Press and hold the left, right and main pedal of the wireless foot control shown below simultaneously for approximately 10 seconds.

6. Pairing is successful when a beep sounds and "OK" is displayed.
Make sure is displayed, and the device is connected.

7. Press twice to end setting mode.

⚠ CAUTION

- Pairing may not be possible if something is obstructing the signal between the Wireless foot control and control unit.

NOTICE

- If pairing cannot be accomplished, replace the batteries with new ones.
- Pairing can be done within 10 minutes after inserting the battery.
- Pairing is not possible after about 10 minutes. Remove the battery and then reinsert it.
- When the wired foot control is connected, the wireless foot control will not operate.

NOTICE

- If pairing is not completed within 5 minutes, “NG” is displayed. If this happens, turn the main power switch OFF, and ON again to repeat pairing operation again from step 1.

7.3 Programming system operation

Nine programs can be memorized and recalled according to treatment procedures.

Each program number is assigned to each mode as shown in the table below. Modes other than those set for each program number cannot be set.

Program No.	Mode	Burst function*1	Power range	Coolant flow (5 levels)	Lighting brightness	
1	SURG	Available	10 – 150 %	20 - 76 mL/min*2	OFF LOW HIGH	
2		OFF	In multiples of 10 %			
3		(B1: 10Hz, B2: 30Hz, B3: 60Hz)				
4						
5						
6	ENDO		Not available	5 - 100 % Below 50 %, multiples of 5%. 50 % and larger, multiples of 10 %.		10 - 46 mL/min*2
7						
8	PERIO					
9						

*1: Burst function: Vibration levels can be changed at certain intervals. Levels can be selected depending on patient's bone hardness (density).

*2: The amount of water flow may vary slightly depending on the condition of the irrigation tube.

NOTICE

- When the power is turned on again, the program No. that was displayed before the power was turned off is displayed.

7.4 Link function

This product can be connected using Bluetooth to NSK Surgic Pro2.

Operate both VarioSurg 4 and Surgic Pro2 for linking.

■ Procedure

VarioSurg 4

1. Press and hold for 2 seconds or more to switch to setting mode, and "BLE SET" is displayed.
2. Press .
3. Press to select "LK".
4. "LK" and "LINK" flash.
5. Press .
6. "PAIRING" and "LK" are illuminated. and "LINK" flash.

Surgic Pro2

7. Press and hold for 2 seconds or more to switch to setting mode, and "BLE SET" is displayed.
8. Press .
9. Press to select "LK".
10. "LK" and "VARIOSURG" flash.
11. Press .
12. "PAIRING" and "LK" are illuminated. and "VARIOSURG" flash.

VarioSurg 4 / Surgic Pro2

13. Pairing is successful when this product beeps and device displays "OK".

Make sure * is displayed and the device is connected.

VarioSurg 4

Surgic Pro2

14. Press  twice to end setting mode for Surgic Pro2 and VarioSurg 4.

NOTICE

- When link mode is enabled, operation of the foot control on the Surgic Pro2 is prioritized.
 - If pairing is not completed within 5 minutes, "NG" is displayed. If this happens, turn the main power switch OFF, and ON again to repeat pairing operation again from step 1.
-

8 Troubleshooting

8.1 Error code

When an error occurs to check the status of the control unit and to easily identify the cause of the error, an error code is displayed on the LCD along with a notification sound. The error code will blink until the abnormality is resolved.

■ Error check

Depress the main pedal of the foot control or press .

If there is no error, error status is released.

If not connected the foot control, press  and error is released.

Error code	Type of error	Cause of error	Action
E0	System error	Self-check error.	Request repair.
E1	Oscillator system error	Tip is overloaded.	Be careful not to overload the tip. In the event of a situation where the tip is pinched during cutting, step on the main pedal on the foot control and slowly pull out the tip while oscillating it (do not gouge it).
		Tip not mounted, tip not tightened enough	Install the tip and tighten it with the tip replacement wrench until it "clicks" and idles.
		Handpiece cord is not securely connected. Disconnection in the handpiece cord.	Check the handpiece cord connection. If it is not resolved, contact an authorized NSK dealer.
E2	Power supply error	Insufficient contact of the power cord. Failure of the power supply circuit.	Securely re-connect the power cord. If the error cannot be eliminated, request repair.
E4	Control unit interior overheating error	Overheating by extended use under heavy load. Use immediately after the control unit is placed at high temperature (such as inside a car under the blazing sun or in a locker exposed to direct sunlight).	Allow it to cool down before use. To allow heat to be sufficiently radiated, the periphery of the control unit should be well-ventilated wherever possible. If the error cannot be eliminated, request repair.
E7	Irrigation pump error	Irrigation tube caught on pump roller. Irrigation pump failure.	Check the irrigation tube. If the error cannot be eliminated, request repair.
E8	Connection error in the control unit	Failure of Bluetooth communication circuit (LED handpiece system) inside the control unit.	Request repair.

Error code	Type of error	Cause of error	Action
E9	Wireless foot control error	Failure of sensor in wireless foot control. Detection of disconnection of Bluetooth communication during ultrasonic vibration.	Request repair.
E10	Handpiece connection error	Handpiece or tip not connected	Check that the handpiece connector is securely connected or that the tip is securely connected using a tip wrench. If not resolved, contact an authorized NSK dealer.
E11	Operation panel error	Abnormal operation panel initialization	Request repair.
E12	Connection error in the control unit	Failure of the Bluetooth communication circuit (link system with Surgic Pro2) inside the control unit	Request repair.
E13	Wired foot control error	Failure of wired foot control	Make sure the connector of the wired foot control is securely connected. If not resolved, contact an authorized NSK dealer.

8.2 Problems and solutions

Every 3 months perform periodical maintenance checks, referring to the check sheet below. If any abnormalities are found, contact an authorized NSK dealer.

Problem	Cause of problem	Action
Screen is blank.	Plug improperly inserted to power source	Check to make sure the plug is in.
	Malfunction of power switch	Request repair.
The power level and " - - - " are alternately displayed on LCD.	Power turned on while pressing the foot control	Remove foot from foot control.
	Foot control malfunction	If the problem is solved by removing the batteries, it is a foot control malfunction. If not, it is a control unit malfunction. Please contact an authorized NSK dealer.
	Control unit malfunction	

8 Troubleshooting

Problem	Cause of problem	Action
After the power is ON, "CONNECT" or "PUSH FC" is displayed and the screen is not switched.	Wireless foot control is not connected.	Keep pressing the main pedal of the wireless foot control while "PUSH FC" is displayed on the LCD. If the screen does not change after 5 minutes, press . The normal screen is displayed. Replace the wireless foot control batteries to check if it can be connected. => "3.2 Installing batteries (Only for wireless foot control)". If it cannot be connected after replacing the batteries, pairing the control unit and wireless foot control. => "7.2 Pairing with wireless foot control". If this does not resolve the issue, please contact an authorized NSK dealer.
Will not vibrate even if the Wireless foot control is pressed.	Wireless foot control is not connected.	Perform the following after confirming that nothing is blocking the signal between the wireless foot control and control unit. • Turn the main power switch OFF, and ON again to connect to the control unit. • Replace the wireless foot control batteries with new ones. => "3.2 Installing batteries (Only for wireless foot control)". • Perform pairing the control unit and wireless foot control. => "7.2 Pairing with wireless foot control". If this does not resolve the issue, please contact an authorized NSK dealer.
	Wireless foot control pedal and button response is poor due to not being used for a long period of time.	Calibrate the foot control. => "6.2 Foot control calibration".

NOTICE

- If there is no response when any key is pressed, the control unit is automatically adjusting the sensitivity. Wait for a while and the sensitivity will automatically return.

9 Spare Parts and Option List

Model	REF	Remarks
AC Power Cord	U260414	
Coolant Solution Hanger Post	U370152	Bottle prop
FC-86	ZF16040001	Wireless foot control
FC-78	Z1102003	Wired foot control
Battery Cover	Z1401068	For FC-86
Tube Holder	Y900767	Pack of 7
Handpiece Stand	Z1029201	
Irrigation Tube	Y900113	Pack of 5
VS4-LED-HPSC	EA14170001	LED handpiece with cord
VA-SG-TIP HOLDER	20001327	
Tips Holder	Z221086	Only for Europe.
CR-30	10000977	Tip Wrench
Connector Cap	UA14260145	
Link Stand3	ZA16230001	For installation of Surgic Pro2 and this product.
Carrying Case	Y1500783	For components and accessories of this product.
iCart L	Z299	For installation of Surgic Pro2 and this product and accessories.
iCart Duo	S9090	For installation of Surgic Pro2 and this product and accessories. Only for Europe.
E TIP WRENCH	Z217399	For V10-S
CR-40	Z305350	VS Tip Wrench. For SG1A, SG8A, Socket Lift Tip.
VA-SG-CASE	Z313102	Sterilization Cassette. Only for Europe.
VA-SG-CASE	20001326	Sterilization Cassette.

10 Specification

10.1 System specifications

■ Control unit

Model	VarioSurg 4
Power Supply Voltage	AC100 - 240 V
Frequency	50/60 Hz
Drive frequency	28 - 32 kHz
Maximum output	25 W
Power consumption	51 - 110 VA
Max. Pump Output	76 mL / min
Dimensions	W 245 x D 235 x H 90 mm
Weight	2.2 kg

■ LED handpiece with handpiece cord

Model	VS4-LED-HPSC
Type of vibration	Piezo type
Optic	White LED
Consumption current (LED)	0.17 A (3.5 V)
Cord length	2 m
Dimensions	Ø20 x 129 mm (with cord)
Weight	170 g (with cord)

■ Wireless foot Control

Model	FC-86
Electric power source	Three AAA/Micro/LR03 1.5 V alkaline batteries 3 pcs.
Dimensions	W 260 x D 185 x H 65 mm
Weight	1.1 kg (with hanger)

■ Wired foot control (option)

Model	FC-78
Cord length	2 m
Dimensions	W 268 x D 230 x H 103 mm (without cord)
Weight	1.1 kg

■ Use, transportation and storage environment

	Temperature	Humidity	Atmospheric pressure
Use Environment			
Transportation and Storage Environment			

No moisture condensation in the control unit.

Use outside of these limits may cause malfunction.

10.2 Bluetooth specifications

Frequency band	2.4 GHz ISM band (2.402-2.480 GHz)
Transmitting power	Class3 1 mW
Modulation	GFSK
Channels	40 channel 2 MHz spacing
Compatibility	EN 300 328, EN 300 489-1, EN 301 489-17, EN 62479:2010

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by NAKANISHI INC. is under license. Other trademarks and trade names are those of their respective owners.

10.3 Classification of equipment

- Type of protection against electric shock:
 - Class I equipment
- Degree of protection against electric shock:
 - Type BF applied part ☒ (Applied part: tip, handpiece)
- Method of maintenance (cleaning) and sterilization recommended by the manufacturer:
 - => “5 Post-use Maintenance”
- Degree of protection against ingress of water as detailed in IEC 60529 (EN 60529):
 - Foot control: IPX8 (Protected against the effects of continuous immersion in water)
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:
 - Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Mode of operation:
 - Non-continuous operation

10.4 Operation principle

A sinusoidal electrical signal, at ultrasonic frequency ($f > 20$ kHz), is delivered by the generator. This signal is applied to the ‘piezoelectric ceramic’ located inside the transducer. Piezoelectric ceramic converts this signal into mechanical vibrations. These vibrations are at the same ultrasonic frequency as the electrical signal. The mechanical vibrations are propagated towards the distal end of the transducer. The “TIP” insert, which is attached at the distal end of the transducer, vibrates at ultrasonic frequencies and makes it possible to achieve the aimed purpose.

10.5 Factory settings

No.	Name	Factory setting
1	Power setting	10%
2	Program selection	1
3	Burst	OFF
4	Mode switching	SURG
5	Coolant solution flow level selection	5
6	Light level selection	HIGH

11 EMC Information (Electromagnetic Compatibility Information)

Guidance and manufacturer's declaration - Electromagnetic Emissions.

The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that is used in such an environment. Guidance and manufacturer's declaration electromagnetic emissions.

Emissions test	Compliance	Electromagnetic environment - guidance
Radiated emissions CISPR 11/EN 55011	Group 1 Class B	The product uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emissions CISPR 11/EN 55011	Group 1 Class B	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2/EN 61000-3-2	Class A (except 120-240 V)	
Voltage fluctuations/flicker emissions IEC 61000-3-3/EN 61000-3-3	Complies (except 120-240 V)	

Guidance and manufacturer's declaration - Electromagnetic Immunity

The product is intended for use in the electromagnetic environment specified below.

The customer or the user of the product should assure that it is used in such an environment.

Immunity test	IEC/EN 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2/EN 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4/EN 61000-4-4	Input a.c. power PORT ± 2 kV Signal input/output parts PORT ± 1 kV	Input a.c. power PORT ± 2 kV Signal input/output parts PORT ± 1 kV
Surge IEC 61000-4-5/EN 61000-4-5	Line to line ± 0.5 kV, ± 1 kV Line to ground ± 0.5 kV, ± 1 kV, ± 2 kV	Line to line ± 0.5 kV, ± 1 kV Line to ground ± 0.5 kV, ± 1 kV, ± 2 kV
Voltage dips IEC 61000-4-11/EN 61000-4-11	0 % U_T ; 0.5cycle (0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0 % U_T ; 1cycle and 70 % U_T ; 25 (50 Hz)/30 (60 Hz)cycles Single phase ; 0°	0 % U_T ; 0.5cycle (0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0 % U_T ; 1cycle and 70 % U_T ; 25 (50 Hz)/30 (60 Hz)cycles Single phase ; 0°

11 EMC Information (Electromagnetic Compatibility Information)

Voltage interruptions IEC 61000-4-11/EN 61000-4-11	0 % U_T ; 250 (50 Hz) /300 (60 Hz) cycle	0 % U_T ; 250 (50 Hz) /300 (60 Hz) cycle
Power frequency (50/60Hz) magnetic field IEC 61000-4-8/EN 61000-4-8	30 A/m	30 A/m
Proximity magnetic field IEC 61000-4-39	134.2 kHz (Pulse modulation) 2.1 kHz, 50 % Duty, 13.56 MHz (Pulse modulation) 50 kHz 50 % Duty	134.2 kHz (Pulse modulation) 2.1 kHz, 50 % Duty, 13.56 MHz (Pulse modulation) 50 kHz 50 % Duty

NOTE U_t is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - Electromagnetic Immunity			
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.			
Immunity test	IEC/EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6/EN 61000-4-6	3 V rms 0.15 MHz - 80 MHz 6 V rms ISM bands between 0.15 MHz and 80 MHz	3 V rms 6 V rms	WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the product, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3/EN 61000-4-3	3 V/m 80 MHz - 2.7 GHz	3 V/m	
	9 V/m in ISM bands 710 MHz, 745 MHz, 780 MHz, 5.24 GHz, 5.5 GHz, 5.785 GHz	9 V/m	
	27 V/m in ISM band 385 MHz	27 V/m	
	28 V/m in ISM bands 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1.72 GHz, 1.845 GHz, 1.97 GHz, 2.45 GHz	28 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

11 EMC Information (Electromagnetic Compatibility Information)

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product series is used exceeds the applicable RF compliance level stated above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product series.
- b Over the 150 kHz to 80 MHz frequency range, the field strength should be less than 3V/m.

Cables and accessories	Maximum length	Complies with
Handpiece Cord	2.0 m (Unshielded)	RF emissions, CISPR 11: Class B/ Group 1
Wired Foot Control Cord	2.0 m (Unshielded)	Electrostatic discharge (ESD): EN/IEC61000-4-2
		Electrical fast transient/burst: EN/IEC61000-4-4
		Surge: EN/IEC61000-4-5
		Voltage dips, short interruptions and voltage variations on power supply input lines: EN/IEC61000-4-11
AC Power Cord	2.0 m (Unshielded)	Power frequency (50/60 Hz) magnetic field: EN/IEC61000-4-8
		Conducted RF: EN/IEC 61000-4-6
		Radiated RF: EN/IEC61000-4-3

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NSK Europe GmbH  
Elly-Beinhorn-Str. 8, 65760 Eschborn, Germany

Visit our
website

We have corrected the quantity of the irrigation tubes included.

Original

Corrected

EN	P12		
	No.	Part name	Quantity
	11	Irrigation tube	5
DE	P64		
	Nr.	Bauteilbezeichnung	Anzahl
	11	Kühlmittelschlauch	5
FR	P120		
	N°	Nom de la pièce	Quantité
	11	Tube d'irrigation	5
ES	P176		
	N.º	Nombre de la pieza	Cantidad
	11	Tubo de irrigación	5
IT	P232		
	N.	Nome parte	Quantità
	11	Tubo di irrigazione	5
PT	P288		
	N.º	Nome da peça	Quantidade
	11	Tubo de irrigação	5

EN	P12		
	No.	Part name	Quantity
	11	Irrigation tube	3
DE	P64		
	Nr.	Bauteilbezeichnung	Anzahl
	11	Kühlmittelschlauch	3
FR	P120		
	N°	Nom de la pièce	Quantité
	11	Tube d'irrigation	3
ES	P176		
	N.º	Nombre de la pieza	Cantidad
	11	Tubo de irrigación	3
IT	P232		
	N.	Nome parte	Quantità
	11	Tubo di irrigazione	3
PT	P288		
	N.º	Nome da peça	Quantidade
	11	Tubo de irrigação	3

Additional information

EN NOTICE

Substances which are carcinogenic, mutagenic or toxic to reproduction ("CMR"), of category 1A or 1B, in accordance with Part 3 of Annex VI to Regulation (EC) No 1272/2008 of the European Parliament and of the Council, or substances having endocrine-disrupting properties for which there is scientific evidence of probable serious effects to human health and which are identified either in accordance with the procedure set out in Article 59 of Regulation (EC) No 1907/2006 of the European Parliament and of the Council or, once a delegated act has been adopted by the Commission pursuant to the first subparagraph of Article 5 of Regulation (EU) No 528/2012 of the European Parliament and the Council, in accordance with the criteria that are relevant to human health amongst the criteria established therein: lead (Pb), cobalt (Co).

This product can be safely used in the treatment of children, or other patient groups considered particularly vulnerable to such substances.

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