Joline®

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Mixer

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Contents

- 1. Device description
 - 1.1 General device description, clinical benefit and mode of operation
 - 1.2 Performance features/ technical specification
 - 1.3 Compatibility with other devices
- 2. Indication
- 3. Contraindication
- 4. Complications
- 5. Patient group
- 6. User profile
- 7. Environment for use
- 8. Warnings and precautions
- 9. Packaging, shelf life, storage and disposal
- 10. Summary of safety and clinical performance
- 11. Reporting obligations for severe incidents
- 12. Use

1. Device description

1.1 General device description, clinical benefit and mode of operation

The Mixer is a device for mixing and transferring bone cement to the Vertebra Filler Nozzles of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System that is used to inject bone cement during kyphoplasty. The individual components of the Mixer are the mixing system (comprising the mixing cylinder and lid), one-way stopcock, spindle unit, funnel and stand. The Joline Mixer is intended for single use only. The period of use is usually 10 minutes.

The Mixer is intended to enable the user to uniformly mix and transfer the bone cement to the Vertebra Filler Nozzles of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System. The closed construction is intended to minimise any unpleasant odours that the user is exposed to due to the liquid monomer components of PMMA bone cements.

Using the Mixer the bone cement is mixed and transferred to the Vertebra Filler Nozzles of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System through a Luer lock connector. Mixing and transferring the bone cement with the Mixer is carried out as follows:

The bone cement is mixed manually using a mixing blade that is fixed to the piston rod and can be moved up and down the length of and rotated in the mixing cylinder. For the mixing of the two individual components of the bone cement, they are homogeneously combined, which starts the polymerisation of the bone cement. The mixing process is stopped as soon as the desired viscosity of the bone cement is reached (see IFU of the bone cement).

For the transfer, the bone cement is transferred to several Vertebra Filler Nozzles using the Luer lock connector. To do so, the spindle unit is connected via the piston rod with the mixing system and the one-way stopcock is opened. Using manual rotating movements of the threaded spindle, the piston seal is moved forwards axially, which transfers doses of the bone cement dough into the Vertebra Filler Nozzles.

1.2 Performance features/ technical specification

- Closed system
- Filling volume: 1 packaging unit of bone cement (approx. 25 g polymer powder + approx. 10 ml monomer liquid)
- One-way stopcock with sealing option
- Luer lock connector to connect the Vertebra Filler Nozzle of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System
- Spindle unit for simple transfer of the bone cement

1.3 Compatibility with other devices

The Mixer is designed to be used in combination with the following devices:

- Radiopaque bone cement based on polymethylmethacrylate (PMMA) or calcium phosphate
- Vertebra Filler Nozzles of the following Vertebra Filler Devices of the Allevo Kyphoplasty System:
 - Vertebra Filler Device (frontal opening)
 - Vertebra Filler Device (side opening)
 - Vertebra Filler Device flexPlunger (frontal opening)

2. Indication

The Mixer is intended for mixing and transferring bone cement during kyphoplasty.

3. Contraindication

There are no contraindications known for use of the Mixer.

4. Complications

The following complications are associated with the Mixer:

- a) Infection, including deep or superficial wound infection
- b) Inhomogeneity of the bone cement with associated consequences for the patient

5. Patient group

Provided that the manufacturer of the bone cement does not specify any limitations on the patient group in the intended use, there are also no limitations for the Mixer because it does not make direct contact with the patient.

6. User profile

The Mixer is designed to be operated by OR (operating room) personnel. It should only be used by or under the supervision of an experienced user or in accordance with his or her instructions.

7. Environment for use

The Mixer is used in a sterile environment such as an operating room during kyphoplasty. The Mixer is used under sterile conditions.

Due to the functional principle and properties of the bone cements, the Mixer is only suitable for use within the application time specified for the bone cement because otherwise the function of the Mixer may be impaired.

For application of the bone cement used, the specification of the bone cement manufacturer must also be complied with.

8. Warnings and precautions

- 8.1 See the 'Explanation of symbols' to understand the symbols used on the packaging. The contents specified on the label are translated under 'Translation of the label'.
- 8.2 The medical techniques and procedures described here make NO claim to completeness. They do not illustrate all the medically possible procedures and, most of all, do not replace the doctor's own experience regarding the use and the treatment of certain patients.
- 8.3 The device may only be used by an experienced user, under his or her supervision or in accordance with his or her instructions.
- 8.4 The device has been sterilised with ethylene oxide.
- 8.5 The device is only sterile and pyrogen-free if the sterile packaging is unopened and undamaged. Do not use the device if the packaging has already been opened or is damaged.
- 8.6 The device must be used before the expiry date stated on the packaging ('Use-by' date). The device must no longer be used after the expiry date has passed and must be disposed of.
- 8.7 Always ensure a sterile surgical field and sterile procedures when handling or using the device.
- 8.8 Do not use the device if it appears to be damaged or defective.
- 8.9 **FOR SINGLE USE ONLY!** The device is designed for single use. It must not be reused, resterilised or reprocessed. Reuse, resterilisation or reprocessing can negatively impact the performance of the device and may result in injury to the patient if this device is used. Reuse, resterilisation or reprocessing may also increase the risk of contamination of the device and/or infection of the patient, including but not limited to transmission of infectious diseases from one patient to another.
- 8.10 The device is designed for use with the Vertebra Filler Nozzle of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System and with the PMMA bone cements distributed by Joline.
- 8.11 In all cases, follow the instructions for use of the bone cement manufacturer.
- 8.12 Do not exceed the maximum recommended capacity of the Mixer of one packaging unit of the bone cement (approx. 25 g polymer powder and approx. 10 ml monomer liquid) to ensure that the bone cement is mixed homogeneously.
- 8.13 Small glass fragments may result from breaking the neck of the ampoule. To prevent them from getting into the bone cement and consequently into the vertebral body, do not open the ampoule directly over the mixing cylinder but open it instead to the side.
- 8.14 Only attach the funnel after the monomer liquid has been added because otherwise the polymer powder will stick to the sides when poured through the wet funnel.
- 8.15 Close the one-way stopcock so that bone cement does not escape during mixing.
- 8.16 When mixing, hold the mixing cylinder upright with the lid facing upwards to prevent air inclusions and thus inhomogeneity of the bone cement.
- 8.17 Mix the components of the bone cement thoroughly with the mixing blades to prevent inhomogeneity of the bone cement.
- 8.18 Holding the mixing cylinder in the palm of your hand could heat up the bone cement, which may significantly shorten the processing times for the bone cement.
- 8.19 The two lateral fixtures on the spindle unit must clearly and audibly latch into place so that the bone cement can subsequently be transferred safely.
- 8.20 Fully vent the Mixer to prevent air inclusions and thus inhomogeneity of the bone cement. When doing so, hold the Mixer away from your face in case dust from any residual polymer powder escapes.

8.21 Transfer the bone cement after mixing as quickly as possible for as long as the increasing viscosity of the bone cement allows transfer.

9. Packaging, shelf life, storage and disposal

The device is supplied sterile and pyrogen-free. It is individually packaged. Sterility is guaranteed as long as the packaging is not opened or damaged.

Note the expiry date on the corresponding packaging. The device must no longer be used after the expiry date has passed and must be disposed of.

Store and transport the device in a cool, dark and dry place until use. Keep away from sunlight and heat.

Do not use the device if the packaging has already been opened or is damaged. If damage of any sort is identified, the device must no longer be used.

The packaging material used, comprising paper, cardboard and plastic, can be recycled and safely disposed of through existing disposal channels. It is considered non-problematic waste.

There are no hazards associated with the device before its use. If the device is inadvertently opened or is damaged, it can be disposed of according to hospital hygiene guidelines.

Residual bone cement should be allowed to cure fully in the Mixer before it can be disposed of with hospital waste in accordance with general and hospital-specific precautions.

10. Summary of safety and clinical performance

The summary of safety and clinical performance pursuant to Article 32 of Regulation (EU) 2017/745 is only to be drawn up for implantable devices and class III devices that are not custom-made devices or investigational devices. Due to the fact that the Mixer does not fall under this definition, a summary of safety and clinical performance is not applicable and thus also was not drawn up.

11. Reporting obligations for severe incidents

The user must immediately report all severe incidents occurring in causal relationship with the device (death of a patient, user or other person; temporary or permanent severe deterioration in the health of a patient, user or other person; serious risk to public health) to the manufacturer as well as the competent authority of the country in which the user is established. The unique device identifier (UDI) required for this purpose as well as other device-related information can be found on the packaging label.

The affected device that is causally linked to a severe incident must be sent immediately to the manufacturer for the necessary inspections and risk assessment by the manufacturer in accordance with Article 89 of Regulation (EU) 2017/745.

12. Use

The Mixer is used in accordance with Figures 1–11 of the supplementary sheet included in the packaging.

- I Filling
- 12.1

Pull the piston rod out to the stop. Place the stand on an even surface. Place the Mixer, without the lid, in the stand with the opening facing upwards.
- 12.2

Pour the liquid monomer component of the bone cement into the mixing cylinder. Attach the funnel to the mixing cylinder. Fill the powder polymer component of the bone cement into the mixing cylinder through the funnel.

WARNING: Small glass fragments may result from breaking the neck of the ampoule. To prevent them from getting into the bone cement and consequently into the vertebral body, do not open the ampoule directly over the mixing cylinder but open it instead to the side.

WARNING: Only attach the funnel after the monomer liquid has been added because otherwise the polymer powder will stick to the sides when poured through the wet funnel.

WARNING: Do not exceed the maximum recommended capacity of the Mixer of one packaging unit of the bone cement (approx. 25 g polymer powder and approx. 10 ml monomer liquid) to ensure that the bone cement is mixed homogeneously.

II Preparing

- 12.3

Remove the funnel from the mixing cylinder. Screw the lid with the closed one-way stopcock onto the mixing cylinder. Remove the mixing system from the stand. The funnel and stand are no longer required and can be disposed of.

CAUTION: Close the one-way stopcock so that bone cement does not escape during mixing.

12.4

Remove the piston rod lock by pulling it sideways away from the piston rod. The lock is no longer required and can be disposed of.

III Mixing

- 12.5

Hold the mixing system in one hand by the grip handles. Mix the components of the bone cement by using the other hand to push the piston rod back and forth along its entire length while simultaneously performing oscillating rotations.

WARNING: When mixing, hold the mixing cylinder upright with the lid facing upwards to prevent air inclusions and thus inhomogeneity of the bone cement.

WARNING: Mix the components of the bone cement thoroughly with the mixing blades to prevent inhomogeneity of the bone cement.

CAUTION: Holding the mixing cylinder in the palm of your hand could heat up the bone cement, which may significantly shorten the processing times for the bone cement.

CAUTION: In all cases, follow the instructions for use of the bone cement manufacturer.

12.6

Once the recommended mixing time has elapsed (as specified in the instructions for use of the bone cement manufacturer), pull the piston rod out to the stop. Hold the piston rod with one hand and use the other hand to forcefully break off the handle of the piston rod and dispose of it.

CAUTION: In all cases, follow the instructions for use of the bone cement manufacturer.
- IV Venting
- 12.7

Screw the threaded spindle of the spindle unit anticlockwise out to the stop. Push the spindle unit over the piston rod, which has been pulled out to the stop, and connect to the mixing system by pressing hard. Two clicks indicate that the two fixtures on the spindle unit are securely latched into place on the mixing cylinder.

CAUTION: The two lateral fixtures on the spindle unit must clearly and audibly latch into place so that the bone cement can subsequently be transferred safely.

12.8

To vent, open the one-way stopcock by a one-quarter turn. Hold the mixing system vertically with the lid facing upwards as described in step 5. Push the piston forwards in a clockwise direction by turning the threaded spindle until bone cement exits the one-way stopcock.

WARNING: Fully vent the Mixer to prevent air inclusions and thus inhomogeneity of the bone cement. When doing so, hold the Mixer away from your face in case dust from any residual polymer powder escapes.

12.9

To avoid the bone cement escaping, turn the threaded spindle back about one-half anticlockwise rotation.

V Transferring

12.10

Connect the Vertebra Filler Nozzle to the Luer lock connector of the one-way stopcock on the lid. Transfer the bone cement to the filler nozzle by turning the threaded spindle clockwise. One complete rotation of the threaded spindle is sufficient to transfer about 1.5 ml of bone cement equivalent to the capacity of one filler nozzle.

CAUTION: Transfer the bone cement after mixing as quickly as possible for as long as the increasing viscosity of the bone cement allows transfer.

CAUTION: In all cases, follow the instructions for use of the bone cement manufacturer.

12.11

Once the bone cement emerges from the distal end of the Vertebra Filler Nozzle, turn the threaded spindle back about one-half anticlockwise rotation to prevent the bone cement escaping. Remove the Vertebra Filler Nozzle from the one-way stopcock.

12.12

Repeat steps 10 and 11 with additional Vertebra Filler Nozzles until the required quantity of bone cement has been transferred.

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Joline GmbH & Co. KG