Joline[®]

de **Gebrauchsanweisung
Dilatator**

en **Instructions for use
Dilator**

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Contents

1. Device description
 - 1.1 Intended purpose
 - 1.2 General device description, clinical benefit and mode of operation
 - 1.3 Performance features/ technical specification
 - 1.4 Compatibility with other devices
2. Indication
3. Contraindication
4. Undesireable side-effects
5. Patient group
6. User profile
7. Environment for use
8. Warnings and precautions
9. Packaging, shelf life, storage and disposal
10. Summary of safety and clinical performance
11. Reporting obligations for severe incidents
12. Use

1. Device description

1.1 Intended purpose

The dilator is used to temporarily dilate the puncture tract in preparation for the implantation of a dialysis catheter.

1.2 General device description, clinical benefit and mode of operation

General device description

The dilator is a flexible, sterile and non-pyrogenic vascular dilator for single use.

It consists of a distally tapered plastic tube (shaft) made of polyethylene (PE) with an atraumatic tip. A colour-coded Luer (Lock) connector is located at the proximal end. The dilator has an inner lumen through which it is inserted into the tissue or vessel via a guidewire.

Clinical benefit

The dilator is required for implanting a Joline dialysis catheter which is performed using the Seldinger technique.

Mode of operation

The dilator is used to widen the percutaneous opening and the punctured vessel, thus preparing it for insertion of the dialysis catheter. The plastic tube is hollow internally and tapers distally to an atraumatic tip so that the dilator can be inserted safely into the vessel via the guidewire.

The proximal end can be connected to a syringe via the Luer (-Lock) connection to flush the lumen with sterile physiological saline solution (0.9% NaCl).

1.3 Performance features/ technical specification

The performance features and technical specifications are shown in the detailed description (see page 9).

1.4 Compatibility with other devices

When using the dilator, the following combinations with other medical devices may occur:

- Luer(-Lock) connection with Luer(-Lock) disposable syringes
- Guidewire (Seldinger Technique)

2. Indication

The dilator is used to prepare percutaneous vascular access (Seldinger technique) in patients who require vascular access for extracorporeal blood purification procedures.

3. Contraindication

In cases of known allergy or contact allergy to the dilator material, the dilator must not be used.

4. Undesireable side-effects

The following undesirable side-effects are associated with the use of the dilators:

- | | |
|--|--|
| a) Vessel injury | g) Product-related infections |
| b) Allergic reaction to product material | h) Thrombosis |
| c) Cardiac arrhythmias | i) Death |
| d) Bleeding / Hemorrhage | j) Injury of adjacent organs or structures |
| e) Hemothorax | k) Injury of heart chambers |
| f) Pericardial tamponade | |

5. Patient group

The patient group includes newborns, infants, children, adolescents and adults. In stable patients who can undergo the operation, age, weight or gender do not matter. Depending on the patient's anatomical and physiological conditions, which are usually reflected in the patient's weight, the treating physician decides which dilator size to use.

6. User profile

The dilators are intended to be used by qualified physicians.

7. Environment for use

The dilator is used either in a hospital intensive care unit or in a sterile environment such as an operating room. It should be performed in a moderate climate, typical for operating rooms, i.e. room temperature and moderate humidity.

8. Warnings and precautions

- 8.1 See the 'Explanation of symbols' to understand the symbols used on the packaging. The contents specified on the label are translated under 'Translation of the label'.
- 8.2 The medical techniques and procedures described here make NO claim to completeness. They do not illustrate all the medically possible procedures and, most of all, do not replace the doctor's own experience regarding the use and the treatment of certain patients.
- 8.3 The device may only be used by an experienced user, under his or her supervision or in accordance with his or her instructions.
- 8.4 The device has been sterilised with ethylene oxide.
- 8.5 The device is only sterile and non-pyrogenic if the sterile packaging is unopened and undamaged. Do not use the device if the packaging has already been opened or is damaged.
- 8.6 The device must be used before the expiry date stated on the packaging ('Use-by' date). The device must no longer be used after the expiry date has passed and must be disposed of.
- 8.7 Always ensure a sterile surgical field and sterile procedures when handling or using the device.
- 8.8 Do not use the device if it appears to be damaged or defective.
- 8.9 **FOR SINGLE USE ONLY!** The device is designed for single use. It must not be reused, resterilised or reprocessed. Reuse, resterilisation or reprocessing can negatively impact the performance of the device and may result in injury to the patient if this device is used. Reuse, resterilisation or reprocessing may also increase the risk of contamination of the device and/or infection of the patient, including but not limited to transmission of infectious diseases from one patient to another.
- 8.10 When selecting the dilator, consider the patient's anatomical conditions.
- 8.11 When selecting the dilator, ensure the compatibility with the selected guidewire.
- 8.12 Flush the dilator with sufficient sterile physiological saline solution (NaCl 0.9 %) before use to reduce the risk of embolisms.
- 8.13 To avoid vascular injury, do not insert the dilator deeper into the vessel than necessary or against any resistance.
- 8.14 To avoid vascular injury, do not under any circumstances insert the dilator beyond the tip of the guide wire.

9. Packaging, shelf life, storage and disposal

The device is supplied sterile and non-pyrogenic. It is individually packaged. Sterility is guaranteed as long as the packaging is not opened or damaged.

Note the expiry date on the corresponding packaging. The device must no longer be used after the expiry date has passed and must be disposed of.

Store and transport the device in a cool, dark and dry place until use. Keep away from sunlight and heat.

Do not use the device if the packaging has already been opened or is damaged. If damage of any sort is identified, the device must no longer be used.

The packaging material used, comprising paper, cardboard and plastic, can be recycled and safely disposed of through existing disposal channels. It is considered non-problematic waste.

The device poses no biological hazards before its use. If the device is unintentionally opened or is damaged, it can be disposed of in accordance with hospital hygiene guidelines. In order to prevent any injuries from cuts, dispose of the device in its original protective packaging.

The device poses potential biological hazards (e.g. microorganisms, viruses, allergens) after its use. It is therefore considered to be potentially infectious waste. If the device comes into contact with pathogens which are subject to mandatory reporting during use, it is considered to be hazardous waste and should be marked with the 'biohazard' symbol upon disposal. It must be handled and disposed of in accordance with recognised medical procedures and in line with local, state and federal laws and regulations.

Once it has been used and therefore contaminated, immediately place the device into tear-resistant, moisture-resistant and sealed containers at the waste-generation site and take it to a central collection point in sealed containers without being transferred or sorted.

10. Summary of safety and clinical performance

The summary of safety and clinical performance pursuant to Article 32 of Regulation (EU) 2017/745 is only to be drawn up for implantable devices and class III devices that are not custom-made devices or investigational devices. Due to the fact that the Dilator does not fall under this definition, a summary of safety and clinical performance is not applicable and thus also was not drawn up.

11. Reporting obligations for severe incidents

The user must immediately report all severe incidents occurring in causal relationship with the device (death of a patient, user or other person; temporary or permanent severe deterioration in the health of a patient, user or other person; serious risk to public health) to the manufacturer as well as the competent authority of the country in which the user is established. The unique device identifier (UDI) required for this purpose as well as other device-related information can be found on the packaging label.

The affected device that is causally linked to a severe incident must be sent immediately to the manufacturer for the necessary inspections and risk assessment by the manufacturer in accordance with Article 89 of Regulation (EU) 2017/745.

12. Use

NOTE: Follow the instructions for use of the dialysis catheter to be implanted and the Seldinger technique described therein.

WARNING: Always ensure a sterile surgical field and sterile procedures when handling or using the device.

CAUTION: When selecting the dilator, ensure the compatibility with the selected guidewire.

WARNING. When selecting the dilator, consider the patient's anatomical conditions.

WARNING: Flush the dilator with sufficient sterile physiological saline solution (NaCl 0.9 %) before use to reduce the risk of embolisms.

- a) Insert the dilator into the vein via the guide wire and make sure that the guide wire does not advance any further into the vein.

DANGER: To avoid vascular injury, do not insert the dilator deeper into the vessel than necessary or against any resistance.

DANGER: To avoid vascular injury, do not under any circumstances insert the dilator beyond the tip of the guide wire.

- b) Pull out the dilator. If necessary, compress the vessel manually to prevent increased blood loss.
- c) Dispose of the dilator after use in accordance with chapter 9.

Joline GmbH & Co. KG
Neue Rottenburger Str. 50
72379 Hechingen
Germany
phone: +49 (0) 7471 9881-0
fax: +49 (0) 7471 9881-111
www.joline.de
e-mail: info@joline.de

Art. no.: 11826
Revision: 00
Date of issue: 2026-02