Joline®

- de** Gebrauchsanweisung
Biopsiezange Knipsa
- en** Instructions for use
Biopsy Forceps Knipsa
- es** Instrucciones de uso
Pinzas de Biopsia Knipsa
- sv** Bruksanvisning
Biopsitång Knipsa

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Contents

1. Device description
 - 1.1 Intended purpose
 - 1.2 General device description, clinical benefit and mode of operation
 - 1.3 Performance features/ technical specification
 - 1.4 Compatibility with other devices
2. Indication
3. Contraindication
4. Undesireable side-effects
5. Patient group
6. User profile
7. Environment for use
8. Warnings and precautions
9. Packaging, shelf life, storage and disposal
10. Summary of safety and clinical performance
11. Reporting obligations for severe incidents
12. Use
 - 12.1 Device selection and preparation
 - 12.2 Biopsy sampling

1. Device description

1.1 Intended purpose

The Biopsy Forceps Knipsa is used for tissue sampling from the endomyocardium (right ventricle). It is intended for single use.

1.2 General device description, clinical benefit and mode of operation

General device description

The Biopsy Forceps Knipsa (referred to in the following as the 'device') is a radiopaque endomyocardial biopsy forceps (also called a biopptome) made of stainless steel, an FEP sheathing and medical silicone oil as a lubricant. The device is supplied sterile and non-pyrogenic. The grip can be freely rotated. A bend indicator is provided on the grip. This allows the user to more easily guide the shaft, which has already been pre-bent.

This is a single-use device which is intended for use in a single intervention. The intervention usually involves approx. 4-6 tissue samples being taken from different regions of the right ventricle.

The device consists of the main components/functional elements shown in Figure 1 (see page 8).

Clinical benefit

The tissue samples taken from the endomyocardium are used for pathological diagnostics to determine suitable therapeutic approaches. They are used, among other things, for follow-up monitoring of heart transplants and for diagnosis in suspected cases of heart muscle disease, storage disease of the heart (e.g. amyloidosis), cardiomyopathy or unclear cases of inflammation of the heart muscle (myocarditis).

Functionality

Using the grip [1], the user can open and close the cutting jaws [3] of the biopsy forceps with one hand by exerting force on the thumb ring and the Bowden cable inside the shaft [2] (Figures 2.1 and 2.2, see page 8).

1.3 Performance features/ technical specification

The performance features and technical specifications are shown in the detailed description (of the individual instruments (see page 8).

1.4 Compatibility with other devices

The Biopsy Forceps Knipsa is introduced into the vascular system through a suitable introducer sheath with an inner diameter of ≥ 6 Fr. The introducer sheath must be provided by the user.

In case of a femoral introduction an introducer sheath with an appropriate length to position the tip at least in the right atrium and before the tricuspid valve is recommended.

2. Indication

Usage of the Biopsy Forceps Knipsa is indicated in cases in which samples of the (right ventricular) heart muscle tissue are needed for further diagnostic purposes.

This includes, among other things:

- Follow-up care after heart transplantations
- Heart muscle inflammation (myocarditis)
- Cardiac storage diseases (e.g. cardiac amyloidosis)
- Cardiomyopathy

The KNIPSA500 is preferably intended for the introduction via the right internal jugular vein, while the KNIPSA1100 is intended for the introduction via one of the femoral veins.

3. Contraindication

Absolute contraindications:

- | | |
|--------------------------|---|
| a) Intracardiac thrombus | c) Severe tricuspid or pulmonary stenosis |
| b) Ventricular aneurysm | d) Tricuspid mechanical prosthesis |

Relative contraindications:

- | | |
|--|---|
| a) Active bleeding | f) Uncontrolled hypertension |
| b) Infection and fever | g) Thin ventricular wall (for the biopsy of the myocardium) |
| c) Infective endocarditis | h) Coagulopathy |
| d) Pregnancy | i) Contrast media hypersensitivity |
| e) Recent cerebrovascular accident/transient ischaemic attack (TIA) (<1 month) | j) Uncooperative patient |

4. Undesireable side-effects

The following undesirable side-effects may occur during an endomyocardial biopsy:

- | | |
|--|--|
| a) Cardiac arrhythmia (e.g. atrioventricular conduction abnormalities, complete heart block, tachycardia such as e.g. atrial fibrillation) | l) Embolism (thromboembolism/air embolism) |
| b) Vasovagal reaction (hypotonia) | m) Death |
| c) Perforation (ventricular, vascular, septum) | n) Pain |
| d) Pericardial tamponade | o) Haemorrhaging/haematoma |
| e) Injury to the tricuspid valve | p) Infection/sepsis |
| f) Pericardial effusion | q) Stroke |
| g) Coronary fistula | r) Neurological events |
| h) Arterial puncture | s) Allergic reaction |
| i) Pneumothorax | t) Shifting of a pacemaker electrode |
| j) Pseudoaneurysm | u) Cardiac infarction |
| k) Haemothorax | v) Pericardial fibrosis |

5. Patient group

The Biopsy Forceps Knipsa is intended for use on HTX patients as follow-up monitoring after heart transplants, for the diagnosis of suspected heart muscle inflammation (myocarditis), suspected cardiac storage disease (e.g. amyloidosis), or if there is existing cardiomyopathy. The patient group contains adults. Gender and weight are irrelevant for the selection of patients.

6. User profile

The Biopsy Forceps Knipsa for the myocardial biopsy is only intended to be used by experienced doctors who are trained in interventional techniques and in dealing with possible complications.

7. Environment for use

The Biopsy Forceps Knipsa is usually used for the endomyocardial biopsy in a sterile environment, such as an operating theatre. Other requirements for carrying out the procedure include monitoring the patient's vital signs and using suitable diagnostic imaging, e.g. a C-arm for diagnostic radiography and/or ultrasound. The procedure is carried out under sterile conditions.

8. Warnings and precautions

- 8.1 See the 'Explanation of symbols' to understand the symbols used on the packaging.
- 8.2 The medical techniques and procedures described here make NO claim to completeness. They do not illustrate all the medically possible procedures and, most of all, do not replace the doctor's own experience regarding the use and the treatment of certain patients.

- 8.3 The device may only be used by an experienced user, under his or her supervision or in accordance with his or her instructions.
- 8.4 The device has been sterilised with ethylene oxide.
- 8.5 The device is only sterile and non-pyrogenic if the sterile packaging is unopened and undamaged. Do not use the device if the packaging has already been opened or is damaged.
- 8.6 The device must be used before the expiry date stated on the packaging ('Use-by' date). The device must no longer be used after the expiry date has passed and must be disposed of.
- 8.7 Always ensure a sterile surgical field and sterile procedures when handling or using the device.
- 8.8 Do not use the device if it appears to be damaged or defective.
- 8.9 **FOR SINGLE USE ONLY!** The device is designed for single use. It must not be reused, resterilised or reprocessed. Reuse, resterilisation or reprocessing can negatively impact the performance of the device and may result in injury to the patient if this device is used. Reuse, resterilisation or reprocessing may also increase the risk of contamination of the device and/or infection of the patient, including but not limited to transmission of infectious diseases from one patient to another.
- 8.10 The performance features of the device must be suitable for use in the intended intervention.
- 8.11 Do not bend the pre-bent shaft excessively during manual modification. This may cause damage to the core wire and the device to lose function as a result.
- 8.12 During the intervention, the heart should be monitored by means of an electrocardiogram (ECG) and other vital signs.
- 8.13 During the intervention, suitable diagnostic imaging, e.g. a C-arm for diagnostic radiography and/or ultrasound, should be used at all times to check whether the device is in the correct position.
- 8.14 If unusual resistance is encountered during the intervention, the intervention must be discontinued. Before the intervention can be continued, the cause of the resistance must be identified. If the cause of the resistance cannot be identified, the device must be removed and the intervention must be discontinued.
- 8.15 Particular caution should be taken when passing aneurysmal bulging in the vascular system and the heart in the area of application.
- 8.16 Always operate the device with caution. The device may be damaged if excessive force is exerted on the grip. It is not necessary to exert excessive force to cut through the tissue. The cutting jaws are sharp and cut with relatively little pressure.
- 8.17 Unintentionally opening the cutting jaws after taking the tissue samples could lead to the tissue sample being lost and result in an embolism.
- 8.18 Clean the device thoroughly in a sterile physiological saline solution (NaCl 0.9%) during the intervention after each tissue sample is taken in order to remove residual tissue and to prevent a possible subsequent embolism.

9. Packaging, shelf life, storage and disposal

The device is supplied sterile and non-pyrogenic. It is individually packaged. Sterility is guaranteed as long as the packaging is not opened or damaged.

Note the expiry date on the corresponding packaging. The device must no longer be used after the expiry date has passed and must be disposed of.

Store and transport the device in a cool, dark and dry place until use. Keep away from sunlight and heat.

Do not use the device if the packaging has already been opened or is damaged. If damage of any sort is identified, the device must no longer be used.

The packaging material used, comprising paper, cardboard and plastic, can be recycled and safely disposed of through existing disposal channels. It is considered non-problematic waste.

The product poses no biological hazards before its use. If the device is unintentionally opened or is damaged, it can be disposed of in accordance with hospital hygiene guidelines. In order to prevent any injuries from cuts, dispose of the device in its original protective packaging.

The product poses potential biological hazards (e.g. microorganisms, viruses, allergens) after its use. It is therefore considered to be potentially infectious waste. If the device comes into contact with pathogens which are subject to mandatory reporting during use, it is considered to be hazardous waste and should be marked with the 'biohazard' symbol upon disposal. It must be handled and disposed of in accordance with recognised medical procedures and in line with local, state and federal laws and regulations.

Once it has been used and therefore contaminated, immediately place the device into tear-resistant, moisture-resistant and sealed containers at the waste-generation site and take it to a central collection point in sealed containers without being transferred or sorted.

10. Summary of safety and clinical performance

The summary of safety and clinical performance pursuant to Article 32 of Regulation (EU) 2017/745 can be retrieved from the European database on medical devices (EUDAMED) using the Basic UDI-DI 42502039 8020000 54 as soon as it becomes available (<https://ec.europa.eu/tools/eudamed>). Alternatively, the summary of safety and clinical performance can be requested using the manufacturer's contact details on the last page of these instructions for use.

11. Reporting obligations for severe incidents

The user must immediately report all severe incidents occurring in causal relationship with the device (death of a patient, user or other person; temporary or permanent severe deterioration in the health of a patient, user or other person; serious risk to public health) to the manufacturer as well as the competent authority of the country in which the user is established. The unique device identifier (UDI) required for this purpose as well as other device-related information can be found on the packaging label.

The affected device that is causally linked to a severe incident must be sent immediately to the manufacturer for the necessary inspections and risk assessment by the manufacturer in accordance with Article 89 of Regulation (EU) 2017/745.

12. Use

12.1. Device selection and preparation

- The physician uses his or her knowledge and experience to decide at his or her own discretion, via which access site the product is used.
- Select the correct length of the product for the selected access site (internal jugular vein or femoral veins).
- Open the packaging and place the device on a sterile operating table.

 **WARNING: Always ensure a sterile surgical field and sterile procedures when handling or using the device.**

- Check that the device is working properly by opening and closing the cutting jaws (Figures 2.1 and 2.2, see page 8).

 **WARNING: Do not use the device if it appears to be damaged or defective.**

- Depending on the selected access site the already pre-bent shaft of the device can be further modified manually if necessary. To do this, hold the distal end between your thumb and index finger (Figure 3.1). Run your fingers along the shaft towards the tip and bend it to the desired angle while doing so (Figures 3.2 and 3.3). Check that the device is working properly after the modification.

 **WARNING: Do not bend the pre-bent shaft excessively during manual modification. This may cause damage to the core wire and the device to lose function as a result.**


Figure 3.1


Figure 3.2


Figure 3.3

- Before inserting the device for the first time, wipe it with a compress moistened with sterile physiological saline solution (NaCl 0.9%).

12.2. Biopsy sampling

- Insert the device through the introducer sheath, which has already been positioned.

 **DANGER: During the intervention, the heart should be monitored by means of an electrocardiogram (ECG) and other vital signs.**

 **DANGER: During the intervention, suitable diagnostic imaging, e.g. a C-arm for diagnostic radiography and/or ultrasound, should be used at all times to check whether the device is in the correct position.**

 **WARNING: If unusual resistance is encountered during the intervention, the intervention must be discontinued. Before the intervention can be continued, the cause of the resistance must be identified. If the cause of the resistance cannot be identified, the device must be removed and the intervention must be discontinued.**

 **WARNING: Particular caution should be taken when passing aneurysmal bulging in the vascular system and the heart in the area of application.**

- In addition to visually observing the device, the course of the bend of the pre-bent shaft can also be identified on the basis of the moulded arrow on the grip, which acts as a bend indicator (Figure 1, ref. 4, see page 8). The pre-bent shaft is oriented by rotating the grip. In this process, the shaft of the device is held between the fingers.
- When the tip of the device is positioned in the ventricle, gently push the closed cutting jaws forward until you feel resistance. Confirm that the device is correctly positioned and directed using suitable diagnostic imaging, e.g. a C-arm for diagnostic radiography and/or ultrasound.
- Pull the device back by approx. 1 cm and open the cutting jaws. To do this, gently pull back the thumb ring of the grip and slide away from the index finger and middle finger (Figure 2.1, see page 8).
- Push the opened cutting jaws forward as far as the planned sampling point. Release the tension on the thumb ring, gently exert pressure on the thumb ring, close the cutting jaws (Figure 2.2, see page 8), and thereby extract the tissue sample.

 **WARNING: Always operate the device with caution. The device may be damaged if excessive force is exerted on the grip. It is not necessary to exert excessive force to cut through the tissue. The cutting jaws are sharp and cut with relatively little pressure.**

- With the cutting jaws closed, pull the device back through the introducer sheath and pull out of the introducer sheath.

 **WARNING: Unintentionally opening the cutting jaws could lead to the tissue sample being lost and result in an embolism.**

- Open the cutting jaws to remove the tissue sample. To do this, gently pull back the thumb ring of the grip (Figure 2.1, see page 8).
- After removing the tissue sample and before the device is reinserted, clean the device in a sterile physiological saline solution (NaCl 0.9%) by repeatedly opening and closing the cutting jaws.

 **WARNING: Clean the device thoroughly in a sterile physiological saline solution (NaCl 0.9%) during the intervention after each tissue sample is taken in order to remove residual tissue and to prevent a possible subsequent embolism.**

- Once you have finished using the device, dispose of it in accordance with Section 9.

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Joline GmbH & Co. KG