Joline®

QuattroMix®

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Contents

1. Device description
 - 1.1 General device description, clinical use and functional principles
 - 1.2 Performance characteristics / technical specifications
 - 1.3 Compatibility with other devices
2. Indication
3. Contraindication
4. Complications
5. Patient group
6. User profile
7. Use environment
8. Warnings and precautions
9. Packaging, shelf life, storage and disposal
10. Summary of safety and clinical performance
11. Reporting obligations for serious incidents
12. Use

1. Device description

1.1 General device description, clinical benefit and mode of operation

General device description

The QuattroMix® is a device for mixing and transferring bone cement to the Vertebra Filler Nozzles of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System that is used to inject bone cement during kyphoplasty. The individual components of the QuattroMix® are the mixing cylinder, mixing unit, union spindle, funnel and stand. The QuattroMix® is intended for single use only. The period of use is usually 10 minutes.

Clinical benefit

The QuattroMix® is intended to enable the user to uniformly mix and transfer the bone cement to the Vertebra Filler Nozzles of the Vertebra Filler Devices in the Joline ALLEVO Kyphoplasty System. The closed construction is intended to minimise any unpleasant odours that the user is exposed to due to the liquid monomer components of PMMA bone cements.

Functionality

The QuattroMix® is used to mix and transfer the bone cement to the Vertebra Filler Nozzles. The bone cement is mixed and transferred using the QuattroMix® as follows:

Mixing is carried out manually using the mixing unit. The unit is placed on the mixing cylinder and the mixing propeller can be moved axially and rotationally in the mixing cylinder. When mixing the two individual components of the bone cement, they are homogeneously combined, which starts the polymerisation of the bone cement. The mixing process is complete once the desired viscosity of the bone cement is reached (see Instructions for use of the bone cement).

For the transfer, the bone cement is transferred to several Vertebra Filler Nozzles using the Luer lock connector. To do so, the union spindle is placed over the mixing unit and the one-way stopcock is opened. Using manual rotating movements of the union spindle, the unit is moved forwards axially, which transfers doses of the pasty bone cement into the Vertebra Filler Nozzles.

1.2 Performance features/ technical specification

- Closed system
- Filling volume: 1 packaging unit bone cement (see the instructions for use for the bone cement)
- One-way stopcock with sealing option
- Luer lock connector to connect the Vertebra Filler Nozzle of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System
- Union spindle for simple transfer of the bone cement

1.3 Compatibility with other devices

The QuattroMix® is designed to be used in combination with the following devices:

- Radiopaque bone cement based on polymethylmethacrylate (PMMA)
- Vertebra Filler Nozzles of the following Vertebra Filler Devices of the Allevo Kyphoplasty System:
 - Vertebra Filler Device (frontal opening)

- Vertebra Filler Device (side opening)
- Vertebra Filler Device flexPlunger (frontal opening)

2. Indication

The QuattroMix® is intended for mixing and transferring bone cement during kyphoplasty.

3. Contraindication

There are no contraindications known for use of the QuattroMix®.

4. Complications

The following complications are associated with the QuattroMix®:

- Infection, including deep or superficial wound infection
- Inhomogeneity of the bone cement with associated consequences for the patient

5. Patient group

Provided that the manufacturer of the bone cement does not specify any limitations on the patient group in the intended use, there are also no limitations for the QuattroMix® because it does not make direct contact with the patient.

6. User profile

The QuattroMix® is intended to be operated by OR personnel. It should only be used by or under the supervision of an experienced user or in accordance with his or her instructions.

7. Use environment

The QuattroMix® is used in a sterile environment such as an operating theatre during kyphoplasty. The QuattroMix® is used under sterile conditions.

Due to the functional principle and properties of the bone cements, the QuattroMix® is only suitable for use within the application time specified for the bone cement because otherwise the function of the QuattroMix® may be impaired.

For application of the bone cement used, the specification of the bone cement manufacturer must also be complied with.

8. Warnings and precautions

- See the 'Explanation of symbols' to understand the symbols used on the packaging. The contents specified on the label are translated under 'Translation of the label'.
- The medical techniques and procedures described here make NO claim to completeness. They do not illustrate all the medically possible procedures and, most of all, do not replace the doctor's own experience regarding the use and the treatment of certain patients.
- The device may only be used by an experienced user, under his or her supervision or in accordance with his or her instructions.
- The device has been sterilised with ethylene oxide.
- The device is only sterile and pyrogen-free if the sterile packaging is unopened and undamaged. Do not use the device if the packaging has already been opened or is damaged.
- The device must be used before the expiry date stated on the packaging ('Use-by' date). The device must no longer be used after the expiry date has passed and must be disposed of.
- Always ensure a sterile surgical field and sterile procedures when handling or using the device.
- Do not use the device if it appears to be damaged or defective.
- FOR SINGLE USE ONLY!** The device is designed for single use. It must not be reused, resterilised or reprocessed. Reuse, resterilisation or reprocessing can negatively impact the performance of the device and may result in injury to the patient if this device is used. Reuse, resterilisation or reprocessing may also increase the risk of contamination of the device and/or infection of the patient, including but not limited to transmission of infectious diseases from one patient to another.
- The device is designed for use with the Vertebra Filler Nozzles of the Vertebra Filler Devices in the Joline Allevo Kyphoplasty System and with the PMMA bone cements distributed by Joline.
- Always follow the instructions for use of the bone cement manufacturer.
- Do not exceed the maximum recommended capacity of the QuattroMix® of one packaging unit of the bone cement (approx. 25 g polymer powder and approx. 10 ml monomer liquid) to ensure that the bone cement is mixed homogeneously.
- Small glass fragments may result from breaking the neck of the ampoule. To prevent them from getting into the bone cement and consequently into the vertebral body, do not open the ampoule directly over the mixing cylinder but open it instead to the side.
- Only attach the funnel after the monomer liquid has been added because otherwise the polymer powder will stick to the sides when poured through the wet funnel.
- Close the one-way stopcock (horizontal position) so that no bone cement escapes during mixing.
- When mixing, hold the mixing cylinder upright with the lid facing upwards to prevent air inclusions and thus inhomogeneity of the bone cement.
- Mix the components of the bone cement thoroughly with the mixing blades to prevent inhomogeneity of the bone cement.
- Holding the mixing cylinder in the palm of your hand could heat up the bone cement, which may significantly shorten the processing times for the bone cement.

- 8.19 The two lateral fixtures on the mixing unit must clearly and audibly latch into place so that the bone cement can subsequently be transferred safely.
- 8.20 Fully vent the QuattroMix® to prevent air inclusions and thus inhomogeneity of the bone cement. When doing so, hold the QuattroMix® away from your face in case dust from any residual polymer powder escapes.
- 8.21 After mixing, transfer the bone cement as quickly as possible for such time that the increasing viscosity of the bone cement allows transfer.
- 8.22 Do not add the liquid monomer component via the funnel because otherwise the polymer powder will stick to the sides of the wet funnel. Instead, add the liquid monomer component directly to the mixing cylinder and then place the funnel on the cylinder.
- 8.23 In order to ensure that the mixing unit is firmly secured to the mixing cylinder, apply pressure when placing it on the mixing cylinder until it reaches the stop. Then turn it in the indicated direction until the bayonet mount audibly clicks into place.
- 8.24 Ensure that the one-way stopcock is closed before mixing.
- 8.25 The one-way stopcock must be opened before filling the Vertebra Filler Nozzles in order to prevent damage to the device due to any pressure that may build up when turning the union spindle.

9. Packaging, shelf life, storage and disposal

The device is supplied sterile and pyrogen-free. It is individually packaged. Sterility is guaranteed as long as the packaging is not opened or damaged.

Note the expiry date on the corresponding packaging. The device must no longer be used after the expiry date has passed and must be disposed of.

Store and transport the device in a cool, dark and dry place until use. Keep away from sunlight and heat.

Do not use the device if the packaging has already been opened or is damaged. If damage of any sort is identified, the device is no longer permitted to be used.

The packaging material used, comprising paper, cardboard and plastic, can be recycled and safely disposed of through existing disposal channels. It is considered non-problematic waste.

There are no hazards associated with the device before its use. If the device is unintentionally opened or is damaged, it can be disposed of in accordance with hospital hygiene guidelines.

Residual bone cement should be allowed to cure fully in the QuattroMix® before it can be disposed of with hospital waste in accordance with general and hospital-specific precautions.

10. Summary of safety and clinical performance

The summary of safety and clinical performance pursuant to Article 32 of Regulation (EU) 2017/745 is only to be drawn up for implantable devices and class III devices that are not custom-made devices or investigational devices. Since the QuattroMix® does not fall under this definition, a summary of safety and clinical performance is not applicable and was therefore not compiled.

11. Reporting obligations for serious incidents

The user must immediately report all serious incidents occurring in causal relationship with the device (death of a patient, user or other person; temporary or permanent serious deterioration in the health of a patient, user or other person; serious risk to public health) to the manufacturer as well as the competent authority of the country in which the user is established. The unique device identifier (UDI) required for this purpose as well as other device-related information can be found on the packaging label.

The affected device that is causally linked to a serious incident must be sent immediately to the manufacturer for the necessary inspections and risk assessment by the manufacturer in accordance with Article 89 of Regulation (EU) 2017/745.

12. Use

The QuattroMix® is used in accordance with the figures in the Quick Guide included in the package as follows:

I Filling

- 12.1 Place the stand on a level and slip-resistant surface. Place the mixing cylinder in the stand with the one-way stopcock closed and with the opening facing upwards.

 **CAUTION: Close the one-way stopcock (horizontal position) so that no monomer component escapes during mixing.**

- 12.2 Pour the liquid monomer component of the bone cement into the mixing system. Place the funnel on the mixing cylinder. Fill the powder polymer component of the bone cement into the mixing cylinder through the funnel.

 **WARNING: Small glass fragments may result from breaking the neck of the ampoule. To prevent them from getting into the bone cement and consequently into the vertebral body, do not open the ampoule directly over the mixing cylinder but open it instead to the side.**

 **CAUTION: Do not add the liquid monomer component via the funnel because otherwise the polymer powder will stick to the sides of the wet funnel. Instead, add the liquid monomer component directly to the mixing cylinder and then place the funnel on the cylinder.**

 **WARNING: Do not exceed the maximum recommended capacity of the QuattroMix® of one packaging unit of the bone cement (see the instructions for use for the bone cement) to ensure that the bone cement is mixed homogeneously.**

 **CAUTION: Check to ensure that the entire quantity of polymer powder has been added to the mixing cylinder and that no residues are stuck to the funnel.**

Note: The QuattroMix® is exclusively intended for use with radiopaque bone cement based on PMMA.

II Preparing

- 12.3 Remove the funnel from the mixing cylinder. Place the mixing cylinder on the mixing unit and turn it in the direction indicated by an arrow until it is firmly locked in place. Remove the mixing cylinder from the stand.

 **CAUTION:** In order to ensure that the mixing unit is firmly secured to the mixing cylinder, apply pressure when placing it on the mixing cylinder until it reaches the stop. Then turn it in the indicated direction until the bayonet mount audibly clicks into place.

 **CAUTION:** The two lateral fixtures on the mixing unit must clearly and audibly click into place so that the bone cement can subsequently be transferred safely.

Note: If the mixing unit is not firmly placed on the mixing cylinder, it may become detached due to pressure from the mixing cylinder.

III Mixing

- 12.4 Hold the mixing cylinder in one hand. Mix the components of the bone cement by using the other hand to push the mixing unit back and forth along its entire length while simultaneously performing oscillating rotations. Ensure that the piston is also moved all the way forward here in order to move the entire bone cement mass.

 **CAUTION:** When mixing, hold the mixing cylinder upright with the one-way stopcock facing upwards to prevent air inclusions and thus inhomogeneity of the bone cement.

 **WARNING:** Mix the components of the bone cement thoroughly with the mixing blades to prevent inhomogeneity of the bone cement.

 **CAUTION:** Holding the mixing cylinder in the palm of your hand could heat up the bone cement, which may significantly shorten the processing times for the bone cement.

 **CAUTION:** Ensure that the one-way stopcock is closed before mixing.

- 12.5 Once the recommended mixing time has elapsed (as specified in the instructions for use of the bone cement manufacturer), pull back on the piston on the mixing unit until it reaches its stop point.

 **CAUTION:** Always follow the instructions for use of the bone cement manufacturer.

IV Venting and transferring

- 12.6 Slide the union spindle over the mixing unit that has been pulled back to the stop point. Apply heavy pressure to connect it to the mixing unit. A click will indicate that the union spindle is connected to the mixing unit.

 **CAUTION:** The one-way stopcock must be opened before filling the Vertebra Filler Nozzles in order to prevent damage to the device due to any pressure that may build up when turning the union spindle.

- 12.7 Remove the plunger from the Joline Vertebra Filler Device. Connect the Vertebra Filler Nozzle to the Luer lock connector of the one-way stopcock to vent. Hold the mixing cylinder upright as described in step 12.4. Push the piston forwards by turning the union spindle in a clockwise direction until bone cement exits the Vertebra Filler Nozzle.

Note: The capacity of 1.5 ml is reached once bone cement exits the Vertebra Filler Nozzle.

 **WARNING:** Fully vent the QuattroMix® to prevent air inclusions and thus inhomogeneity of the bone cement. When doing so, hold the QuattroMix® away from your face in case dust from any residual polymer powder escapes.

- 12.8 Close the one-way stopcock in order to prevent any excess bone cement from exiting.

- 12.9 Remove the Vertebra Filler Nozzle from the one-way stopcock.

Note: Wipe the excess bone cement from the Luer lock connector before connecting any further Vertebra Filler Nozzles.

- 12.10 Repeat steps 12.7 through 12.9 with additional Vertebra Filler Devices until the required quantity of bone cement has been transferred.

 **CAUTION:** Transfer the bone cement after mixing as quickly as possible for as long as the increasing viscosity of the bone cement allows transfer.

 **CAUTION:** Always follow the instructions for use of the bone cement manufacturer.

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Joline GmbH & Co. KG