

Name	Maria Musterfrau	Scope of analysis	285 Allergens
Age	37 Years	Analysis procedure	285 ALEX 2 - Allergy Explorer
Sample-ID	Test_ALEX2_Mit		IgE-multiplex-diagnostics
Date of analysis	28. February 2020	Analysis laboratory	Macro Array Diagnostics, Vienna

Thank you very much for choosing igevia for your allergy testing! The laboratory analysis of your blood for the presence of immunoglobulin E (short IgE) antibodies determines the Total IgE and the Allergen-specific IgE.

1. Your Total IgE

62 kU/L

The Total IgE-value shows your personal allergy tendency.

Less than 20 kU/L <i>Allergy unlikely</i>	20 - 100 kU/L <i>Allergy possible</i>	More than 100 kU/L <i>Allergy likely</i>
--	--	---

2. Your Specific IgE at a glance

The overview shows the total number of tested allergens and their assignment to the RAST classes 0-4 based on your laboratory results. RAST classes show the level of the Specific IgE-level.

	0 Negative	1 Low	2 Moderate	3 High	4 Very high
Inhalative allergen sources	0	1	2	3	4
Pollen					
27 Tree Pollen	22	3	2		
12 Grass Pollen	3	2	3	4	
19 Weed Pollen	19				
Mites & Cockroaches					
19 Storage- & House dust mites	17	1	1		
7 Cockroach	7				
Animals					
18 Pets	16		2		
7 Farm Animals	7				
Microorganisms & Spores					
11 Fungal Spores & Yeast	10	1			
Other allergen sources	0	1	2	3	4
6 Latex	6				
3 Others	3				
Ingested allergen sources	0	1	2	3	4
Plant Food					
11 Vegetables	10	1			
15 Cereals	15				
7 Spices	7				
15 Legumes	15				
29 Nuts & Seeds	28	1			
19 Fruits	16	3			
Animal-Derived Food					
7 Egg	7				
31 Fish & Seafood	31				
12 Meat	12				
8 Milk	8				

You can find more details about these allergen sources on the next page in section 3.

3.

Your detailed results for the Specific IgE

Tested allergen	E / C ¹⁾	Value ²⁾	Allergen family / cross-allergy ³⁾	Options for action				Technical term for research ⁸⁾
				Cat ⁴⁾	Avoidance ⁵⁾	Cook ⁶⁾	Immuno-therapy ⁷⁾	
4 Allergens with high IgE-levels (5 - 15 kU_A/L)								
Cultivated rye, Pollen	E	6.28		🔔			👉	Secale cereale
Perennial Ryegrass (Lol p 1)	C	6.24	Beta-Expansin	🔔			👉	Lolium perenne
Timothy grass (Phl p 1)	C	14.33	Beta-Expansin	🔔			👉	Phleum pratense
Timothy grass (Phl p 5.0101)	C	9.9	Grass Group 5/6	🔔			👉	Phleum pratense
8 Allergens with moderate IgE-levels (1 - 5 kU_A/L)								
Bahia grass	E	1.1		🔔				Paspalum notatum
Bermuda grass (Cyn d 1)	C	3.67	Beta-Expansin	🔔				Cynodon dactylon
Timothy grass (Phl p 6)	C	3.93	Grass Group 5/6	🔔			👉	Phleum pratense
Silver birch (Bet v 1)	C	2.35	PR-10	🔔			👉	Betula verrucosa
Hazel pollen (Cor a 1.0103)	C	1.72	PR-10	🔔			👉	Corylus avellana
Rabbit, epithel (Ory c 2)	C	1.2	Lipophilin	🔔	🚫			Oryctolagus cuniculus
Rabbit, epithel (Ory c 3)	C	2.64	Uteroglobin	🔔	🚫			Oryctolagus cuniculus
American house dust mite (Der f 2)	C	1.27	NPC2 Family	🔔			👉	Dermatophagoides farinae
11 Allergens with low IgE-levels (0.3 - 1 kU_A/L)								
Common reed	E	0.94		🔔				Phragmites communis
Timothy grass (Phl p 12)	C	0.65	Profilin	🔔			👉	Phleum pratense
Hazel pollen	E	0.65		🔔			👉	Corylus avellana
Silver birch (Bet v 2)	C	0.35	Profilin	🔔			👉	Betula verrucosa
Date palm (Pho d 2)	C	0.61	Profilin	🔔				Phoenix dactylifera
European house dust mite (Der p 2)	C	0.92	NPC2 Family	🔔			👉	Dermatophagoides pteronyssinus
Melon (Cuc m 2)	C	0.53	Profilin	🔔	🚫	🔥		Cucumis melo
Strawberry (Fra a 1+3)	C	0.82	PR-10+LTP	🔔	🚫			Fragaria ananassa
Apple (Mal d 1)	C	0.31	PR-10	🔔	🚫	🔥		Malus domestica

Name:

Celery (Api g 1)	C	0.91	PR-10	☞	⊘	🔥		Apium graveolens
Hazelnut (Cor a 1.0401)	C	0.85	PR-10	☞	⊘	🔥		Corylus avellana
262 Allergens with negative IgE-levels (less than 0.3 kU _A /L)								
For these allergens no sensitization was found in your blood.								

4. Your indicator for cross allergies




The laboratory analysis shows that you react positively to 2 or more allergens from the same allergen family
This is a clear indication of a cross-allergy. Please bear this in mind when choosing further options for action or when drawing up a therapy plan with your allergy specialist.

Additional information and explanations for your detailed results




- 1) E = Extract: An extract is a natural substance that is prepared as a whole for testing.
K = Component: A component is an allergenic molecule and thus part of an extract. It provides more specific results than an extract and is an important aid in determining cross-reactions between allergens.
- 2) Value = result of the blood analysis for the IgE-level, measured in kilo units per litre (kU_A/L)
- 3) Allergen family / cross-allergy: Allergens within an allergen family are very similar. The probability is high that one reacts to several allergens within an allergen family of the same name. This is called a cross-allergy or cross-reaction.
- 4) Cat (category): ☞ = Inhalative allergens (via air), ☞ = Ingested allergens (via food), ☞ = Other allergens
- 5) Avoidance ⊘: With these allergens, there is the possibility of alleviating symptoms by taking a leave of absence from allergens. This means avoiding contact with the allergenic substance. Important: Never change your diet fundamentally without professional guidance.
- 6) Cook 🔥: These allergens can be rendered harmless by cooking (heating above 100 degrees Celsius) and then consumed.
- 7) Immunotherapy ☞: For these allergens there is the possibility of immunotherapy / hyposensitization. Please talk to your allergy specialist about it.
- 8) Technical term for research: These technical terms simplify further research on the single allergen sources.

A list of all tested allergens can be found on the next page.

List of all tested allergens with the number of tested extracts (E) and components (C)

 Inhalative allergen sources	E	C	 Inhalative allergen sources	E	C	 Ingested allergen sources	E	C
Pollen			Mites & Cockroaches			Plant Food		
Tree Pollen	11	16	Storage- & House dust mites	2	17	Vegetables	5	5
Acacia	1		Acarus siro	1		Carrot	1	1
Arizona-Cypress		1	American house dust mite		2	Potato	1	
Silver birch		3	Blomia tropicalis		3	Garlic	1	
Beech		1	European house dust mite		9	Celery		3
Date palm		1	Glycyphagus domesticus		1	Tomato	1	1
Alder		2	Lepidoglyphus destructor		1	Onion	1	
Ash	1	1	Tyrophagus putrescentiae	1	1	Cereals	10	5
Tree of Heaven	1		Cockroach	1	6	Common buckwheat	1	1
Hazel	1	1	American Cockroach	1	1	Spelt	1	
Sugi		1	German Cockroach		5	Barley	1	
Mulberry tree	1		Animals			Oat	1	
Olive		2	Pets	1	16	Millet	1	
Paper mulberry	1		Djungarian hamster		1	Lupine seed	1	
Cottonwood	1		Dog		6	Sweet corn	1	1
Plane tree		3	Dog Urine	1		Quinoa	1	
Elm	1		Rabbit		3	Rice	1	
Walnut	1		Cat		4	Cultivated rye	1	
Mountain cedar	1		House mouse		1	Wheat		3
Cypress	1		Guinea Pig		1	Spices	6	1
Grass Pollen	4	8	Farm Animals	4	4	Anise	1	
Bahia grass	1		Horse		3	Caraway	1	
Bermuda grass	1	1	Rat	1		Oregano	1	
Timothy grass		6	Cattle		1	Peppers	1	
Cultivated rye	1		Sheep	1		Parsley	1	
Common reed	1		Pig	1		Mustard	1	1
Ryegrass		1	Goat	1		Legumes	4	11
Weed Pollen	9	10	Microorganisms & Spores			Pea	1	
Common Pigweed	1		Fungal Spores & Yeast	3	10	Peanut		7
Mugwort	1	2	Alternaria alternata		2	Green bean	1	
Annual mercury		1	Aspergillus fumigatus		4	Chickpea	1	
Nettle	1		Yeast	1		Lentil	1	
Glass herb	1	1	Cladosporium herbarum	1	1	Soya		4
Hemp		1	Malassezia sympodialis		3			
Hemp (CBD)	1		Penicillium chrysogenum	1				
Russian thistle	1	1						
Ribwort	1	1						
Ragweed	1	2						
Lamb's quarter	1	1						

List of all tested allergens with the number of tested extracts (E) and components (C)

 Ingested allergen sources	E	C	 Ingested allergen sources	E	C	 Other allergen sources	E	C
Plant Food			Animal-Derived Food					
Nuts & Seeds	5	17	Egg	2	5	Latex		6
Fenugreek seeds	1		Chicken egg yolk	1	1	Latex		6
Cashew	1	2	Egg white	1	4	Others	1	2
Hazelnut		5	Fish & Seafood	16	15	Birch fig	1	
Pumpkin seed	1		Anisakis simplex		2	Hom s Lactoferrin		1
Macadamia	1	1	Atlantic cod	1	2	Grape tick		1
Almond	1		Atlantic herring	1	1			
Poppy seed	1	1	Oyster	1				
Brazil nut	1	1	Black Tiger Shrimp		4			
Pecan nut	1		Northern prawn	1				
Pistachio		3	Lobster	1				
Sesame	1	1	Scallop	1				
Sunflower seed	1		Carp		1			
Walnut		5	Crab	1				
Fruits	9	11	Salmon	1	1			
Apple		3	Mackerel	1	1			
Avocado	1		Common mussel	1				
Banana	1		Thornback ray	1	1			
Pear	1		Brown shrimp		1			
Strawberry		1	Swordfish		1			
Fig	1		Shrimp mix	1				
Blueberry	1		Tuna	2				
Cherry	1		Squid	1				
Kiwi		4	Clam	1				
Mango	1		Meat	10	2			
Melon		1	Cricket	1				
Orange	1		Chicken	1				
Papaya	1		Rabbit	1				
Peach		1	Lamb	1				
Grape		1	Mealworm	1				
			Horse	1				
			Cattle	1	1			
			Pig	1	1			
			Turkey	1				
			Migratory locust	1				
			Milk	5	3			
			Camel milk	1				
			Cow milk	1	3			
			Sheep milk	1				
			Mare's milk	1				
			Goat milk	1				

Important information for igevia-Customers at home

As part of the laboratory analysis, your blood was examined for the presence of Immunoglobulin E (short IgE) antibodies . They are responsible for the occurrence of allergic reactions and therefore play a decisive role in the detection of allergies.

Important note for you

Your personal allergen-profile is not a substitute, but rather the basis for a profound expert advice.
Please always consult your allergy specialist if you have any medical questions.

This is an allergy

An allergy is a hypersensitive reaction of the immune system to normally harmless, exogenous substances (e.g. pollen, food). After contact with these substances (= allergens), the body starts to form Specific IgE-antibodies as a "defence reaction" and is thus "sensitised". In case of further contact with the allergenic substance, the body may (but does not have to) react with allergic symptoms (e.g. watering eyes, sneezing fits, itching).

Only if Allergen-specific IgE-antibodies occur together with allergic symptoms it is called a clinically manifested allergy . That is why it is important not only to analyse the blood, but also to record the symptoms precisely. This can be done directly during a conversation with a qualified person or by means of a symptom questionnaire. You can use the symptom questionnaire in your personal login area under www.igevia.com/login.

The analysis procedure

For the analysis of your blood sample, our partner laboratory used the testing method "ALEX 2 - Allergy Explorer" from the company Macro Array Diagnostics. This is the first multiplex allergy test that can simultaneously test Total and Specific IgE for a wide range of allergen sources (extracts and components). This allows to create an almost complete allergen-profile for each person tested. As a result tailor-made dietary recommendations can be made and allergen-specific immunotherapy (AIT) can be implemented to target your specific allergy.

The next steps with your allergen-profile

With this allergen-profile and a positive symptom questionnaire you have an important basis to get your allergy under control.

- Step 1:** Try to avoid allergens or reduce them as much as possible. You will find numerous tips and recommendations in your login area on the website www.igevia.com/login.
- Step 2:** Talk to your allergy specialist or pharmacist about the options for immunotherapy or symptomatic treatment (e.g. antihistamines).

Are you happy with igevia?

We would love you to recommend us to friends or family.

Please visit our [Facebook-Page](#) or spread it via word of mouth.
If you have any questions regarding your results, please feel free to contact us at any time:

Telephone +43 5 99 07 99
Email contact@igevia.com

Your igevia-Team