

**FACTOR
VIII**

Pre-analytical
heat treatment
step, to inactivate
endogenous
Factor VIII

Ready
to use
components

Excellent
repeatability and
reproducibility

Reduce
inter-laboratory
variation



NEW CRYOcheck™ Factor VIII Inhibitor Kit

For laboratories wishing to determine and quantify the presence of functional FVIII inhibitors, now there is a ready-to-use sample preparation kit to facilitate the performance of a modified Nijmegen-Bethesda assay.

With standardised components and a validated protocol, it brings you excellent repeatability and reproducibility.

α alpha laboratories
supplying quality to science

CE IVD


PrecisionBioLogic

Increase Confidence in Your Results

The **CRYOcheck™** FVIII Inhibitor Kit, delivers excellent repeatability and reproducibility to reduce intra and inter laboratory variation.

Repeatability

- Inhibitor-negative and inhibitor positive plasma from patients with haemophilia A were combined to yield a panel of test plasmas at four different levels of inhibitor.
- A **3 lot x 20 day x 2 run x 2 repeat** precision study was performed using the **CRYOcheck™** FVIII Inhibitor Testing Kit to determine the repeatability and within laboratory imprecision.
- Below, the mean measured titre is shown as a solid line, and the 2 standard deviations limit as a dotted line.



Sample	N	Mean Value (BU/mL)	Repeatability		Between-Run		Between-Day		Within Lab	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
Kit Negative Control	80	0	0.0	–	0.0	–	0.0	–	0.1	–
Kit Positive Control	80	1.6	0.1	8.8%	0.1	3.5%	0.1	5.5%	0.2	11%
Below Cut-off Test Plasma	80	0.3	0.1	–	0.0	–	0.0	–	0.1	–
Low Titre Test Plasma	80	1.2	0.1	6.4%	0.0	3.8%	0.0	2.3%	0.1	7.8%
Mid Titre Test Plasma	80	5.3	0.3	5.9%	0.2	3.9%	0.1	1.3%	0.4	7.2%
High Titre Test Plasma	80	8.6	0.6	7.6%	0.4	4.6%	0	0%	0.7	8.9%

The History of Haemophilia

1803

A bleeding disorder that mainly affects men is recognised

1828

This bleeding disorder is given the name haemophilia

1840

The first blood transfusion to treat haemophilia

1941

The existence of FVIII inhibitors are reported

1947

It's reported that there is more than one type of haemophilia

1950

The Haemophilia Society is established as a charity

1952

The first identification of 'Christmas disease', now known as haemophilia B

1963


The World Federation of Haemophilia is established

1964

The first successful use of cryoprecipitate to treat patients with haemophilia


Reproducibility

- The same panel of inhibitor positive plasmas were subjected to a **3 site x 5 day x 2 run x 3 repeat** precision study using the **CRYOcheck™** FVIII Inhibitor Testing Kit.
- The testing was performed on three different analysers with three different operators and a single lot of chromogenic FVIII kit.
- Below, the mean measured titre is shown as a solid line, and the 2 standard deviations limit as a dotted line.



Sample	N	Mean Value (BU/mL)	Repeatability		Between-Run		Between-Day		Between Site		Reproducibility	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
Kit Negative Control	90	0	0.1	–	0	–	0	–	0	–	0.1	–
Kit Positive Control	90	1.8	0.2	8.8%	0.1	3.4%	0.1	3.4%	0.2	12.4%	0.3	15.9%
Below Cut-off Test Plasma	90	0.3	0.1	–	0	–	0.0	–	0.0	–	0.1	–
Low Titre Test Plasma	90	1.4	0.1	10.1%	0	0%	0.0	1.9%	0.0	1.1%	0.1	10.3%
Mid Titre Test Plasma	90	5.5	0.4	7.8%	0.1	2.2%	0.1	2.0%	0	0%	0.5	8.3%
High Titre Test Plasma	89	9.6	0.7	7.8%	0.4	4.1%	0.3	3.3%	0.9	9.2%	1.3	13.1%

Data source: Douglas CD, et al. ¹



Factor VIII Inhibitor Kit

- For use in performing a modified Nijmegen-Bethesda assay
- Standardised components with a validated method, working to reduce inter laboratory variation
- Excellent repeatability and reproducibility
- Platform independent - FVIII activity can be measured using a 1-stage clot based assay or a chromogenic assay on any coagulation analyser
- Ready to use components, simply thaw in a water bath in a few minutes - saving time, improving workflow and eliminating reconstitution errors
- Includes a pre-analytical heat treatment step to inactivate endogenous FVIII
- Bovine Serum Albumin (BSA) has been validated as a practical alternative to FVIII deficient plasma, improving stability and reducing inter laboratory variation from different sources of FVIII deficient plasma²
- Includes an additional inhibitor positive and negative control, to increase confidence in your results



Ordering Information

Product Code	Description	Pack Size
CCIK08	Factor VIII Inhibitor Kit	Contains 5 vial sets, enough to prepare five high titer patient samples. Each vial set comprises: 2 x 1.5ml Imidazole Buffered Pooled Normal Plasma (●) 2 x 1.5ml Imidazole Buffered Bovine Serum Albumin (●) 1 x 0.5ml Negative FVIII Inhibitor Control (●) 1 x 0.5ml Positive FVIII Inhibitor Control (●)

References

1. Douglas CD, et al. A Standardised Kit for a Chromogenic Modified Nijmegen-Bethesda Assay: Repeatability, Reproducibility, and Analytical Sensitivity. Presented at 60th ASH Annual Meeting, December 1-4, 2018 San Diego, California, USA.
2. Miller CH, et al. Validation of Nijmegen-Bethesda assay modifications to allow inhibitor measurement during replacement therapy and facilitate inhibitor surveillance. Journal of Thrombosis and Haemostasis. 2012; 10: 1055-61



T: +44 (0) 2380 483000
E: sales@alphalabs.co.uk
W: www.alphalabs.co.uk