

# cobas® 8100

## automated workflow series

### Input station

Mandatory, maximum 1 station

#### Input Buffer (IPB)

##### Basic specs

- 50, 100 position generic trays
- Either 13 mm or 13 & 16 mm flex-trays
- Optional sample loading with 5 pos. racks and 100 pos. trays
- Loading area for up to 400 samples
- Tube type identification (tube dimensions & cap color)
  - used to cross-check the tube type with the tests/ specimen type
- Throughput 800 samples/hour
- Up to 2 modules per **cobas** 8100 to increase throughput to 1,100 samples/hour

##### Dimension (W x D x H)

1,050 x 1,100 x 1,500 mm

##### Weight

330 kg

##### Heat emission

482 kcal/h  
0,56 kW

##### Power rating

1,0 kVA



#### Automatic Centrifuge Unit (ACU)

##### Basic specs

- 36 tubes per run (9 tubes per bucket)
- Prioritization of STAT samples
- Individual parameter setting per ACU
- Prioritized unloading of STAT samples
- Up to 3,000 g centrifugation force
- Throughput of 300 samples/hour @ 5 min. spin time
- Up to 3 modules per **cobas** 8100

##### Dimension (W x D x H)

900 x 1,100 x 1,500 mm

##### Weight

250 kg

##### Heat emission

1496 kcal/h  
1,74 kW

##### Power rating

2,0 kVA



---

## Sample Check Module (SCM)

---

### Basic specs

- **Sample volume detection / LLD**
  - Accurate determination after separation
  - Detectable through labels
  - Compares the volume in the tube to the required volume
- **Sample quality assessment functionality**
  - Qualitative measurement of Serum Indices without direct access to the sample
  - Window of 6.5mm to see the sample required
  - Up to 1 module per **cobas** 8100

---

<b>Dimension (W x D x H)</b>	300 x 1,100 x 1,500 mm
------------------------------	------------------------

---

<b>Weight</b>	160 kg
---------------	--------

---

<b>Heat emission</b>	260 kcal/h 0,32 kW
----------------------	-----------------------

---

<b>Power rating</b>	0,5 kVA
---------------------	---------

---



---

## Destopper (DSP)

---

### Basic specs

- Decapping of rubber, hemogard, Sarstedt screw caps (including new FBT) and flex caps in any combination
- Throughput: 800 samples/hour
- Up to 1 module per **cobas** 8100

---

<b>Dimension (W x D x H)</b>	300 x 1,100 x 1,500 mm
------------------------------	------------------------

---

<b>Weight</b>	150 kg
---------------	--------

---

<b>Heat emission</b>	241 kcal/h 0,28 kW
----------------------	-----------------------

---

<b>Power rating</b>	0,5 kVA
---------------------	---------

---



---

COBAS is a trademark of Roche.  
© 2023 Roche

### Published by

Roche Diagnostics International Ltd  
CH-6343 Rotkreuz  
Switzerland

[diagnostics.roche.com](https://diagnostics.roche.com)

---

# cobas® 8100

## automated workflow series

Aliquot station

Optional, maximum 1 station

### Barcode Labeler / Tube Feeder (BCL)

#### Basic specs

- Thermal transfer printer
- Compatible with all standard bar codes (Codabar (NW7), interleaved 2 of 5 (ITF), Code 39, Code 128)
- 3,000 aliquot tubes onboard
- Bulk loading of aliquot tubes
  - Optical sensor
- Throughput: 800 samples/hour
- Up to 2 modules per **cobas** 8100 allowing for feeding of two different aliquot tubes (e.g. standard and false bottom tube)

#### Dimension (W x D x H)

600 x 1,100 x 1,500 mm

#### Weight

250 kg

#### Heat emission

464 kcal/h  
0,54 kW

#### Power rating

1,0 kVA



### Aliquot Module (AQM)

#### Basic specs

- Aliquoting into labeled aliquot tubes
- LLD, clot detection
- Aliquot tip volume: 2,000 µL
- Automatic tip feeding (1,000 tips onboard)
- Waste capacity: up to 1,200 tips
- Up to 10 aliquots per sample
- Throughput: 800 samples/hour (max. 1,600 aliquots/hour)
- Up to 1 module per cobas 8100

#### Dimension (W x D x H)

900 x 1,100 x 1,500 mm

#### Weight

250 kg

#### Heat emission

585 kcal/h  
0,68 kW

#### Power rating

1,0 kVA



---

COBAS is a trademark of Roche.  
© 2023 Roche

**Published by**

Roche Diagnostics International Ltd  
CH-6343 Rotkreuz  
Switzerland

[diagnostics.roche.com](https://diagnostics.roche.com)

---

# cobas® 8100

## automated workflow series

### Output station

Mandatory station, maximum 2 stations

#### Restopper Flex-cap (RSF) / Restopper Screw Cap (RSS)

<b>Basic specs</b>	<ul style="list-style-type: none"> <li>Storage capacity at least 2,500 caps               <ul style="list-style-type: none"> <li>– Optical sensor</li> </ul> </li> <li>RSF: flex caps (Sarstedt) for 13 and 16 mm tubes</li> <li>RSS: screw caps</li> <li>Throughput: 600 samples/hour</li> <li>Up to 2 modules per output station</li> </ul>
--------------------	---

<b>Dimension (W x D x H)</b>	300 x 1,100 x 1,500 mm
------------------------------	------------------------

<b>Weight</b>	150 kg
---------------	--------

<b>Heat emission</b>	378 kcal/h 0,44 kW
----------------------	-----------------------

<b>Power rating</b>	0,5 kVA
---------------------	---------



#### Add-on / Output Buffer (AOB)

<b>Basic specs</b>	<ul style="list-style-type: none"> <li>Generic trays with 100 pos.</li> <li>Samples can be retrieved automatically for add-on, rerun, reflex/repeat tests as requested</li> <li>Storage area for up to 1000 samples (10 trays)</li> <li>Manual retrieval of samples possible</li> <li>Multiple buffering targets available per module</li> <li>Throughput: 800 samples/hour</li> <li>Up to 3 modules per output station</li> </ul>
--------------------	--

<b>Dimension (W x D x H)</b>	900 x 1,100 x 1,500 mm
------------------------------	------------------------

<b>Weight</b>	320 kg
---------------	--------

<b>Heat emission</b>	568 kcal/h 0,66 kW
----------------------	-----------------------

<b>Power rating</b>	1,0 kVA
---------------------	---------



---

## Output Buffer / Sorter (OBS)

---

### Basic specs

- Sample sorting into 25, 50, 100 position generic trays
- Either 13 mm or 13 & 16 mm flex-trays
  - Group multiple targets into logical targets
  - Split sorting targets into smaller logical targets
  - Sysmex racks, Hitachi 5 pos. racks
- Sorting area for up to 400 samples
- Default sorting area for tubes with problems  
e.g. BC not readable, no test requests, aliquot error
- Throughput: 800 samples/hour
- Up to 2 modules per output station

---

### Dimension (W x D x H)

600 x 1,100 x 1,500 mm

---

### Weight

230 kg

---

### Heat emission

292 kcal/h  
0,34 kW

---

### Power rating

0,5 kVA

---



---

COBAS is a trademark of Roche.  
© 2023 Roche

### Published by

Roche Diagnostics International Ltd  
CH-6343 Rotkreuz  
Switzerland

[diagnostics.roche.com](https://diagnostics.roche.com)

---

# cobas® 8100

## automated workflow series

### Installation specifications

#### Environmental conditions

**The following environmental specifications should be followed in order to ensure correct operation of the system.**

<b>1</b>	Indoor use only
<b>2</b>	Avoid exposure to dusts and provide adequate ventilation
<b>3</b>	No direct sunlight
<b>4</b>	The floor is level (angle: less than 1/200)
<b>5</b>	The floor is strong enough to hold the weight of the instrument (at least 500 kg/m <sup>2</sup> )
<b>6</b>	Room temperature: 15 to 32°C
<b>7</b>	Any temperature changes should be not greater than ±2°C when using the system
<b>8</b>	Room humidity: 30 to 80 % (non condensing)
<b>9</b>	Ambient air pressure: 106 kPa to 80 kPa
<b>10</b>	No perceptible vibration
<b>11</b>	No significant power supply fluctuation (see table power source requirements for details)
<b>12</b>	No machines discharging ultrahigh frequencies (e.g. electric discharger)
<b>13</b>	No equipment generating electromagnetic waves in near vicinity (e.g. cell phones, transceivers, cordless phones, etc.)
<b>14</b>	A three-wire power supply cord with a safety ground (ground resistance of less than 10Ω)
<b>15</b>	Maximum installation site altitude is 2,000 m
<b>Noise level</b>	Operation: <65 dB (A) Stand-by: <55 dB (A)
<b>Waste</b>	Solid waste in DSP and AQM
<b>Interface</b>	USB-printer, standard Ethernet

---

COBAS is a trademark of Roche.  
© 2023 Roche

**Published by**

Roche Diagnostics International Ltd  
CH-6343 Rotkreuz  
Switzerland

[diagnostics.roche.com](https://diagnostics.roche.com)

---