





FormthoticsTM
Foot Science

Heel Pain: Clinical Insights and Effective Mechanical Management



Heel pain remains one of the most frequent presentations in lower-limb clinical practice. While the condition may appear routine, its underlying pathology, diagnostic complexity, and therapeutic response vary widely – making a thoughtful, evidence-based approach essential.

This overview synthesises current best practice, particularly for practitioners integrating **Formthotics** into their management protocols. As a medically credible, adaptable orthotic system, Formthotics supports rigorous clinical standards while enabling efficient delivery of customised care.





Terminology Matters: Fasciitis, Fasciosis, or Something Else?

Historically referred to as “policeman’s heel” or “heel spur syndrome,” the modern diagnostic landscape has moved towards **plantar fasciitis** – though this term may mislead. There is now broad consensus that chronic heel pain typically involves **degenerative change (fasciosis)** rather than inflammation alone. This shift is critical in guiding treatment.

Diagnostic accuracy directly affects therapeutic outcomes. Failed interventions often reflect diagnostic oversight rather than mechanical mismanagement. Correct identification of the pain generator – be it plantar fasciopathy, fat pad atrophy, or stress reaction – is the prerequisite for success.



Red Flags and Key Questions

Clinical success begins with comprehensive assessment. Ask:

- Duration: >6 months suggests transition to chronicity
- Age: older patients may exhibit different pathology
- Medical history: systemic conditions (e.g., diabetes, inflammatory arthritis)
- Load profile: changes in activity, footwear, work demands
- Previous treatments: what's failed — and why?

Imaging, particularly ultrasound, may reveal **fascia thickening, partial tearing, or fibroma** – all consistent with degenerative disease. In such cases, regenerative modalities (e.g., **ESWT, LLLT, PRP, dry needling**) may be appropriate. Orthoses remain useful for **load redistribution**, but they should be viewed as adjunctive within a broader treatment plan.





The Mechanical Model: Tension and Traction

The prevailing model for mechanical heel pain implicates excessive tension within the plantar fascia at its calcaneal insertion. This tension may result from:

- Intrinsic biomechanical contributors (e.g., equinus, hypermobility, STJ alignment)
- Extrinsic loading factors (e.g., occupational, athletic demands)
- Gait-related stressors (e.g., abrupt footwear transition, training changes)

Formthotics provides a reliable platform for mechanical intervention. Modifications can be implemented based on clinical findings, such as:

- Rearfoot inversion posting for STJ control
- Medial arch fill or first ray cut-out to address foot posture
- Lateral forefoot wedging to reduce fascial tension (supported by emerging data)
- Heel lifts for equinus or Achilles-fascia coupling

The mechanical interface must reflect patient-specific biomechanics – not generic templates. Formthotics enables that level of individualisation.



Considerations for the Older Patient

In older populations, heel pain etiology often shifts:


- Fat pad atrophy increases impact sensitivity
- Plantar calcaneal spurs may become symptomatic due to cushioning loss
- Stress fractures of the calcaneus must be excluded
- Subtalar or inflammatory arthropathies may present atypically

In these cases, cushioning and shock attenuation take precedence over motion control. Models such as Formthotics ShockStop, paired with appropriately designed footwear, can meaningfully reduce loading stress. Orthotic design should reflect both tissue integrity and loading thresholds.



The Role of Footwear: Not Just an Afterthought

Footwear innovation – particularly in **high-stack PEBA** midsoles and **rocker sole geometry** – offers a valuable adjunct to orthotic care. These features reduce sagittal plane loading, assist forward propulsion, and **decrease plantar fascia strain** by limiting MTPJ dorsiflexion. Educating patients about appropriate footwear selection enhances outcomes and compliance.



When Orthotics Underperform



Orthotic “failure” should prompt re-evaluation:

- Has the diagnosis changed or evolved?
- Are mechanical corrections insufficient or misapplied?
- Is the device appropriately matched to patient size/activity?
- Are lifestyle factors (e.g., footwear, hours worn, occupation) impeding success?

Even in cases where regenerative therapies are pursued, Formthotics often retains value – stabilising biomechanics and addressing persistent mechanical loading. As the primary inciting factor may remain, orthotic therapy continues to play a preventative and supportive role post-intervention.





Delivering Tailored Care

– Without Custom Orthotics

Formthotics offers a compelling alternative to full custom devices – enabling precise, **evidence-aligned modifications** based on clinical findings. Rearfoot/forefoot position, sagittal plane mobility, leg length asymmetry, and bodyweight considerations can all be addressed in-practice, without compromising on performance or credibility.

Correct density selection is crucial – a 100kg rugby player will clearly require different material properties than a recreational walker. With a range of densities and configurations, Formthotics supports **scalable personalisation** at the point of care.



Closing Thoughts



Heel pain is multifactorial – and sometimes persistent – but rarely unsolvable. With accurate assessment, thoughtful biomechanical planning, and appropriate device selection, most cases respond well to conservative care.

As a trusted tool in your clinical arsenal, Formthotics reinforces your ability to deliver precise, patient-centred care – supporting both your standards and your outcomes.



References

1. **Whittaker et al.** Foot orthoses for plantar heel pain: a systematic review and meta-analysis. British Journal of Sports Medicine. March 2018
2. **Rasenberg, et al.** Efficacy of foot orthoses for the treatment of plantar heel pain: a systematic review and meta-analysis. British Journal of Sports Medicine. Vol 52, 2018
3. **Schuitema et al.** Effectiveness of Mechanical Treatment for Plantar Fasciitis: A Systematic Review. Journal of Sport Rehabilitation. October 2019.
4. **Caglar, O., Aydin, A.** Comparison of extracorporeal shock wave therapy with custom foot orthotics in plantar fasciitis treatment: A prospective randomized one-year follow-up study. Journal of Musculoskeletal Neuronal Interaction. 2019 Jun 1;19(2):178-186.
5. **Heide et al.** Is radial extracorporeal shock wave therapy (rESWT), sham-rESWT or a standardised exercise programme in combination with advice plus customised foot orthoses more effective than advice plus customised foot orthoses alone in the treatment of plantar fasciopathy? A double-blind, randomised, sham-controlled trial. British Journal of Sports Medicine. July 2024
6. **Pfeffer et al.** Comparison of Custom and Prefabricated Foot Orthoses in the Initial Treatment of Plantar Fasciitis. Foot and Ankle International. April, 1999.
7. **Landorf, A M Keenan.** Effectiveness of Foot Orthoses to Treat Plantar Fasciitis-A Randomized Trial. Archives of Internal Medicine. June, 2006.
8. **Whittaker et al.** Foot orthoses for plantar heel pain: a systematic review and meta-analysis. British Journal of Sports Medicine. March 2018.
9. **Landorf, A M Keenan.** Efficacy of foot orthoses: What does the literature tell us? Journal of the American Podiatric Medical Association. Vol 32, 1998.
10. **Seligman DAR, Dawson D, Streiner DL, Seligman DJ, Davis A.** Treating Heel Pain in Adults: A Randomized Controlled Trial of Hard vs Modified Soft Custom Orthotics and Heel Pads. Arch Phys Med Rehabilitation. 2021 Mar;102(3):363-370.



26 Dakota Crescent, Wigram,
Christchurch, 8042, New Zealand

www.formthotics.com