

# Rapid whole-genome sequencing, *de novo* assembly, and characterisation of bacterial isolates

The nanopore-only microbial isolate sequencing solution (NO-MISS), combined with the EPI2ME™ wf-bacterial-genomes workflow, enables rapid sequencing and automated analysis of bacterial isolates, including efficient *de novo* assembly, identification, sequence typing, and predicted antimicrobial resistance profiling

More information can be found at: [epi2me.nanoporetech.com/epi2me-docs/workflows/wf-bacterial-genomes](https://epi2me.nanoporetech.com/epi2me-docs/workflows/wf-bacterial-genomes)




Fig. 1 Overview of the end-to-end process for characterisation of bacterial isolates.

## Bacterial isolates characterisation

Whole-genome sequencing and characterisation of microbial isolates is a valuable approach across many settings, including in public health, food safety, biopharma, and research sectors. In particular, predicting the antimicrobial resistance status or other phenotypic properties of strains is important, especially during outbreaks. Oxford Nanopore sequencing overcomes limitations that exist for short-read bacterial sequencing, and enables high-quality, contiguous genome assemblies that resolve complex genomic structures, including repetitive elements and plasmids. Here we demonstrate the nanopore-only microbial isolate sequencing solution (NO-MISS) approach, combined with analysis using the wf-bacterial-genomes workflow, on 32 foodborne pathogens from an FDA proficiency panel.



Fig. 3 a) The EPI2ME platform was used to run the wf-bacterial-genomes workflow and generate a report. The report includes b) quality metrics, and c) *de novo* assembly statistics with further characterisation results.

## Analysis with wf-bacterial-genomes via EPI2ME

wf-bacterial-genomes can be run either on the command line or via the EPI2ME Desktop Application, which provides a simple graphical interface for setting parameters and running the workflow. It accepts either FASTQ or BAM files as input. In this example, FASTQ files from MinKNOW were used with default parameters, and with additional isolate characterisation analysis selected. The report generated includes quality control metrics and results for each sample, indicating that the sequencing was completed successfully and with high accuracy and coverage.




Fig. 2 a) PCR-free library preparation with the Rapid Barcoding Kit. A transposase is used to fragment and attach barcodes to DNA before adding a sequencing adapter. b) Sequencing is carried out on a benchtop GridION™ device using two MinION™ Flow Cells, each run for 72 hours.

## Library preparation and sequencing

Following DNA extraction with Maxwell RSC Cultured Cells DNA Kit, libraries were prepared using the Rapid Barcoding Kit (SQK-RBK114.24). NO-MISS allows for up to 24 samples to be multiplexed together, so the 32 samples were split into two sets of 16 samples and loaded onto two MinION Flow Cells. They were sequenced for 72 hours to reach the minimum target coverage of 50x. Basecalling was performed using high accuracy (HAC) mode on MinkNOW™.




Fig. 4 Multilocus sequence typing (MLST), serotyping, and antimicrobial resistance prediction report.

## Further analyses and antimicrobial resistance characterisation

An additional report was generated for each sample showing further analyses. The example in Fig. 4 was identified by MLST as belonging to the senterica\_achtman\_2 scheme and therefore further characterised by SeqSero2, which identified the isolate as belonging to the Heidelberg serotype. Additionally, ResFinder was used to identify genes and SNVs associated with antimicrobial resistance, and in this case, identified eight features associated with resistance. These results demonstrate how NO-MISS can be combined with wf-bacterial-genomes to enable efficient, high-resolution characterisation of bacterial isolates.