
Flowpression Boots Ultra Red Light

FLOWLIFE

Brochure

Flowpression Boots Ultra Red Light

Infrared Compression Boots

Table of Content

4 —	Manifest
6 —	Flowpression Boots Ultra Red Light
7 —	Accessories
8 —	Flowpression Boots Ultra Red Light: A Comprehensive Guide
10 —	Starting your Journey
16 —	The Four Modes Flow · Pulse · Progressive · Zone
26 —	For Muscle Soreness
30 —	For General Recovery
34 —	Before, During, and After your Workout
39 —	Contact Information
40 —	Manual
42 —	• English • Svenska • Dansk • Norsk • Suomi • Deutsch • Français • Nederlands • Español • Italiano

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.

For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowpression Boots Ultra Red Light?

Revolutionize leg recovery with boots that compress, flush, and revive your muscles — and go anywhere you do.

Advanced Gradient Compression Technology – Flowpression Boots Ultra Red Light applies precisely calibrated sequential pressure across five independent chambers from foot to hip, simulating the natural action of muscles pumping blood and lymph toward the heart. With 15 adjustable pressure levels between 80 and 220 mmHg, each session is tailored to the exact stimulus your body needs.


Four Precision Recovery Modes – Not all recovery is the same. Flowpression Boots Ultra Red Light offers four distinct compression programmes — Flow, Pulse, Progressive, and Zone — each engineered for a specific moment in your training or daily recovery cycle. Select the right mode for the right moment, and let the boots do the work.

Integrated Infrared Light Therapy — Flowpression Boots Ultra Red Light combines pneumatic compression with infrared light therapy across 108 LED chips per boot — 18 operating at 660 nm (red light) and 90 at 850 nm (near-infrared). Together, these two scientifically established wavelengths support cellular repair, reduce local inflammation, and enhance circulation alongside any compression mode.

Fully Cordless with Removable Batteries — Each boot operates independently with a removable 2600 mAh Li-ion battery and USB-C charging. No hub, no cables, no external controller. Designed for unrestricted use at home, at the track, or while travelling.


Accessories

Flowpression Boots Ultra Red Light comes with everything you need for a complete, cable-free compression therapy session — whether at home or on the move.


Flowpression Boots Ultra Red Light Compression Boots

Two integrated boot sleeves with five independent compression chambers each, 108 built-in LED chips per boot (18 × 660 nm red light + 90 × 850 nm near-infrared), and an onboard LCD control unit. Available in two sizes: S/M (860 mm) and L/XL (938 mm).


Rechargeable Li-ion Battery Packs (x2)

One removable battery per boot. Each 2600 mAh pack snaps into the control unit with a click and charges via USB-C. Provides up to 150 minutes of runtime per charge with compression only, or approximately 80 minutes when the infrared light therapy function is active.


USB-C Charging Cables (x2)

One cable per boot. Compatible with any 5V / 2A USB-A power source, power bank, or phone charger. Full charge in approximately 3 hours.


Storage and Carry Bag

Structured carry bag designed to hold both boots, cables, and battery packs. Compact enough for a training bag or suitcase.

Flowpression Boots Ultra Red Light: A Comprehensive Guide

This essential guide is designed to enhance leg health and maximise the recovery benefits and performance of your Flowpression Boots Ultra Red Light by Flowlife. The boots feature advanced sequential air compression technology with integrated infrared light therapy — 660 nm red light spectrum and 850 nm near-infrared light spectrum — 15 adjustable pressure levels, and four precision recovery modes, each built around a specific phase of training and daily recovery.

Starting your
Journey

Prepare Yourself

Ensure your legs are clean and dry before starting a session with Flowpression Boots Ultra Red Light. This essential step guarantees a hygienic recovery experience and prolongs the product's lifetime. Wearing slim or no clothing under the boots will ensure full contact pressure.

First Charge

Before your first session, remove the battery packs and charge them fully via the included USB-C cables. Allow approximately 3 hours for a complete charge. Do not use the device without a charged battery installed.

Installing the Batteries

Align the bottom slide of the battery pack with the base slot on the control unit. Gently push the battery inside until it is flush with the control unit. When you hear a click, the battery is installed in place. Ensure it is installed correctly for the device to work properly.

Wearing the Boots

The boots are not marked left or right. To identify orientation, position each boot so the control unit and battery pack sit on the outside of the leg — this tells you which boot belongs on which leg. Pull the zipper fully open from foot to hip, place your leg inside, and zip completely closed before powering on. Always ensure the zipper is fully zipped before use.

Starting Your Session

Power On:

Press the power button briefly when the device is off to turn on the power. The display screen will light up, indicating the device has started up.

Select Your Mode:

When the device is in standby mode, press the M/IR button briefly to select the mode. Each press changes the mode once, cycling between the four programmes.

Selecting Your Compression Mode

Choose from four tailored programmes to address your legs' specific needs:

Flow

Classic sequential compression that builds progressively from foot to hip. The foundational compression therapy mode. Primary recommendation for post-training recovery and muscle soreness relief.

Pulse

Wide two-chamber pulsating compression covering large sections of the leg at once. Broad, even pressure — ideal for pre-training warm-up and comfortable massage-style use without high pinpoint intensity.

Progressive

Slow ascending mode that completes four full cycles through the lower leg before reaching the thigh. More focus on feet and calves, gentler on the thigh — ideal for sensitive users or lower-leg targeted sessions.

Zone

Single-chamber sustained hold, one section at a time. Works through the leg like a methodical massage — one zone fully treated before moving to the next.

Adjusting Compression Intensity

Use the Intensity Up and Intensity Down buttons to adjust compression across 15 levels between 80 and 220 mmHg. Intensity can be adjusted while the device is running.

Low Intensity (80–110 mmHg)

Begin with a low-intensity setting, which provides gentle pressure. Ideal for first-time users, sensitive legs, or a subtle daily recovery session.

Medium Intensity (115–160 mmHg)

The standard working range — balanced pressure suitable for regular use, post-training recovery, and moderate soreness relief.

High Intensity (165–220 mmHg)

For significant muscle soreness or deep tissue work. Delivers deep and sustained pressure for significant tension relief in stiff or heavily worked areas. Increase gradually and stop if discomfort arises.

Setting the Session Duration

Press the Time button to set session length between 15 and 60 minutes in 5-minute steps.

15 Minutes

Perfect for a quick recovery session or warm-up.

30 Minutes

Ideal for a standard recovery session offering thorough muscle relaxation.

60 Minutes

Best for an intensive recovery session, providing extended relief and rejuvenation.

Infrared Light Therapy

During any active session, press the M/IR button briefly to activate the integrated infrared light therapy function. Each boot contains 108 LED chips — 18 emitting at 660 nm (red light) and 90 emitting at 850 nm (near-infrared) — working simultaneously across the full sleeve. The corresponding icon on the LCD screen will illuminate. Press M/IR again to deactivate. The infrared function can be activated independently of the current compression mode and is available across all four modes.

Syncing Both Boots

The boots are factory-paired and connect automatically on startup. If they do not pair, press and hold the M/IR button on both units simultaneously for more than 5 seconds until you hear a beep. Three beeps confirm a successful pairing. The sync icon will illuminate on both screens.

Battery Indicators

Low Battery: The indicator will alert you to recharge when the battery is running low. When the first indicator light flashes, the battery has dropped below 10% and must be charged before the next session.


Charging: When the battery is charging, the battery indicator lights will pulse in sequence.

Full Charge: All battery indicator lights will stop pulsing and remain fully illuminated once charging is complete.


Charging

Connect the output end of the USB-C cable to the charging port on the battery pack, then connect the input end to a 5V / 2A power source — a power bank, phone charger, or wall adapter. Charging time is approximately 3 hours. You can charge the battery at any level at any time. It is recommended to charge every 3 months if the product is not in regular use.


Concluding Your Session


1. Power Off: Press and hold the power button for 2 seconds to begin the automatic air evacuation cycle. The boots will deflate fully over approximately 150 seconds before powering off automatically.

2. Skip Evacuation: If you need to power off immediately, press and hold the power button for 2 seconds again during the evacuation cycle to shut down immediately.


3. Remove the Boots: Unzip and remove the sleeves gently.

4. Store Properly: Deflate fully, then fold and place into the carry bag. Follow the illustrated folding guide.

Emergency Pressure Relief

If you experience any discomfort or excessive pressure during a session, push the Emergency Pressure Relief slider on the control unit to the ON position. This instantly releases all pressure in the sleeves. After deflating, push the slider back to the OFF position for normal operation.

The Four Modes

Flow · Pulse · Progressive · Zone

Each mode is engineered around a specific phase of training, recovery, and daily wellbeing. Select the right programme for the right moment, and let the boots do the work.

Flow

Classic

Compression Therapy

Flow is the most fundamental expression of what compression boot therapy was originally designed to do. It applies progressive sequential pressure from the foot upward, systematically squeezing blood and lymphatic fluid out of the lower limbs and back toward the heart and core — the same mechanism the body relies on naturally, amplified and made consistent.

Every time the chambers build to their full-leg peak, the boots are doing exactly what decades of compression therapy research was built on: pushing deoxygenated blood and metabolic waste products up and out, creating space for fresh, oxygen-rich blood to flow back in. Flow is effective both after training and after any long day where blood has been pooling in the lower legs against gravity.

When to Use


After training, Flow is the single most effective tool for your legs in the first hour. The moment your workout ends, blood has been pushed into the periphery — pooling in the feet, calves, and thighs under the effects of exertion and gravity.

Flow's cumulative chamber build actively reverses this, squeezing blood from the foot upward in a sustained five-stage wave and driving it back to the core where the body can process and recirculate it. This is the mechanism that initiates the recovery cascade — clearing metabolic debris, restoring fresh oxygenated blood, and reducing the delayed-onset muscle soreness that sets in over the following hours.

The earlier you use it after training, the more work it saves your body from doing on its own. Flow is equally effective outside of training — use it after long days of standing, walking, or travel whenever the lower limbs need a full circulatory reset. It is the most versatile of the four modes and the natural first choice whenever you reach for the boots.

Recommended Settings

Duration: 20–40 minutes. Pressure: 90–150 mmHg. Post-training, aim for the upper half of the pressure range to maximise circulatory stimulus. For daily fatigue, the mid-range is sufficient and comfortable for longer sessions.

Pulse

Wide Coverage Compression

Pulse works across large sections of the leg at once. Two chambers inflate simultaneously in a pulsating rhythm, creating a broad, even pressure that covers a wide area rather than focusing on a single concentrated point. The result is a full, generous compression that feels more like a wide-handed massage than a targeted press — comfortable at any intensity level.


Because the pressure is distributed across two overlapping zones at once, Pulse delivers a strong stimulus without the intensity of single-chamber modes. This makes it one of the most versatile modes in the range — effective before training, between sessions, or whenever the legs need movement without high-pressure point work.

When to Use

Use Pulse as a pre-training warm-up to activate circulation and raise tissue temperature before physical output. Use it between workout efforts to keep the legs fresh. It is also the natural first choice for users who prefer a broad, comfortable massage sensation over deep concentrated pressure — effective without ever needing to push the intensity high.

Recommended Settings

Duration: 15–30 minutes for warm-up. 20–40 minutes for recovery or massage use. Pressure: 90–150 mmHg. The wide coverage means you feel the compression clearly without needing to increase intensity — start in the mid-range and adjust to comfort.

Progressive

Lower Leg to Full Body

Progressive is a slow, deliberate ascending mode that starts at the foot and works its way upward — completing four full cycles through the lower leg before reaching the thigh. This means the feet, ankles, and calves receive sustained, repeated compression long before the session ever moves to the upper leg.


This makes Progressive the right choice for anyone who is sensitive around the thighs or who simply wants more time and attention on the lower legs. Where Flow builds across the full leg in a single pass, Progressive stays low, repeats, and earns its way upward — making it ideal for users who find full-leg modes too intense at the top, or who carry their tension primarily from the knee down.

When to Use

Use Progressive after long days of standing, walking, or travel when the feet and calves need the most attention. Use it as an alternative to Flow for users who are sensitive around the thigh area — the mode delivers all the lower-leg benefit of a full session before it ever reaches the upper leg. It is also well suited as an evening recovery session or a pre-sleep wind-down.

Recommended Settings

Duration: 30–60 minutes. Pressure: 80–140 mmHg. Progressive works well at moderate pressure — the repeated lower-leg cycles deliver significant stimulus without needing high intensity. Start low and increase gradually.

Zone

Section by Section

Zone works through the leg one section at a time. Each of the five chambers inflates to full pressure, holds, deflates completely, and then the sequence moves to the next zone — from foot to hip. There is no overlap, no travelling wave, no simultaneous action. Just one section of the leg receiving the full selected pressure, one at a time.


This is the mode that most closely resembles a hands-on massage. The sustained single-zone hold gives each area of the leg extended time under pressure — longer than any other mode provides — making Zone the right choice when you want to work through the leg methodically and feel each section separately.

When to Use

Use Zone on rest days, on low-activity days, or in the morning after prolonged sitting or sleeping. Use it when the legs feel stiff, sluggish, or unresponsive to faster modes. It is also the recommended mode for users who prefer a massage-like experience — deliberate, unhurried, and focused.

Recommended Settings

Duration: 30–60 minutes. Pressure: 130–220 mmHg. Zone benefits from higher pressure — the single-zone hold is most effective when the pressure is sufficient to create a clear compressive sensation in each section. Increase gradually and find the level that feels like firm, comfortable pressure on each zone.

For Muscle
Soreness

Muscle soreness — whether from a training session that pushed you to the limit, a competition, or simply a day that demanded more than expected — responds directly to targeted compression therapy. Here is how to get the most from Flowpression Boots Ultra Red Light when your muscles need it most.

Select Mode and Intensity

Flow Mode

For overall leg soreness after training or competition, Flow is the primary recommendation. The classic sequential compression builds progressively from foot to hip, systematically driving metabolic waste and lactic acid out of the tissue and pulling fresh, oxygenated blood back in — directly addressing the root cause of post-exercise soreness. For best results, use within one to two hours of finishing your session.

Zone Mode

If soreness is deep, localised, or recurring in a specific area, Zone delivers the concentrated pressure needed to address it. The single-chamber sustained hold penetrates more deeply than travelling wave modes, making it the right choice for significant or chronic tension in specific zones.

Adjust Compression Intensity

Low Intensity

(80–110 mmHg)

Start here if soreness is mild, if using immediately after very intense activity, or if using Flowpression Boots Ultra Red Light for the first time on sore muscles.

Medium Intensity

(115–160 mmHg)

The standard working range for moderate post-training soreness. Provides sufficient stimulus to accelerate recovery without discomfort.

High Intensity

(165–220 mmHg)

For significant muscle soreness or deep tension in stiff areas. Delivers deep and sustained pressure for thorough relief. Increase gradually and stop if discomfort arises.

Set the Optimal Duration

15 Minutes

Quick relief for mild soreness or a targeted flush between sessions.

30 Minutes

Thorough treatment for moderate soreness after a hard training day.

60 Minutes

Best for intensive recovery — significant soreness, accumulated fatigue, or multiple days of consecutive training.

For General
Recovery

LOWLIFE

Flowpression Boots Ultra Red Light provides a comprehensive recovery solution for legs that have endured the demands of training, competition, or daily activity. Regular sessions deliver benefits that extend well beyond immediate relief — supporting the long-term health and performance of your legs.

Enhanced Circulation

Flowpression Boots Ultra Red Light enhances circulation through advanced sequential air compression technology. The boots actively drive blood and lymphatic fluid from the lower limbs back toward the heart, increasing the delivery of oxygen and nutrients to fatigued muscle tissue while simultaneously removing metabolic waste. By boosting blood flow in this way, the boots facilitate faster recovery and support overall leg health.

Muscle Relaxation

The four compression programmes effectively alleviate muscle tension across the full length of the leg — from the feet to the thighs — promoting relaxation and reducing the risk of overuse injuries. Regular sessions break the cycle of accumulated tightness, ease tension from repetitive activity or prolonged sitting, and help the legs maintain their natural range of motion.

Flexibility and Tissue Health

Regular use of Flowpression Boots Ultra Red Light is associated with maintained and improved leg flexibility. By preventing the accumulation of excess fluid, keeping circulation active between training sessions, and reducing chronic tightness in the lower limbs, consistent compression therapy supports the flexibility and long-term resilience of the muscles and connective tissue in the legs — keeping them in optimal condition for whatever comes next.

Infrared Light Therapy

The boots integrate 108 LED chips per sleeve — 18 emitting at 660 nm (visible red light) and 90 at 850 nm (near-infrared) — distributed across the full length of the boot to deliver light energy directly through the skin and into the underlying tissue during any compression session.

660 nm — Red Light

Red light at 660 nm penetrates the superficial layers of skin and tissue. At this wavelength, light energy is absorbed by cells near the surface, where it supports cellular metabolism, promotes local microcirculation, and helps reduce surface-level inflammation. You can see the red glow when the function is active.

850 nm — Near-Infrared


Near-infrared at 850 nm is invisible to the eye but penetrates significantly deeper — reaching muscle tissue, tendons, and joints. At this depth, the light energy supports mitochondrial activity, promotes tissue repair, and helps reduce deeper inflammation and recovery time after physical load.

Calibrated Irradiance for Longer Sessions

The irradiance is set at 75 W/m^2 — intentionally calibrated for extended wear rather than peak intensity. High-irradiance devices deliver a large light dose quickly but are only used for short sessions of a few minutes. At 75 W/m^2 , the boots are designed to be worn for a full 30–40 minute compression session, delivering a therapeutic joule dosage accumulated gradually over time. The result is a meaningful dose of light energy delivered comfortably alongside the compression therapy — without the need to stop, adjust, or switch between devices.

Combined Effect

Using infrared alongside compression creates a compounding benefit: compression drives blood and lymphatic fluid through the tissue while infrared light simultaneously supports cellular repair and reduces inflammation in that same tissue. The two therapies are fully compatible and can be used together across all four compression modes.

A black and white photograph of a person's legs and feet. The person is standing on a light-colored surface, with their right leg in the foreground and their left leg slightly behind it. The skin is tanned and shows some hair. The muscles of the legs are well-defined, particularly the quadriceps and calves. The feet are bare, and the toes are visible. The background is a plain, light-colored wall.

Before, During,
and After your
Workout

Flowpression Boots Ultra Red Light is designed to support your training at every stage — not just recovery. Here is how to integrate the boots into your exercise routine for maximum benefit at each phase.

Before Your Workout Warm-Up Enhancement

Using Flowpression Boots Ultra Red Light before training activates circulation in the lower limbs and increases tissue temperature ahead of physical load. The gentle but progressive compression increases blood flow and oxygen delivery, which helps reduce the risk of injury and prepares the legs for activity.

How to Use:

1. Put on the boots and zip both sleeves fully closed.
2. Select Flow mode — the foundational sequential compression that builds blood flow progressively from foot to hip, priming the working muscles before load.
3. Set a low to medium intensity (90–135 mmHg).
4. Use for 15–20 minutes to effectively warm up and prepare the legs.

During Your Workout Active Recovery Between Efforts

Incorporating Flowpression Boots Ultra Red Light during breaks in your workout aids active recovery between efforts. It reduces muscle fatigue and keeps the legs fresh throughout longer sessions — allowing you to perform at your best across multiple efforts or rounds.

How to Use:

5. Put on the boots during breaks of 10–15 minutes.
6. Select Flow for a fast circulatory flush, or Progressive if time allows for a more thorough ascending treatment.
7. Use a low to medium intensity (90–135 mmHg) to keep the legs active without adding fatigue.

After Your Workout

Post-training is where Flowpression Boots Ultra Red Light delivers its most significant benefit. Using the boots after exercise drives metabolic waste and lactic acid out of the muscle tissue and accelerates the return of fresh, oxygenated blood — the essential first step in recovery.

How to Use:

8. After completing your workout, sit or lie comfortably and put on the boots.
9. Select Flow mode — the most effective post-training programme for driving metabolic waste out of the legs and bringing fresh blood back in.
10. Set medium to high intensity (130–180 mmHg) to maximise the circulatory stimulus.
11. Use for 20–45 minutes to reduce muscle stiffness and promote blood flow.

Benefits

- **Alleviates Muscle Fatigue** - Eases discomfort and reduces sensations of stiffness and soreness — both immediately after training and in the hours and days that follow.
- **Promotes Circulation** - Enhances blood flow, flushes metabolic waste from the muscles, and brings essential oxygen and nutrients to recovering tissue.
- **Supports Flexibility** - Regular sessions help maintain and improve leg flexibility and tissue pliability, keeping you ready for the next session.


Manufacturer:
Flowlife AB
Runebergsgatan 7
114 30 Stockholm
Sweden

Contact:
hello@flowlife.com
www.flowlife.com

Manual

Flowpression Boots Ultra Red Light

ENGLISH

INTRODUCTION

Welcome to your Flowpression Boots Ultra Red Light. This guide will help you use the boots safely and effectively. Read all instructions before first use and retain this manual for future reference.

In this document, Flowpression Boots Ultra Red Light is referred to as the product. The product must be used in accordance with the following safety instructions to reduce the risk of fire, electric shock, property damage, and personal injury.

Flowpression Boots Ultra Red Light is a personal-use sequential pneumatic compression device with integrated infrared light therapy, intended for use on the lower limbs. It is designed to support general wellness, circulation, and recovery in healthy adult users. It is not a medical device and is not intended to diagnose, treat, cure, or prevent any disease or medical condition.

SAFETY INSTRUCTIONS

SAFETY HAZARDS — REDUCE THE RISK OF ELECTRIC SHOCK

- Do not use the product near water or in wet environments.
- Never immerse the boots, control unit, or battery packs in water or any other liquid.
- Use only the original Flowlife USB-C cable and approved battery packs supplied with the product.
- Disconnect the charging cable if the control unit or battery pack appears damaged or wet.
- Do not handle the charging port or battery contacts with wet hands.
- Do not modify, repair, or attempt to open the control unit or battery packs.

DANGER — REDUCE THE RISK OF INJURY

- Do not walk or stand while wearing the boots. Use only while seated or lying down.
- Always zip the sleeves fully before powering on. Never unzip the sleeves during an active session.
- Stop use immediately if you experience pain, numbness, tingling, or any unusual sensation.
- Do not puncture or otherwise damage the boot sleeves — this may cause the device to operate incorrectly.
- Do not insert foreign objects into any openings on the boots or control unit.
- The boots are designed for single-person use only. Do not share between users.
- Use by unconscious or incapacitated persons may be dangerous without supervision.
- Keep out of reach of children under 13. Users aged 13–15 must be supervised by an adult.
- Ensure the Power button is always accessible during use so the device can be switched off if needed.
- If the device emits an unusual smell or appears to overheat, switch it off and contact Flowlife support.

SKIN PRECAUTIONS

- Do not use on broken skin, open wounds, active infections, or areas of recent burns or injury.
- Do not use if you are experiencing inflammation, pain of unknown origin, or bleeding at or near the site of application.
- Ensure legs are clean and dry before each session. Remove any topical creams or oils.
- Start at lower intensity and increase gradually. Stop use if redness, bruising, or persistent discomfort occurs.
- Do not use on areas of reduced skin sensation — you may not detect excessive pressure.
- Consult your doctor before use on sensitive skin.

CONSULT A DOCTOR BEFORE USE IF YOU HAVE

- Suspected or confirmed deep vein thrombosis (DVT) or thrombophlebitis
- Episodes of pulmonary embolism or acute pulmonary oedema
- Acute congestive cardiac failure, heart disease, or labile hypertension
- Acute infections of the limbs
- Wounds, lesions, or tumours at or near the site of application
- Bone fractures or dislocations at or near the site of application
- Any condition where increased venous and lymphatic return is undesirable
- Pacemaker, ICD, or any active implanted electronic device
- Pregnancy
- Active cancer or ongoing cancer treatment in the lower limbs
- Recent surgery to the lower limbs (within 6 months)
- Any other condition — when in doubt, consult your doctor before starting

INSTRUCTIONS FOR USE

WHAT IS IN THE BOX

- 2 × Flowpression Boots Ultra Red Light
- 2 × USB-C charging cables
- 1 × User manual
- 2 × Removable Li-ion battery packs (2600 mAh each)
- 1 × Storage and carry bag

BEFORE YOUR SESSION

- Read this manual before first use.
- Charge both battery packs fully via the included USB-C cables. Allow approximately 3 hours for a full charge.
- Ensure legs are clean and dry. Remove thick socks or bulky clothing under the boots.
- Inspect the boots, control units, and cables for any damage before each use. Do not use if damaged.

INSTALLING THE BATTERIES

- Align the base of the battery pack with the slot on the underside of the control unit.
- Slide the battery in firmly until you hear a click. Ensure it sits flush.
- Repeat for the second boot.

WEARING THE BOOTS

- The boots are not marked left or right. To determine orientation, position each boot so the control unit and battery pack sit on the outside of the leg — this is how you identify which boot goes on which leg.
- Pull the zip fully open from foot to hip before putting on.
- Place your leg inside and zip completely closed before powering on.
- Always ensure the zip is fully closed before starting a session.

BUTTONS AND CONTROLS

Power button Press briefly to power on. Press and hold for 2 seconds to begin the air evacuation and shutdown cycle.

M/IR button Press briefly to cycle through the four compression modes. During an active session, press briefly to toggle the infrared light therapy function on or off.

Intensity Up button Increases compression intensity by one level. Can be adjusted during an active session.

Intensity Down button Decreases compression intensity by one level. Can be adjusted during an active session.

Time button Sets session duration between 15 and 60 minutes in 5-minute steps.

Chamber selection Hold the Intensity Up button to remove active chambers starting from the thigh downwards, allowing you to focus compression on the lower leg or any specific zone.

COMPRESSION MODES

M1 — Flow Cumulative sequential compression building from foot to hip. Classic compression therapy.

M2 — Pulse Pulsating dual-chamber wave from foot to hip. The rhythmic wave action creates a relaxing massage effect that promotes circulation and eases tension throughout the leg.

M3 — Progressive Five-stage ascending cycle from lower leg to full body. Lower leg focus and travel recovery.

M4 — Zone Single-chamber sustained hold, one zone at a time. Deep targeted compression.

ADJUSTING COMPRESSION INTENSITY

- Use the Intensity Up and Intensity Down buttons to select from 15 levels between 80 and 220 mmHg.
- Intensity can be adjusted at any time during an active session.
- New users should begin at low intensity (80–110 mmHg) and increase gradually over several sessions.
- Stop use and reduce intensity immediately if any pain, numbness, or discomfort arises.

INFRARED LIGHT THERAPY

- During any active session, press the M/IR button briefly to activate infrared light therapy.
- The corresponding icon on the LCD screen will illuminate to confirm the function is active.
- Press M/IR again to deactivate. Infrared can be used independently of any compression mode.
- Each boot contains 108 LED chips: 18 × 660 nm red light + 90 × 850 nm near-infrared.

SYNCING BOTH BOOTS

- The boots are factory-paired and connect automatically on startup.
- If the boots do not sync, press and hold the M/IR button on both units simultaneously for more than 5 seconds.
- Three beeps confirm a successful pairing. The sync icon will illuminate on both screens.

DURING YOUR SESSION

- Remain seated or lying down for the duration of the session. Do not walk or stand.
- Intensity and mode can be adjusted at any time during the session.
- If you experience any discomfort or excessive pressure, use the Emergency Pressure Relief slider on the control unit to release all pressure instantly.
- After deflating, push the slider back to the OFF position to resume normal operation.

CONCLUDING YOUR SESSION

- Press and hold the Power button for 2 seconds. The boots will enter the automatic air evacuation cycle and deflate fully over approximately 150 seconds before powering off.
- To force an immediate shutdown during evacuation, press and hold the Power button again for 2 seconds.
- Remove the boots by unzipping fully and lifting out each leg carefully.
- Check the skin for any redness, marks, or irritation after each session.

SESSION LENGTH AND FREQUENCY

- Minimum session: 15 minutes
- Recommended frequency: 3–6 sessions per week
- Maximum session: 60 minutes

Daily use is safe and well-tolerated at moderate intensity settings. Do not exceed 60 minutes per session. Allow at least 60 minutes between sessions on the same area.

CLEANING & STORAGE

- Power off and remove the battery pack before cleaning.
- Wipe the outer surface of the sleeves with a soft, dry cloth after each session.
- For deeper cleaning, use a slightly damp cloth and allow to dry fully before storage.
- Do not immerse the boots, control unit, or battery pack in water or any other liquid.
- Do not use alcohol, solvents, bleach, or abrasive cleaners on any component.
- Do not machine wash the boots.
- Allow the boots to deflate fully before folding and storing.
- Fold and store as follows: deflate and lay flat, fold the two ends toward each other, fold in half again, then place inside the carry bag.
- Store in a cool, dry place at room temperature, away from direct sunlight and heat sources.
- Do not store at temperatures below -10°C or above 45°C .
- It is recommended to charge the battery packs every 3 months during extended periods of non-use.
- Inspect cables and battery contacts before each use. Stop use if any damage is observed.

RECYCLING & DISPOSAL

- Do not dispose of this product as unsorted household waste.
- Flowpression Boots Ultra Red Light contains electronic components and Li-ion batteries subject to separate collection under WEEE Directive 2012/19/EU and applicable local legislation.
- Return the product to a designated collection point for electrical and electronic equipment recycling.
- Remove battery packs before disposal and take them to an approved battery collection point.
- Follow all local regulations for the disposal of electronic and electrical equipment.

TECHNICAL SPECIFICATIONS

Pressure range: 80–220 mmHg (15 levels)

Compression chambers: 5 per boot

Operation modes: 4 — Flow (M1), Pulse (M2), Progressive (M3), Zone (M4)

Session timer: 15–60 minutes (5-minute steps)

Battery capacity: 2 × 2600 mAh Li-ion (removable, one per boot)

Charge input: 5V / 2A via USB-C

Charge time: Approximately 3 hours

Runtime: Up to 150 minutes (compression only) · Approx. 80 minutes (infrared light therapy active)

Infrared light therapy: 108 LED chips per boot — 18 × 660 nm red light + 90 × 850 nm near-infrared

Irradiance: 75 W/m²

Rated power: 25 W

Boot sizes: S/M — 860 mm boot length | L/XL — 938 mm boot length

Product weight: S/M: 3.2 kg | L/XL: 3.5 kg (batteries included)

Operating temperature: 10°C – 40°C

Storage temperature: -10°C – 45°C

Certifications: FCC / CE / ISO 9001 / ISO 13485 / RoHS

Manufacturer: Flowlife AB, Runebergsgatan 7, 114 29 Stockholm, Sweden.

hello@flowlife.com — www.flowlife.com

WARRANTY

Flowpression Boots Ultra Red Light is covered by a 2-year warranty from the date of purchase, covering material and manufacturing defects under normal conditions of use as described in this manual. The warranty applies to purchases made directly from Flowlife. Purchases made from unauthorised third-party retailers are not covered by this warranty.

The warranty does not cover damage caused by misuse, improper storage, water exposure, dropping or impact, unauthorised repair or modification, use of non-approved accessories, or use outside the instructions in this manual.

Proof of purchase is required for all warranty requests. For support, contact hello@flowlife.com.

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Produttore: Flowlife AB, Runebergsgatan 7, 114 29 Stoccolma, Svezia.