Flowcontrast Spot Go

FLOWLIFE

Brochure

Flowcontrast Spot Go

Portable Contrast Therapy Device

Table of Content

- 4 — Manifest
- 6 — Why Flowcontrast Spot Go
- 7 — Accessories
- 8 — Flowcontrast Spot Go - A Comprehensive Guide
- 12 — Starting your Journey
- 16 — For General Recovery
- 18 — For Muscle Soreness
- 18 — Before, During, or After Your Workout
- 25 — Contact Information
- 26 — Manual
 - English
 - Swedish / Svenska
 - Danish / Dansk
 - Norwegian / Norsk
 - German / Deutsch
 - French / Français
 - Dutch / Nederlands
 - Italian
 - Spanish / Espanol

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.


We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.

For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowcontrast Spot Go?

- **Portable contrast therapy for targeted recovery** Flowcontrast Spot Go delivers precise, alternating hot–cold therapy designed to reset your recovery response at the neuromuscular level. It uses time-based contrast cycles to support circulation, reduce pain, and improve mobility—all without water, ice, or cords.
- **Targeted contrast therapy, anywhere**
Alternates cold and heat in a pre-programmed cycle to support tissue repair, circulation, and inflammation management—where you need it most.
- **Adjustable comfort, clinical precision**
Choose between three intensity levels for both heat and cold. The default protocol (4 min cold, 7 min heat, 4 min cold) mirrors contrast techniques used in sports medicine.
- **Wireless. Wearable. Intuitive**
Designed for real life. Strap it to your knee, shoulder, or back with no setup time, no melting ice, and no waiting for cooldowns.


Accessories


What's Included

Frame Strap

helps keep device flat on wide areas

USB-C Charging Cable

3 Velcro Bands (S/M/L)

for customized fit on various body parts

Manual

Flowcontrast Spot Go: A Comprehensive Guide

How It Works

This device uses semiconductor thermal regulation to deliver focused hot and cold therapy without fluids or refrigeration. Instead of passive heat or disposable ice packs, Flowcontrast Spot Go creates a responsive temperature gradient that promotes circulation, reduces inflammation, and relaxes muscle tone. The system runs on pre-set minute sessions and automatically shuts down after completion.

Mode 1: Contrast Therapy Mode (4 min cold – 7 min heat – 4 min cold)

Contrast therapy works by alternating vasoconstriction (cold) and vasodilation (heat), creating a pumping action that stimulates blood flow, flushes metabolic waste, and helps relieve deep-tissue pain.

- Cold Phase
- Level 1 (15°C): Gentle cooling (surface-level calming)
 - Level 2 (10°C): Moderate cold (supports inflammation control)
 - Level 3 (5°C): Deep cold (maximizes constriction, ideal for acute pain or swelling)

- Heat Phase
- Level 1 (38°C): Mild warmth to support relaxation and comfort
 - Level 2 (42°C): Therapeutic heat for soft tissue pliability and mobility
 - Level 3 (45°C): Deep heat for fascial release and chronic tension relief

The default mode when starting is set on level 2 for the cold and heat phases. You can only adjust intensity during the active element (cold or heat) to preserve sequence integrity

Mode 2: Manual Cold

Contrast therapy works by alternating vasoconstriction (cold) and vasodilation (heat), creating a pumping action that stimulates blood flow, flushes metabolic waste, and helps relieve deep-tissue pain.

Activates cold-only therapy immediately at Level 1.

Tap the cold button to cycle through:

- Level 1 (15°C): Surface cooling for recovery after light strain
 - Level 2 (10°C): Moderate cold for post-training soreness or circulation reset
 - Level 3 (5°C): Max cooling for acute inflammation, joint pain, or swelling
- Best used within 24 hours of intense training, strain, or inflammation flare-up.

Mode 3: Manual Heat

Activates heat-only therapy immediately at Level 1.

Tap the heat button to cycle through:

- Level 1 (38°C): Soft warmth for nervous system calming (good for pre-sleep use)
- Level 2 (42°C): Improves range of motion and muscle elasticity
- Level 3 (45°C): Deep-tissue heat for recovery from chronic tightness or mobility restrictions

Use before stretching, mobility drills, or soft tissue therapy for optimal effect.

When to Use Contrast vs Manual

Therapy Goal	Recommended Mode	Why
Inflammation control after workout	Contrast or Cold Only	Flushes metabolites, reduces swelling
Morning stiffness or restricted movement	Heat Only	Prime muscles, increase elasticity
Chronic joint discomfort	Contrast	Alternating input stimulates tissue repair
Nervous system reset or calming	Heat Only, Level 1-2	Triggers the parasympathetic state
Acute injury (within 48 hrs)	Cold Only, Level 3	Slows blood flow, reduces inflammation
Pre-training warmup	Heat Only, Level 2-3	Flushes metabolites, Increases circulation and readiness

Starting Your Journey

Set it up. Strap it on—press start.

Flowcontrast Spot Go is built for real-life recovery. Whether you're dealing with soreness, inflammation, or just trying to stay mobile, getting started takes less than a minute.


Here's how to begin your first session:

1. Place the device on your body

Choose the muscle or joint you want to treat—commonly knees, shoulders, back, or calves.

2. Secure with the right strap

- Use one of the three Velcro bands to hold the device in place

- Small: for wrists, ankles, arms or elbows
- Medium: for knees, calves, or thighs
- Large: for back, quads, glutes or shoulders

- If you want the device to sit flat on a wide surface, use the frame strap to prevent bending or tilting

3. Start your mode

- Hold both cold + heat buttons → Start contrast mode
- Hold just cold → Cold only
- Hold just heat → Heat only

LED lights will tell you what's active:

- Blue: Cooling
- Red: Heating

4. Adjust intensity (optional)

- Tap the active button to switch between Level 1, 2, or 3
- Only the active function (cold or heat) can be adjusted

5. Let it run

Each mode has a pre-set timer, heat & cold mode runs for 10 minutes, contrast mode runs for 24 minutes, then shuts off automatically. You can relax, read, work, or simply sit still while it works.

6. After your session

- Unstrap and let the device cool before storing
- Wipe it down with a soft, dry cloth if needed
- Recharge via USB-C when the battery runs low

For General Recovery

How Flowcontrast Therapy Supports General Recovery. After long hours at your desk, standing on your feet, or pushing through demanding days, your body needs more than rest. It needs restoration. Recovery is what keeps you mobile, focused, and ready to go again tomorrow.

Flowcontrast Spot Go supports circulation—the foundation of recovery—by using alternating cold and heat to reset the body's internal rhythm. This contrast effect boosts blood flow, reduces low-grade tension, and clears out fatigue built up in overused muscles.

It's your daily reset when your day doesn't stop.

How to Use Flowcontrast Spot Go for General Recovery

- Use contrast mode on common tension zones like quads, calves, or the lower back
- Start at cold level 2 and heat level 2 for a balanced circulatory reset
- Best used in the evening after long days of standing, sitting, or commuting
- For calm recovery before sleep, use heat-only mode at level 1–2
- For joint stiffness, apply contrast before movement or daily heat-only at level 2

Why Contrast Therapy Works

When your day doesn't stop, neither does the load on your body. Prolonged tension, poor posture, or constant low-grade stress can slow circulation and fatigue your system. Contrast therapy reactivates healthy blood flow, flushes built-up waste, and restores a relaxed state—without needing intense effort or long downtime.

For Muscle Soreness

Training was the easy part. Recovery is where progress happens. Sore muscles need more than rest—they need circulation. Cold reduces inflammation. Heat brings in oxygen and nutrients. The contrast effect helps flush waste, ease tension, and accelerate tissue repair.

It's how your body resets after effort.

How to Use Flowcontrast Spot Go for Muscle Soreness

- Use contrast mode within 24 hrs of exercise
- Cold level 2–3, heat level 2–3 based on tolerance
- Target sore zones (quads, glutes, shoulders, knees, calves)
- For daily stiffness: use heat-only at level 2
- For swelling: cold-only at level 3 after activity

Why Contrast Therapy Works

Soreness signals inflammation and micro-damage. Contrast therapy supports the healing response: cold limits excess inflammation while heat increases nutrient delivery and waste clearance. Together, they accelerate tissue recovery without pushing your body further.

Before, During, or After Your Workout

Support every rep—from warm-up to recovery. Every workout has three parts: how you prepare, how you perform, and how you recover. Flowcontrast Spot Go supports all three. Use heat to prime muscles before training. Apply cold to limit fatigue mid-session. Use contrast to accelerate recovery when you're done. Whether it's a run, a heavy lift, or a long session on the field—this tool helps you stay ready, stay mobile, and stay in control.

How to Use Flowcontrast Spot Go Around Your Workout

- Use heat-only mode (level 2–3) for 10–15 minutes before activity to warm up
 - Apply to hips, calves, or shoulders to improve range and reduce restriction
- During breaks, cold-only mode (level 1–2) helps reduce fatigue or joint stress
- After training, use contrast mode to support faster recovery
- Focus on high-strain zones like quads, knees, or lower back post-workout

Movement begins where blood flows freely.

Healthy blood flow is essential to recovery. It delivers oxygen, nutrients, and warmth to your tissues—while removing waste like lactic acid and inflammatory byproducts. But when you're injured, inflamed, or overtrained, circulation can slow down or become blocked, delaying healing. That's where contrast therapy works differently.

By alternating cold and heat, Flowcontrast Spot Go creates a powerful “vascular pump” effect. Cold constricts blood vessels. Heat opens them back up. This rhythm boosts blood flow, clears inflammation, and triggers your body's natural repair signals—all in one cycle.

What it helps with:

- Muscle tightness after training
- Joint stiffness in the morning
- Swelling from minor injuries
- Recovery from long runs or lifting sessions
- Circulatory issues in high-use areas (knees, shoulders, calves)

When to use:

- After hard workouts
- For joint stiffness (knees, elbows, shoulders)
- To manage flare-ups from tendon strain or minor injury
- When dealing with long-standing pain that worsens with movement

Why Spot Go makes it easier:

Traditional contrast therapy needs two tubs or a lot of prep.

Spot Go gives you the same physiological effect—without water, ice, or guesswork.


Just strap it on, press start, and let the cycle do its work.

hello@flowlife.com
www.flowlife.com
Luntmakargatan 90
113 51 Stockholm

Manual

Flowcontrast Spot Go

ENGLISH

INTRODUCTION

Welcome to your Flowcontrast Spot Go. This guide will help you get started and use the device safely and effectively for contrast therapy, recovery, and mobility. Flowcontrast Spot Go delivers targeted hot and cold therapy in one compact, wireless unit. With adjustable intensity levels and a time-based contrast cycle, it supports circulation, relieves pain, and restores mobility—anytime, anywhere.

SAFETY INSTRUCTIONS

Safety Hazards – Reduce the Risk of Electric Shock

- Do not use or charge the device near water or in wet environments.
- Never immerse the device in water or clean it with running water.
- Use only the provided USB-C charging cable.
- Do not charge the device while in use.
- Disconnect immediately if the device or cable becomes damaged or wet.

Danger – Reduce the Risk of Burns, Fire, or Injury

- Do not leave the device unattended during use.
- Only use accessories supplied or approved by Flowlife.
- Do not block air vents or insert foreign objects into openings.
- Do not place heavy objects on the device or drop it.
- Do not disassemble, crush, or attempt to modify the device.

Heat & Cold Precautions

- Test skin sensitivity before extended use.
- Do not apply to broken, irritated, or numb skin.
- Avoid prolonged use on the same area (max 30 minutes per session).

- Stop use if redness, discomfort, or numbness occurs.
- Keep away from extreme cold or heat during storage.

Risk Groups – Consult a Doctor Before Use

- If pregnant or nursing
- If you have cardiovascular disease, cold intolerance, or Raynaud's
- If you have diabetes, nerve damage, or circulatory disorders
- If using medical implants or pacemakers
- If under 15 years old (unless supervised)

INSTRUCTIONS FOR USE

Using the Velcro Straps and Frame Strap

The device includes three Velcro straps to help you secure it comfortably to different body areas:

- Small (45×3 cm): wrists, ankles, or narrow joints
- Medium (75×3 cm): thighs, calves, or knee
- Large (110×3 cm): upper back, lower back, or shoulders

The strap should fit snugly to hold the device in place, but not so tightly that it restricts circulation. Adjust for a firm yet comfortable fit.

Use the frame strap when you want the device to remain flat—ideal for straight surfaces like the back or upper leg.

Remove the frame strap to allow the device to flex and curve, which helps it adapt better to rounded areas like the knee, shoulder, or calf.

Powering On & Off

- Hold the heat and cold buttons together for 2 seconds to enter Contrast Mode.
- Hold cold button (2 sec) for Cold Mode.
- Hold heat button (2 sec) for Heat Mode.
- Tap the active element button to cycle between Levels 1–3.

Mode Overview

- Contrast Mode: 26 min timer.
4 min cold → 7 min heat → 4 min cold --> 7 min heat --> 4 min cold
- Cold Mode: 10 min timer.
Direct cold, Level 1–3
 - Level 1: 15°C
 - Level 2: 10°C
 - Level 3: 5°C
- Heat Mode: 10 min timer.
Direct heat, Level 1–3
 - Level 1: 38°C
 - Level 2: 42°C
 - Level 3: 45°C

Note: You can only adjust the temperature of the active mode.

Auto Shut-Off

RECOMMENDED USAGE

- General Recovery: Use contrast mode in the evening or after long workdays.
- Muscle Soreness: Use contrast within 24 hrs of training. Cold-only helps inflammation.
- Joint Mobility: Use heat-only on stiff joints. Use contrast before stretching or movement.
- Max use per area: 30 minutes
- Let the temperature do the work—no pressure needed.

BATTERY & CHARGING

- Use the included USB-C cable
- Charging time: approx. 2.5 hours
- Battery capacity: 3200 mAh (7.2V)
- Power input: DC 5V = 2A
- Power output: 5V = 2A

BATTERY RUNTIME (APPROXIMATE)

- Cooling, Level 3: ~1 hour
- Heating, Level 3: ~5.5 hours
- Contrast Mode (Cooling + Heating), Level 3: ~1 hour

BATTERY INDICATOR

- Red: 0–33%
- Green: 34–66%
- White: 67–100%
- LED turns off 3 seconds after full charge

CLEANING & STORAGE

- Wipe the metal surface with a soft, damp cloth
- Do not use harsh chemicals or submerge the device
- Store in a cool, dry place
- Charge monthly if unused for extended periods

TECHNICAL SPECIFICATIONS

- Product Name: Flowcontrast Spot Go
- Size: 150 × 80 × 50 mm
- Weight: 384 g
- Temperature Range:
 - Cold: 5°C, 10°C, 15°C
 - Heat: 38°C, 42°C, 45°C
- Battery: 3200 mAh, 7.2V
- Input: 5V USB-C
- Output: 5V = 2A
- Materials: Aluminum contact plate, ABS housing
- Velcro Strap Sizes:
 - Small: 45×3 cm
 - Medium: 75×3 cm
 - Large: 110×3 cm
- Included Accessories: USB-C cable, manual, frame strap, 3 velcro straps (S/M/L)

WARRANTY & SUPPORT

Flowcontrast Spot Go includes a 2-year warranty covering manufacturing defects.

Warranty does not cover:

- Damage from misuse or accidents
- Modifications or unauthorized repairs
- Improper charging or liquid exposure
- Use not in accordance with this manual
- Missing proof of purchase

For questions or service, contact: support@flowlife.com

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