
Flowlight Wearable Knee 60 Two Waves

FLOWLIFE

Brochure

Flowlight Wearable Knee 60 Two Waves

660nm + 850nm NIR | 60 LEDs | Wearable | Wireless | Rechargeable

Table of Content

- 4 — Manifest
- 6 — Flowlight Wearable Knee 60 Two Waves
- 7 — Accessories
- 8 — Flowlight Wearable Knee: A Comprehensive Guide
- 12 — Starting your Journey
- 14 — Muscle Soreness and Recovery
- 18 — Joint Comfort and Mobility
- 22 — Active Recovery and Tendon Support
- 26 — Priming for Performance
- 30 — Structured Rest and Recovery Support
- 36 — Manual
 - English • Deutsch • Italiano
 - Svenska • Français • Dutch
 - Norsk • Espanol • Suomi
 - Dansk

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.


For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowlight Wearable Knee 60 Two Waves?

The Flowlight Wearable Knee 60 Two Waves delivers 660nm red and 850nm near-infrared light directly to the knee area — designed to support post-activity recovery, everyday comfort, and joint wellbeing. With dual-wavelength LEDs in direct skin contact and a fully wireless, wearable design, it fits naturally into your daily routine — at home, at the gym, or on the move.

- **Dual-wavelength light, worn on the body** - 60 dual-wavelength LEDs deliver 660nm red light and 850nm near-infrared across the full knee area. The flexible wrap ensures even coverage over the joint — used comfortably during rest or gentle movement.
- **Built for life in motion** - Whether you're winding down after a workout, supporting everyday knee comfort, or preparing for your next session — the Flowlight Wearable Knee 60 Two Waves wraps on in seconds. No clinic visit. No cables. No waiting room.
- **Wireless. Wearable. Rechargeable** - A removable battery pack powers full sessions with no cords. Charge independently, snap in, and go — or wear during rest, work, or gentle movement while the light does its work.

Accessories

What's Included

Flowlight Wearable Knee 60 Two Waves wrap — 60 triple-LED chips (180 total emitters), performance fabric with adjustable velcro straps

User manual

Two extra straps is included: — adjust for leg size and lock the wrap in place during movement

Works With

- Flowlight Wearable Battery pack — powers all Flowlight wearables, snap-in design for quick swapping between sessions

Flowlight Wearable Knee
60 Two Waves:
A Comprehensive Guide

The Flowlight Wearable Knee 60 Two Waves combines two wavelengths of light — 660nm red and 850nm near-infrared — designed to reach both the surface and deeper tissue around the knee in a single session.

Why it matters:

- Red and near-infrared light is absorbed by tissue — supporting the body's natural recovery processes during rest
- Dual-wavelength design delivers light across both the surface and deeper tissue around the knee in a single session
- Wearable wrap format allows consistent daily use at home — no clinic visits or appointments needed
- Wireless design means sessions fit around your schedule — not the other way around

660nm — Surface and Connective Tissue

Reaches the surface tissue around the knee, including the skin and outer connective layer. Supports local circulation and the body's natural recovery response in the tissues closest to the surface.

850nm — Deep Joint Penetration

Reaches deeper into the knee area — muscle and surrounding joint tissue. Supports the body's natural recovery processes at depth, helping the knee feel more comfortable after physical activity or prolonged use.

Combined Mode


The recommended starting point. Both wavelengths work together in a single session — covering the surface and deeper knee tissue for broad-spectrum post-activity recovery and everyday comfort.

Pulse Mode (Variable Frequency)

Cycles the LEDs on and off at a set frequency. Different frequencies may suit different preferences and use cases. Use + / - on the controller to adjust between 1 and 1000Hz and find what feels best for you.

Starting
Your Journey

1. Charge your battery pack before first use. Connect the Type-C cable to the battery pack. Full charge takes approximately 2.5 hours – the LED indicator shows 100% when complete.
2. Prepare your knee – clean and dry the skin. Remove any lotion, oils, or thick clothing from the joint area. Direct skin contact is essential for light to penetrate effectively.
3. Wrap and secure – place the LED device over the knee and fasten the velcro straps snugly. Firm contact across the joint, not restrictive.
4. Select your mode – Combined for general post-activity recovery, 660nm for surface-focused use, 850nm for deeper coverage, or Pulse to adjust light delivery frequency. Start with 10 minutes for your first session.
5. Press power to begin – the session timer counts down automatically. The device powers off when your session is complete.
6. After your session – remove the wrap, wipe clean with a soft dry cloth. Recharge the battery if below 30% before your next session.

A black and white photograph of a person's legs from the knees down, positioned in a dynamic, athletic pose. The legs are extended forward and slightly to the right, with the feet flat on the ground. The skin is smooth and the lighting creates soft shadows, emphasizing the contours of the muscles. The background is a plain, light-colored surface.

Muscle Soreness and Recovery

The body is built to repair. Help it do it faster.

Intense training creates microscopic tears in muscle fibers and the connective tissue surrounding the knee. Soreness peaks 24–72 hours after exercise, reducing performance in the sessions that follow and raising injury risk during the recovery window.

Accelerating Cellular Repair

Red and near-infrared light is absorbed by muscle and connective tissue around the knee, supporting the body's natural recovery process after physical exertion. With regular use, post-activity soreness duration may be reduced and recovery between sessions improved.

Why it matters:

Consistent use during recovery periods supports better readiness for your next session and helps manage cumulative fatigue across training cycles.

Clearing Inflammation and Restoring Circulation

660nm red light supports local circulation in the tissue around the knee after training. Improved circulation helps the area feel more comfortable and less heavy in the hours after physical activity.

Result:

Less stiffness. Reduced morning soreness. Faster return to full performance.

How to use

- Apply within 2 hours of training for best effect
- Wrap snugly with LEDs in full contact over the joint — no gaps
- Use Combined mode (660nm + 850nm) for post-training recovery
- Run a 15–20 minute session seated or lying down
- Pair with gentle elevation for additional swelling reduction

RECOMMENDED USE

Positioning	Wrap snugly over the knee with LEDs in full skin contact
Mode	Combined (660nm + 850nm)
Duration	15–20 minutes
Frequency	After every training session
Best time	Post-session — within 2 hours of exercise

Joint Comfort and Mobility

Keep moving. Feel better every day.

Everyday activity, age, and prolonged periods on your feet can leave the knee feeling stiff, uncomfortable, and less mobile. Without consistent support, daily movement becomes harder and recovery between active days slower.

Reaching Deeper Knee Tissue

850nm near-infrared reaches deeper into the knee area — warming the surrounding tissue and supporting the body's natural recovery process at depth. Regular use contributes to a feeling of greater comfort and ease of movement.

Why it matters:

Daily use supports long-term joint comfort and keeps the knee feeling its best through active periods and rest days alike.

Supporting Surface Comfort

660nm red light supports local circulation and surface tissue comfort around the knee. Improved circulation brings fresh oxygen and nutrients to the area — helping the knee recover from daily demands.

Result:

Reduced stiffness. Improved sense of mobility. Ongoing joint comfort with daily use.

How to use

- Use daily for cumulative wellbeing support — results build over weeks of consistent use
- Ensure full knee coverage with the LED device in firm contact across the joint
- Morning sessions help ease stiffness before daily movement begins
- Evening sessions support recovery from the demands of the day
- Try Pulse mode alongside Combined for variety in your routine

RECOMMENDED USE

Positioning	Full knee coverage — LED device in firm contact across the entire joint
Mode	Combined (660nm + 850nm), or add Pulse for a different delivery pattern
Duration	15–20 minutes
Frequency	Daily
Best time	Morning to ease into the day / Evening to wind down and recover

Active Recovery and Tendon Support

Strength begins with the tissue that holds it all together.

Tendons and connective tissue around the knee can take time to recover from intense activity, overuse, or sudden strain. Stiffness, soreness, and reduced comfort often linger longer than expected — slowing the return to full movement.

Supporting the Tissue Around the Joint

660nm red light supports local circulation and tissue health in the connective structures surrounding the knee. Consistent use during rest sessions helps maintain the area in better condition through the demands of an active lifestyle.

Why it matters:

Keeping connective tissue supported through regular use — not just when soreness strikes — contributes to better resilience over time.

Deeper Tissue Comfort

850nm NIR reaches deeper into the area around tendons and ligaments — supporting circulation and the body's natural recovery response in tissue that may take longer to bounce back after heavy or sustained use.

Result:

Reduced post-activity soreness. Better tissue comfort through the recovery period. More consistent readiness for the next session.

How to use

- Position the LED device over the area of the knee you want to focus on
- For the front of the knee: align LEDs below the kneecap
- For the inner or outer knee: use full wrap coverage across the joint line
- Use daily for consistent support — do not skip sessions
- For more intensive recovery periods: morning and evening sessions work well

RECOMMENDED USE

Positioning	LED device in firm contact over the area you want to focus on
Mode	Combined (660nm + 850nm)
Duration	20 minutes
Frequency	Daily during active recovery periods
Best time	Morning and evening during intensive recovery / Once daily for general maintenance

Priming for
Performance

The best sessions start before the first rep.

Cold, stiff joints and unprepared connective tissue reduce mechanical efficiency at the start of physical activity. Traditional warm-up protocols raise heart rate and temperature — but adding a light session beforehand may support how ready the knee feels before the first movement.

Warming the Joint from Within

850nm near-infrared supports local circulation and tissue warmth around the knee — a useful addition to a pre-activity routine that helps the joint feel ready before load is placed on it.

Why it matters:

A joint that is warm and well-circulated responds better to physical demand. Pre-activity use helps bridge the gap between rest and full performance.

Elevating Cellular Readiness

Applying red and near-infrared light before physical activity may help you feel more prepared from the first movement. Used as part of a pre-session routine, it complements other warm-up steps and supports tissue readiness.

Result:

Better performance from the first movement. Reduced warm-up time needed. Lower risk of early-session strain.

How to use

- Apply 15–30 minutes before activity begins
- Use while completing other preparation — taping, kit, mental preparation
- Follow the light session immediately with dynamic movement
- Pair with your existing warm-up for maximum benefit
- Use consistently before demanding sessions to build cumulative protection

RECOMMENDED USE

Positioning	Wrap over the knee 15–30 minutes before activity
Mode	Combined (660nm + 850nm)
Duration	10–15 minutes
Frequency	Before demanding training sessions
Best time	Pre-session — 15–30 minutes before activity

Structured Rest and Recovery Support

Recovery is not passive. It is a daily investment.

Extended recovery periods — whether after medical procedures, significant physical events, or demanding training phases — require consistent support. Keeping the knee comfortable, maintaining local circulation, and supporting the body's own processes through a structured rest period can make a meaningful difference to how quickly you feel ready to return to full activity.

Supporting Circulation and Surface Comfort

660nm red light supports local circulation and surface tissue comfort — making it well-suited for use during rest-based recovery periods when the knee needs gentle, consistent support without physical stress.

Why it matters:

Maintaining good local circulation through a recovery period supports the body's natural processes and helps keep the knee feeling as comfortable as possible day to day.

Deeper Tissue Warmth and Comfort

850nm NIR reaches deeper into the knee area — supporting local warmth, circulation, and the body's recovery response in tissue that may take longer to feel normal again after an intensive period. Use during rest sessions, in a comfortable seated or lying position.

Result:

Better comfort during the recovery period. Consistent daily support for the knee. A structured routine that complements the rest of your recovery plan.

How to use

- Always follow the guidance of your healthcare professional regarding any recovery period
- Apply to clean, intact skin only — do not use over any area that is not fully healed
- Use while resting with the leg supported and comfortable
- Start with shorter sessions (10 minutes) and increase gradually as comfort allows
- Progress to Combined mode once you are ready for broader coverage

RECOMMENDED USE


Positioning	LED coverage over the knee — use on intact, healthy skin only
Mode	660nm only initially; progress to Combined (660nm + 850nm) as comfort allows
Duration	10 minutes initially — increase to 20 minutes over time
Frequency	Daily, 5–7 days per week
Best time	Morning or evening — whenever rest and device use can be combined comfortably

hello@flowlife.com
www.flowlife.com
Flowlife AB
Runebergsgatan 7
114 30 Stockholm
Sweden

Manual

Flowlight Wearable Knee 60 Two Waves

1. magnetic clip

2. Velcro

3. Infrared light chips

ENGLISH

INTRODUCTION

Welcome to your Flowlight Wearable Knee 60 Two Waves. This guide will help you use your device safely and effectively. Read all instructions before first use and retain this manual for future reference.

SAFETY INSTRUCTIONS

Safety Hazards – Reduce the Risk of Electric Shock

- Do not use or charge near water or in wet environments
- Never immerse the device in water or clean with running water
- Use only the original Flowlife power adapter and Type-C cable
- Do not charge the device while in use
- Disconnect immediately if the device, cable, or battery pack appears damaged or wet
- Never store the device or cable on hot surfaces
- Unplug the charger when not in use
- Do not leave the charger unattended for extended periods in unsafe environments

Danger – Reduce the Risk of Burns, Fire, or Injury

- Do not leave the device unattended during use
- Do not cover the device or battery pack during use or charging – risk of overheating
- Do not disassemble, crush, or attempt to modify the device or battery pack
- Never insert objects into any openings on the device
- Stop use immediately if you detect unusual heat, a burning smell, or visible damage
- Never use the device where aerosol products are used or where oxygen is administered
- Never attempt to repair the device yourself – contact support@flowlife.com
- Never use if there is visible or suspected damage to the device, cable, or battery pack
- Keep out of reach of children under 15
- Children aged 15 and above may use the product only when properly informed about safe usage

Light Safety

- Never look directly at the LEDs during operation
- Do not direct the device toward the face or eyes

Skin & Application Precautions

- Do not use on open wounds, broken skin, active infections, or areas of recent burns
- Do not apply for more than 30 minutes per session to the same area
- Allow at least 6 hours between sessions on the same site
- Clean and dry skin before each session — remove all lotion or oils
- Stop use if redness, irritation, blistering, or persistent discomfort occurs
- Start with shorter sessions and increase gradually

Risk Groups – Consult a Doctor Before Use

- Pregnant or nursing — avoid use on the abdomen
- Active or suspected cancer at or near the treatment area — theoretical risk of increased cell activity
- Epilepsy or history of light-triggered seizures
- Cardiovascular disease, pacemaker, ICD, or any active implanted electronic device
- Hyperthyroidism — avoid applying to the neck and upper chest
- Fever or active infection — red and near-infrared light use may worsen hyperthermia
- Systemic lupus erythematosus (SLE) or other photosensitivity conditions
- Severe bleeding disorders or recent blood loss
- Photosensitising medications — consult your doctor before and after every session
- Peripheral neuropathy or reduced sensation in the treatment area — check skin carefully before and after every session
- Melasma or hyperpigmentation concerns — start slowly
- Under 15 years old — use only under supervision of a healthcare professional
- Any other health concern not listed — consult your doctor before starting

INSTRUCTIONS FOR USE

Before Your Session

1. Snap the battery pack in until seated
2. Clean and dry the skin — no lotions or oils

Applying the Device

1. Place the LED device over the knee with LEDs in full contact with skin
2. Fasten the velcro straps snugly – firm, not restrictive
3. Use the additional velcro straps if needed for a more secure fit or larger leg size

Using the Controller

1. Hold the power button to turn on
2. Press Mode – use + / – to select: Red / NIR / Combined
3. Press Mode again – use + / – to set intensity
4. Press Mode again – use + / – to set Pulse.
5. From the home screen, use + / – to set session timer (1–30 min)
6. Press Start to begin the session

After Your Session

1. Device powers off automatically when timer ends
2. Check skin for any redness after every session
3. Wipe the LED device with a soft, dry cloth

BATTERY & CHARGING

- Use the Type-C USB cable (included with the Flowlight Wearable Battery pack) to charge the battery pack
- Charging time: approx. 2.5 hours
- Battery type: removable lithium-ion pack – 6000mAh
- LED indicator shows 100% when fully charged
- Do not use the device while charging
- Dispose of used battery packs at a certified recycling point – do not place in household waste

CLEANING & STORAGE

- Wipe the LED device and neoprene with a soft, dry cloth after each session
- Do not use water, alcohol, or chemical cleaners on the LED device
- The neoprene may be wiped with a slightly damp cloth – dry fully before storage
- Do not machine wash
- Store in a cool, dry place away from direct sunlight
- Do not store at temperatures below –10°C or above 45°C
- Charge monthly if unused for extended periods

TECHNICAL SPECIFICATIONS

Specification	Value
Product Name	Flowlight Wearable Knee 60 Two Waves
Spectrum	660 nm (Red) + 850 nm (Near-Infrared)
LED Count	60 triple-LED chips — 180 total emitters (60×660nm + 120×850nm)
Irradiance	approx. 120 mW/cm ² at direct contact
Light Intensity	(10–100%)
Modes	Red / NIR / Combined / Pulse (1–1000Hz)
Pulse Frequency	1–1000 Hz
Session Timer	1–30 min
Power Source	Flowlight Wearable Battery pack
Shell Material	Neoprene
Dimensions	765 × 260 × 10 mm (L×W×H) ±2mm
Weight (excl. battery)	185g ±5g
Rated Power	60W
Working Current	Input & Output: 5V / 0.8A
Certifications	CE, UKCA
Warranty	2 years

WARRANTY & SUPPORT

Flowlight Wearable Knee 60 Two Waves includes a 2-year warranty covering manufacturing defects from the date of purchase.

Warranty does not cover:

- Physical damage from misuse, drops, or accidents
- Water damage or use in wet environments
- Missing proof of purchase
- Use with non-original accessories
- Unauthorized modifications or repair attempts

For questions or service, contact: support@flowlife.com | www.flowlife.com

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