
Flowlight Wearable Foot 60 Two Waves

FLOWLIFE

Brochure

Flowlight Wearable Foot 60 Two Waves

660nm + 850nm NIR | 60 LEDs | Wearable | Wireless | Rechargeable

Table of Content

- 4 — Manifest
- 6 — Flowlight Wearable Foot 60 Two Waves
- 7 — Accessories
- 8 — Flowlight Wearable Foot: A Comprehensive Guide
- 12 — Starting your Journey
- 14 — Heel and Arch Recovery
- 18 — Achilles and Heel Recovery
- 22 — Foot Circulation and Daily Comfort
- 26 — Ankle Recovery and Support
- 30 — Post-Activity Recovery
- 35 — Contact Information
- 37 — Manual
 - English • Deutsch • Dutch
 - Svenska • Français • Suomi
 - Norsk • Espanol
 - Dansk • Italiano

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.


For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowlight Wearable Foot 60 Two Waves?

The Flowlight Wearable Foot 60 Two Waves delivers 660nm red and 850nm near-infrared light to the sole, heel, and ankle — the areas that carry every load and rarely get the recovery they deserve. With dual-wavelength LEDs in direct skin contact, it supports post-activity recovery, everyday foot comfort, and general wellbeing — without clinic visits or extended downtime.

- **60 LEDs — full plantar and ankle coverage** - The flexible neoprene wrap conforms to the natural shape of the foot, delivering both wavelengths across the plantar fascia, heel, arch, and ankle simultaneously.
- **Designed for the foot's most demanding conditions** - From plantar fasciitis and Achilles tendinopathy to everyday foot discomfort — the Flowlight Wearable Foot 60 Two Waves targets the conditions that limit daily movement and performance.
- **Wireless. Wearable. Grounded.** - A removable battery pack powers full sessions with no cords. Charge independently, strap on, and recover — wherever your day takes you.

Accessories

What's Included

Flowlight Wearable Foot 60 Two Waves wrap — 60 triple-LED chips (180 total emitters), performance fabric with adjustable velcro straps

User manual

Works With

- Flowlight Wearable Battery pack — powers all Flowlight wearables, snap-in design for quick swapping between sessions

Flowlight Wearable Foot
60 Two Waves:
A Comprehensive Guide

The Flowlight Wearable Foot 60 Two Waves combines two wavelengths of light — 660nm red and 850nm near-infrared — designed to reach both the surface and deeper tissue of the foot in a single session.

Why it matters:

- The Flowlight Wearable Foot 60 Two Waves combines two wavelengths of light — 660nm red and 850nm near-infrared — designed to reach both the surface and deeper tissue of the foot in a single session.
- Red and near-infrared light is absorbed by tissue — supporting the body's natural recovery processes in the foot and ankle during rest
- Dual-wavelength design delivers light across both the surface and deeper tissue of the foot in a single session
- Wearable wrap format allows consistent daily use at home — no clinic visits or appointments needed
- Wireless design means sessions fit around your schedule — not the other way around

660nm — Surface and Plantar Tissue

Reaches the surface tissue across the plantar sole and heel. Supports local circulation and the body's natural recovery response in the tissue closest to the surface — including the arch and surrounding skin.

850nm — Deep Tissue and Nerve Pathways

Reaches deeper into the heel, Achilles area, and surrounding tissue. Supports local circulation and the body's natural recovery processes at depth — helping the foot and ankle feel more comfortable after sustained activity or prolonged use.

Combined Mode


The recommended starting point. Both wavelengths work together in a single session, covering the surface and deeper foot tissue for broad-spectrum post-activity recovery and everyday comfort.

Pulse Mode (Variable Frequency)

Pulse mode cycles the LEDs on and off at a set frequency. Different frequencies suit different preferences. Use + / – on the controller to adjust between 1 and 1000Hz and find what works best for you.

Starting
Your Journey

1. Charge your battery pack before first use. Connect the Type-C cable to the battery pack. Full charge takes approximately 2.5 hours — the LED indicator shows 100% when complete.
2. Prepare your foot — clean and dry the skin. Remove any lotion or oils. Direct skin contact is essential.
3. Wrap and secure — place the LED device flat against the sole of the foot and fasten the straps snugly. Firm contact, not restrictive.
4. Select your mode — Combined for general post-activity recovery, 660nm for surface coverage, 850nm for deeper tissue use, or Pulse to adjust light delivery frequency. Start with 10–15 minutes for your first session.
5. Press power to begin — the session timer counts down automatically. The device powers off when your session is complete.
6. After your session — remove the wrap, wipe clean with a soft dry cloth. Recharge the battery if below 30% before your next session.

A black and white photograph of a person's legs and feet, showing the heel and arch of the foot. The person is sitting on the floor, with their legs extended forward. The feet are flat on the ground, and the arches are visible. The lighting is soft, highlighting the contours of the legs and feet. The background is a plain, light-colored surface.

Heel and Arch Recovery

Every step begins here.

Heel and arch soreness — especially the kind that's worst with the first steps of the day — is one of the most common complaints among active people, those who spend hours on their feet, and anyone whose feet absorb daily impact without enough recovery time.

Reaching the Heel and Arch

850nm NIR reaches into the heel and arch area — supporting local circulation in tissue that absorbs a lot of daily load. 660nm supports the plantar surface tissue closest to the skin, helping the underside of the foot feel more comfortable during rest.

Why it matters:

Daily use at rest supports the foot's natural recovery, especially for those whose feet rarely get a real break.

Supporting Comfort Through Consistent Use

Red and near-infrared light together support circulation and the body's natural recovery response across the plantar surface. Used consistently — particularly in the morning before the first steps — the foot may feel more comfortable and mobile through the day.

Result:

Reduced morning stiffness. Greater comfort with first steps. Gradual improvement in foot comfort with regular use.

How to use

- Place the wrap with LEDs flat against the plantar surface — full heel and arch coverage
- Sit or lie during sessions — no weight-bearing
- Use Combined mode for most sessions
- Morning sessions help ease stiffness before first steps / Evening sessions support overnight recovery

RECOMMENDED USE

Positioning	LED device flat against the plantar surface — full heel and arch coverage
Mode	Combined (660nm + 850nm)coverage
Duration	15–20 minutes
Frequency	Daily
Best time	Morning — before first steps / Evening — to support overnight comfort

Achilles and Heel Recovery

What connects you to the ground deserves attention.

The Achilles tendon and heel area take significant load with every step. After intense activity, prolonged running, or periods of overuse, stiffness and soreness in this area can linger — slowing recovery and making the first movements of the day noticeably uncomfortable.

Reaching the Achilles and Heel Area

850nm NIR reaches deeper into the heel and Achilles area — supporting circulation in tissue that is under regular load and slow to recover from intense activity. 660nm supports the surrounding tissue and surface circulation around the heel.

Why it matters:

Consistent rest-based use helps the Achilles area stay in better condition through the demands of an active lifestyle — especially during high-mileage or high-load periods.

Supporting Recovery Between Sessions

Red and near-infrared light used consistently during rest may help the heel and Achilles area feel more comfortable and better prepared for the next bout of activity. Used morning and evening during intensive periods, it fits into a structured recovery routine

Result:

Reduced morning stiffness in the heel. Greater comfort during the recovery period. Better readiness for the next activity.

How to use

- Extend LED coverage up over the heel and lower Achilles area
- Ensure the ankle strap holds the wrap in firm contact with the back of the foot
- Use Combined mode throughout the recovery period
- Apply in the morning before activity and again in the evening for overnight recovery

RECOMMENDED USE

Positioning	LED coverage up over the heel and lower Achilles — ankle strap firm against the back of the foot
Mode	Combined (660nm + 850nm)
Duration	20 minutes
Frequency	Daily throughout recovery
Best time	Morning before activity / Evening to support overnight recovery

Foot Circulation and Daily Comfort

Keep your feet feeling their best every day.

For people who experience tingling, heaviness, sensitivity, or reduced circulation in the feet — whether from extended standing, physical demands, or underlying health factors — consistent foot care and support are essential to maintaining comfort and quality of daily movement.

Supporting Circulation in the Foot

850nm NIR reaches deeper into the foot area, supporting local circulation and the body's natural recovery response. 660nm supports surface tissue and local circulation across the plantar sole — beneficial for feet that experience heaviness or reduced warmth.

Why it matters:

Improved local circulation supports general foot comfort and may help the foot feel more vital and comfortable through the day.

Pulse Mode for Personalised Comfort

Pulse mode delivers light in on/off cycles at an adjustable frequency. Some users find this delivery pattern particularly comfortable for feet that feel sensitive or tingly. Use + / – on the controller to adjust the frequency and find what suits you.

Result:

Improved foot comfort with regular use. Better circulation support. A simple, non-invasive addition to a daily foot care routine.

How to use

- Apply to clean, intact skin only — check the skin before and after every session
- If you have an underlying health condition affecting your feet, consult your doctor before starting
- Use 850nm or combined
- Evening sessions work well when feet feel most in need of recovery after a full day

RECOMMENDED USE

Positioning	Apply to clean, intact skin only – inspect skin before and after every session
Mode	850nm or Combined
Duration	15–20 minutes
Frequency	Daily
Best time	Evening – after a full day on your feet

A black and white photograph of a person's lower leg and foot, showing the ankle and calf muscles. The leg is positioned vertically, with the foot flat on the ground. The skin is light-colored, and the muscles are well-defined. The background is a plain, light-colored surface.

Ankle Recovery and Support

A stable ankle is the foundation of everything above it.

Ankle soreness, stiffness, and swelling after physical activity or a misstep are common complaints — and without consistent recovery support, the ankle can take much longer than expected to feel stable, comfortable, and ready for full movement again

Supporting the Ankle Area

660nm red light supports local circulation and surface tissue comfort around the ankle. 850nm NIR reaches deeper into the ankle area — supporting circulation and the body's natural recovery response in the surrounding tissue.

Why it matters:

Consistent daily use keeps the ankle area supported through the recovery period — contributing to better comfort and gradual return to normal movement.

Reducing Post-Activity Swelling and Discomfort

660nm light supports local circulation around the ankle, which may help the area feel less heavy and more comfortable after activity or a period of swelling.

Result:

Reduced ankle soreness. Improved comfort during the recovery period. Greater readiness to return to full movement.

How to use

- Full wrap coverage over the ankle and surrounding area
- Combined
- Elevate the foot during sessions where possible
- Daily use throughout recovery — morning and evening works well for active recovery periods

RECOMMENDED USE

Positioning	Full wrap coverage with LED contact over the ankle and surrounding area
Mode	Combined (660nm + 850nm)
Duration	15–20 minutes
Frequency	Daily throughout recovery
Best time	Morning and evening during intensive recovery / Once daily for general maintenance

Post-Activity Recovery

After every mile, give your feet better than rest.

After prolonged standing, running, or high-impact activity, the foot can feel heavy, tight, and in need of real recovery. Without active support, this fatigue compounds across training weeks — leaving the foot less ready for the demands ahead.

Supporting Post-Activity Recovery

850nm NIR supports local circulation in fatigued plantar tissue, helping the foot feel less heavy after activity. 660nm supports the surface tissue across the sole — contributing to an overall feeling of lightness and reduced tightness after long sessions.

Why it matters:

Faster cellular restoration means more capacity available for the next session — and less cumulative tissue stress over weeks and months of training.

Improving Foot Circulation

Consistent post-activity red and near-infrared light use improves the microcirculatory response in the foot. Improved local circulation means better oxygen and nutrient delivery to plantar structures, and better removal of inflammatory byproducts that accumulate during high-load activity.

Result:


Less foot fatigue at day's end. Reduced plantar tension across training weeks. Lower cumulative injury risk with consistent use.

How to use

- Apply immediately after activity or within 2 hours
- Wrap snugly with full plantar contact
- Use Combined mode for comprehensive recovery
- Follow with gentle foot elevation for additional circulatory benefit

RECOMMENDED USE


Positioning	Apply immediately after activity or within 2 hours — wrap snugly with full plantar contact
Mode	Combined (660nm + 850nm)
Duration	15–20 minutes
Frequency	Every training or high-activity day
Best time	Post-activity — within 2 hours of exercise or prolonged standing

hello@flowlife.com
www.flowlife.com
Flowlife AB
Runebergsgatan 7
114 30 Stockholm
Sweden

Manual

Flowlight Wearable Foot 60 Two Waves

1. Velcro
2. magnetic clip
3. Infrared light chips

ENGLISH

INTRODUCTION

Welcome to your Flowlight Wearable Foot 60 Two Waves. This guide will help you use your device safely and effectively. Read all instructions before first use and retain this manual for future reference.

SAFETY INSTRUCTIONS

Safety Hazards – Reduce the Risk of Electric Shock

- Do not use or charge near water or in wet environments
- Never immerse the device in water or clean with running water
- Use only the original Flowlife power adapter and Type-C cable
- Do not charge the device while in use
- Disconnect immediately if the device, cable, or battery pack appears damaged or wet
- Never store the device or cable on hot surfaces
- Unplug the charger when not in use
- Do not leave the charger unattended for extended periods in unsafe environments

Danger – Reduce the Risk of Burns, Fire, or Injury

- Do not leave the device unattended during use
- Do not cover the device or battery pack during use or charging – risk of overheating
- Do not disassemble, crush, or attempt to modify the device or battery pack
- Never insert objects into any openings on the device
- Stop use immediately if you detect unusual heat, a burning smell, or visible damage
- Never use the device where aerosol products are used or where oxygen is administered
- Never attempt to repair the device yourself – contact support@flowlife.com
- Never use if there is visible or suspected damage to the device, cable, or battery pack
- Keep out of reach of children under 15
- Children aged 15 and above may use the product only when properly informed about safe usage

Light Safety

- Never look directly at the LEDs during operation
- Do not direct the device toward the face or eyes

Skin & Application Precautions

- Do not use on open wounds, broken skin, active infections, or areas of recent burns
- Do not apply for more than 30 minutes per session to the same area
- Allow at least 6 hours between sessions on the same site
- Clean and dry skin before each session — remove all lotion or oils
- Stop use if redness, irritation, blistering, or persistent discomfort occurs
- Start with shorter sessions and increase gradually

Risk Groups – Consult a Doctor Before Use

- Pregnant or nursing — avoid use on the abdomen
- Active or suspected cancer at or near the treatment area — theoretical risk of increased cell activity
- Epilepsy or history of light-triggered seizures
- Cardiovascular disease, pacemaker, ICD, or any active implanted electronic device
- Hyperthyroidism — avoid applying to the neck and upper chest
- Fever or active infection — red and near-infrared light use may worsen hyperthermia
- Systemic lupus erythematosus (SLE) or other photosensitivity conditions
- Severe bleeding disorders or recent blood loss
- Photosensitising medications — consult your doctor before and after every session
- Peripheral neuropathy or reduced sensation in the treatment area — check skin carefully before and after every session
- Melasma or hyperpigmentation concerns — start slowly
- Under 15 years old — use only under supervision of a healthcare professional
- Any other health concern not listed — consult your doctor before starting

INSTRUCTIONS FOR USE

Before Your Session

1. Snap the battery pack in until seated
2. Clean and dry the skin — no lotions or oils

Applying the Device

1. Place the LED device flat against the sole of the foot — full heel and arch coverage
2. For heel or Achilles focus, align LEDs over the heel area
3. Fasten the straps snugly — firm, not restrictive

Using the Controller

1. Hold the power button to turn on
2. Press Mode — use + / – to select: Red / NIR / Combined
3. Press Mode again — use + / – to set intensity
4. Press Mode again — use + / – to set Pulse.
5. From the home screen, use + / – to set session timer (1–30 min)
6. Press Start to begin the session

After Your Session

1. Device powers off automatically when timer ends
2. Check skin for any redness after every session
3. Wipe the LED device with a soft, dry cloth

BATTERY & CHARGING

- Use the Type-C USB cable (included with the Flowlight Wearable Battery pack) to charge the battery pack
- Charging time: approx. 2.5 hours
- Battery type: removable lithium-ion pack — 6000mAh
- LED indicator shows 100% when fully charged
- Do not use the device while charging
- Dispose of used battery packs at a certified recycling point — do not place in household waste

CLEANING & STORAGE

- Wipe the LED device and neoprene with a soft, dry cloth after each session
- Do not use water, alcohol, or chemical cleaners on the LED device
- The neoprene may be wiped with a slightly damp cloth — dry fully before storage
- Do not machine wash
- Store in a cool, dry place away from direct sunlight
- Do not store at temperatures below –10°C or above 45°C
- Charge monthly if unused for extended periods

TECHNICAL SPECIFICATIONS

Specification	Value
Product Name	Flowlight Wearable Foot 60 Two Waves
Spectrum	660 nm (Red) + 850 nm (Near-Infrared)
LED Count	60 triple-LED chips – 180 total emitters (60×660nm + 120×850nm)
Irradiance	approx. 120 mW/cm ² at direct contact
Light Intensity	10–100%
Modes	Red / NIR / Combined / Pulse (1–1000Hz)
Pulse Frequency	1–1000 Hz
Session Timer	1–30 min
Power Source	Flowlight Wearable Battery pack
Shell Material	Neoprene
Dimensions	470 × 30 mm
Weight (excl. battery)	278g
Rated Power	60W
Working Current	Input & Output: 5V / 0.8A
Certifications	CE, UKCA
Warranty	2 years

WARRANTY & SUPPORT

Flowlight Wearable Foot 60 Two Waves includes a 2-year warranty covering manufacturing defects from the date of purchase.

- Physical damage from misuse, drops, or accidents
- Water damage or use in wet environments
- Missing proof of purchase
- Use with non-original accessories
- Unauthorized modifications or repair attempts

For questions or service, contact: support@flowlife.com | www.flowlife.com

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