
Flowpemb IR Mat Pro

FLOWLIFE

Brochure

Flowpemf IR Mat Pro

Pulsed Electromagnetic Field Mat

Table of Content

4 —	Manifest
6 —	Flowpemf IR Mat Pro
7 —	Accessories
8 —	Eight Coils, Full-Body Coverage
10 —	Far-Infrared Heat, Head To Feet
12 —	The Science Inside The Mat
16 —	A Comprehensive Guide
18 —	Starting Your Journey
24 —	Temperature Guide
26 —	Your First Sessions
30 —	Building A Recovery Ritual
37 —	Contact Information
38 —	Manual
40 —	• English • Svenska • Norsk • Dansk • Deutsch • Français • Español • Italiano • Dutch • Suomi

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.

For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowpemf IR Mat Pro?

Recovery is no longer something you fit around training. It is part of training. Flowpemf IR Mat Pro turns recovery into a daily ritual you can keep — predictable, repeatable, and built into the rhythm of an active life.

Eight copper PEMF coils — Eight copper PEMF coils running the full length of the mat — head to feet — producing a wide, continuous electromagnetic field across the entire body. No zones. No drop-off. Wherever you lie, you're in the field.

True far-infrared heat — not surface warmth. Flowpemf IR Mat Pro emits far-infrared radiation in the 5–14 μm range. It penetrates past the skin into connective tissue and muscle, dilating blood vessels at depth and relaxing the body from the inside out. You feel warm.

Three guided programs — Sleep, Circulation, Recovery. P1, P2, P3 on the controller. One for winding down at night. One for clearing the static of a sedentary day. One for the deep reset after training. No guesswork, no overthinking.

Full manual control when you want it. Dial PEMF from 1 to 30 Hz, set heat from 30 to 70°C, run sessions from one to twelve hours. The presets cover most days. Manual mode is there when they don't.



Heated mineral therapy at the core. Natural amethyst and black tourmaline are layered into the mat above the heating element. Warmed by it, they hold and redistribute the heat across the surface — and re-radiate it as far-infrared in the 5–14 μm range, the wavelength the body absorbs most readily.

Low-EMF heating system. Heated mats often introduce more electromagnetic noise than people realise. Flowpemf IR Mat Pro uses a silicone-resistance heating system that measures just 0.2–0.8 milligauss — a low, near-passive reading, well below typical heated mats and most household appliances. The deliberate field — the PEMF — only engages when you choose to activate it.

A serious piece of recovery equipment. Full-body footprint. Substantial weight. Built for a dedicated space in your home — a recovery corner, a studio, a quiet room. Folds down for storage and travels in its carry bag when you need to move it.

Accessories

Flowpemf IR Mat Pro comes with everything you need to build a recovery practice at home.

Flowpemf IR Mat Pro Recovery Mat

A full-body 17-layer recovery surface, built around eight copper PEMF coils, far-infrared heat, and heated mineral therapy. Non-toxic synthetic PU leather on the surface, breathable polyester mesh, integrated graphene and far-infrared composite layers inside, and an even thermal field from edge to edge.

Flowpemf IR Mat Pro Controller


A seven-button controller that puts every setting one press away. Power, decrease, increase, Temp, Time, PEMF, and Mode — each function on its own dedicated key. Press Mode to load one of the three preset programs (P1 Sleep, P2 Circulation, P3 Recovery), or use the dedicated keys to set heat between 30 and 70°C, PEMF between 1 and 30 Hz, and the timer between 1 and 12 hours. Built-in safety with automatic overheat shutoff and a high-stability temperature sensor.

PEMF Tester

A small plastic cylinder with free-moving magnetic balls inside. Hold it close to the mat surface while a program is running and the balls visibly move and click in response to the magnetic field — a simple, hands-on way to confirm the PEMF you can't feel is actually working. No batteries, no electronics.

Carry and Storage

The mat folds along its built-in fold lines and stows inside the included carry bag. Designed primarily for a permanent place at home — and rated to handle the occasional move between rooms or trip away.

EIGHT COILS,
FULL-BODY
COVERAGE

The thing inside the mat that nobody sees — and the difference you feel.

Flowpemf IR Mat Pro is built around eight copper PEMF coils running the full length of the mat — from head to legs. A wide, deliberately engineered array that produces a continuous, even electromagnetic field across the entire body. Wherever you lie, you're in the field.

Coils are the heart of any PEMF system. The more coils, and the better their distribution, the more uniform the stimulus and the more of the body each session reaches. Eight is a deliberate choice — wide enough to cover head to feet, evenly enough to keep the field consistent from one end of the mat to the other.

Wide coverage

Eight coils running head-to-feet — the whole body in the field, every session.

Even distribution

A continuous electromagnetic field from head to lower back to legs, with no drop-off between zones.

Copper construction

Copper is the standard for high-precision electromagnetic coils — efficient, stable, and consistent over years of use.

A black and white photograph of a woman performing a backbend (Urdhva Dhanurasana) on a dark, textured mat. She is wearing a dark, strappy backless top. Her head is tucked down towards her knees, and her arms are extended forward, with her hands resting on the mat. The background is a plain, light color. The text "FAR-INFRARED HEAT, HEAD TO FEET" is overlaid in white, sans-serif capital letters on the left side of the image.

FAR-INFRARED
HEAT, HEAD TO
FEET

Most heated mats warm your skin. Flowpempf IR Mat Pro warms the body.

There is heat and there is far-infrared heat — and they are not the same thing. Conventional heating mats raise the temperature of the surface you're lying on. Flowpempf IR Mat Pro emits far-infrared radiation in the 5 to 14 μm range — the same wavelength the human body emits, and the same range the body absorbs most efficiently. It's sometimes called the "vital ray." Instead of stopping at the skin, it penetrates past the surface and into connective tissue and muscle. The result is a deeper, more even warmth that dilates blood vessels at depth, supports circulation, and helps muscles release without an aggressive surface temperature. You feel warm. You don't feel scorched.

5–14 μm wavelength

The body's optimal absorption range. The wavelength the body resonates with naturally — and the reason far-infrared feels different from surface heat.

Multiple sources, one field

Multiple sources, one field. A non-EMF silicone-resistance heating wire, a dedicated far-infrared composite layer, and the natural emission of heated amethyst and black tourmaline — every layer in the mat works toward the same far-infrared output.

Reflective build

An aluminium foil layer reflects far-infrared back toward the body rather than letting it escape downward. More of the heat reaches you, less of it is lost to the floor.

30 to 70°C, your choice

From a low, sleep-friendly warmth to a full far-infrared sauna effect — adjustable across the full range, with safety cutoff and a high-stability temperature sensor as standard.

THE SCIENCE INSIDE THE MAT

Three technologies, working together. None of them new. All of them held to the same standard: do they actually help the body recover?

Pulsed Electromagnetic Fields (PEMF)

PEMF stands for pulsed electromagnetic fields — low-frequency, low-intensity magnetic pulses delivered through coils embedded in the mat. Unlike heat, you don't feel them. The body does.

Flowpemf IR Mat Pro uses eight copper coils to generate this field, running the full length of the mat. The result is a wide, continuous electromagnetic field across the body — and a consistent stimulus regardless of where you lie.

The 1–30 Hz range used in Flowpemf IR Mat Pro is sometimes called the body's "window" — frequencies that early bioelectromagnetic research (Adey, Bawin) identified as the range cells respond to most readily. Within this window, PEMF has been studied as a way to support circulation, reduce perceived muscle guarding, and help the nervous system transition out of a stress-dominant state.

Three of those frequencies are particularly well-mapped, and they form the basis of Flowpemf IR Mat Pro's three preset programs.

FREQUENCY MAP

5 Hz — Sleep

Calming and grounding. Used to support sleep, downregulate the nervous system, and ease tension after long mental or emotional days.

15 Hz — Circulation

The frequency most studied for supporting blood flow and easing tissue stiffness. Built for sedentary days, long flights, and the dull fatigue that comes from not moving enough.

30 Hz — Recovery

A higher-frequency stimulus used to support post-training recovery — accelerating circulation, easing muscle soreness, and shortening the gap between sessions.

Far-Infrared Heat

Far-infrared heat is experienced as a deep, even warmth rather than surface heat. As temperature rises, blood vessels dilate, circulation improves, and muscles and connective tissue settle. The body reads it as a clear physical signal — it is safe to rest.

Flowpemf IR Mat Pro generates far-infrared heat through its internal heating system, an integrated far-infrared composite layer, and the natural emission of heated minerals. The temperature range is 30 to 70°C. Most users will find their sweet spot somewhere between 35 and 55°C; higher is not better.

Heated Mineral Therapy

Two natural minerals sit inside the mat — not as decoration, but as functional layers.

Amethyst

When heated, amethyst emits far-infrared radiation in the 5–14 μm range — the wavelength range the human body absorbs most efficiently. It also contributes to the mat's steady, even thermal field and has long been associated with calm, focus, and emotional regulation.

Black Tourmaline


Black tourmaline also emits far-infrared when heated and is traditionally associated with grounding and detoxification. Functionally, it helps the mat deliver a consistent heat distribution and adds to the depth of the thermal experience.

17 Functional Layers

A recovery mat is only as good as the experience you can repeat. Flowpemp IR Mat Pro's 17-layer construction exists to make every session feel the same — stable, comfortable, predictable. Every layer is functional.

1. Non-toxic synthetic PU leather — durable, easy to clean, designed for daily contact.
2. 100% polyester high-density mesh — breathable and heat-resistant.
3. Natural amethyst and tourmaline — far-infrared emission and steady heat distribution.
4. Non-woven fabric layer — structural separation between mineral and ion layers.
5. Negative ion composite layer — emits anions into the immediate environment.
6. Aluminium foil layer — disperses heat evenly across the full surface, shields against stray electromagnetic waves, and reflects far-infrared back into the body.
7. Non-EMF / silicone-resistance heating system — high-precision warmth without conventional EMF spikes.
8. Non-woven fabric — fixes the heating elements in place for long-cycle stability.
9. High-stability temperature sensor — continuously reads the mat's thermal state.
10. Overheat protector — automatic shutoff if the temperature ever exceeds the safe threshold.
11. Far-infrared composite layer — emits infrared radiation for deeper, more even penetration.
12. Graphene composite layer — moisture absorption and circulatory support.
13. Yellow mud mineral layer — traditionally associated with detoxification.
14. Nano charcoal lamination — adds far-infrared emission and supports surface hygiene.
15. Eight-coil PEMF system — generates the pulsed magnetic field across the full body.
16. Eco-friendly compressed cotton / sponge — thermal protection and body comfort.
17. Waterproof bottom material — durability and base protection.

The result is not one magic feature. It is the sum of clean engineering choices that make sessions feel consistent — every time.

A Comprehensive Guide

Flowpempf IR Mat Pro is built for one thing — to make recovery something you actually do, not something you mean to do. This guide covers setup, the three preset programs, manual mode, and how to build a recovery ritual around training, work, and sleep.

The principle is simple: lie down, choose a setting, breathe. The mat does the work. There are no complicated routines, no equipment to manage, no choreography to remember. Just consistency.

Recovery is cumulative

Three to six moderate sessions a week will move you further than two intense ones. Lower the bar to use it. Raise the floor of how you feel.

FLOWLIFE

Starting Your
Journey

Prepare Your Body and Setup

Choose a calm, stable environment. The mat works best when the body can settle without distraction.

Lay the mat on a flat surface — a bed, the floor, or a yoga mat. Wear light clothing, or use a thin towel between you and the surface if preferred. Ensure the controller cable is free of bends or compression.

Silence your phone. Lie down in a comfortable position and let your breathing slow before starting the session.

Powering On

Connect the mat to the controller and the power source. Power on. Select your preferred preset, or switch to manual mode and set heat and PEMF yourself. The mat warms gradually; PEMF begins automatically with the chosen program.

During Your Session

Stay still, or move gently if you like. Many users combine sessions with slow breathing, meditation, or light stretching. P2 Circulation and P3 Recovery run for 30 minutes each. P1 Sleep runs for 8 hours, with PEMF active during the first 30 minutes. In manual mode you can set the timer anywhere from 1 to 12 hours. Comfort and consistency matter more than intensity.


Ending Your Session

Turn off the controller and unplug the mat. Let it cool before folding or storing. Wipe the surface with a soft, dry — or slightly damp — cloth as needed. Store flat or loosely folded in a dry place.

SAFETY FIRST

PEMF therapy is generally well-tolerated, but Flowpemf IR Mat Pro is not intended for use during pregnancy, by anyone with an implanted electronic device (pacemaker, defibrillator, insulin pump, cochlear implant), or by anyone with active malignancies, severe cardiac arrhythmia, or active bleeding. If you are uncertain, speak to your doctor before starting. Full safety information is included in the Manual section of this booklet.

USING THE CONTROLLER

Seven buttons. Three preset programs. Full manual control. Every setting is one press away — heat, time, PEMF, and program selection each sit on a dedicated key.

1 · Power

Press the power button to turn the controller and the product's functionality on or off.

2 & 3 · Decrease ▼ and Increase ▲

The shared adjustment keys for every value on the controller. Once you press Temp, Time, or PEMF the setting name starts to blink — use ▼ and ▲ to step the value down or up. Changes apply immediately and start when the setting name stops blinking and glows steady.

4 · Temp

Toggles the heating function. Press once, then use ▼ and ▲ to set the temperature between 30 and 70°C. The value flashes while you adjust, then returns to the current actual temperature.

5 · Time

Toggles the session timer. Press to activate, then use ▼ and ▲ to set the session length between 1 and 12 hours. The session ends automatically when the timer reaches zero.

6 · PEMF

Toggles the pulsed electromagnetic field on or off. Once active, use ▼ and ▲ to dial the frequency between 1 and 30 Hz in single-Hz steps. PEMF works in both preset and manual modes. Pressing PEMF puts the field in standby — it activates only once heat and time are running (started via the Power button). Use the PEMF tester to confirm the field is active during a session.

7 · Mode

Cycles through the three preset programs — P1 Sleep, P2 Circulation, P3 Recovery. Each press selects the next preset; the active program is shown on the display. The mat loads that program's frequency, temperature, and timer in a single step.

CONTROLLER AT A GLANCE

P1 — Sleep	5 Hz · 35°C · PEMF 30 min · heat 8 h
P2 — Circulation	15 Hz · 40°C · 30 min · daytime reset
P3 — Recovery	30 Hz · 45°C · 30 min · post-training
Manual mode	PEMF 1–30 Hz · Heat 30–70°C · Timer 1–12 h

P1 — SLEEP 5 HZ · 35°C · PEMF 30 MIN · HEAT 8 H

Sleep is the overnight mode. A short PEMF wind-down at the start, then a steady, gentle warmth that runs alongside you through the night.

Designed for nervous system downregulation and deep sleep support, Sleep pairs a 5 Hz PEMF signal — a frequency long associated with deep relaxation and the transition into rest — with a gentle far-infrared warmth of 35°C, matching the body's natural night-time core temperature. PEMF runs for the first 30 minutes to help the body downshift, then stops automatically. The 35°C heat continues for the full 8-hour cycle, providing a clear, low-grade signal to the nervous system that it is safe to rest. The mat shuts off automatically at the end of the cycle.

Sleep is the mode for evenings, for long mental days, and for the night itself. It is not stimulating. It is not meant to be. The work of Sleep is the work of letting go.

When to Use

At bedtime. Place the mat on or under your sheet, set the program, and lie down. The PEMF burst eases you out of the day in the first half hour; the heat keeps the conditions for rest in place for the rest of the night.

Heat: 35°C, fixed, runs for 8 hours. PEMF: 5 Hz, runs for the first 30 minutes then stops automatically. The mat shuts off when the 8-hour cycle completes.

How It Works

5 Hz sits within the theta brainwave band — the frequency range associated with deep relaxation, creative states, and the transition into sleep. Running PEMF only at the start of the cycle gives the nervous system a downregulation signal where it matters most — the half hour before sleep — without continuing the stimulus through the night. 35°C matches the body's natural night-time core temperature, so the heat blends into rest rather than disturbing it.

P2 — CIRCULATION 15 HZ · 40°C · 30 MIN

Circulation is the everyday session — the one you'll come back to most. The answer to long flights, full days at a desk, the dull fatigue that comes from not moving enough.

Circulation runs at 15 Hz — the PEMF frequency with the deepest research base for supporting blood flow and lymphatic movement — paired with mid-range far-infrared heat at 40°C. The combination creates a passive recovery environment that helps the body do what it would otherwise do on its own — just more efficiently, and without the requirement to move.

Use it after travel, after long meetings, after a day on your feet, or anywhere along the fatigue curve. The shorter the session, the more flexible the placement — fit one in over lunch, after work, before dinner.

When to Use

Mid-afternoon, evening, after travel, or as part of a structured recovery week. Three to five sessions per week is a good baseline; daily during high-workload periods or long-haul travel.

Default duration: 30 minutes. Heat: 40°C, fixed. PEMF: 15 Hz, runs for the full session.

How It Works

15 Hz has been studied for its effects on red blood cell behaviour and nitric oxide release — both factors that influence how efficiently oxygen and nutrients reach soft tissue. The far-infrared heat does the parallel work mechanically

P3 — RECOVERY 30 HZ · 45°C · 30 MIN

Recovery is the post-training mode — the answer to sore muscles, heavy legs, and the deep fatigue that comes after a hard session. A higher frequency, a deeper heat, and a clear job to do.

At 30 Hz, the PEMF stimulus sits at the upper edge of the body's response window — a frequency associated with circulation, tissue repair, and physical readiness. Combined with far-infrared heat at 45°C, Recovery accelerates the body's natural post-load processes: clearing metabolic byproducts, easing muscle soreness, and shortening the gap between training sessions.

Use Recovery within one to two hours of training, after a competition, or in the evening after a physically demanding day. The earlier you use it after physical load, the more it helps.

When to Use

Within 1–2 hours of training, in the evening after a physically demanding day, or as part of a structured recovery week. Three to five sessions per week during heavy training blocks.

Default duration: 30 minutes. Heat: 45°C, fixed. PEMF: 30 Hz, runs for the full session.

How It Works

30 Hz is widely used in long-wave therapy for its association with circulation and tissue repair. The 45°C heat dilates peripheral blood vessels and reduces residual muscle tone. The two together create a strong but passive recovery stimulus — the kind of session that lets you skip the cool-down and head straight to dinner.

MANUAL MODE YOUR SETTINGS · 1–30 HZ · 30–70°C · 1–12 H

The presets cover most days. Manual mode is there when they don't.

Manual mode opens up the full operating range of Flowpempf IR Mat Pro. PEMF: dial anywhere from 1 to 30 Hz in single-Hz steps. Heat: set anywhere from 30 to 70°C. Timer: choose a session length from one to twelve hours.

When to Use Manual

When you know exactly what you want — a specific frequency, a precise temperature, a non-standard session length. If you are using the mat as part of a coach- or therapist-guided protocol, manual mode is where those settings live.

Frequency Notes

The full 1–30 Hz range corresponds to the body's response window. Generally: 1–4 Hz for deep rest and sleep preparation; 4–7 Hz for cool-down and emotional regulation; 8–12 Hz for daytime calm and balance; 13–20 Hz for activation and circulation; 20–30 Hz for post-training and pre-movement. The controller offers all 30 single-Hz settings — 30 PEMF programs in addition to the three named presets.

Heat Notes

30–40°C is sleep-friendly and suitable for long sessions. 40–55°C is the recovery and circulation range — comfortable for 20–45-minute sessions. 55–70°C is the deep heat range; use shorter sessions and start lower until you know how your body responds.

Build your own.

Manual mode is the mode that turns Flowpempf IR Mat Pro from a recovery tool into a recovery practice. Find the combination of heat and frequency that works for your body — and run it consistently.

A black and white photograph of a woman's neck and shoulder. The woman is wearing a small hoop earring and a dark, thin strap over her shoulder. The text "TEMPERATURE GUIDE" is overlaid in white, sans-serif capital letters. The background is a plain, light color.

TEMPERATURE GUIDE

Heat is a tool. The right temperature depends on what you want the session to do — and how long you plan to stay on the mat.

The four ranges below cover the full 30–70°C spectrum of Flowpempf IR Mat Pro. Start at the lower end of any range and increase gradually to comfort. Higher is not better.

35–40°C Extended Sessions and Overnight

The sleep-friendly range. Ideal for long sessions, evening wind-downs, and overnight use. Supports deeper sleep and full-body relaxation. Recommended for users new to heated mats or anyone with heat sensitivity.

45–55°C Sessions Up to 2 Hours

The comfort range. Helps temporarily relieve joint pain and stiffness, promote muscle relaxation, and give a short-term boost to local circulation. Where most people will find their natural session temperature.

55–65°C Sessions of 30–90 Minutes


A higher temperature for shorter, more intensive sessions. A noticeable step up from comfortable warmth — bring water. Use no more than twice a day at this range.

65–70°C Sessions of 30–60 Minutes (Max Once a Day)

The far-infrared sauna range. Retains heat for a sweating session that mimics the cardiovascular effects of a sauna. Wrap yourself in a towel to keep the warmth in. Use once a day, maximum. Always hydrate before, during, and after.

Heating time.

Allow up to 40 minutes for the mat surface to reach its highest temperature. For shorter sessions, preheat the mat for 15 to 20 minutes before lying down.

A black and white photograph of a person's legs in a dynamic, athletic pose. The legs are bent and crossed, with one foot flat on the ground and the other leg raised and bent. The skin is smooth and the lighting creates soft shadows, emphasizing the contours of the legs. The background is a plain, light-colored surface.

YOUR FIRST SESSIONS

What to expect when you're new to PEMF — and how to set up the first week so it sticks.

What PEMF Feels Like

Most people don't feel the electromagnetic field itself. Some describe a faint sensation of warmth, a gentle vibration, or a feeling of relaxation that settles in as the session progresses. There's no sharp signal — it's subtle, by design. Use the included PEMF tester — the small cylinder with magnetic balls inside — to see the field at work.

If You're Sensitive

A small number of people are highly sensitive to electromagnetic fields and may feel discomfort at higher frequencies. If that's you, start with shorter sessions at 2 Hz (in manual mode), build tolerance gradually, and step up only when the body settles into the lower frequencies.

What People Notice in the First Week

PEMF reaches the body at a cellular level, so some of the effects are subtle and some are surprisingly

- A faint tingling sensation during or just after a session — normal, particularly in the first few sessions.
- More frequent trips to the bathroom in the first week. This is the body's metabolic rate adjusting; it settles.
- A clearer sense of relaxation in the hour or two after a session.
- Slightly faster nail and hair growth over a few weeks of consistent use.
- Better sleep, especially on nights where you've used a 2 Hz or 8 Hz session before bed.

Setting Up Your First Week


Run three to four sessions in the first seven days. Start with P2 Circulation (15 Hz, 40°C, 30 min) once a day, midday or after work. Add P1 Sleep (5 Hz, 35°C, PEMF 30 min, heat 8 h) overnight on a few of those days. Don't go above 50°C in week one. Drink water before and after each daytime session.

If You're Treating a Specific Area

For lower-back or spinal discomfort, lie on the mat with your knees bent or supported by a pillow. This places the spine in a neutral position and brings it into full contact with the mat. For other areas, position yourself so the area you're targeting sits directly above one of the PEMF coils — they run the full length of the mat.

Alternate frequencies

Once you're comfortable with the presets, try alternating between Sleep, Circulation, and Recovery across your week. The body responds well to varied stimulus — different sessions train different recovery pathways.

A black and white photograph showing the back and shoulder of a person. A hand is placed on their waist. The person is wearing a dark, strappy garment. The text "BUILDING A RECOVERY RITUAL" is overlaid in white, sans-serif capital letters.

BUILDING A RECOVERY RITUAL

Three to six daytime sessions a week, thirty minutes each. Add P1 Sleep overnight as often as you like. That's the only rule that matters. Everything else is personal.

The patterns below are starting points — not prescriptions. Use them to find your own rhythm.

Muscle Soreness and Post-Training

Muscle soreness after training is a natural response to micro-stress in muscle fibres. It's necessary for adaptation. It's also the thing that gets in the way of the next session.

Far-infrared heat at 45°C increases blood flow and oxygen delivery to stressed tissue. PEMF at 30 Hz adds a circulation-supporting electromagnetic stimulus on top. Together they help muscles relax, ease perceived stiffness, and reduce the time you spend feeling broken.

Pattern:

P3 Recovery, 30 minutes, within two hours of training. Three to five times per week during heavy training blocks.

Circulation and Tissue Recovery

Reduced circulation is a quiet contributor to stiffness, fatigue, and slow recovery — especially after prolonged sitting, long-haul travel, or repetitive movement.

Far-infrared heat at 40°C increases local blood flow through vasodilation. PEMF at 15 Hz supports microcirculation and fluid movement in tissue. The combination creates a passive recovery environment that doesn't require active movement to work.

Pattern:

P2 Circulation, 30 minutes, mid-afternoon or evening. Three to five times per week — daily during high training or workload periods.

Stress and Nervous System Regulation

Chronic stress keeps the nervous system in a low-grade state of alert. Over time it shows up as muscle tone you can't release, shallow breathing, poor recovery, and broken sleep.

Sleep mode is built to interrupt that pattern. The first 30 minutes deliver a 5 Hz PEMF downregulation cue, then the gentle 35°C warmth continues for the rest of the cycle, providing a steady, low-grade signal that it is safe to rest.

Pattern:	
	P1 Sleep — start in the evening and leave running overnight. The PEMF eases the first half hour; the heat carries you through the night.

Deep Sleep

Poor sleep is often not about being tired enough — it's about not being able to transition into rest. A body that's still activated can't fall asleep, no matter how hard it tries.

Used overnight, Sleep mode delivers a 30-minute PEMF downregulation cue at 5 Hz at the start of the cycle, then maintains a gentle 35°C far-infrared warmth — matched to the body's natural night-time temperature — for the remaining 8 hours. The mat shuts off automatically when the cycle ends.

Pattern:	
	P1 Sleep, overnight, on or under the sheet. Daily, or as needed during disrupted sleep periods.

FREQUENTLY ASKED QUESTIONS

Where can I use the mat?

On any flat, stable surface — a bed, a massage table, the floor, or on top of a yoga mat. Make sure the mat lies completely flat with no folds, and that the controller cable isn't pinched or compressed.

How often can I use it?

Daily use is safe and well-tolerated at low and moderate heat settings. Most users find a rhythm of three to six sessions per week. At higher temperatures (above 55°C), use no more than once a day.

Can I lie directly on the mat?

Yes. The mat is designed for direct skin contact through light clothing or a thin towel. For high-heat sessions or extended use, place a clean cotton towel or sheet over the surface to absorb sweat and protect the mat.

What should I wear?

Light, breathable clothing made of natural fibres — cotton, bamboo, linen. Remove jewelry, watches, and metal accessories before a session.

Do I need to drink water?

Yes — especially at higher temperatures or for longer sessions. Drink a glass of water before lying down and another after. Far-infrared heat encourages mild sweating, and hydration is part of the recovery.

Can I use the mat if I have a pacemaker or other implanted device?

No. PEMF therapy is contraindicated for anyone with an implanted electronic device, including pacemakers, defibrillators, insulin pumps, and cochlear implants. The full list of contraindications is in the Manual section of this booklet.

How do I clean the mat?

Power off and unplug, then wipe the surface with a soft, dry — or slightly damp — cloth. For routine surface drying after a sweat session, run a high-heat setting for 15–20 minutes to evaporate moisture, then power down. Never use solvents, alcohol, or harsh cleaners.

If I have a medical condition, can I use the mat?

If you have a chronic condition, heart disease, hypertension, diabetes, or any serious medical concern, speak with your doctor before using Flowpemp IR Mat Pro. The Manual section lists the full set of conditions that require medical consultation before use.


Manufacturer:
Flowlife AB
Runebergsgatan 7
114 30 Stockholm
Sweden

Contact:
hello@flowlife.com
www.flowlife.com

Manual

Flowpemp IR Mat Pro

1. Flowpemp IR Mat Pro
2. Controller
3. Connection plug

ENGLISH

INTRODUCTION

Welcome to your Flowpemp IR Mat Pro. This guide will help you use your mat safely and effectively. Read all instructions before first use and retain this manual for future reference.

In this document, Flowpemp IR Mat Pro is referred to as the product. The product must be used in accordance with the following safety instructions to reduce the risk of fire, electric shock, property damage, and personal injury. Read all safety instructions before using Flowpemp IR Mat Pro.

SAFETY INSTRUCTIONS

SAFETY HAZARDS — REDUCE THE RISK OF ELECTRIC SHOCK

- Do not use the product near water or in wet environments
- Never immerse the mat, controller, cable, or plug in water or any other liquid
- Use only the original Flowlife power cable supplied with the product
- Disconnect immediately if the mat, controller, cable, or plug appears damaged or wet
- Never store the product or cable on hot surfaces
- Unplug the product when not in use, before cleaning, and before folding or moving
- Do not modify, repair, or open the mat, controller, or cable
- Do not handle the plug with wet hands

DANGER — REDUCE THE RISK OF BURNS, FIRE, OR INJURY

- Do not leave the product unattended while powered on
- Do not cover the mat with thick blankets, duvets, or heavy objects during use
- Do not fold, bunch, or crease the mat while it is powered on
- Do not insert pins, needles, or sharp objects into the mat or controller
- Do not stand, jump, or place heavy items on the mat
- Stop use immediately if you detect unusual heat, a burning smell, or visible damage
- Never attempt to repair the product yourself — contact support@flowlife.com
- Never use where aerosol products are used or where oxygen is administered
- Allow the mat to cool fully before folding, moving, or storing
- Keep out of reach of children under 15
- Children aged 15 and above may use the product only when properly informed about safe usage and under adult supervision

HEAT SAFETY

- Always start at a lower heat setting and increase gradually to comfort
- Do not use the mat if you are unable to feel heat or pain in the area of contact
- Stop use if any part of the surface feels hotter than the rest, or if the heat becomes uncomfortable

SKIN & CONTACT PRECAUTIONS

- Do not use on broken skin, open wounds, active infections, or areas of recent burns
- Clean and dry the skin before each session — remove lotion, oil, or topical medication
- Do not use the mat with creams or ointments containing heat-producing ingredients
- Do not use on skin that is sensitive to heat, such as the face
- Start with shorter sessions and increase gradually as the body adapts
- Stop use if redness, rash, blistering, or persistent discomfort occurs
- Allow at least 60 minutes between sessions on the same area

RISK GROUPS — CONSULT A DOCTOR BEFORE USE

- Pregnancy, or trying to become pregnant
- Pacemaker, ICD, insulin pump, cochlear implant, or any active implanted electronic device
- Metal or other body implants in the area of contact with the mat
- Active or suspected cancer, or undergoing cancer treatment
- Cardiovascular disease, severe arrhythmia, or any unstable heart condition
- High or low blood pressure
- Diabetes or impaired blood circulation
- Epilepsy or any history of seizures
- History of stroke or transient ischemic attack
- Severe bleeding disorders or recent blood loss
- Fever or active infection — heat may worsen hyperthermia
- Surgery in the past six months
- Recent fractures or open wounds in the area of contact
- Skin conditions, including eczema, psoriasis, or open lesions in the area of contact
- Peripheral neuropathy or reduced sensation in the area of contact — check skin carefully before and after every session
- Hemophilia or anhidrosis

- Mental illness, or under the influence of alcohol or substance use
- Photosensitising medications, or medication affecting circulation, blood pressure, or pain perception
- Under 15 years old — use only under supervision of a healthcare professional
- Any other health concern not listed — consult your doctor before starting

INSTRUCTIONS FOR USE

WHAT IS IN THE BOX

- 1 × Flowpemp IR Mat Pro recovery mat
- 1 × Flowmat controller + 1 × Power cable
- 1 × Carry bag
- 1 × PEMF tester
- 1 × User manual

BEFORE YOUR SESSION

- Review this manual before first use
- Remove jewelry, watches, and metal accessories
- Wear clean, thin, comfortable clothing made of natural fibres such as cotton
- Drink a glass of water
- Inspect the mat, controller, and cable for any wear, fraying, or damage — do not use if damaged

SETTING UP THE MAT

- Place the mat on a flat, stable, heat-resistant surface — a bed, the floor, or a yoga mat
- Make sure the mat lies flat with no folds or creases
- Plug the controller cable into the connector jack on the mat
- Plug the power cable into a standard grounded outlet (220–240 V)
- Place the controller on a heat-resistant surface beside the mat — never on top of the mat

USING THE CONTROLLER

- Press the Power button to turn on the heat and time settings
- Press Mode to cycle through the three preset programs — P1 Sleep, P2 Circulation, P3 Recovery
- Or set heat, PEMF, and timer manually:
- Press Temp, then use ▼ / ▲ to set the temperature (30 to 70°C)
- Press Time, then use ▼ / ▲ to set the session length (1 to 12 hours)
- Press PEMF, then use ▼ / ▲ to set the frequency (1 to 30 Hz)
- Changes apply immediately. The session starts automatically when all setting names stop blinking and glow steady
- Lie down on the mat — the mat warms gradually
- Allow up to 40 minutes for the mat to reach the highest temperature settings

PRESET PROGRAMS

- P1 Sleep — 5 Hz · 35°C · PEMF 30 min, heat 8 h — overnight wind-down and sleep support
- P2 Circulation — 15 Hz · 40°C · 30 min — daytime reset and circulation support
- P3 Recovery — 30 Hz · 45°C · 30 min — post-training recovery

MANUAL MODE

- PEMF: 1 to 30 Hz, adjustable in single-Hz steps
- Heat: 30 to 70°C, adjustable
- Timer: 1 to 12 hours, adjustable

TEMPERATURE GUIDELINES

- 35–40°C — overnight or extended sessions, ideal for users with heat sensitivity
- 45–55°C — sessions up to 2 hours, the comfortable working range for most users
- 55–65°C — sessions of 30–90 minutes, no more than twice a day
- 65–70°C — sessions of 30–60 minutes maximum, once per day, hydrate before, during, and after

DURING YOUR SESSION

- Stay still, or move gently if you wish
- Settings can be adjusted at any time — changes take effect immediately
- Drink water if using high heat or long sessions
- Hold the PEMF tester close to the mat to see the magnetic balls move — confirmation the field is active

AFTER YOUR SESSION

- The mat shuts off automatically when the timer ends
- Press the Power button to turn the controller off if required
- Drink water
- Allow the body to cool for at least 15 minutes before showering or exercising
- Allow the mat to cool fully before folding or moving
- Unplug the power cable from the mains and the controller cable from the mat
- Check the skin for any redness or irritation
- Wipe the mat surface with a soft, dry cloth

SESSION LENGTH AND FREQUENCY

- P2 Circulation and P3 Recovery: 30 minutes each
- P1 Sleep: PEMF active for the first 30 minutes, heat continues for 8 hours, automatic shutoff at the end of the cycle
- Manual mode: 1 to 12 hours per session
- Recommended use: 3 to 6 sessions per week
- Daily use is suitable for most users at low to moderate settings
- At settings above 55°C, do not use more than once per day

CLEANING & STORAGE

- Power off the mat and unplug the power cable before cleaning
- Allow the mat to cool fully before cleaning, folding, or storing
- Wipe the surface with a soft, dry cloth after each session
- For deeper cleaning, wipe with a slightly damp cloth — dry fully before use
- Do not immerse the mat or controller in water or any other liquid
- Do not use alcohol, solvents, bleach, or abrasive cleaners on the mat or controller
- Do not machine wash the mat
- Allow the mat to cool fully before folding
- Fold the mat only along its built-in fold lines
- Store the mat flat, or loosely folded inside the carry bag
- Store in a cool, dry place at room temperature, away from direct sunlight
- Do not store at temperatures below -10°C or above 45°C
- Do not place heavy objects on top of the mat during storage
- Inspect the cable before each use for cuts, fraying, or exposed wires — do not use if the cable is damaged.

RECYCLING & DISPOSAL

- Do not dispose of this product as unsorted household waste
- Return the product to a designated collection point for electronic equipment recycling
- Follow local regulations for the disposal of electronic and electrical equipment
- Contact your local authority for recycling locations and procedures

TECHNICAL SPECIFICATIONS

PRODUCT NAME	Flowpemf IR Mat Pro	P1 SLEEP	5 Hz · 35°C · PEMF 30 min, 8 h
COLOUR	Black	P2 CIRCULATION	15 Hz · 40°C · 30 min
DIMENSIONS	80 × 190 cm	P3 RECOVERY	30 Hz · 45°C · 30 min
WEIGHT	19 kg	CONTROLLER	Seven-button interface (Power · ▼ · ▲ · Temp · Time · PEMF · Mode)
CONSTRUCTION	17 functional layers	SAFETY	High-stability temperature sensor with automatic overheat shutoff
SURFACE	Non-toxic synthetic PU leather	CABLE LENGTH	180 cm
MINERALS	Natural amethyst and black tourmaline	INPUT POWER	220–240 V · 50 Hz
PEMF Coils	8 copper coils, full-body coverage	RATED POWER	280 W
PEMF Frequency	1–30 Hz, adjustable in single-Hz steps	OUTPUT CURRENT	1.3 A
HEATING SYSTEM	Non-EMF silicone-resistance heating wire	OPERATING TEMPERATURE	5°C — 35°C ambient
HEAT RANGE	30–70°C, adjustable	STORAGE TEMPERATURE	-10°C — 45°C
FAR-INFRARED WAVELENGTH	5–14 µm	HUMIDITY	10% — 85% non-condensing
TIMER	1–12 h, adjustable	CERTIFICATIONS	CE · FCC · RoHS · UKCA · FDA
PROGRAMS	30 PEMF settings + 3 preset programs	CE LVD Number	B-S210335561
		CE EMC Number	B-E210335562

WARRANTY

Flowpemf IR Mat Pro is covered by a 2-year warranty from the date of purchase, covering material and manufacturing defects. The warranty does not cover damage caused by misuse, improper storage, water exposure, unauthorised repair or modification, or use outside the instructions in this manual. Proof of purchase is required for all warranty requests. For support, contact support@flowlife.com.

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