
Flowroller Ball Go

FLOWLIFE

Brochure

Flowroller Ball Go

Vibration Therapy Ball

Table of Content

- 4 – Manifest
- 6 – Why Flowroller Ball Go
- 7 – Accessories
- 8 – Flowroller Ball Go - A Comprehensive Guide
- 10 – Starting your Journey
- 14 – For Muscle Soreness
- 16 – For General Recovery
- 18 – During Your Workout – Maintain Mobility
- 20 – After Your Workout – Recover and Repair
- 22 – Contact Information
- 24 – Manual
 - English
 - Swedish / Svenska
 - Danish / Dansk
 - Norwegian / Norsk
 - French / Français
 - German / Deutsch
 - Spanish / Español
 - Finnish / Suomi
 - Dutch / Nederlands

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.

For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.


Why Flowroller Ball Go?

Compact vibration therapy for targeted relief, deep muscle activation, and total mobility control.

- **Powerful vibration in a palm-sized tool** - Flowroller Ball Go delivers up to 3600 RPM of deep tissue vibration in an 80 mm format. It helps unlock tight muscles and stimulate circulation in small and hard-to-reach areas. With four intensity levels, you can adapt the pressure to suit your recovery needs. The compact shape makes it ideal for the feet, glutes, calves, and scapula.
- **Focused muscle activation and myofascial release** - Unlike traditional massage balls, Flowroller Ball Go combines vibration with static pressure to stimulate muscle response and support trigger point release. It helps increase range of motion, ease deep-seated tension, and accelerate warm-up and recovery.
- **A stall force that holds up under pressure** - With a stall force of 100 kilograms, the motor continues to vibrate even when compressed against a wall or floor. You achieve consistent output throughout the entire session, without dropouts or interruptions.
- **Smart timer and durable construction** - The built-in 10-minute shutoff helps prevent overstimulation, keeping treatments safe and efficient. The structure is impact-resistant and designed for long-term use. Portable size, lasting performance - Flowroller Ball Go weighs just 281 grams and fits in your hand, gym bag, or desk drawer. The battery lasts up to 10 hours, depending on intensity, and recharges via USB-C in under three hours.
- **A precision tool for conscious recovery** - Flowroller Ball Go is designed for users who want targeted recovery without compromise. It gives you complete control over pressure, depth, and location, wherever you are.


Accessories

Enhance your Flowroller Ball Go experience with essential accessories designed for maximum convenience and care.


Functional Storage Pouch

Flowroller Ball Go comes with a soft, drawstring storage pouch designed for both protection and function. The dual-thread system allows you to hold or suspend the device during use, which can be especially helpful when treating hard-to-reach areas. For users who find floor or wall pressure too intense, the pouch allows for manual positioning and control of the ball—without requiring a hard surface.

Silicone Charging Port Protector


An integrated silicone cover protects the USB-C charging port from moisture, dust, and debris. This helps ensure long-term safety and performance during charging. The port is easily accessible, and the cover remains in place throughout everyday use.

Type C Charging Cable

A compact USB-C cable is included for fast, efficient charging. The battery reaches full charge in approximately three hours and supports up to 10 hours of use, depending on vibration level. The universal port makes charging simple, whether you are at home, in the gym, or on the move.

All accessories are designed to support flexible, portable, and consistent recovery without compromise.

Flowroller Ball Go: A Comprehensive Guide

This guide will help you unlock the full benefits of your Flowroller Ball Go. With four calibrated vibration levels, an 80 millimeter diameter, and a stall force of up to 100kg, it is designed for high-performance recovery in a compact format. Whether you are treating trigger points, preparing for training, or managing daily mobility, this tool gives you clinical precision and total control.

Preparing for Your Session

Flowroller Ball Go is suitable for pre-activity warm-ups, post-training recovery, or as part of a structured mobility routine. It can be used directly on skin or through clothing, depending on preference and sensitivity. For best results, use the device on a stable surface such as a floor, yoga mat, or wall. If you prefer a softer or more controlled application, the included drawstring pouch allows you to position and move the ball without applying bodyweight.

Before use, ensure that the silicone charging port protector is securely in place to keep dust, sweat, and moisture away from the charging contact.

Powering On and Starting Your Session

Press and hold the power button for two seconds. The motor will activate at Level 1. Tap the power button to cycle through the four intensity levels. On the fifth press, the device will shut off completely. The motor will also turn off automatically after ten minutes of continuous use.

Positioning the Device

Place the Flowroller Ball Go under the target muscle or fascia point. Use your body weight to apply pressure as needed. Adjust your position slowly to maintain consistent contact for areas such as the glutes, hamstrings, or upper back. The drawstring pouch allows manual guidance without excessive pressure. Always start with a shorter session and gradually increase duration and intensity based on your tolerance.

Adjusting Vibration Intensity

Flowroller Ball Go features four speed levels, each with a defined RPM range. These levels are engineered to support different stages of recovery, from warm-up to deep tissue release.

Level 1 1500 RPM

Low intensity for warm-up, circulation, and use on sensitive areas, such as the feet, shoulders, or neck. Recommended for first-time users.

Level 2 2200 RPM

Moderate intensity for general muscle relaxation and improved blood flow. Suitable for daily maintenance and full-body mobility routines.

Level 3 2900 RPM

High intensity for deep muscle activation and myofascial release. Effective for hamstrings, glutes, calves, and postural tension around the spine.

Level 4 3600 RPM

Maximum intensity for experienced users and advanced recovery needs. Ideal for treating stubborn trigger points or delayed onset muscle soreness (DOMS) after intense physical activity.

Recommended Usage and Session Duration

5 to 10 minutes

Recommended for pre-workout activation or quick recovery, targeting specific areas.

10 to 15 minutes

Ideal for post-training recovery, full-body mobility work, or general tension release.

15 minutes
(maximum per area)

Extended deep tissue sessions should be used with care. Avoid overuse of the same area within a single session.

The built-in ten-minute auto shut-off function ensures safe usage and helps regulate treatment duration.

Ending Your Session

To power off the device, press and hold the button until all lights turn off and the motor stops. If you reach Level 4, one additional press will also power it down. Allow the device to cool for a few minutes before storing. Wipe clean with a dry or slightly damp cloth if used on bare skin or in a high-sweat setting.

Battery and Charging Instructions

Charging the Device

Insert the included USB-C cable into the charging port by gently lifting the silicone cover. After charging, always reseal the port by pressing the silicone cover firmly into place. This ensures protection from dust, moisture, and wear over time. Click the button with a fast click to get a indication of the battery level with a constant light.

One light

0 to 25
percent charge

Two lights

25 to 50
percent charge

Three lights

50 to 75
percent charge

Four lights


75 to 100
percent charge

To maintain battery health, recharge the device at least once a month if it is not in regular use. The device cannot be operated while charging to preserve battery integrity.

Storage and Maintenance

Store the Flowroller Ball Go in its protective box or fabric pouch, away from direct sunlight and heat. Do not expose the device to wet environments. Avoid dropping the unit or applying excessive force to the housing. Clean the surface with a non-abrasive cloth. Always allow the device to cool completely before storing it after extended use.

For Muscle Soreness

What happens in the muscles.

After high-intensity exercise, muscle fibers develop microscopic tears as part of the body's natural repair process. This triggers inflammation and swelling, which restricts mobility and creates discomfort known as delayed onset muscle soreness. Blood flow increases to deliver oxygen and nutrients to damaged tissue, but recovery can be slowed by localized tension, pressure buildup, and poor circulation.

Without intervention, soreness can limit joint movement, alter posture, and cause compensatory strain in nearby muscle groups. This not only prolongs recovery but also increases the risk of imbalance and injury over time.

Why Flowroller Ball Go helps

Flowroller Ball Go uses concentrated vibration therapy to accelerate muscle recovery in small, high-load areas. By combining high RPM vibration with focused pressure, it supports recovery through several physiological pathways.

- Increases local blood flow, improving oxygen and nutrient delivery to healing tissue.
- Relieves tightness by relaxing overactive muscle fibers and reducing protective tension.
- Flushes out lactic acid and other waste products through improved lymphatic drainage.
- Stimulates sensory nerves to override pain signals and promote neuromuscular relaxation.

How to use Flowroller Ball Go for muscle soreness

Start by placing the device under or against the sore muscle group. For larger zones such as the glutes or hamstrings, begin with slow, controlled pressure to allow the vibration to take effect. Move the ball gradually, no more than a few centimeters per second.

Pause on tender or tight spots for 10 to 20 seconds. This gives the tissue time to respond and release tension. Use the pouch if full bodyweight feels too intense, or when targeting harder-to-reach areas.

Begin with Level 2 or 3 and increase to Level 4 only if soreness allows. A session of 10 to 15 minutes is typically sufficient to reduce stiffness and restore movement.

Flowroller Ball Go helps reduce soreness faster by treating the areas that matter most. The compact design allows you to manage recovery on your own terms, restoring mobility and preparing you for your next session.

For General Recovery

What happens in the body

Even outside of training, muscles accumulate low-level tension from everyday stressors, such as prolonged sitting, poor posture, and repetitive movements. Over time, this leads to reduced circulation, shortened muscle fibers, and inefficient movement patterns. These changes often go unnoticed until they begin to affect performance, posture, or mobility.

Muscles that are not regularly activated or released tend to lose elasticity. This limits their ability to generate force, adapt to movement, or recover from fatigue. In the long run, this contributes to stiffness, joint discomfort, and an increased risk of strain or imbalance.

Why Flowroller Ball Go helps

By integrating Flowroller Ball Go into your daily routine, you can support long-term muscle health and prevent passive tension from accumulating.


- Stimulates blood flow in high-load and sedentary zones, keeping tissue oxygenated and reducing fatigue.
- Releases low-grade muscular tension that builds from static posture and repetitive movement.
- Maintains soft tissue pliability and helps preserve joint range of motion.
- Supports movement quality and postural balance by restoring muscle responsiveness and improving overall joint function.

How to use the Flowroller Ball Go for general recovery

Place the device on or against the muscle groups most affected by daily activity—commonly the glutes, lower back, calves, feet, or upper traps. Use slow, steady contact with light to moderate pressure, allowing the vibration to work passively through the tissue.

Select Level 2 or 3 for a calming, restorative effect. Spend 5 to 10 minutes per session, moving between areas as needed. If preferred, use the drawstring pouch to adjust position without loading the tissue directly.

Making Flowroller Ball Go part of your daily recovery practice helps reduce background tension, improves circulation, and supports long-term mobility and muscular function.

During Your Workout – Maintain Mobility

What happens in the body

As muscles contract repeatedly during training, they build up local tension and metabolic waste. This can limit movement, restrict joint range, and reduce force output. In high-volume or technique-driven sessions, even mild tightness can impact form and control.

Why Flowroller Ball Go helps

Brief use between sets helps prevent mobility loss and keep movement efficient.

- Relieves accumulating tightness to preserve the range of motion
- Stimulates local circulation to reduce early fatigue
- Releases surface tension that may disrupt technique or posture

How to use Flowroller Ball Go mid-workout

Apply the device to the active muscle group for 10 to 15 seconds between sets if stiffness appears. Use Level 1 or 2 to maintain circulation without causing muscle inhibition. Keep the contact light and dynamic. This is not a recovery session but a mobility reset.

Using Flowroller Ball Go mid-workout helps maintain quality in your movement, so you can train longer and more effectively.

After Your Workout — Recover and Repair

What happens in the body

Following training, muscles enter a recovery phase to repair microtears and restore energy. This process involves inflammation, swelling, and temporary stiffness as the body mobilizes resources to damaged tissue. Without active recovery, these responses can limit flexibility, prolong soreness, and reduce performance in subsequent sessions.

Why Flowroller Ball Go helps

Post-workout application of vibration therapy helps accelerate the recovery process and reduce discomfort.

- Improves blood flow to damaged tissue, supporting nutrient delivery and cellular repair
- Flushes metabolic waste like lactic acid to reduce soreness and swelling
- Relaxes overactive muscle tone and reduces residual tension after exertion

How to use Flowroller Ball Go after a workout

Use the device on fatigued or sore muscle groups within 30 minutes of completing your session. Start at Level 2 or 3 for circulation and tissue response. For areas with deeper tension, increase to Level 4 if comfortable. Target each area for 1 to 2 minutes, with a total session time between 10 and 15 minutes.

Position the device on the floor or wall to apply steady pressure, or use the drawstring pouch for a more adjustable and lighter approach. Allow the vibration to penetrate without rushing the process.


This method supports faster recovery, reduces post-exercise discomfort, and prepares your body for the next training session.

hello@flowlife.com
www.flowlife.com
Luntmakargatan 90
113 51 Stockholm

Manual

Flowroller Ball Go

ENGLISH

INTRODUCTION

Welcome to your Flowroller Ball Go. This guide will help you get started and ensure you use the device safely and effectively for recovery, mobility, and performance.

Flowroller Ball Go delivers compact, targeted vibration therapy in an 80 mm format, with up to 3600 RPM of deep-tissue stimulation and a stall force of 100 kg. Whether you're warming up, treating sore muscles, or maintaining mobility, Flowroller Ball Go gives you precision recovery on your own terms.

SAFETY INSTRUCTIONS

Safety Hazards – Reduce the Risk of Electric Shock

- Never use the product while it is plugged into an electrical outlet
- Do not expose the product to water or use near bathtubs, showers, or pools.
- Store in a dry place, away from liquids and high humidity.
- If the product comes into contact with water, disconnect immediately.

Danger – Reduce the Risk of Burns, Fire, and Injury

- Do not leave the product unattended while in operation.
- Only use accessories provided or approved by Flowlife.
- Discontinue use if the product, charger, or cable is damaged.
- Do not insert foreign objects into openings
- Do not use the device near aerosol sprays or oxygen therapy equipment.
- Not recommended for individuals under 15 years without supervision.
- Do not apply excessive pressure on sensitive areas (head, front of neck).
- Stop use immediately if discomfort, dizziness, or irritation occurs.

Risk Groups – Consult a Doctor Before Use

- If you have a fever, infection, or inflammation.
- If you suffer from osteoporosis, epilepsy, or neurological conditions.
- If you are pregnant, have a pacemaker, heart disease, or diabetes complications.

INSTRUCTIONS FOR USE

Powering On & Off

- Press and hold the power button for 2 seconds to turn on/off.
- The LED indicators show battery status at startup.

Adjusting Speed Levels

- Tap the power button to cycle through 4 vibration levels:
 - Level 1: 1500 RPM – Gentle warm-up & sensitive areas.
 - Level 2: 2200 RPM – General relaxation & mobility.
 - Level 3: 2900 RPM – Deep muscle activation & release.
 - Level 4: 3600 RPM – Maximum intensity for experienced users.
- On the fifth press, the device shuts off.

Auto Shut-Off

- Flowroller Ball Go turns off automatically after 10 minutes of continuous use for safety.

RECOMMENDED USAGE

- 5–10 minutes: Quick warm-up or recovery on specific areas.
- 10–15 minutes: Post-training recovery or mobility routine.
- Max 15 minutes per area: Avoid overstimulation.

Let the vibrations do the work—avoid applying unnecessary pressure.

BATTERY & CHARGING

- Use the included USB-C cable.
- Charging time: ~3 hours.
- Up to 10 hours of use depending on intensity.
- LED battery indicators:
 - 1 light = 0–25%
 - 2 lights = 25–50%
 - 3 lights = 50–75%
 - 4 lights = 75–100%

Always reseal the silicone charging port protector after charging.

TECHNICAL SPECIFICATIONS

- Product Size: Diameter 80 mm
- Product Weight: 281 g ± 5 g
- Package Size: 115 × 110 × 90 mm
- Power Input: 100–240V ~ 50/60 Hz
- Power Output: DC 5V=2A
- Battery Capacity: 2400 mAh
- Charging Port: USB-C
- Charging Time: 2.5–3 hours
- Auto Shut-Off: 10 minutes
- Battery Life (Use Time):
 - Level 1: 10 hours
 - Level 2: 8 hours
 - Level 3: 6 hours
 - Level 4: 4 hours
- Noise Level: ≤ 55 dB
- Stall Force: 100 kg
- Power Consumption: ~10 W
- Speed Levels: 4 adjustable settings
- Speed Range: 1500–3600 RPM
 - Level 1: 1500 RPM
 - Level 2: 2200 RPM
 - Level 3: 2900 RPM
 - Level 4: 3600 RPM
- Accessories: USB-C charging cable, storage pouch, storage box

CLEANING & STORAGE

- Wipe with a soft, damp cloth after use.
- Do not submerge in water or use harsh cleaners.
- Store in a cool, dry place, away from direct sunlight.
- Use the included pouch or storage box for protection.
- Charge at least once per month if not used regularly.

WARRANTY & SUPPORT

Flowroller Ball Go comes with a 2-year warranty covering manufacturing and material defects. Exclusions:

- Damage due to misuse, accidents, or improper handling.
- Unauthorized repairs or modifications.
- Storage in unsuitable environments.
- Use not in accordance with this manual.
- Lack of proof of purchase.

For support, contact Flowlife customer service.

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