# IMAGE XPLOERER

## UPUTA ZA KORIŠTENJE

### SADRŽAJ

|        |                                                           |   |
|--------|-----------------------------------------------------------|---|
| I.     | ODRICANJE.....                                            | 3 |
| II.    | ODRICANJE ODGOVORNOSTI U VEZI S JEZIKOM.....              | 3 |
| III.   | IZJAVA O ODGOVORNOSTI .....                               | 3 |
| IV.    | POJMOVI I DEFINICIJE .....                                | 4 |
| V.     | PREDVIĐENA NAMJENA IMAGEXPLOERER .....                    | 5 |
| VI.    | PROIZVOĐAČ.....                                           | 5 |
| VII.   | PODACI O IZVEDBI .....                                    | 5 |
| VII.1  | KALIBRACIJA TESTA .....                                   | 5 |
| VII.2  | MJERNI RASPON .....                                       | 5 |
| VII.3  | KONTROLA KVALITETE .....                                  | 5 |
| VII.4  | ANALIZA PODATAKA .....                                    | 6 |
| VII.5  | REZULTATI .....                                           | 6 |
| VII.6  | OGRANIČENJE POSTUPKA.....                                 | 6 |
| VII.7  | OČEKIVANE VRIJEDNOSTI .....                               | 6 |
| VII.8  | KARAKTERISTIKE IZVEDBE .....                              | 6 |
| VIII.  | PRINCIP POSTUPKA.....                                     | 6 |
| VIII.1 | ALEX <sup>2</sup> , ALEX <sup>3</sup> PRINCIP TESTA ..... | 6 |
| VIII.2 | FOX PRINCIP TESTA.....                                    | 7 |
| IX.    | SERVIS .....                                              | 7 |
| X.     | JAMSTVO .....                                             | 7 |
| XI.    | INFORMACIJE O NARUČIVANJU .....                           | 8 |
| XII.   | SIGURNO RUKOVANJE .....                                   | 8 |
| XII.1  | KVALIFIKACIJA OPERATERA .....                             | 8 |
| XII.2  | RADNA SIGURNOST .....                                     | 9 |
| XII.3  | DEKONTAMINACIJA.....                                      | 9 |
| XIII.  | RJEČNIK SIMBOLA .....                                     | 9 |

|          |                                                                                                 |    |
|----------|-------------------------------------------------------------------------------------------------|----|
| XIV.     | SIGURNOSNE PORUKE .....                                                                         | 10 |
| XV.      | PRAVNI ZAHTJEVI .....                                                                           | 10 |
| XV.1     | MEĐUNARODNI STANDARDI .....                                                                     | 10 |
| XV.2     | CE SUKLADNOST .....                                                                             | 10 |
| XV.3     | ELEKTROMAGNETSKA KOMPATIBILNOST (EMC), SUZBIJANJE RADIO<br>SMETNJI I OTPORNOST NA SMETNJE ..... | 11 |
| XVI.     | ŽIVOTNI CIKLUS .....                                                                            | 11 |
| XVI.1    | DOSTAVA.....                                                                                    | 11 |
| XVI.1.1  | ŠTEĆENJA TIJEKOM TRANSPORTA.....                                                                | 11 |
| XVI.1.2  | OPSEG ISPORUKE.....                                                                             | 11 |
| XVI.2    | RASPOLAGANJE .....                                                                              | 11 |
| XVI.3    | SIGURNOSNA KOPIJA PODATAKA .....                                                                | 12 |
| XVII.    | OPIS.....                                                                                       | 13 |
| XVII.1   | POVEZIVANJE IMAGEXPLORERA S RAČUNALOM.....                                                      | 14 |
| XVII.2   | RAPTOR SERVER POSTAVLJANJE SOFTVERA .....                                                       | 15 |
| XVII.3   | PREUZIMANJE SOFTVERA IMAGEXPLORER AGENT I IMAGEXPLORER<br>DEFINICIJE.....                       | 16 |
| XVII.4   | PODEŠAVANJE POSTAVKI IMAGEXPLORERA .....                                                        | 19 |
| XVII.5   | SNIMANJE I ANALIZA MJERENJA .....                                                               | 25 |
| XVII.5.1 | UMETANJE SPREMNIKA U IMAGEXPLORER .....                                                         | 25 |
| XVII.5.2 | SNIMANJE SLIKE, ČITANJE QR KODA I PRONALAZENJE MREŽE.....                                       | 26 |
| XVII.6   | INTERNA KONTROLA KVALITETE.....                                                                 | 27 |
| XVII.7   | TEHNIČKA PODRŠKA .....                                                                          | 28 |
| XVII.8   | TEHNIČKI PODACI I ZAHTJEVI.....                                                                 | 28 |
| XVII.9   | ODRŽAVANJE.....                                                                                 | 29 |
| XVII.10  | RJEŠAVANJE PROBLEMA.....                                                                        | 30 |
| XVIII.   | PROMIJENI POVIJEST .....                                                                        | 32 |

## I. ODRICANJE

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**ImageXplorers**

Use this to configure ImageXplorers, or add new ones. Every ImageXplorer needs to be known to Raptor Server, otherwise it won't be able to connect.

[Manage ImageXplorers](#)

---

**RAPTOR** [New measurements](#) [Approved measurements](#) [Start measurement](#) [Automats](#)

Add new ImageXplorer

General information

Name:

Key:

This key is required when installing an Agent.

ImageXplorer SerialNumber:

Serialnumber of connected ImageXplorer.

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

In total there are 17 new measurements in the selected time period.

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[New measurements](#)
[Approved measurements](#)
[Start measurement with ImageXplorer](#)

The ImageXplorer IX training cannot be used, because it requires a ConfigXplorer scan.

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**RAPTOR**
[New measurements](#)
[Approved measurements](#)
[Start measurement with ImageXplorer](#)

## ConfigXplorer run report for ImageXplorer

[← Back to ImageXplorer](#)
[Download Report](#)
[Apply detected settings](#)

Please check the automatically detected values for plausibility before applying them to your device.

### Results

|                                        |                         |
|----------------------------------------|-------------------------|
| Device                                 | 22124137                |
| QR-Code                                | 30BAC11E24122024;53173; |
| Exposure                               | 1546                    |
| Offset X                               | 402                     |
| Offset Y                               | 370                     |
| Width                                  | 618                     |
| Height                                 | 288                     |
| Linearity Test [ $R^2$ ] (> 0.95)      | 0.999479013663686       |
| Homogeneity Test [CV] (< 5)            | 1.3903116236964623      |
| Duration                               | 00:00:42.1146637        |
| Completed                              | now                     |
| Acceptance Criteria ( $R^2$ and CV ok) | Passed                  |

Redovito pokretanje ConfigXplorer Scan će osigurati da ImageXplorer koristi optimalnu konfiguraciju. Stoga će se nakon 60 dana pojaviti poruka koja vas poziva da ponovite ImageXplorer ConfigXplorer skeniranje unutar 30 dana.

**Raptor**
[New measurements](#)
[Approved measurements](#)
[Start measurement](#)
[Automats](#)
Adliance

The last testarray scan for this ImageXplorer is more than 68 days ago. Please make sure to complete another ImageXplorer testarray scan in the next 22 days.

You requested a new measurement from ImageXplorer "Hannes' ImageXplorer". Please wait until the measurement is completed ...

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## Configure "Neuer IX "

## General information

|                                   |                                                                                                                           |
|-----------------------------------|---------------------------------------------------------------------------------------------------------------------------|
| Name:                             | Neuer IX                                                                                                                  |
| Key:                              | hOfOxvDXlxF1xrp4p9ttjmNfH9qait                                                                                            |
|                                   | <small>This key is required when installing an Agent.</small>                                                             |
| ImageXplorer Serialnumber:        | 23167844                                                                                                                  |
|                                   | <small>Serialnumber of connected ImageXplorer.</small>                                                                    |
| Notes:                            |                                                                                                                           |
| Internal Notes:                   |                                                                                                                           |
|                                   | <small>For internal use only.</small>                                                                                     |
| Disabled:                         | Not disabled                                                                                                              |
| Ignore ConfigXplorer scan checks: | <input checked="" type="checkbox"/> Enforce to ignore ConfigXplorer scan checks (required e.g. > 60 days) on this device. |

## QR-Code image configuration

|                   |      |
|-------------------|------|
| QR-Code Exposure: | 3000 |
|-------------------|------|

## Calibration configuration

|       |                        |
|-------|------------------------|
| Mode: | Use automatic exposure |
|-------|------------------------|

## Array image configuration

|                 |     |
|-----------------|-----|
| Array Offset X: | 362 |
| Array Offset Y: | 326 |
| Array Width:    | 634 |
| Array Height:   | 298 |
| Array Exposure: | 834 |

Konfiguracija kalibracije je prema zadanim postavkama postavljena na "Use automatic exposure (Koristi automatsku ekspoziciju)", koja ostaje za ekspoziciju od 2500. Međutim, može se postaviti na "Use manual exposure (Koristi ručnu ekspoziciju)". U ovom načinu rada, ekspozicija se može korigirati gore ili dolje.

## Add new ImageXplorer

### General information

**Name:**

**Key:**   
This key is required when installing an Agent.

**ImageXplorer Serialnumber:**   
Serialnumber of connected ImageXplorer.

**Notes:**

**Internal Notes:**   
For internal use only.

**Disabled:**

### QR-Code image configuration

**QR-Code Exposure:**

### Calibration configuration

**Mode:**

### Array image configuration

**Array Offset X:**

**Array Offset Y:**

**Array Width:**

**Array Height:**

**Array Exposure:**

**RAPTOR** [New measurements](#) [Approved measurements](#) [Start measurement](#)

## Add new ImageXplorer

### General information

**Name:** ImageXplorer

**Key:** U7gPsUGvNeP00306Kd9CupIKFYub4o  
This key is required when installing an Agent.

**ImageXplorer Serialnumber:**   
Serialnumber of connected ImageXplorer.

**Notes:**   
For internal use only.

**Internal Notes:**   
For internal use only.

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



General information

Name: IX IFU

Key: WQZQZif8xd7WBamER1MvTXtpRKv6y5  
This key is required when installing an Agent.

ImageXplorer Serialnumber: 40739900  
Serialnumber of connected ImageXplorer.

Notes:

Internal Notes:

Disabled: Not disabled

Require ImageXplorer ConfigXplorer scan: Warn after 60 days, enforce after 90 days.

Ignore ConfigXplorer scan checks: ☐ Enforce to ignore ConfigXplorer scan checks (required e.g. > 60 days) on this device.

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

Settings

Please enter the settings for your ImageXplorer. You can find your ImageXplorer key when you log in to Raptor. Please contact support if you're not sure what the proper settings are for your ImageXplorer.

Raptor URL:

ImageXplorer Key:

Continue

### Početak mjerenja nije moguć:

**RAPTOR** New measurements Approved measurements Start measurement with ImageXplorer

The ImageXplorer IX training cannot be used, because it requires a ConfigXplorer scan.

Go to ImageXplorer configuration

### Moguća rješenja:

- Provjerite je li ImageXplorer povezan
- izvršeno posljednje potrebno ConfigXplorer skeniranje

Ako i dalje imate problema ili imate bilo kakvih drugih pitanja o proizvodima i uslugama MacroArray Diagnostics, obratite se svom lokalnom distributeru.

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| Verzija | Opis                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Datum izdavanja | Zamjenjuje |
|---------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|------------|
| 01      | First compilation in this language; adaptation to the English IFU version 16<br>Prva kompilacija na ovom jeziku; prilagodba engleskoj verziji IFU-a 16                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 09-2024         | -          |
| 02      | Prilagodba engleskoj verziji Uputa za uporabu (IFU) 17.1: Dodana je jezična napomena za Odjeljak II; dodan je ALEX <sup>3</sup> i u cijelom dokumentu je uklonjena svaka referenca na ALEX <sup>2</sup> (50 testova); u odjeljku VI je uklonjena oznaka; QualityXplorer je uklonjen iz cijelog dokumenta; u odjeljku XII je dodana daljnja izjava u skladu s normom EN IEC 61326-2-6; uklonjene su informacije o specifičnim točkama smjera; snimke zaslona ažurirane su u cijelom dokumentu; u odjeljku XVIII je dodana cijela povijest promjena; XVII.6 Prilagodba klauzuli o ispitivanju kontrole kvalitete (identično u IX i MAX 9/45k) | 09-2026         | 01         |

©Copyright by MacroArray Diagnostics

MacroArray Diagnostics (MADx)

Lemböckgasse 59, Top 4

1230 Vienna, Austria

+43 (0)1 865 2573

[www.madx.com](http://www.madx.com)

Version number: 11-IFU-01-HR-02