A black and white profile photograph of a woman with her hair pulled back into a bun. She is looking upwards and to the right. The lighting is soft, highlighting the contours of her face and neck. The background is a plain, light color.

Flowlight Laser Mask Ultra+

FLOWLIFE

Brochure

Flowlight Laser Mask Ultra+

Red & Blue Light Therapy Mask

Table of Content

4 —	Manifest
6 —	Why Flowlight Laser Mask Ultra+?
7 —	Accessories
8 —	Flowlight Laser Mask Ultra+ – A Comprehensive Guide
10 —	Getting Started
12 —	Understanding the Modes
13 —	Using the Display Remote
14 —	Intensity Levels
15 —	Voice Control
16 —	For Sleep Quality
18 —	For Stress Relief
20 —	For Energy & Cognitive Clarity
22 —	For Skin Health
25 —	Contact Information
26 —	Manual
	• English
	• Svenska
	• Norsk
	• Dansk
	• Suomi
	• Deutsch
	• Polski
	• Français
	• Italiano
	• Español
	• Nederlands

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.

For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowlight Laser Mask Ultra+?

Four precision light modes. Four wavelengths. A new level of at-home light experience.

Four Targeted Wavelengths – 460nm Blue Light, 665nm Red Laser, 850nm Near Infrared, and 1064nm Near Infrared. Four wavelengths used and explored by people around the world as part of their daily routines.

Four Light Modes – Blue Light for skin clarity, Red Light for overall skin feel and tone, Near Infrared for deep relaxation, All for the full-spectrum experience.

Voice Control – Hands-free operation with the Hi Flowlight wake phrase. Set mode, intensity, and timer, or start and stop your session, without touching a button.

Display Remote Controller – The new wireless display remote shows your current mode, intensity, and timer countdown at a glance. Simple, clean, and precise.

Adjustable Intensity – Three brightness levels (High, Medium and Low) allow you to customise the depth and comfort of every session to your needs.


Flexible Timer – Choose between 5, 10, or 20-minute sessions. Short enough for a daily habit, long enough for a deep treatment.

Rechargeable Battery – 1800mAh built-in lithium-ion battery with USB-C charging. Approximately 40 minutes of runtime per full charge.

236 Laser-Emitting Diodes – 72 pcs 460nm blue Light + 72 pcs 665nm red laser + 72 pcs 850nm near-infrared laser + 20 pcs 1064nm near-infrared laser. Precision-placed across the full face and neck area.


Accessories

Everything you need to begin your light therapy practice on day one.

Flowlight Laser Mask Ultra+

236 laser-emitting diodes across four wavelengths: 72 pcs 460nm blue light, 72 pcs 665nm red laser, 72 pcs 850nm near-infrared laser, and 20 pcs 1064nm near-infrared laser. Ergonomic face-conforming design with adjustable head strap.

Wireless Display Remote Controller

The new rechargeable display remote with a built-in screen showing mode, intensity, and timer. Four dedicated buttons: Mode (top), Intensity (right), Timer (left), On/Off (bottom).

USB-A to USB-C Charging Cable

Compatible with any standard USB power adapter. Charges both the mask and the display remote.


Protective Carry Pouch

Slim, durable, and structured for storage or travel. Keeps the mask and accessories protected without bulk.

User Manual

Comprehensive safety and usage guide available in 11 languages: English, Svenska, Norsk, Dansk, Suomi, Deutsch, Polski, Français, Italiano, Español, and Nederlands.

A high-contrast, black and white close-up photograph of a person's face. The person is wearing a laser mask, which is visible as a dark, textured surface covering the lower half of the face. The skin on the upper half of the face is illuminated, showing fine details and texture. The overall mood is dramatic and technological.

Flowlight Laser Mask Ultra+: A Comprehensive Guide

This guide walks you through every feature of the Flowlight Laser Mask Ultra+, from your first session to advanced voice control. Read it once, then let the experience speak for itself.

Getting Started

Charging the Mask

Connect the included USB-C cable to the charging port on the side of the mask. A full charge gives approximately 40 minutes of runtime.

Charging indicator: Red flashing light = Charging. White light = Fully charged.

Always charge fully before your first session.

Wearing the Mask

Cleanse your face before use. Place the mask over your face and adjust the head strap for a comfortable, secure fit. The mask should rest gently against the skin without pressure.

Remove contact lenses and avoid applying serums or oils immediately before use.

Starting a Session

Hold the white button on the mask to power it on. The mask starts immediately in default mode: Red Light and Near Infrared at High intensity for 20 minutes.

To adjust your session, use the display remote to select mode, intensity, or timer. Or short-press the white button on the mask to open a voice control window and speak a command. The mask powers off automatically when the timer ends.

New to Light Therapy?

If you are just starting out, consider beginning with 5 or 10-minute sessions for the first week before moving up to 20 minutes. Many people find this helps their skin adjust to the higher brightness levels comfortably.

Most people find 3 to 4 sessions per week a good rhythm to settle into.

During Your Session

Relax and disconnect. Close your eyes and breathe slowly. Consider pairing with breathwork or a quiet room for the full experience.

Monitor comfort. If the light feels intense, reduce intensity or duration. Your comfort level is the right level.

After Your Session

The mask powers off automatically when the timer ends.

Clean the mask surface with a dry microfiber cloth after use. Store in the protective carry pouch to keep it clean and safe.

Recharge as needed via USB-C. The indicator light confirms charging status.

Understanding the Modes

The Ultra+ features four light modes, each combining specific wavelengths favoured by the light therapy community. Press the Mode button to cycle through them.

Blue Light

460nm Blue Light + 850nm Near Infrared

460nm blue is a wavelength widely favoured in the light therapy community for its surface-level skin experience. Combined with 850nm near-infrared, many people use this mode as part of their daily skin clarity routine.

Red Light

665nm Red Laser + 850nm Near Infrared + 1064nm Near Infrared

665nm is one of the most popular wavelengths in the light therapy community, widely used in skin and recovery routines. Paired with 850nm and 1064nm near-infrared, this is the mode most people reach for as their everyday go-to.

Near Infrared

850nm + 1064nm Near Infrared (invisible light)


Both near-infrared wavelengths are invisible to the human eye. This mode will appear inactive but is fully operational. 850nm and 1064nm are deeply valued in the light therapy community for their relaxation and recovery experience, often described as the most restorative mode.

All

460nm + 665nm + 850nm + 1064nm, full spectrum

All four wavelengths simultaneously. The full-spectrum experience: everything the Ultra+ has to offer in one session. Popular among people who want to cover every wavelength in a single routine.

Using the Display Remote

The wireless display remote has four dedicated buttons and a circular display showing your current mode, intensity, and timer at a glance.

Mode Button (M)

Top button. Cycles through the four light modes: Blue Light → Red Light → Near Infrared → All.

Intensity Button

Right button. Cycles through brightness levels: High → Medium → Low. Default on mode activation is always High.

Timer Button

Left button. Cycles through session durations: 5 minutes → 10 minutes → 20 minutes. Default on mode activation is 20 minutes. The display shows remaining time during a session.

On / Off Button

Bottom button. Hold to wake the remote and start your session. Hold again to pause the session and turn off the remote. The mask powers off automatically when the timer reaches zero.

Intensity Levels

Three brightness levels let you adjust the session to your comfort and preference.

High

Maximum output. The setting most people prefer once they have found their comfort level.

Medium

Balanced output. A good starting point for daily use.

Low

Gentle output. A comfortable starting point for those new to light therapy or who prefer a softer experience.

Note: when any light mode is activated (via remote or voice) the device defaults to 20 minutes and High intensity. Adjust with the Intensity button or voice command if needed.

If your skin feels sensitive to higher brightness, starting on Medium or Low is a perfectly valid approach. Most people move up to High naturally over time as they get used to the experience.

Timer Settings

Select a session length to fit your schedule and routine. The mask powers off automatically when the timer ends.

5 min

A quick daily maintenance session. Good for mornings when time is limited.

10 min

The standard session length for most use cases

20 min

Full treatment session. Default on mode activation. The length most people settle into once they have found their rhythm.

A general guide used by many in the light therapy community: 3 to 4 sessions per week tends to feel like a sustainable and effective routine. Daily use is fine too. Listen to your skin and adjust as needed.

Voice Control

Control the mask entirely hands-free. Short-press the white button on the mask to open a 1-minute voice window, then speak any command below.

Wake Phrase

"Hi Flowlight"

Opens the voice control session. Speak this first after pressing the white button. The mask responds: Hi.

Light Modes

"Red light on" / "Blue light on" / "Infrared on" / "All light on"

Activates the named mode at default settings (20 min, High). The mask confirms with: Red on / Blue on / Infrared on / All on.

Intensity

"High intensity" / "Medium intensity" / "Low intensity"

Sets the brightness level. The mask confirms: High set / Medium set / Low set.

Timer


Say: "Five" / "Ten" / "Twenty"

Sets the session timer. The mask confirms: Five minutes / Ten minutes / Twenty minutes.

Session Control

Say: "Start" / "Stop" / "Goodbye"

Start begins the session. If no mode was set first, the mask defaults to Red Light (665nm + 850nm + 1064nm). Stop ends it early. Goodbye exits voice control mode. The mask responds: Starting / Turning off.

For Sleep Quality

Many people in the light therapy community report that near-infrared sessions in the evening have become a favourite part of their wind-down routine. Consistent use is where people say they feel the difference.

The Session


Use the mask 20–30 minutes before bed. Lie down in a dimly lit or dark room. Breathe slowly and let the session become part of your wind-down ritual. The Near Infrared mode operates invisibly. The absence of visible light makes it ideal for pre-sleep use.

Consistency

People who get the most out of the mask tend to use it regularly. Daily sessions, even short ones, add up. Build it into your evening the same way you would skincare.

Recommended Settings

Mode: Near Infrared
Intensity: Low or Medium
Timer: 20 minutes

A high-contrast, black and white close-up photograph of a person's face, focusing on the right eye. The eye is light-colored and looking directly at the camera. The skin is textured, and the lighting is dramatic, with deep shadows and bright highlights. The text "For Stress Relief" is overlaid in white, sans-serif font across the lower half of the image.

For Stress Relief

A dedicated moment to yourself, with warm light and no distractions. People in the light therapy community consistently describe their sessions as one of the simplest ways to reset during or after a demanding day.

Create the Space

Even a 10-minute session is effective when approached intentionally. Sit or lie comfortably, close your eyes, and let the session serve as a dedicated pause.

Timing

Stress relief sessions work well at any point in the day: mid-afternoon, post-work, or before any demanding task. Evening sessions double as preparation for sleep.

Recommended Settings

Mode: Red Light
Intensity: Medium
Timer: 10–20 minutes

A black and white photograph showing the back and shoulder of a person. The skin is smooth and the lighting creates soft shadows, highlighting the contours of the body. The text is overlaid in the center-left area.

For Energy &
Cognitive Clarity

Among the most commonly shared experiences in the light therapy community: morning sessions with red light described as starting the day with a clearer head and a stronger sense of readiness. Many have made it the first thing they reach for.

Morning Routine

Use the mask as part of your morning ritual, before coffee or directly after waking. Even a 5–10 minute session is enough to feel the difference in how you approach the day.

Stacking

Pair it with movement, hydration, and natural light for the full morning stack. The mask fits neatly into any routine without adding much time.

Recommended Settings

Mode: Red Light or All
Intensity: High
Timer: 10–20 minutes

For Skin Health

From daily skin maintenance to targeted routines, the Ultra+ covers the full wavelength range that the light therapy community favours for skin. Blue for surface clarity, red for overall skin feel and tone, near-infrared when you want to go deeper.

Skin Preparation

Always start with a clean, dry face. Remove makeup and apply no active serums or retinoids immediately before a session. Light therapy is most effective on bare skin.

After the Session

Follow with your usual skincare routine. Many people find their skin feels more receptive immediately after a session, a good moment for serums or moisturisers.

Recommended Settings

Skin clarity: Blue Light mode, High intensity, 20 minutes
Skin feel and tone: Red Light mode, High intensity, 20 minutes
Deep relaxation and recovery: Near Infrared or All mode, Medium–High, 20 minutes


Manufacturer:
Flowlife AB
Runebergsgatan 7
114 30 Stockholm
Sweden

Contact:
hello@flowlife.com
www.flowlife.com

Manual

Flowlight Laser Mask Ultra+

1. Adjustable head strap
2. Mask
3. On/Off
4. USB Charging Port
5. Breathing Vents

6. Mode Selection
7. Timer
8. Light Intensity
9. On/Off

ENGLISH

INTRODUCTION

- Read this manual carefully before use.
- Keep this manual for future reference.
- This product is intended for personal, non-medical use only.

SAFETY INSTRUCTIONS

DANGER — ELECTRIC SHOCK

- Do not immerse the mask or remote in water or any liquid.
- Do not use in wet environments such as bathrooms or near swimming pools.
- Do not attempt to disassemble or repair the device yourself.
- Keep away from children under 14 years of age.
- Inspect the charging cable before each use. Do not use if damaged.

DANGER — LASER / LIGHT INJURY

- Do not look directly into the diodes while the mask is active.
- Always close your eyes fully during a session. Do not use without eye protection if you are photosensitive.
- Near Infrared mode (850nm + 1064nm) is invisible to the human eye. The mask may appear inactive but is emitting light. Do not remove the mask mid-session without first pressing Stop.
- Blue light (460nm) is visible and bright. Avoid prolonged direct eye exposure.
- Do not use on broken, irritated, or sunburned skin.
- Keep the mask flat and stationary during use. Do not drop or impact the diode panel.
- Do not cover the mask with fabric or other materials during a session.
- Do not leave the mask operating unattended.
- Do not use near oxygen or aerosol products.
- If you experience discomfort, pain, or excessive heat, stop immediately and remove the mask.

RISK GROUPS

Do not use this product if you belong to any of the following groups without first consulting a medical professional:

- Pregnancy
- Heart disease or pacemaker use
- Active cancer or light-sensitive conditions
- Epilepsy or migraines triggered by light
- Recent facial surgery or open wounds
- Active skin infections or severe rosacea
- Medication that causes photosensitivity
- Diabetes with complications such as neuropathy
- Any other concerning medical conditions

Medication note: Certain medications, including blood thinners, anti-inflammatory drugs, and medications that affect circulation or sensation, may increase the risk of adverse reactions from light therapy. If you take regular medication, consult your doctor or pharmacist before use.

RECYCLING AND SAFETY

- This product must not be disposed of with regular household waste.
- Dispose of at a certified electronic waste recycling facility.
- The lithium-ion battery must be removed before disposal. Contact Flowlife for guidance.
- Do not incinerate or puncture the battery.
- Store the product in a cool, dry location away from direct sunlight.

INSTRUCTIONS FOR USE, MAINTENANCE AND STORAGE

- Charge the mask fully before first use using the included USB-C cable. Full charge takes approximately 3 hours.
- Charging indicator: Red flashing light = Charging. White light = Fully charged.
- Always start with clean, dry skin. Remove all makeup before use.
- Select your mode, intensity, and timer with either the remote control or voice commands.
- For voice control: short-press the white button on the mask to open a voice window, then speak your command.
- Default session settings on mode activation: 20 minutes, High intensity.
- If you are new to light therapy, consider starting with 5 or 10-minute sessions for the first week before moving up to 20 minutes. This gives your skin time to adjust at its own pace.
- A general guideline: 3 to 4 sessions per week. Daily use is also fine. Adjust based on how your skin feels.
- Sessions of up to 20 minutes per day are generally what most users stay within. If your skin feels sensitive, reduce duration or intensity.
- If high brightness feels intense on your skin, try Medium or Low. Many people find these levels equally enjoyable.
- Clean the mask surface with a soft, dry or slightly damp cloth. Do not use alcohol or chemical cleaners on the diode panel.
- Store in the original packaging or a protective pouch when not in use.
- Keep away from extreme heat, cold, or humidity during storage.
- The mask is for single-user use only.

WARRANTY

If your product breaks within two years from the date of purchase due to material or manufacturing defects, Flowlife will repair or replace the necessary components or the product free of charge.

Flowlife AB is not responsible for warranty or returns of products purchased through other retailers.

The warranty does not apply if:

- the product is damaged due to misuse, improper handling, accident, or transport
- the product has been repaired by an unauthorized service provider
- the product has not been used according to Flowlife's instructions
- the damage exceeds the value of the product
- the product has deteriorated due to improper storage
- a dated proof of purchase cannot be provided

PACKAGE CONTENTS

- Flowlight Laser Mask Ultra+
- Wireless Display Remote Controller
- USB-A to USB-C Charging Cable
- Protective Carry Pouch
- User Manual (11 languages)

TECHNICAL SPECIFICATIONS

Model	Flowlight Laser Mask Ultra+
Wavelengths	460nm (Blue Light), 665nm (Red Laser), 850nm (Near Infrared Laser), 1064nm (Near Infrared Laser)
Light Modes	Blue Light / Red Light / Near Infrared / All
Diodes	236 total -- 72 pcs 460nm + 72 pcs 665nm + 72 pcs 850nm + 20 pcs 1064nm
Intensity	High / Medium / Low
Timer	5 / 10 / 20 minutes
Output Power	6W
Battery	1800mAh lithium-ion, built-in
Charging	USB-C, 5V/2A -- approx. 3 hours to full charge
Operating time	Approx. 40 minutes per full charge
Voice control	Yes -- wake phrase: Hi Flowlight
Remote	Wireless display remote, rechargeable via USB-C
Wireless range	Up to 10 metres
Dimensions	640 × 220 × 8.8mm
Weight	370g
Certifications	CE, RoHS, FCC, CB, IE
Beam Angle	20°
Irradiance	Approx. 90 mW/cm ² at skin surface
Charging Input	AC 100-240V / DC 5V-1A (standard USB adapter)

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Flowlife AB is niet verantwoordelijk voor de garantie of retouren van producten die bij andere retailers zijn gekocht.