A black and white photograph of a person's midsection, showing the waist and hips. The person is wearing a dark, form-fitting athletic top and shorts. The lighting is soft, highlighting the contours of the body. The text is overlaid on the image.

Flowplate Vibe Pro

FLOWLIFE

Brochure

Flowplate Vibe Pro

Whole-Body Vibration Platform by Flowlife

Table of Content

4	—	Manifest
6	—	Why Flowplate Vibe Pro
7	—	In the Box
8	—	Flowplate Vibe Pro: A Comprehensive Guide
10	—	Getting Started
12	—	Speed Settings
14	—	Vibration Modes & Timer
16	—	For Muscle Recovery
18	—	For General Recovery
20	—	Before Your Workout
22	—	Position Is Treatment
24	—	Manual
		• English • Svenska • Norsk • Dansk • Suomi • Deutsch • Français • Italiano • Español • Nederlands

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.

For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowplate Vibe Pro?

Unlock whole-body recovery with a vibration platform engineered for precision, adaptability, and therapeutic depth. The Flowplate Vibe Pro is built for those who take recovery seriously — delivering oscillating vibration therapy that concentrates treatment exactly where your body needs it most.

- **30-Level Precision Control.** Fine-grained speed adjustment across 30 levels — from the gentlest activation to intense therapeutic vibration. Use the included remote control to change intensity without stepping off the platform mid-session.
- **Oscillating Vibration Platform.** The pivotal oscillation mechanism delivers alternating stimulation to each side of the body, mimicking natural walking mechanics and engaging muscles from feet through the core.
- **Position Is Treatment.** Standing concentrates therapy in the legs, glutes, and core. Placing hands or forearms on the surface redirects treatment to the shoulders, arms, and upper back. One platform — the entire body.
- **Wireless Remote Control.** Four-function remote: Power, Speed+, Speed-, Mode. Adjust intensity and cycle vibration modes without stepping off. A single beep confirms each press.
- **Four Vibration Modes.** Constant, Dynamic Peaks, Burst, and Wave — four distinct rhythm patterns from steady therapeutic vibration to slow rolling waves. Each mode targets a different physiological response.
- **10-Minute Auto-Shutoff.** The Flowplate Vibe Pro runs a fixed 10-minute session and powers down automatically — no timer button to press, no session to manage. Just step on and go.
- **Non-Slip Platform Surface.** A textured standing surface ensures secure footing across all 30 speed levels, including the highest intensities.
- **Quiet 200W Motor.** The high-torque motor operates below 50 dB — quiet enough for home use at any hour, and powerful enough to deliver effective therapy at every speed level.

In the Box

Flowplate Vibe Pro comes complete with everything required for immediate use.

Flowplate Vibe Pro Vibration Platform

A full-body oscillating vibration platform with 30-level precision speed control and a textured non-slip surface. A beep confirms each speed change.

Remote Control


Compact wireless remote with four buttons: Power, Speed+, Speed-, Mode. Requires 2 x AAA batteries (not included).


Power Cable

Standard IEC cable. Connect to the side port and flip the red toggle switch to power on.

User Manual

Multilingual guide in 10 languages: EN, SV, NO, DA, FI, DE, FR, IT, ES, NL.

Flowplate Vibe Pro: A Comprehensive Guide

Explore each speed level and session type to discover the combinations that work best for your body — the right approach makes a meaningful difference to how quickly and completely you recover.

This essential guide will help you maximize the recovery benefits and therapeutic performance of your Flowplate Vibe Pro. Take the time to explore each speed level and session type, and discover the combinations that work best for your body — the right approach can make a meaningful difference to how quickly and completely you recover.

Getting Started

Setting Up

Place Flowplate Vibe Pro on a flat, stable, non-slip surface with adequate space on all sides. A rubber mat beneath the platform reduces vibration transfer to hard floors.

Connect the IEC power cable to the port on the side of the unit. Flip the red toggle switch to the ON position. The unit is ready when it beeps.

Starting Your Session

Press the Power button on the remote control. The platform starts vibrating at the last-used speed level and a beep confirms the start.

Use Speed + / Speed – to adjust intensity — each press is confirmed by a beep. The unit runs for 10 minutes and shuts down automatically. Press Power at any time to stop early.

Standing Position

Stand with feet shoulder-width apart, knees slightly bent. This stance engages the leg muscles and allows oscillation to travel from feet through the lower body and into the core.

Avoid locking your knees. Soft, relaxed joints transmit the vibration more effectively into the muscle tissue.

Upper Body Position

Sit in a chair in front of the platform and place both hands, forearms, or wrists flat on the surface. The vibration transmits directly into the arms, wrists, elbows, and shoulders.

Vary the position — palms flat, knuckles down, or forearms resting — to target different areas of the upper limb.

Speed Settings

Control vibration intensity across 30 precision levels, ranging from 2 to 13 Hz. Start low and increase gradually to find the settings that produce the best therapeutic effect for your body.

- | | |
|--------------------|--|
| Level 1–5 | Minimal vibration — ideal for warming up joints, activating the lymphatic system, and first-time use. Suitable for those returning from injury or with high sensitivity. |
| Level 6–15 | Gentle to moderate vibration — promotes circulation, reduces morning stiffness, and provides a light warm-up. The recommended daily maintenance range for most users. |
| Level 16–25 | Medium to firm vibration — active recovery and deep muscle engagement. Effective for post-workout recovery, reducing DOMS, and stimulating the neuromuscular system. |
| Level 26–30 | Maximum intensity — for experienced users. Delivers the strongest whole-body stimulation at 13 Hz for advanced recovery and deep muscle treatment. |

Vibration Modes

Cycle through four modes using the remote. Each mode delivers a different rhythm — from steady constant vibration to slow, rolling waves. A beep confirms each mode change. Choose the pattern that best matches your body's needs in the moment.

Mode 1 — Constant

Steady, uninterrupted vibration at your chosen intensity level. The most direct and predictable therapeutic input. Ideal for focused recovery, first-time users, and any session where consistent stimulation matters most.

Mode 2 — Dynamic Peaks

Short rhythmic cycles between low and high intensity — the vibration rises to a peak, drops, and climbs again. Prevents the body from adapting to a fixed rhythm and keeps the muscles actively responding. Good for warm-up and general daily recovery.

Mode 3 — Burst


Rapid, high-intensity pulses with quick cycling between peaks. More aggressive stimulation pattern that challenges the neuromuscular system. Best for experienced users, post-workout recovery, and sessions where deep muscle activation is the goal.

Mode 4 — Wave

Long, gradual rises and falls — a slow, rolling wave pattern. The gentlest dynamic mode. Promotes relaxation, supports parasympathetic recovery, and works well for cool-down sessions and evening use.

Session Timer

The unit runs a fixed **10-minute session** and shuts down automatically with a beep. No timer button to press — the session length is built in. Press Power at any time for an immediate stop.

For Muscle Recovery

Muscle soreness following training is the result of micro-damage and inflammation in the tissue. The most effective response is to increase local circulation, promote lymphatic drainage, and reduce neuromuscular tension. Vibration therapy does all three simultaneously.

Standing Recovery

After leg day, running, cycling, or any lower-body training, stand on the platform at medium-to-high speed. The oscillating motion stimulates the quadriceps, hamstrings, calves, glutes, and surrounding connective tissue — accelerating clearance of metabolic by-products and reducing next-day stiffness.


Upper Body Recovery

After upper body training or manual work, sit and place forearms or hands on the platform. The vibration travels directly into the wrists, forearms, biceps, and shoulders.

Recommended Settings

Speed: Level 16–25 — medium-firm vibration reduces inflammation and breaks down lactic acid build-up.


Timer: 10 minutes — let the auto-shutoff end the session naturally for maximum therapeutic effect.

For General Recovery

Recovery does not only happen after intense training. The cumulative effect of standing, walking, and everyday movement places continuous load on the body. Used as a daily routine, the Flowplate Vibe Pro actively reduces accumulated tension and supports healthy circulation.

Enhanced Circulation	Gentle vibration combined with oscillating motion stimulates blood flow throughout the lower body, facilitating oxygen and nutrient delivery essential for tissue repair and long-term musculoskeletal health.
Lymphatic Activation	The rhythmic muscle contractions produced by vibration therapy act as a mechanical pump for the lymphatic system — supporting waste clearance, reducing fluid retention, and enhancing the body's natural recovery processes.
Recommended Settings	Speed: Level 6–15 — gentle, maintenance-level vibration without inducing fatigue. Timer: 10 minutes — a daily habit that compounds meaningful results over time.

Before Your Workout

A pre-workout vibration session raises tissue temperature, lubricates the joints with synovial fluid, and activates the neuromuscular system — preparing the body to perform at a higher level from the first rep.

Neuromuscular Activation

Vibration at medium intensity activates motor units throughout the lower body simultaneously — a level of muscle recruitment that static stretching alone cannot achieve. The result is a body ready to contract powerfully and safely from the first movement.

Joint Preparation

The oscillating motion stimulates synovial fluid production and distribution in the knee, hip, and ankle joints, reducing friction and the risk of impact-related injury during high-load training.

Recommended Settings

Speed: Level 10–20 — activating without fatiguing. Increase to Level 20+ for high-intensity sessions.

Timer: The first 5–10 minutes of the 10-minute session is the ideal warm-up window. Step off once you feel activated.

Position Is Treatment

The Flowplate Vibe Pro is not a fixed-function device. The area of the body in contact with the platform determines where the therapeutic effects are concentrated. One platform — the entire body.

Standing — Lower Body

Both feet on the platform, knees soft. Vibration travels into the feet, calves, hamstrings, quadriceps, glutes, hip flexors, and lower back.

Best for: leg recovery, lower back relief, general warm-up, circulation, lymphatic drainage.

Single-Leg — Targeted

One foot on the platform, the other on the floor. Concentrates vibration for targeted recovery or tension on a specific side.

Best for: unilateral recovery, calf and hamstring isolation.

Arms / Hands — Upper Body

Sit and place both hands or forearms flat on the surface. Vibration travels through the wrists, forearms, elbows, biceps, triceps, and shoulders.

Best for: upper body recovery, desk tension, forearm flexibility.


Calves Elevated — Passive

Sit and rest both calves on top of the platform. Targets the lower legs and promotes drainage from the feet upward.

Best for: end-of-day foot and calf recovery, reducing swelling.

Manual

Flowplate Vibe Pro

ENGLISH

Flowplate Vibe Pro USER MANUAL ENGLISH

- In this document, Flowplate Vibe Pro is referred to as the product or Flowplate Vibe Pro.
- The product must be used in accordance with the following safety instructions.
- This reduces the risk of fire, electric shock, property damage, and personal injury.
- Read all safety instructions before using Flowplate Vibe Pro.

SAFETY INSTRUCTIONS

DANGER – REDUCE THE RISK OF ELECTRIC SHOCK

- Never keep the product connected to a power outlet or other power source while cleaning it.
- Never use the product near water, such as in a bathtub, shower, pool, or similar environment.
- Never place the product anywhere it could fall into water.
- Always store the product in a dry place, away from contact with water or liquids.
- Immediately disconnect the product from power if it comes into contact with water.

DANGER – REDUCE THE RISK OF BURNS, FIRE, ELECTRIC SHOCK AND INJURY

- Never leave the product unattended when it is connected to power.
- The product must never be covered with fabric or similar material during use. This may cause overheating, fire, electric shock, or injury.
- Provided they have been instructed on safe use and the related risks, children from 14 years of age, individuals with functional limitations, and inexperienced users may use the product.
- Never use the product for any purpose other than that described in these safety instructions. The product is intended solely for whole-body vibration exercise and therapy. Do not use accessories not supplied by the manufacturer.
- Never use the product if it is visibly damaged, or if the cable or power source is damaged.
- If the product has been in contact with water, dropped, or otherwise damaged, contact customer support for inspection.
- Do not use the product until you have read all customer support information and recommendations.
- Never use the cord as a carrying handle.
- Never store the product or the cable on hot surfaces.
- Do not use the product in areas where aerosols are used or where oxygen is administered.
- Always stand securely on the platform during use. Keep knees slightly bent. Do not sit or lie on the vibrating platform.
- Children must not use, clean, or maintain the product without adult supervision.
- Never attempt to repair the product yourself. Do not open or disassemble it.
- In the event of obvious or suspected damage, always contact customer support.
- If the product becomes excessively warm, turn it off and allow it to cool before further use.
- Do not exceed the maximum user weight of 150 kg.

RISK GROUPS

- If you belong to any of the following risk groups, ALWAYS consult your doctor before using Flowplate Vibe Pro.
- Use the product only after approval from a healthcare professional:

Pregnancy, heart disease, medication*, diabetes with complications such as neuropathy or retinal damage, pacemaker use, recent surgery or injury, migraines or epilepsy, disc herniation, vertebral displacement, spondylosis, recent joint replacements, intrauterine devices, metal implants, osteoporosis, acute thrombosis, severe varicose veins, or any other concerning medical conditions.

- If unsure, consult a doctor.

Medication note: Certain medications, including blood thinners, anti-inflammatory drugs, and medications that affect circulation or sensation, may increase the risk of adverse reactions from vibration therapy. If you take regular medication, consult your doctor or pharmacist before use.

RECYCLING AND SAFETY

- Connect and disconnect the cable carefully.
- Flowplate Vibe Pro must be disposed of at an approved electronic waste facility.
- If Flowplate Vibe Pro is exposed to extreme conditions or handled improperly, it may pose a fire hazard.
- Do not expose Flowplate Vibe Pro to external temperatures above 45°C.
- Never disassemble the product. Children must not operate the product unsupervised.

INSTRUCTIONS FOR USE, MAINTENANCE AND STORAGE

- Place Flowplate Vibe Pro on a flat, stable, non-slip surface.
- Connect the IEC power cable to the power input port on the side of the unit and plug into a standard wall outlet. Flip the red toggle power switch to the ON position. The unit will beep once when ready.
- Press the Power button on the remote control to start a session. The unit beeps once to confirm start.
- Use the Speed + / Speed – buttons on the remote control to adjust vibration intensity (Level 1–30). Each press is confirmed by a single beep.
- The unit runs a 10-minute session and shuts down automatically when complete.
- Stand with feet shoulder-width apart, knees slightly bent. For upper body treatment, sit and place forearms or hands on the platform surface.
- Press the Power button at any time for an immediate shutdown.
- Wipe the platform surface with a damp cloth after each use. Do not use abrasive cleaners or solvents.
- Allow the unit to cool completely before storing. Store in a dry place, away from direct sunlight.
- Do not submerge the product in water. The product is not waterproof.

WARRANTY

- If your product breaks within two years from the date of purchase due to material or manufacturing defects, Flowlife will repair or replace the necessary components or the product free of charge.
- Flowlife AB is not responsible for warranty or returns of products purchased through other retailers.

The warranty does not apply if:

- the product is damaged due to misuse, improper handling, accident, or transport
- the product has been repaired by an unauthorized service provider
- the product has not been used according to Flowlife's instructions
- the damage exceeds the value of the product

- the product has deteriorated due to improper storage
- a dated proof of purchase cannot be provided

PACKAGE CONTENTS

- 1 x Flowplate Vibe Pro Vibration Platform
- 1 x Remote Control (requires 2 x AAA batteries, not included)
- 1 x IEC Power Cable
- 1 x User Manual

TECHNICAL SPECIFICATIONS

Model Flowplate Vibe Pro

Vibration Type Oscillating (Pivotal)

Motor Power 200W

Speed Levels 1–30

Frequency Range 2–13 Hz

Amplitude 9 mm

Vibration Modes 4 (Constant / Dynamic Peaks / Burst / Wave)

Session Timer 10 min (auto-shutoff)

Max User Weight 150 kg

Power Connection IEC 3-pin cable, 110V / 220V, 50/60 Hz

Power Switch Red toggle switch (on unit)

User Feedback Audible beep (confirms each button press)

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