



LYPHOCHEK INTERNAL QUALITY CONTROL

Lyphochek® Allergen IgE

Panel A

Lot 22680

Release date of expected values: 07.01.2026

DESCRIPTION:

The Lyphochek Allergen IgE, Panel A (Lot 22680) by Bio-Rad Laboratories, Inc. can be used as an internal quality control material to monitor the test performance of ALEX², ALEX Air, ALEX Food or ALEX³ in your lab. The control should be processed like a standard patient serum sample according to the Instructions for Use of the respective product.

Please find expected values and acceptance ranges in the table below:

Allergen	Code	Target Value (kU _A /l)	Expected Min (kU _A /l)	Expected Max (kU _A /l)
Bet v 1	t215	2,78	1,20	5,30
Fel d 1	e94	9,90	4,9	16,9
Pen m 2	f545	4,44	1,50	8,50
Pru p 3	f420	2,04	0,30	4,00

IMPORTANT NOTE:

Lyphochek Allergen IgE, Panel A (Lot 22680) shows more positive values than listed in this document. Allergens from represented allergen sources/families on ALEX², ALEX Air, ALEX Food or ALEX³ were chosen.