

Use only for ALEX²!

Lyphochek® Allergen IgE

Panel A

Lot 22650

Release date of expected values: June 17th, 2022

DESCRIPTION:

The Lyphochek Allergen sIgE, Panel A (Lot 22650) by Bio-Rad Laboratories, Inc can be used as an internal quality control material to monitor the Allergy Explorer² (ALEX²) test performance in your lab.

The control should be processed like a standard patient serum sample according to the ALEX²

Instructions for Use. Please find expected values and acceptance ranges in the table below:

Allergen	Code	target Value (kU _A /l)	Expected Min (kU _A /l)	Expected Max (kU _A /l)
Alt a 1	m229	9.20	5.00	14.38
Art v 1	w231	0.00	0.00	0.30
Asp f 3	m220	1.44	0.30	2.73
Bet v 1	t215	2.84	1.02	5.17
Bos d 2	e100	0.00	0.00	0.30
Bos d 8	f78	1.62	0.50	2.70
Can f 1	e101	5.22	2.65	8.78
Can f 2	e102	4.24	1.72	7.41
Can f 3	e221	7.99	3.79	12.38
Der p 1	d202	4.12	1.58	6.76
Der p 2	d203	9.01	4.64	14.87
Fel d 1	e94	6.60	3.07	10.99
Fel d 2	e220	2.90	1.19	4.94
Phl p 1	g205	5.55	2.30	10.56
Phl p 5.0101	g215	5.40	2.65	8.65
Phl p 7	g210	4.72	2.47	7.87
Pru p 3	f420	1.03	0.30	1.96