



LYPHOCHEK INTERNAL QUALITY CONTROL

Lyphochek® Allergen IgE

Panel A

Lot 22670

Release date of expected values: December 18th, 2024

DESCRIPTION:

The Lyphochek Allergen sIgE, Panel A (Lot 22670) by Bio-Rad Laboratories, Inc. can be used as an internal quality control material to monitor the test performance of ALEX², ALEX Air or ALEX³ in your lab. The control should be processed like a standard patient serum sample according to the Instructions for Use of the respective product used. This product is not compatible with ALEX Food.

Please find expected values and acceptance ranges in the table below:

Allergen	Code	Target Value (kU _A /l)	Expected Min (kU _A /l)	Expected Max (kU _A /l)
Bet v 1	t215	2,22	0,30	4,30
Der p 1	d202	5,12	2,20	9,20
Fel d 1	e94	7,74	2,70	14,50
Phl p 1	g205	4,65	1,90	10,00

IMPORTANT NOTE:

Lyphochek Allergen sIgE, Panel A (Lot 22660) shows more positive values than listed in this document. Allergens from the four most prevalent allergen sources were chosen.