



macroarraydx.com

MAD
MACRO ARRAY DIAGNOSTICS

FOX
EXPLORE FOOD
INTOLERANCE THE
COMPREHENSIVE WAY

FOX TEST FOR 286 FOOD ANTIGENS

- ALL-IN-ONE ARRAYS
- MEASUREMENT OF sIgG
- KIT WITH 50 ARRAYS

BENEFITS

- Covers all relevant food groups in one assay run
- Result interpretation supported by RAVEN
- Individualised dietary recommendations included

KEY FEATURES

ALLERGEN PANEL

- ELISA-based IgG multiplex assay
- 286 food antigens from 13 food groups
- Staple & novel foods included

TECHNOLOGY

- Nano-bead technology
- Detection of IgG antibodies (IgG subclasses 1-4)
- Semi-quantitative results
- Design your food antigen panel via on-demand multiplex feature



The connection between food intake, elevated IgG levels and chronic disorders has been described in peer reviewed publications and case studies. Nonetheless, this connection is still debated in the scientific community and a consensus has not been reached thus far.



MADx tests are processed and analysed by our proprietary hardware (MAX 45k, MAX 9k, ImageXplorer) and software systems (RAPTOR Server Analysis Software, RAVEN Interpretation Software).

EXPLORE MORE

80-MKT-05-EN-01