

Use only for ALEX²!

Lyphochek® Allergen IgE

Panel A

Lot 22640

Release date of expected values: February 16th, 2021

DESCRIPTION:

The Lyphochek Allergen sIgE, Panel A (Lot 22640) by Bio-Rad Laboratories, Inc can be used as an internal quality control material to monitor the Allergy Explorer² (ALEX²) test performance in your lab.

The control should be processed like a standard patient serum sample according to the ALEX²

Instructions for Use. Please find expected values and acceptance ranges in the table below:

Allergen	Code	target Value (kU _A /l)	Expected Min (kU _A /l)	Expected Max (kU _A /l)
Alt a 1	m229	11.89	5.21	25.31
Ara h 2	f423	0.85	0.30	1.97
Ara h 9	f427	1.08	0.30	1.78
Art v 1	w231	0.00	0.00	0.30
Asp f 3	m220	2.27	1.03	3.89
Bet v 1	t215	2.08	0.50	3.49
Bos d 2	e100	0.00	0.00	0.30
Bos d 8	f78	1.71	0.50	2.88
Can f 1	e101	7.03	3.17	11.84
Can f 2	e102	7.54	3.30	12.71
Can f 3	e221	9.94	4.77	15.65
Der p 1	d202	5.42	2.15	9.81
Der p 2	d203	8.98	4.19	14.38
Fel d 1	e94	8.60	4.27	21.50
Fel d 2	e220	4.49	1.97	7.87
Phl p 1	g205	5.23	2.11	10.83
Phl p 5.0101	g215	3.81	1.97	8.94
Phl p 7	g210	7.95	4.07	13.02
Pru p 3	f420	1.13	0.30	2.41