

Coloured soft contact lenses

Colouring of soft contact lenses is done only on ML-Soft designs made of Hydrogel material (ML 65 or ML75). Si-H material is not coloured.

COLOURING

Type	Design	Specification	
Colouring of the iris <i>Please always include a specification regarding if you would like an open or closed pupil. If open is selected, then please also specify the size.</i>	Standard colour. One of the following colours and intensity is selected:	Iris: Ø 8.0 - 16.0 mm (1.0 step) and Ø 10.5, 11.5*, 12.5 and 13.5 mm Pupil open: Ø 2.5 - 6.0 mm (0.5 step) including 7.0 mm (Ø 4.5*mm)	
	Blue <i>Light, Dark</i>		
	Green <i>Handling, Light, Dark</i>		
	Aqua <i>Handling, Light, Dark</i>		
	Brown <i>Handling, Light, Medium, Dark, Max**</i>		
	Grey <i>Handling, Light, Medium, Dark, Max**</i>		
	Filter colour / Not standar colour		
	Filter colour ML 450		
	Filter colour ML 550		
	Filter Colour ML Ruby		
	Mixed iris colour		
	Black Iris***		
	Iris opaque		
	Iris Black Opaque***		

*Standard diameter. If no diameter is specified on the order, this diameter will be used. NB! The difference between the color diameter and the total diameter of the lens must be min. 2 mm.

** Max refers to double colouring of Dark. Please note: The price for Max is the same as for Filter Colour / Non-Standard Colour.

*** Black lenses ordered without power are produced as standard with a power of -1.5 D to make them easier to handle.

Type	Design	Specification	
Advanced colouring according to photo	An iris print is colored according to your request (reference photo). Specify: Open/closed pupil and with or without white back surface. Contact us for advice when ordering for the time.	Ø 11.3 and 12.3 mm Spherical powers: -10.0 / +15.0	

Type	Design	Specification	
Colouring of pupil	Pupil, standard colour. One of the following colours is selected:	Ø 3.0 - 5.0 mm in 0.5 mm step. Ø 6.0 - 9.0 and 1.0 mm step (Standard diameter: 4.5 mm)	
	Brown <i>Handling, Light, Medium, Dark, Max</i>		
	Grey <i>Handling, Light, Medium, Dark, Max</i>		
	Black		
	Pupil, ML filter colour. One of the following filter colours is selected:		
	Filter colour ML 450		
	Filter colour ML 550		
	Filter Colour ML Ruby		
	Pupil, black opaque		

*Standard diameter. If no diameter is specified on the order, this diameter will be used.


** Lenses with a black pupil ordered without power are produced as standard with a power of -1.5 D to make them easier to handle.

±5% colour deviation may occur

Add on	
Limbal ring	
Pupilring	
Coloured white back surface	
White back surface x 2	
Black dot (used as a marker)	
Decentration (1,0 mm steps) incl. prism ballast	

Exchange/Return Guarantee	Article	Artnr
	Exchange / Return Guarantee Fee Standard coloured lens	N30404

For iris colour/ print based on photo, please attach images according to the following when placing your order



Picture instructions for ordering Iris colour according to sample
 1. Take a picture from straight ahead, in daylight (without flash)
 2. Crop the image as above and send with the order

Please contact us for advice regarding colouring.