

## Coloured soft contact lenses

Colouring of soft contact lenses is done only on ML-Soft designs made of Hydrogel material (ML 65 or ML75). Si-H material is not coloured.

### COLOURING

| Type                                                                                                                                                                                          | Design                                                                        | Specification                                               |  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------|-------------------------------------------------------------|--|
| <b>Colouring of the iris</b><br><br><i>Please always include a specification regarding if you would like an open or closed pupil. If open is selected, then please also specify the size.</i> | <b>Iris, standard colour. One of the following colours and % is selected:</b> | Iris: Ø 8.0 - 14.0 mm (1.0 step) and Ø 11.5* and 12.5 mm    |  |
|                                                                                                                                                                                               | Blue (5; 10; 20; 30%)                                                         |                                                             |  |
|                                                                                                                                                                                               | Green (5; 10; 20; 30%)                                                        |                                                             |  |
|                                                                                                                                                                                               | Aqua (5; 10; 20; 30%)                                                         |                                                             |  |
|                                                                                                                                                                                               | Brown 5%, and 10 - 100% (10 steps)                                            |                                                             |  |
|                                                                                                                                                                                               | Grey 5%, and 10 - 100% (10 steps)                                             |                                                             |  |
|                                                                                                                                                                                               | <b>Iris, not standard / filter colour</b>                                     | Pupil open: Ø 1.5 - 5.0 mm (0.5 step) and 7.0 mm (Ø 4.5*mm) |  |
|                                                                                                                                                                                               | Mixed iris colour                                                             |                                                             |  |
|                                                                                                                                                                                               | Black Iris                                                                    |                                                             |  |
|                                                                                                                                                                                               | Double colouring of standard colour (high intensity)                          |                                                             |  |
|                                                                                                                                                                                               | Filter colour ML 450                                                          |                                                             |  |
|                                                                                                                                                                                               | Filter colour ML 550                                                          |                                                             |  |
|                                                                                                                                                                                               | <b>Iris opaque</b>                                                            |                                                             |  |
|                                                                                                                                                                                               | Iris, black x2 (opaque)                                                       |                                                             |  |

\*Standard diameter. If no diameter is specified on the order, this diameter will be used. NB! The difference between the color diameter and the total diameter of the lens must be min. 2 mm.

| Type                                         | Design                                                                                                                                                                                     | Specification                                         |  |
|----------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------|--|
| <b>Advanced colouring according to photo</b> | An iris print is colored according to your request (reference photo). Specify: Open/closed pupil and with or without white back surface. Contact us for advice when ordering for the time. | Ø 11.3 and 12.3 mm<br>Spherical powers: -10.0 / +15.0 |  |


| Type                      | Design                                                                           | Specification                                      |  |
|---------------------------|----------------------------------------------------------------------------------|----------------------------------------------------|--|
| <b>Colouring of pupil</b> | <b>Pupil, standard colour. One of the following colours is selected:</b>         | Ø 3.0 - 5.0 mm in 0.5 mm step.<br>Ø 6.0 and 7.0 mm |  |
|                           | Brown 5%, and 10 - 100% (10 steps)                                               |                                                    |  |
|                           | Grey 5%, and 10 - 100% (10 steps)                                                |                                                    |  |
|                           | Black                                                                            | (Standard diameter: 4.5 mm)                        |  |
|                           | <b>Pupil, ML filter colour. One of the following filter colours is selected:</b> |                                                    |  |
|                           | Filter colour ML 450                                                             |                                                    |  |
|                           | Filter colour ML 550                                                             |                                                    |  |
|                           | Pupil, black x2 (opaque)                                                         |                                                    |  |

±5% colour deviation may occur

| Add on                                          |  |
|-------------------------------------------------|--|
| Limbal ring                                     |  |
| Black pupil ring                                |  |
| Coloured white back surface                     |  |
| White back surface x 2                          |  |
| Black dot (used as a marker)                    |  |
| Decentration (1,0 mm steps) incl. prism ballast |  |

| Warranty exchange and right of return | Article                                                       | Artnr  |  |
|---------------------------------------|---------------------------------------------------------------|--------|--|
|                                       | Warranty Exchange & Right of return standard/ML filter colour | N30404 |  |

For iris colour/ print based on photo, please attach images according to the following when placing your order



Picture instructions for ordering Iris colour according to sample  
 1. Take a picture from straight ahead, in daylight (without flash)  
 2. Crop the image as above and send with the order

Please contact us for advice regarding colouring.