


# BeFlared

FEVAR




# BeFlared

## FEVAR



### Elevating the standard in FEVAR procedures

Bentley continues its mission of making on-label solutions available for the treatment of patients. With the BeFlared FEVAR Stent Graft System Bentley wants to move even beyond by implementing innovative features dedicated to bridging stents in FEVAR procedures. As a 2-in-1 device BeFlared thrives to set a new standard in FEVAR.




The image shows the BeFlared FEVAR Stent Graft System. It consists of a long, blue, flexible catheter that is coiled in a large loop on the left side of the frame. The catheter extends to the right, where it terminates in a stent graft assembly. This assembly includes a blue catheter with a yellow tip, a white, mesh-like stent graft, and a blue catheter with a yellow tip. A small blue dot on the white stent graft is connected by a thin blue line to the text 'Precise & standardised positioning'.

**First indicated & dedicated platform**  
for bridging stents in FEVAR

**2-in-1 design**  
that combines inflation & flaring in one step

**Precise & standardised positioning**  
due to a dedicated fenestration marker

# THE FLARED ONE



## Product specifications

Stent Material	CoCr (L-605)
Graft Material	One layer of micro-porous ePTFE tubing (203 ± 25 µm)
Stent Graft Design	Single stent scaffold Open-cell design Patented ePTFE cover fixation
Delivery System	0.035" OTW
Catheter Shaft Length	120 cm
Expanded Stent Graft Diameter Distal/Proximal	Ø 5/8 and 6/8 mm <sup>1</sup> Ø 7/10 and 8/10 mm <sup>2</sup> Ø 9/11 mm <sup>3</sup>
Nominal Stent Graft Length	22, 27, 37 mm <sup>1</sup> 22*, 27, 37 mm <sup>2</sup> 27, 37 mm <sup>3</sup>
Introducer Sheath Compatibility	6 F for Ø 5/8 and 6/8 mm 7 F for Ø 7/10 to 9/11 mm
Shelf Life	3 years

<sup>1</sup>SV= Small Vessel Design, <sup>2</sup>MV= Medium Vessel Design,

<sup>3</sup>LV= Large Vessel Design, \*not available in Ø 8/10 mm

Three different stent designs ensure ideal stent performance.

FEVAR Stent Graft System

## Indication

The BeFlared FEVAR Stent Graft System is indicated for use as a bridging stent graft in patients requiring endovascular abdominal aortic aneurysm repair with a fenestrated aortic endovascular graft to keep the visceral arteries perfused.



## Ordering information

Expanded Stent Graft Diameter Distal/Proximal	Nominal Stent Graft Length	Introducer Sheath Size	Catalogue Number for Catheter Length 120 cm
5/8 mm	22 mm	6 F	BF0522_2
	27 mm		BF0527_2
	37 mm		BF0537_2
6/8 mm	22 mm	6 F	BF0622_2
	27 mm		BF0627_2
	37 mm		BF0637_2
7/10 mm	22 mm	7 F	BF0722_2
	27 mm		BF0727_2
	37 mm		BF0737_2
8/10 mm	27 mm	7 F	BF0827_2
	37 mm		BF0837_2
9/11 mm	27 mm	7 F	BF0927_2
	37 mm		BF0937_2

### Be innovative

We focus on the development of unique products and strive to find creative answers to questions posed by the healthcare industry

### Be dedicated

A marked customer focus and compliance to international standards help us to uphold unsurpassed quality in products and services

### Be responsible

We care for better health, our region and the environment. Every member of our team is motivated to make a vital contribution to the quality of life

### Be effective

We strive for time- and cost efficient solutions to satisfy the versatile interests of our customers and patients worldwide

All information provided in this brochure is not intended to be professional medical advice. The Instructions for Use (IFU) for the product must therefore be consulted before use of any product. Please be aware that the product and/or an indication might not be approved in your country, e.g. US. Please contact your Bentley representative for further information.

All data on file at Bentley InnoMed GmbH.



Bentley InnoMed GmbH  
Lotzenäcker 3, 72379 Hechingen, Germany  
tel +49 7471 984 995 10  
mail [info@bentley.global](mailto:info@bentley.global)

CE 0123

© 2025 Bentley InnoMed GmbH | Germany  
All Rights Reserved | Art.-Nr. 406319-001/004