

# Knipsa

## Endomyocardial Biopsy Forceps




# Knipsa – Endomyocardial Biopsy Forceps

■ Manufactured in Germany

■ Various lengths available for different access points





## ■ Further advantages

### Excellent maneuverability

- Bowden cable is easy to move

### Quick and efficient

- Ambulant and stationary

## Order information

REF	Configuration	Insertion length	Outer diameter shaft	Cutting jaws volume	Recommended sheath size	Packaging
KNIPSA500	pre-bent	50 cm	1.90 mm	2.42 mm <sup>3</sup>	≥ 6 Fr (2.00 mm)	sterile, single packed
KNIPSA1100	pre-bent	109 cm	1.90 mm	2.42 mm <sup>3</sup>	≥ 6 Fr (2.00 mm)	sterile, single packed

The Biopsy Forceps Knipsa is for single use only. It is available in standard sizes and used in cardiology.

# Joline<sup>®</sup> Medical Technology

Joline GmbH & Co. KG  
Neue Rottenburger Str. 50  
72379 Hechingen  
Germany

Phone: +49 (0)7471 9881-0  
Fax: +49 (0)7471 9881-222  
E-mail: [info@joline.de](mailto:info@joline.de)  
[www.joline.de](http://www.joline.de)



Quality first: Joline is EN ISO 13485 certified.  
All Joline products are CE-marked.  
Subject to alterations: 51499/25/02\_R05

CE0124