



 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 1 / 16


1. AUDIT TRAIL

Version	Change Description
Version 2.5.2	2.5.2 is an update for 2.5.1 with optimizations to the report generation engine and bug fixes

2. REMARKS

3. SIGNATURES

Signature Meaning	Name	Role	Date	Signature
Author	Ignas Žakaitis	QA Lead	2025/03/10 17:03:15	 sign-p14-i189230-b20250310160315141
Reviewer	Abbas Dhilawala	Group CTO	2025/03/10 18:02:45	 sign-p14-i189230-b20250310170245922
Reviewer	François Chasseur	Deputy CTO	2025/03/10 17:18:05	 sign-p14-i189230-b20250310161805403
Reviewer	Olena Motrichuk	Product Owner	2025/03/10 16:48:12	sign-p14-i189230-b20250310154812069
Approver	Ann Vankrunkelsven	RAQA Responsible	2025/03/11 10:45:36	 sign-p14-i189230-b20250311094536354

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 2 / 16

4. PURPOSE

Documentation of verification and validation results of Matrix Requirements Release 2.5.2.0-75ff6b3

5. SCOPE

This release with number 2.5.2.0-75ff6b3 is a patch of the 2.5.1 release and contains the changes listed in the "Appendix: Changes" below.

6. REFERENCES AND APPLICABLE DOCUMENTS

- **SIGN-625** Project Kickoff Form
- **SIGN-454** Verification and Validation Plan, Ver. F
- **SIGN-587** Software Development Plan, Ver.H

7. RESPONSIBILITIES

Name	Responsibility
Ignas Žakaitis	creating tests, testing
François Chasseur	coding, creating tests, code review
Siarhei Klimovich	coding, code review, testing
Aloïs Paulus	coding, code review, testing
Daniel Kula	coding, creating tests, code review, testing
Michael Stanton	coding, code review, testing
Sylvain Caillet	coding, code review, testing
Youssef Barnoukh	coding, testing
Nestor Canales	testing

8. VERIFICATION AND VALIDATION SUMMARY


8.0.1. Validation Summary

No validation was required due to the scope of changes.

8.0.2. Verification Summary

- The critical code changes were reviewed and no major issue was found. The minor issues were fixed.

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 3 / 16

- Three remaining non-critical issues found in the test cycle were documented as known issues and will be fixed in one of the upcoming releases.
- The tests were executed using three builds (2.5.2.0-5ab6e45, 2.5.2.0-8af1f7b, 2.5.2.0-75ff6b3) and no new issues were found in the final build

8.0.3. Final Verdict


No new blocking issue was found and the software build 2.5.2.0-75ff6b3 can be released.


9. VERIFICATION DETAILS

All changes have been reviewed and tests have been executed to ensure they are working. A total of 148 manual tests and 40 automated tests have been executed.

The traceability has been reviewed to make sure all changes have been implemented, reviewed and tested as planned.

Three issues found existed already in previous versions and will not be fixed. The other issues found during the testing phase were fixed and the tests were successfully repeated.



 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 4 / 16

9.1. TESTED RELEASE BUILDS

The tests were executed using 3 builds.

Build	Test Cycle	Number of Tests	Passed Tests
2.5.2.0-5ab6e45	First Cycle	117	103
2.5.2.0-8af1f7b	Second Cycle	27	25
2.5.2.0-75ff6b3	Third Cycle	4	4 (2 passed with deviation)
2.5.2.0-75ff6b3	Automated Tests	40	40

9.2. TESTERS

Manual tests were executed by 8 testers:


Tester	Number of Executed Tests
ignas	77
siarhei	10
sylvain	9
mike	21
daniel	9
youssef	13
nestor	8
alois	1

Results of automated tests were reviewed by Ignas.

9.3. TEST ENVIRONMENT

All tests were executed with the latest version of the used browsers

Browser / OS	Number of Tests
Windows 11 / Brave Browser Version 1.75.180 Chromium: 133.0.6943.126 (Official Build) (64-bit)	73
macOS Sonoma 14.2 (23C64) / Chrome 133.0.6943.127	10
Mac OS Sonoma 14.6.1 (23G93) - Brave Version 1.75.180 Chromium: 133.0.6943.126 (Official Build) (x86_64)	9
Mac OS / Chrome Version 133.0.6943.127 (Official Build) (arm64)	18
macOS Sonoma Version 14.5/Chrome Version 133.0.6943.127 (Official Build) (arm64)	9

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 5 / 16

MacOs 15.2 / Chrome Browser Version 1.75.180 Version 132.0.6834.111 (Official Build) (arm64)	21
Chrome Version 133.0.6943.142 (Official Build) (arm64)	3
MacOS 15.3 (24D60) / Version 133.0.6943.142 (Official Build) (arm64)	1
Windows 11 / Brave Browser Version 1.76.73 Chromium: 134.0.6998.45 (Official Build) (64-bit)	4

Notes

- It was verified that above changes were all changes done to the source code (nothing else changed than the planned changes)
- Tests were executed after the last code review
- All FEAT have been tested (see Appendix: Test Results)

10. VALIDATION DETAILS

n/a - see above

11. RISK ANALYSIS DETAILS

0 new risks have been identified with this release.


12. APPENDIX: CHANGES

FEAT-1150 Report Generation Performance Improvement

Report generation takes too long on projects with deep traceability links between items because of a quirk in the implementation of our report generator. We've reimplemented this to vastly reduce the size of intermediate output files and (often) dramatically reduce the amount of time it takes to generate a report.

We have identified both files size / generation performance are especially bad when there is a need of traceability matrix because it generates an intermediate file (j<job_id>_filter.xml) which contains nested downlinks in multiple levels and duplicated many times.

For reports that don't use a custom XSLT transformation file, we've reimplemented the work in Java, such that there is no need to generate the large intermediate filter files.

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 6 / 16

FEAT-1151 Improvement management of database ResultSet objects

The Clouds server makes many database queries as it performs an export. These queries consume resources until they are explicitly closed. An export with history is especially resource intensive. A customer with a large project experienced a server crash trying to perform an export with history.

Rather than waiting until the end of the export process to clean up all the database query objects, we can clean them up right after we've finished with them. This saves on memory, and allows the export to succeed. The report generation process will also benefit from this improved use of resources.

FEAT-1152 Yellow text mark removed after moving to next published item

Yellow text mark removed after moving to next published item.

Yellow text mark for searched text is no longer removed after moving to next published item in live QMS.

FEAT-1153 Risk analysis may not be shown in a DOC if the risk config is missing a key


The RiskStats control will fail to display in a "Risk Analysis" DOC field if there is no "risk_config" project setting. It should be okay that we don't have this setting, for the case where the risk_config is being configured from the field instead of the global setting. With a minor adjustment, this RiskStats control could be made to accept that there is a null/empty risk_config setting.

Before, the user had to workaround this issue by creating an empty risk_config setting.


The software now tolerates an empty/non-existent risk_config setting.

13. APPENDIX: TEST RESULTS

Item	Tested By	Execution Date	Last Modification	Test Result
FEAT-1150 Report Generation Performance Improvement	XTC-5814 Risk Management (TC-6)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5819 Downlink count optimization (TC-401)	2025/02/24	2025/03/10 10:57:32	passed with deviation
	XTC-5970 Controlled Documents (TC-9)	2025/03/05	2025/03/10 16:31:16	passed with deviation
	XTC-5821 Controlled Documents Customization (TC-41)	2025/02/25	2025/03/10 10:57:32	passed

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 7 / 16

Item	Tested By	Execution Date	Last Modification	Test Result
	XTC-5992 Document TOC and Index Test (TC-49)	2025/03/07	2025/03/10 16:31:55	passed with deviation
	XTC-5972 Trace Matrix Section (TC-54)	2025/03/06	2025/03/10 16:34:01	passed
	XTC-5824 Dropdown / (Multi-) Select section (TC-55)	2025/02/26	2025/03/10 10:57:32	passed
	XTC-5825 DOC - List of Items Test (TC-60)	2025/02/26	2025/03/10 10:57:32	passed
	XTC-5973 Table Section Test (TC-63)	2025/03/05	2025/03/10 16:34:21	passed
	XTC-5830 Report Functional Tests (TC-10)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5974 Risk Analysis Section (TC-56)	2025/03/05	2025/03/10 16:34:38	passed
	XTC-5836 Item Table Section (TC-80)	2025/02/25	2025/03/10 10:57:32	passed
	XTC-5841 Folders ID in documents (TC-148)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5975 Items tree for the List of items reference (TC-153)	2025/03/05	2025/03/10 16:34:52	passed
	XTC-5844 Custom column headers in tables (TC-157)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5847 List Of Items Title Numbering (TC-190)	2025/02/27	2025/03/10 10:57:32	passed
	XTC-5848 Risk Table Custom Columns (TC-194)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5849 Advanced Columns In items as table sections (TC-195)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5852 Sub section rendering of packages in SIGN (TC-316)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5853 Risk Pie Chart Labels are clean (TC-291)	2025/02/24	2025/03/10 10:57:32	passed
	XTC-5854 Risk zone rendering (TC-296)	2025/02/26	2025/03/10 10:57:32	passed


 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 8 / 16

Item	Tested By	Execution Date	Last Modification	Test Result
FEAT-1151 Improvement management of database ResultSet objects	XTC-5831 Excel Export Tests (TC-70)	2025/02/24	2025/03/10 10:58:16	passed
	XTC-5958 Branch Item Move Export/Import History Test (TC-400)	2025/03/04	2025/03/10 10:58:16	passed
	XTC-5959 Project Restore Test (TC-61)	2025/03/05	2025/03/10 10:58:16	passed
	XTC-5960 XML item and audit export (TC-11)	2025/03/04	2025/03/10 10:58:16	passed
	XTC-5961 Compose: Export/Import (TC-386)	2025/03/04	2025/03/10 10:58:16	passed
	XTC-5962 Project With Invalid Project Setting Export/Copy (TC-508)	2025/03/04	2025/03/10 10:58:16	passed
FEAT-1152 Yellow text mark removed after moving to next published item	XTC-5953 Highlighting in LiveQMS Searches (TC-595)	2025/02/25	2025/03/10 10:58:31	passed
FEAT-1153 Risk analysis may not be shown in a DOC if the risk config is missing a key	XTC-5872 Empty Risk Config is Tolerated (TC-594)	2025/02/24	2025/03/10 10:58:35	passed

14. APPENDIX: ADDITIONAL REGRESSION TESTS


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5990 Smoke Test (UC-1)	2025/03/07	2025/03/10 09:20:56	passed
XTC-5774 Variations, Labels and Filter (UC-8)	2025/02/24	2025/02/24 14:54:04	passed
XTC-5775 Review With Labels (UC-3)	2025/02/25	2025/02/25 15:02:58	passed
XTC-5776 Items and Field Rendering (UC-5)	2025/02/25	2025/02/25 16:11:16	passed
XTC-5777 Item Tree and Navigation and Keyboard (UC-6)	2025/02/24	2025/03/10 09:18:59	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 9 / 16


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5778 Traceability Test (UC-7)	2025/02/25	2025/02/26 02:33:46	passed
XTC-5779 Testing (UC-9)	2025/02/25	2025/02/25 15:16:27	passed
XTC-5780 Risk Management (UC-10)	2025/02/24	2025/02/26 17:13:53	passed
XTC-5964 Controlled Document Test (UC-11)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5991 Controlled Document Content Test (UC-16)	2025/03/07	2025/03/10 09:21:08	passed
XTC-5783 Audit Log (UC-12)	2025/02/24	2025/03/10 09:18:53	passed
XTC-5784 Searching And Filtering (UC-13)	2025/02/25	2025/02/25 17:00:08	passed
XTC-5785 JIRA tests (UC-14)	2025/02/28	2025/03/01 09:21:22	passed
XTC-5786 Email (UC-15)	2025/02/24	2025/02/24 21:26:30	passed
XTC-5966 Design Reviews (UC-27)	2025/03/06	2025/03/10 09:20:33	passed
XTC-5788 XTC Status overview (UC-29)	2025/02/24	2025/03/10 09:18:45	passed
XTC-5789 Traceability Overview dashboard (UC-31)	2025/02/26	2025/02/26 10:02:46	passed
XTC-5790 Concurrent work (UC-36)	2025/02/25	2025/02/25 17:51:46	passed
XTC-5791 Merging and Branching (UC-38)	2025/02/27	2025/02/27 11:49:56	passed
XTC-5792 Gates (UC-25)	2025/02/26	2025/02/26 13:37:18	passed
XTC-5793 Compose Module (UC-39)	2025/02/26	2025/02/26 09:40:55	passed
XTC-5794 Detailed history view (UC-35)	2025/02/24	2025/03/10 09:18:18	passed
XTC-5795 Cleanup tools (UC-33)	2025/02/26	2025/02/26 10:46:10	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 10 / 16


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5796 Multiple risk categories in one project (UC-32)	2025/02/26	2025/02/26 14:57:07	passed
XTC-5797 FMEA table editing (UC-30)	2025/02/26	2025/02/26 09:32:00	passed
XTC-5798 Frequently Used Functions (UC-37)	2025/02/26	2025/02/26 12:39:53	passed
XTC-5799 Other QMS functions (UC-17)	2025/02/27	2025/02/27 19:44:01	passed
XTC-5800 Install and uninstall a plugin (UC-40)	2025/02/24	2025/02/24 21:21:35	passed
XTC-5803 Audit Log Tests (TC-62)	2025/02/27	2025/02/27 09:05:03	passed
XTC-5804 Item redlining edge cases (TC-555)	2025/02/27	2025/02/27 13:24:13	passed
XTC-5931 Uplinks Fields (TC-47)	2025/02/26	2025/02/26 10:47:14	passed
XTC-5805 Document Packages (TC-167)	2025/02/25	2025/02/25 15:11:02	passed
XTC-5933 Revision Number Column Test (TC-228)	2025/02/24	2025/02/26 16:29:21	passed
XTC-5806 Gate E-Signature Test (TC-230)	2025/02/24	2025/02/26 17:15:18	passed
XTC-5807 Risk Zone Test (TC-283)	2025/02/24	2025/02/26 16:39:52	passed
XTC-5934 smart tag test (TC-325)	2025/02/25	2025/02/26 01:58:08	passed
XTC-5808 Verify image display (TC-340)	2025/02/25	2025/02/25 15:29:21	passed
XTC-5967 Hyperlink field test (TC-415)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5809 Test and Use Cases basics (TC-7)	2025/02/26	2025/02/26 21:35:16	passed
XTC-5810 Test run creation (TC-504)	2025/02/25	2025/02/25 08:50:40	passed
XTC-5968 Risk Management deep traces (TC-506)	2025/03/05	2025/03/10 09:20:33	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2		
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44	
	Internal Document ID: SIGN-634	Page: 11 / 16	


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5937 Risk Customization (TC-28)	2025/02/26	2025/02/26 14:54:51	passed
XTC-5812 Risk Configuration Manual RPN After Control (TC-50)	2025/02/25	2025/02/25 15:49:07	passed
XTC-5813 Risk Readonly Rendering (TC-402)	2025/02/24	2025/02/27 11:38:51	passed
XTC-5873 Risk Mitigation tests (TC-527)	2025/02/27	2025/02/27 10:15:26	passed
XTC-5969 Design Review Label Test (TC-39)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5816 Create and remove new DOC templates based on labels (TC-133)	2025/02/25	2025/02/25 18:31:07	passed
XTC-5817 Creating DOC Template labels creates the correct data structures under labels (TC-134)	2025/02/24	2025/02/25 13:18:25	passed
XTC-5818 Generating Documents with Document Status Labels (TC-135)	2025/02/26	2025/02/26 13:28:33	passed
XTC-5938 Lock Test (TC-57)	2025/02/28	2025/03/01 21:46:37	passed
XTC-5826 Electronic Signatures (TC-16)	2025/02/24	2025/02/27 13:41:24	passed
XTC-5827 Word Template Test (TC-24)	2025/02/27	2025/02/27 09:39:02	passed
XTC-5828 Create DOC from DOCX (TC-58)	2025/02/26	2025/02/26 16:33:25	passed
XTC-5833 Document Redlining (TC-89)	2025/02/26	2025/02/28 16:04:09	passed
XTC-5834 Copy Dashboard to DOC (TC-97)	2025/02/25	2025/02/25 10:05:51	passed
XTC-5835 Documents with item selection based on search expressions (TC-105)	2025/02/24	2025/02/26 10:48:20	passed
XTC-5837 Report Templates (TC-110)	2025/02/24	2025/02/27 11:49:43	passed
XTC-5838 Electronic Signatures - Target folder selection (TC-116)	2025/02/24	2025/02/26 14:00:03	passed
XTC-5839 Label Filter in DOCs (docFilter) (TC-130)	2025/02/25	2025/02/26 02:03:25	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2		
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44	
	Internal Document ID: SIGN-634	Page: 12 / 16	


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5946 Checking the creation of feedback with the "Risks" category (TC-136)	2025/02/25	2025/02/25 18:30:33	passed
XTC-5840 Document Structures (TC-141)	2025/02/28	2025/03/01 21:29:14	passed
XTC-5842 Cache signed SIGN documents (TC-151)	2025/02/24	2025/02/24 21:43:33	passed
XTC-5947 File attachments for DOCs / Signatures for external PDFs/DOCX (TC-158)	2025/02/24	2025/02/27 14:56:39	passed
XTC-5845 PDF Locks (TC-163)	2025/02/26	2025/02/26 10:54:04	passed
XTC-5846 Show checkbox in documents only if checked (TC-166)	2025/02/24	2025/02/24 21:56:22	passed
XTC-5850 Section Headings (TC-258)	2025/02/25	2025/02/26 01:58:16	passed
XTC-5976 Signature Rejection and Signature Hint (TC-261)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5855 Two Signature Table Test (TC-324)	2025/02/25	2025/02/25 12:49:31	passed
XTC-5856 Group name is correctly displayed in locked DOC (TC-407)	2025/02/28	2025/02/28 18:41:45	passed
XTC-5857 Double Gate Test MATRIX-5237 (TC-408)	2025/02/24	2025/02/24 21:19:24	passed
XTC-5858 Page Break Test (TC-587)	2025/02/24	2025/02/26 20:36:09	passed
XTC-5860 QMS Editor test (TC-76)	2025/02/28	2025/03/01 09:11:16	passed
XTC-5861 Document Template Tests (TC-77)	2025/02/28	2025/03/01 21:13:38	passed
XTC-5862 My Documents Test (TC-103)	2025/02/24	2025/02/26 16:19:45	passed
XTC-5863 QMS Revision Table Test (TC-229)	2025/02/28	2025/03/01 22:04:01	passed
XTC-5977 Passwords / User Profile (TC-15)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5864 XSS Tests General (TC-233)	2025/02/28	2025/02/28 21:11:15	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2		
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44	
	Internal Document ID: SIGN-634	Page: 13 / 16	


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5978 XSS Rich Text Editor Test (TC-234)	2025/03/06	2025/03/10 09:20:33	passed
XTC-5866 XSS Tests Plain Text Editors (TC-235)	2025/02/28	2025/03/03 20:18:33	passed
XTC-5993 & escaping test (TC-305)	2025/03/07	2025/03/10 09:21:16	passed with deviation
XTC-5949 DOC History Test (TC-289)	2025/02/24	2025/03/10 09:23:05	passed
XTC-5980 Test conflicting CSS Styling (TC-363)	2025/03/04	2025/03/04 16:36:12	passed
XTC-5981 Ensure that a deleted users full name appears in a signed document (TC-345)	2025/03/06	2025/03/10 09:20:33	passed
XTC-5982 Check docx/pdf generation of risk tables with very little probabilities (TC-350)	2025/03/06	2025/03/10 09:20:33	passed
XTC-5874 Reseting gate with id's of fields (TC-457)	2025/02/25	2025/02/25 18:27:43	passed
XTC-5868 Scaled images on printed paper (TC-458)	2025/02/24	2025/02/26 14:00:55	passed
XTC-5983 Abbreviations in docs (TC-459)	2025/03/06	2025/03/10 09:20:33	passed
XTC-5986 Smart links in custom tables (TC-461)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5984 Items table formating in DOC (TC-470)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5870 Downloading of DOC in old project (TC-475)	2025/02/24	2025/03/10 09:22:59	passed
XTC-5950 Table DOC Sections Configuration Saving (TC-518)	2025/02/27	2025/02/28 14:40:57	passed
XTC-5987 TOC Edge Case with List of Items Folders as Heading (TC-574)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5955 Empty SOP Group in LiveQMS (TC-575)	2025/02/28	2025/03/01 09:04:26	passed
XTC-5951 SIGN item Review Label Race Condition (TC-576)	2025/02/24	2025/02/26 17:13:11	passed
XTC-5988 Verify that URsL are not detected as cross project links (TC-585)	2025/03/05	2025/03/10 09:20:33	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 14 / 16


Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5989 Risk Mitigation Table (TC-593)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5985 Design Review Usability (TC-596)	2025/03/05	2025/03/10 09:20:33	passed
XTC-5956 Compose v2 (ATC-3)	2025/03/07	2025/03/10 10:18:47	passed
XTC-5877 Compose Update (ATC-4)	2025/03/07	2025/03/10 10:18:56	passed
XTC-5878 Compose License Check (ATC-8)	2025/03/07	2025/03/10 10:19:00	passed
XTC-5880 Verify Restricted Access to User Details (ATC-11)	2025/03/07	2025/03/10 10:19:05	passed
XTC-5882 Password Strength settings (ATC-10)	2025/03/07	2025/03/10 10:19:14	passed
XTC-5885 Generate multiple SIGN (ATC-15)	2025/03/07	2025/03/10 10:19:17	passed
XTC-5887 Converting TC to XTC Dead Links (ATC-18)	2025/03/07	2025/03/10 10:19:21	passed
XTC-5888 Help Button (ATC-20)	2025/03/07	2025/03/10 10:19:25	passed
XTC-5889 Project Navigation (ATC-21)	2025/03/07	2025/03/10 10:19:31	passed
XTC-5890 Delete categories in admin (ATC-22)	2025/03/07	2025/03/10 10:19:34	passed
XTC-5891 Project Renaming (ATC-24)	2025/03/07	2025/03/10 10:19:38	passed
XTC-5892 Search text (ATC-25)	2025/03/07	2025/03/10 10:19:42	passed
XTC-5893 XTC creation filter inheritance (ATC-26)	2025/03/07	2025/03/10 10:19:47	passed
XTC-5894 Out of Date Indication (ATC-27)	2025/03/07	2025/03/10 10:19:50	passed
XTC-5896 Project Groups (ATC-30)	2025/03/07	2025/03/10 10:19:54	passed
XTC-5897 MRQL: Review search (ATC-31)	2025/03/07	2025/03/10 10:19:57	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2		
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44	
	Internal Document ID: SIGN-634	Page: 15 / 16	

Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5899 Downlinks in Locked Items (ATC-37)	2025/03/07	2025/03/10 10:20:01	passed
XTC-5900 Copied Category ID and Name (ATC-38)	2025/03/07	2025/03/10 10:20:04	passed
XTC-5901 Re-indexing of items after configuration change (ATC-39)	2025/03/07	2025/03/10 10:20:08	passed
XTC-5902 Published Item Order (ATC-42)	2025/03/07	2025/03/10 10:20:11	passed
XTC-5903 VisibleOption hides field in folder (ATC-43)	2025/03/07	2025/03/10 10:20:15	passed
XTC-5904 Required fields (ATC-44)	2025/03/07	2025/03/10 10:20:19	passed
XTC-5905 Search Tabbing (ATC-46)	2025/03/07	2025/03/10 10:21:18	passed
XTC-5906 Search List Limit (ATC-47)	2025/03/07	2025/03/10 10:21:21	passed
XTC-5907 Risk Factor Options (ATC-48)	2025/03/07	2025/03/10 10:21:24	passed
XTC-5909 Automatic XTC results (ATC-50)	2025/03/07	2025/03/10 10:21:27	passed
XTC-5910 Field Default Values (ATC-51)	2025/03/07	2025/03/10 10:21:31	passed
XTC-5912 Detailed Merge Test (ATC-16)	2025/03/07	2025/03/10 10:21:34	passed
XTC-5913 Folder Merging (ATC-33)	2025/03/07	2025/03/10 10:21:39	passed
XTC-5914 Last Child Searches (ATC-17)	2025/03/07	2025/03/10 10:21:54	passed
XTC-5915 Uplink, Downlink and Folder Search (ATC-9)	2025/03/07	2025/03/10 10:21:58	passed
XTC-5916 MRQL Test Searches (ATC-28)	2025/03/07	2025/03/10 10:22:02	passed
XTC-5917 MRQL risk search (ATC-35)	2025/03/07	2025/03/10 10:22:12	passed
XTC-5918 Author and Editor Searches (ATC-34)	2025/03/07	2025/03/10 10:22:16	passed

All dates and times are in the project's default time zone and formatting.

 matrix requirements	Verification and Validation Report Ver 2.5.2	
	MATRIXSPECS: Validation (VALID)	Date: 2025/03/10 16:38:44
	Internal Document ID: SIGN-634	Page: 16 / 16

Executed Test Case	Execution Date	Last Modification	Test Result
XTC-5919 Dropdown and Label Searches (ATC-19)	2025/03/07	2025/03/10 10:22:23	passed
XTC-5920 Categories Ending in F Searching (ATC-23)	2025/03/07	2025/03/10 10:22:27	passed
XTC-5921 Date Searches (ATC-40)	2025/03/07	2025/03/10 10:22:30	passed
XTC-5925 Calendar of Changes (ATC-14)	2025/03/07	2025/03/10 10:22:34	passed