

#whimpexinc



wh.com/na



# Piezomed Module Promo

The Piezomed module adds a Surgical Piezo to an existing W&H rotary unit that can be used to perform extractions, sinus & bone surgery, in addition to implant surgery. Any W&H glass touch screen system can be combined with the Piezomed module.

Promo Code: PzModTp



#### Modularity

W&H can update the device to adapt & incorporate new clinical requirements of clinicians.



#### Automatic Tip Detection

Piezomed detects the instrument when it is inserted and automatically sets the correct power settings.



#### LED Handpiece & Efficient Cooling

The Piezo handpiece provides excellent illumination & irrigation close to the instrument's area of activity; resulting in efficient of the operating site.



#### Bluetooth Foot Pedal

The wireless foot pedal allows the clinician to choose the optimal location - one less cord in the operating area.



■ REF 90000350 = \$7,734.65 USD  
Piezomed Module set with 1.8m cable

■ REF 90000351 = \$7,941.01 USD  
Piezomed Module set with 3.5m cable

**Promotion includes Piezomed Module & 6 tips - discounted valued over \$1,000 USD**



or



or



or



Does not include tip cassette.



Select 1 tip



Select 4 tips

