

For stripping

Intensiv Swingle Professional Kit Ortho

Efficient and calibrated stripping in Orthodontics

A professional procedure for applying mechanical strips in IPR protocols



2. Space Creation



3. Finishing



1. Opening



20'000 oscillations
1.8 mm stroke

fix

free



- Creation of interproximal space between adjacent teeth in relation to the clinical requirements of IPR (InterProximalReduction) stripping protocols
- Finishing and polishing of interproximal enamel

Clinical pictures: Dr. Francesco Garino, Turin, Italy

Intensiv SA
Via al Molino 107
6926 Montagnola
Switzerland

Tel. + 41 91 986 50 50
Fax + 41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

Intensiv
Swiss Dental Products

Nordic / April 2019

Intensiv Swingle Professional Kit

Order via fax: 091 986 50 59 or via e-mail: order@intensiv.ch

Order now!



Product	Ref.	Price	Quantity
Intensiv Swingle Professional Kit, with light 6 Ortho-Strips included	WG-69 LT Ortho PROF	1'080.00	
Intensiv Swingle Professional Kit, without light 6 Ortho-Strips included	WG-69 A Ortho PROF	750.00	

Choice between different strip grits according to a professional sequence, following clinical needs for the creation of interproximal space.

All packages contain 3 strips



Opening of tight contact areas			
IntensivOrtho-Strips System, 08µm, opener	OS08OP-DS/3	78.00	
Creation of interproximal space			
IntensivOrtho-Strips System, 80µm, extracoarse	OS80XC-DS/3	78.00	
IntensivOrtho-Strips System, Central, 60µm, coarse	OS60C-CEN-DS/3	83.70	
IntensivOrtho-Strips System, Central, 40µm, medium	OS40M-CEN-DS/3	83.70	
Interproximal finishing/polishing			
IntensivOrtho-Strips System, Central, 25µm, fine	OS25F-CEN-DS/3	83.70	
IntensivOrtho-Strips System, Central, 15µm, polishing	OS15POL-CEN-DS/3	83.70	
Space measurement			
Intensiv IPR-DistanceControl Set, 8 elements	IPR-DC Set	68.00	
Tray			
Intensiv Ortho-StripsTray*	OST400	49.00	

* delivered without strips

Orthodontics

Intensiv Swingle Professional Kit, WG-69 LT Ortho PROF, W&H Synea, with light

Contents:

Contra-angle, Ref. WG-69LT*

Ejector, Ref. 053

Spray canal cleaner needle, Ref. 054

+

Intensiv Ortho-Strips System 08µm

Intensiv Ortho-Strips System 80µm

Intensiv Ortho-Strips System, Central, 60µm

Intensiv Ortho-Strips System, Central, 40µm

Intensiv Ortho-Strips System, Central, 25µm

Intensiv Ortho-Strips System, Central, 15µm



**6 Ortho-Strips
for free**

* Contra-angle without light: WG-69 A

Clinic-Stamp:

My dental clinic:

Date

Signature

☐ Please send me the Intensiv Catalogue 2019

