



# Correlation of Irrigation Tubing to W&H Surgical Units (and W&H Private Label Units)

implantMED

SI-95



REF 00993003

implantmed

SI-915



REF 00900102 / 00900103

00900106 / 00900107



REF 16929001

implantmed

SI-1015



REF 30289000



04363600 (6 pcs, 2.2 m)

Disposable (one time use)

irrigation tubing set




REF 90000270 / 271 Classic



04363600 (6 pcs, 2.2 m)


Disposable irrigation tubing set



06338400

Irrigation spike with roller clamp

+



04013900

Pump tubing with nipple

+




00929300 (10 pcs, 2.2 m)

Disposable (one time use) irrigation tubing



Spare pump tubing also available



04014000 (3 pcs, 85 mm)

Spare Pump Tubes

\*\*This is the only tubing compatible with the Implantmed SI-1015 / SI-1023



# Correlation of Irrigation Tubing to W&H Surgical Units (and W&H Private Label Units)

elcomed  
SA-310



REF 15933100 / 15933101  
15933102 / 15933103

piezomed  
SA-320




REF 15933100 / 15933101

elcoMED  
SA-200 / SA-200 C



REF 00994013 / 00994003




04363600 (6 pcs, 2.2 m)

Disposable (one time use) irrigation tubing set



04365300 (6 pcs, 2.2 m)


Disposable (one time use) irrigation tubing set



06338400

Irrigation spike with roller clamp


+



04013900


Pump tubing with nipple

+



00929300 (10 pcs, 2.2 m)


Disposable (one time use) irrigation tubing



06338400

Irrigation spike with roller clamp

+



04035100

Pump tubing with nipple

+



00929300 (10 pcs, 2.2 m)

Disposable (one time use) irrigation tubing


Spare pump tubing also available



04014000 (3 pcs, 85 mm)

Spare Pump Tubes

Spare pump tubing also available



04035200 (3 pcs)

Spare Pump Tubes