

#whdentalwerk



wh.com



# Show casing W&H

Up to  
**60**  
**% OFF**



**Black Friday  
Promotions**

**24<sup>th</sup> to 28<sup>th</sup>  
November 2023**

## Implantmed Plus SI-1023

- › Ultra short, lightweight motor with 6.2 Ncm torque
- › Speed range between 200 – 40,000 rpm
- › Precise torque between 5 – 80 Ncm
- › Automatic thread cutting function
- › Modular system for ultimate upgradeability

# Show casing W&H Up to **60% OFF**



## Osstell Beacon

Implant stability  
measurement ISQ  
Supported by over  
1,200 scientific  
articles

**1**  
Year  
Warranty



Claim  
**FREE**  
Smart  
Pegs\*

**OSSTELL**  
A W&H COMPANY

### Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the colour light indicator, depending on the ISQ value. Reduce treatment times and avoid rejection in high-risk patients.

**£1,952**

\*Register on Osstell Connect to claim your free SmartPegs when you buy an Osstell Beacon

## Surgical Handpieces



- › Straight handpiece 1:1
- › External single spray
- › Dismantleable

### S-11

- › Without light

**£480**

### S-11 L

- › Mini LED+

**£675**

- › Contra-angle 20:1
- › External single spray
- › Dismantleable

### WS-75

- › Without light

**£594**

### WS-75 L

- › Mini LED+

**£649**



### Kit 04 SI-1023 Optic/Wireless

- › Implantmed SI-1023
- › Optic motor
- › Wireless Foot Control

**£2,760**

**2**  
Years  
Warranty

## Turbines & Contra-angles

syneco FUSION



### Syneco Fusion TG-98

- › Without light
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

**£372**



### Syneco Fusion TG-97

- › Without light
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

**£372**



### Syneco Fusion WG-56 A

- › Contra-angle 1:1
- › Without light
- › Single Spray
- › Head Ø 9.5 mm

**£341**



### Syneco Fusion WG-99 A

- › Contra-angle 1:5
- › Without light
- › 4x Spray
- › Head Ø 9.5 mm

**£483**



### Syneco Fusion TG-98 L

- › LED+
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

**£428**



### Syneco Fusion TG-97 L

- › LED+
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

**£428**



### Syneco Fusion WG-56 LT

- › Contra-angle 1:1
- › Optic
- › Single Spray
- › Head Ø 9.5 mm

**£428**



### Syneco Fusion WG-99 LT

- › Contra-angle 1:5
- › Optic
- › 4x Spray
- › Head Ø 9.5 mm

**£541**



LED+ light



Ring LED+



Light



No light



Sprays