


#whdentalwerk



wh.com



Incredible **inside & outside**

Lexa Sterilizer – More Capacity and Efficiency

lexa

High capacity and efficiency

Lexa offers versatile cycle program options for several loads, guaranteed by the closed-door drying process with HEPA filters and the use of fresh water for each cycle.



**18 pounds
sterilized and dried
in 30 minutes!**



Incredible efficiency

Large chamber capacity of 7.4 gallons and short cycle times, including excellent load drying.



Incredible flexibility

Lexa comes with 4 large trays along with a unique rack providing larger loading surface and more loading options.



Incredible protection

Automatic digital cycle recording on USB stick.



Incredible speed

Lexa offers versatile program options for several loads in order to reduce overall cycle duration.



Incredible usability

User-friendly interface and high-end 4.3" color touchscreen for simple navigation to all sterilization option.



Incredible standards

Closed-door vacuum-assisted drying for excellent load drying and single use water management.



W&H Impex Inc
6490 Hawthorne Drive
Windsor, Ontario,
N8T 1J9 Canada
ph +1 519 944 6739
office.ca@wh.com
wh.com

W&H North America
33091 W Jefferson Ave
Brownstown, MI 48173 U.S.A.
ph +1 800 265 6277
office.us@wh.com
wh.com

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.