

#whdentalwerk



wh.com



Show casing W&H

Up to
60
% OFF



13th to 16th
May 2022

Implantmed Plus SI-1023

- › Ultra short, lightweight motor with 6.2 Ncm torque
- › Speed range between 200 – 40,000 rpm
- › Precise torque between 5 – 80 Ncm
- › Automatic thread cutting function
- › Modular system for ultimate upgradeability



Kit 04 SI-1023 Optic/Wireless

- › Implantmed SI-1023
- › Optic motor
- › Wireless Foot Control


£2,503

£6,268 RRP

2
Years
Warranty

Show casing W&H Up to **60** %OFF

Surgical Handpieces



- › Straight handpiece 1:1
- › External single spray
- › Dismantleable

S-11

- › Without light

£435

£1,098 RRP

S-11 L

- › Mini LED+

£612

£1,529 RRP



- › Contra-angle 20:1
- › External single spray
- › Dismantleable

WS-75

- › Without light

£538

£1,345 RRP

WS-75 L

- › Mini LED+

£589

£1,471 RRP



Osstell Beacon

Implant stability
measurement ISQ
Supported by over
1,200 scientific
articles

1
Year
Warranty



Claim
FREE
Smart
Pegs*



Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the colour light indicator, depending on the ISQ value. Reduce treatment times and avoid rejection in high-risk patients.

£1,699

£4,424 RRP

Turbines & Contra-angles

syneo FUSION



Syneo Fusion TG-98

- › Without light
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

£336

£841 RRP



Syneo Fusion TG-97

- › Without light
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

£336

£841 RRP



Syneo Fusion WG-56 A

- › Contra-angle 1:1
- › Without light
- › Single Spray
- › Head Ø 9.5 mm

£309

£773 RRP



Syneo Fusion WG-99 A

- › Contra-angle 1:5
- › Without light
- › 4x Spray
- › Head Ø 9.5 mm

£438

£1,095 RRP



Syneo Fusion TG-98 L

- › LED+
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

£388

£969 RRP



Syneo Fusion TG-97 L

- › LED+
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

£388

£969 RRP



Syneo Fusion WG-56 LT

- › Contra-angle 1:1
- › Optic
- › Single Spray
- › Head Ø 9.5 mm

£388

£969 RRP



Syneo Fusion WG-99 LT

- › Contra-angle 1:5
- › Optic
- › 4x Spray
- › Head Ø 9.5 mm

£490

£1,225 RRP

Lara Vacuum B Sterilizer



Lara B type Sterilizer

- › Intuitive user interface and touchscreen display
- › Cycle traceability available thanks to Wi-fi connection
- › Mobile App for Real Time Remote Monitoring
- › Ergonomic design, with modular feet for flexible positioning
- › Pre-programmable vacuum and helix test function
- › Automatic water fill function
- › Eco Dry technology automatically adapts the drying time to the load
- › Unique optional upgradeability

20* 17l. 22l.

*with fast cycle upgrade

LARA 17L. **£4,121** £5,723 RRP

LARA 22L. **£4,630** £6,431 RRP

Showcasing W&H

Up to
60% OFF



Assistina Twin

- › Two process chambers allowing loading of a 2nd handpiece while the 1st is processed.
- › Superior Maintenance and Lubrication Guaranteed
- › Intuitive and easy to use
- › A flexible range of handpiece adaptors available

£1,340 £2,684 RRP

W&H (UK) Ltd

6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ

t +44 (0) 1727 874990

f +44 (0) 1727 874628

office.uk@wh.com

www.wh.com