

#whdentalwerk



video.wh.com



**Show  
casing  
W&H**

Special offers  
July - December 2025



**New  
Thed**  
Fulfills every wish

Turbines

High speed, low weight

All W&H turbines have a hygienic head' anti-retraction system to prevent cross infection.



3

10.4

122

Without light

3x Spray

Head Ø 10.4 / Ø 12.2 mm

16 / 18 W

390,000 / 330,000 rpm

Roto Quick® Connection

£ 245

RRP £568

The TK-94 L microturbine is perfect for hard to access areas and the mouths of children.



5

9

122

Synea Vision TK-94 L

LED+

5x Spray

MICRO Head Ø 9 mm

14 W/410,000 rpm

Roto Quick® Connection

£ 751

RRP £1,502





4

10

122

Synea Fusion TG-97 / TG-97 L

Without light / LED+

4x Spray

MINI Head Ø 10 mm

18 W/390,000 rpm

Roto Quick® Connection

TG-97

TG-98

£ 503

RRP £1,006

TG-97 L

TG-98 L

£ 579

RRP £1,157

4

11.5

122

Synea Fusion TG-98 / TG-98 L

Without light / LED+

4x Spray

MIDI Head Ø 11.5 mm

21 W/360,000 rpm

Roto Quick® Connection

Special Offer Price

When you buy

3x TG 97/98

£ 1,358 Non-optic

£ 1,562 Optic



Air motors & Couplings



Couplings

Available from

£ 105

RRP £210



Air motors

Available from

£ 247

RRP £448

Restoration & Prosthetics

Contra-angles

High torque and precision

4

10.7

122

Synea Fusion WG-99 LT / WG-99 A

Contra-angle 1:5

Optic / Without light

4x Spray

Head Ø 9.5 mm

£ 695

RRP £1,463

£ 628

RRP £1,308

3

10.7

122

Synea Fusion WG-900 LT / WG-900 A

Contra-angle 1:4

Optic / Without light

3x Spray

Head Ø 10.7 mm

£ 695

RRP £1,489

£ 627

RRP £1,334



Alegra WE-56

Contra-angle 1:1

Without light

Optional External Spray

Head Ø 10.0 mm

£ 135

RRP £356

1

10.7

122

Synea Fusion WG-56 LT / WG-56 A

Contra-angle 1:1

Optic / Without light

Single Spray

Head Ø 9.5 mm

£ 579

RRP £1,157

£ 462

RRP £923

3

10.7

122

Synea Vision WK-900 LT

Contra-angle 1:4

Optic

3x Spray

Head Ø 10.7 mm

£ 875

RRP £1,750



Alegra HE-43

Straight Handpiece 1:1

Without light

Optional External Spray

£ 125

RRP £309

Prophylaxis

Polishing by rotation

Enjoy professional tooth cleaning and give your patients a beaming smile.

Proxeo HP-44 M

Handpiece 4:1

System Doriot

For disposable contra-angles with plastic shank

£ 199

RRP £328





Proxeo WP-64 MU

Prophy contra-angle 4:1

Without light

For universal screw-in and snap on (adaptor) cups and brushes

£ 289

RRP £578

Our turbines are available to fit other manufacturers couplings.

# Implantmed Plus SI-1023

Workflow driven solution

Flexibility to meet your needs!

W&H offer a synergistic solution for oral surgery guiding the clinician through the cycle from implementation to evaluation.

You can enjoy the benefits of utilising the modular Implantmed surgical motor and Piezomed piezosurgery unit alongside the Osstell Connect consultation tool and the Osstell Beacon, which ensures predictable assessment of ideal loading.

Free Surgical WS-75  
Contra-angle Included  
with selected Implantmeds



**Implantmed Plus SI-1023 Kit 01**  
Non-optic/Cable

- Non-optic motor
- Cable Foot Control
- FREE WS-75

£ 3,321  
RRP £7,379



**Implantmed Plus SI-1023 Kit 02**  
Non-optic/Wireless

- Non-optic motor
- Wireless Foot Control
- FREE WS-75

£ 3,634  
RRP £8,076



**Implantmed Plus SI-1023 Kit 03**  
Optic/Cable

- Optic motor
- Cable Foot Control
- FREE WS-75 L

£ 3,836  
RRP £8,524



**Implantmed Plus SI-1023 Kit 04**  
Optic/Wireless

- Optic motor
- Wireless Foot Control
- FREE WS-75 L

£ 4,149  
RRP £9,220

2 Years\*  
Warranty



**Implantmed Plus SI-1023 + Piezomed Module**  
Optic/Wireless

- Optic motor
- Wireless Foot Control
- FREE WS-75 L
- Piezomed Plus Module SA-435
- Extraction Tip Set

£ 6,769  
RRP £13,664

# Implantmed Classic SI-923



**Implantmed Classic SI-923**  
Non-optic/Cable

- Non-optic motor
- Cable Foot Control
- FREE WI-75 E/KM

£ 2,880  
RRP £6,401

# Piezomed

"A real practice builder"



**Piezomed Classic Optic/Cable**

- Piezomed Classic
- Handpiece with LED ring
- Cable Foot Control
- FREE BONE SET Worth £719

£ 5,319  
RRP £9,671



**Piezomed Plus Module**

- Power 24W. Handpiece with light SA-40 L
- Automatic instrument detection
- Instrument changer and instrument library
- FREE EXTRACTION SET INCLUDED

£ 3,411  
RRP £6,201



# Piezomed Instruments

The W&H Piezomed instrument range includes different instruments for bone surgery, periodontology and extraction, for retrograde endodontics and for lateral sinus lift.



Automatic  
instrument  
detection



Bone surgery

£ 741  
RRP £1,348



Extraction

£ 235  
RRP £428



Periodontology

£ 328  
RRP £596



Retrograde  
endodontics

£ 582  
RRP £1,058



Lateral  
sinus lift

£ 532  
RRP £967



Preparation of the implant  
site and crestal sinus lift

£ 1,043  
RRP £1,897



## Surgical Handpieces

High-quality stainless steel with a **scratch-resistant coating** makes W&H surgical handpieces robust and gives them a particularly long service life.



**Self-generating LED Surgical Handpieces available.**



**Handpiece 20:1 Zygoma SZ-75**

- › Without light
- › External single spray
- › Dismantleable

**£ 1,180**

RRP £2,603



**Contra-angle 20:1 WI-75 E/KM**

- › Without light
- › External single spray
- › Not dismantleable

**£ 585**

RRP £1,300



**Contra-angle WS-91 / 1:2.7**

- › Without light
- › 45° head with 3x spray
- › Dismantleable

**£ 681**

RRP £1,513



**Contra-angle 20:1 WS-75**

- › Without light
- › External single spray
- › Dismantleable

**£ 723**

RRP £1,606



**Contra-angle WS-91 L / 1:2.7**

- › Mini LED+
- › 45° head with 3x spray
- › Dismantleable

**£ 753**

RRP £1,674



**Contra-angle 20:1 WS-75 L**

- › Mini LED+
- › External single spray
- › Dismantleable

**£ 791**

RRP £1,757



**Straight handpiece S-11 / 1:1**

- › Without light
- › External single spray
- › Dismantleable

**£ 584**

RRP £1,297



**Straight handpiece S-11 L / 1:1**

- › Mini LED+
- › External single spray
- › Dismantleable

**£ 822**

RRP £1,826



## Osstell Beacon

Measurement of implant stability



### Osstell Classic

Osstell Classic helps you to objectively determine implant stability and measure the degree of osseointegration.

**£ 1,350**

RRP £2,800

### Osstell SmartPegs

The SmartPeg measures the resonance frequency when it is attached to an implant.

**£ 85**

RRP £117




### Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and avoid rejection in high-risk patients.

**£ 1,650**

RRP £5,277

**OSSTELL**  
A W&H COMPANY



### What are the positive outcomes when using the Osstell Implant stability solutions?

- › A positive experience for patients, encouraging trust through understanding
- › Enhanced patient education builds confidence, resulting in increased treatment acceptance
- › Reduction in treatment visits and failures with improved average healing times
- › Improved outcomes for patients with risk factors
- › Enhanced communication between specialists ensuring a more synergistic approach to treatment
- › Improved treatment documenting

Osstell ISQ devices with the Osstell ISQ Scale help implantologists to make informed implant treatment decisions.



**Osstell Connect**  
FREE App for Osstell Beacon

Osstell Connect is designed to document quantified data which ensures traceability, whilst complementing data in your Patient Management System.

Systems with graphical presentation of data and statistics which is great for patient communication. You can easily present treatments in a graphic and totally personalized way, even generating pdf reports, while having access to the world's largest online surgical portal for peace of mind. The patient consultation tool allows you to easily communicate your proposed implant.



Hygiene & Maintenance



Internal cleaning & lubrication

Increase the life of your handpieces and reduce the costs of repairs!



10" quick connect by W&H

Assistina Twin

- Two process chambers allowing loading of a 2nd handpiece while the 1st is processed
- Superior Maintenance & Lubrication
- Intuitive and easy to use
- A flexible range of adaptors available

£ 1,441

RRP £3,202



18" quick connect by W&H

Assistina One

Internal cleaning & lubrication of one handpiece in only 18 seconds, with Assistina One the aerosols stay in the HEPA filter thereby releasing clean air.


£ 1,020

RRP £2,169



Service Oil F1

Ask your dealer about W&H Handpiece Oil



Assistina Twin and Assistina One includes

Care set and Quick RM and ISO adaptors

Ultrasonic cleaner

Cleaning with high-frequency sound waves



New!

U-son 3l

- 3 liter tank capacity
- Stainless steel chamber
- Two transducers
- Programmable temperature
- Programmable timer
- Automatic "degas" function

£ 495

RRP £800

U-son 6l

- 6 liter tank capacity
- Stainless steel chamber
- Three transducers
- Programmable temperature
- Programmable timer
- Automatic "degas" function

£ 755

RRP £1,250

Sterilization



Accessories for sterilization



Osmo

- Supplies multiple sterilizers simultaneously with top quality demineralized water.
- Automatic water filling
  - Most cost-effective solution
  - Increase sterilizer service life
  - Up to 6 devices simultaneously

£ 926

RRP £1,543



3x EliReel Included

Seal2 Starter Pack

- Universal Pouch Sealing
- 12mm hermetic seal
- 3 x EliReel Rolls FREE

£ 747

RRP £1,245

W&H EliReel sealing pouches

From / per roll

£ 14

RRP £22

Pro Service

PremiumCare



Gives you the flexibility to select a plan that best suits the individual needs of your practice.

The standard PremiumCare package offers a flexible plan allowing you to add optional extras as you please or the **PremiumCare Plus or Flex\*\*** plan, which include fully comprehensive service packages and the required pressure vessel inspection and validation testing at no extra cost.

PremiumCare

Maintenance plan for all W&H sterilizer models, with a maximum limit of 4,000 cycles over the 3 year term\*.

PremiumCare Plus & Flex\*\*

A fully comprehensive maintenance plan for all W&H sterilizer models, including PVI certification, unlimited cycles and validation in line with manufacturer's or current HTM 01-05 guidance (model dependent).

Our Technical Service Team are here to support you nationwide!



ProService

Services	Premium Care	Premium Care Plus & Flex**
All routine service visits for three years	✓	✓
All service parts and scheduled sundries	✓	✓
All breakdown calls	✓	✓
All breakdown parts including labour	✓	✓
All engineer travel in the UK Excluding Highlands, Islands and ferries etc,	✓	✓
Free on-site operator training	✓	✓
Free technical software upgrades	✓	✓
Unlimited Cycles - for the duration of the plan	⬜	✓
PVI Certification - in line with Pressure Systems Safety Regulations 2000:128	⬜	✓
Periodic Test Logbook	⬜	✓
Validation in line with manufacturer's or current HTM 01-05 guidance (model dependent)	✗	✓
Installation and Commissioning in line with current guidelines	✗	✓
Helix Test Kit***	✗	✓
Extra Discount	✗	✓

✓ Included ✗ Not Included ⬜ Optional

\*An additional cycle tariff will be levied if the 4,000 cycle limit is exceeded within the 3 year period \*\* PremiumCare Flex Plans are only available for certain models \*\*\* Only applies to B Type Sterilizers



## Thermal washer disinfector devices Thed+ and Thed

The only system with two washing levels, for powerful cleaning, full flexibility, and faster turnaround in even the busiest clinics.

### BePro WDD Starter Pack

£ 65

- › BePro WDD Cleaner TE – 5000 ml
- › BePro WDD Rinse – 1000 ml
- › Canister opener



New!

### Thed

- › Two washing levels
- › Fast cycles with high cleaning performance
- › Active air-drying with HEPA filter
- › Full connectivity options
- › Internal positioning of process chemicals
- › Unique cartridge system (tabletop)
- › Modern design
- › Comprehensive range of accessories
- › User-friendly interface

### Thed Tabletop

with upper and lower rack

£ 5,715

RRP £8,860

### Thed Underbench

with upper and lower rack

£ 6,315

RRP £9,460

### Thed+

Thed features and:

- › EliSense: device status indicator
- › Automatic rack recognition
- › Automatic manifold injector recognition
- › Automatic program recommendation
- › Automatic adjustment of cycle time, water, process chemicals, and power consumption
- › Washing chamber with internal light

### Thed+ Tabletop

with upper and lower rack

£ 6,686

RRP £9,560

### Thed+ Underbench

with upper and lower rack

£ 7,286

RRP £10,160



## Sterilization at its Best!

With a wide range of sterilizers, W&H is able to offer the perfect option to meet individual needs.

### All W&H Sterilizers offer:

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Ergonomic design, with modular feet for flexible positioning
- › Easy to clean surfaces and accessible water
- › Automatic Water Fill function
- › Pre-programmable vacuum and helix (for B Type Sterilizers) test function to run automatically out of surgery hours



S 23'

### Lyla S type Sterilizer New concept sterilizer

- › Mobile App available via upgrade for remote data storage
- › Future-proof upgradability to a type B sterilizer using the new Activation Code concept



Unique upgrade flexibility with **Activation Codes**

- › B Cycle
- › Remote data storage
- › Traceability
- › All in One

£ 3,400

17 l. RRP £4,762

£ 3,800

22 l. RRP £5,008



\*Based on average 2kg load for immediate use



B 20'

New!

### Lara B type Sterilizer

- › Mobile App for Real Time Remote Monitoring
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradability to increase: Performance and Traceability or to add a Fast Cycle



Unique upgrade flexibility with **Activation Codes**

- › Performance\*\*
- › Traceability
- › Fast Cycle
- › All in One
- \*\* Only for Lara 17 and 22 litres

£ 4,650

17 l. RRP £6,857

£ 5,250

22 l. RRP £7,700

£ 5,750

Lara XL 28 l. RRP £8,492

£ 7,200

Lara XXL 38 l. RRP £10,370

\*with fast cycle upgrade

### Lisa Safe Label Printer & QR Barcode Reader

Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Barcode reader

Lisa Safe printer

£ 145

RRP £234

£ 510

RRP £802



### Seal Check

- › Pack of 100 sealing check

£ 31

RRP £62



B 13'

### Lisa B type sterilizer

- › Full extended traceability thanks to EliTrace
- › Mobile App for Real Time Remote Monitoring
- › Eco Dry+ technology automatically adapts the drying time to the mass of the load
- › Elisense user interface generates status, temperature and other cycle information
- › Cycle times from 13 minutes for specific loads
- › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use



B 10'

### NEW Lisa MINI B type sterilizer

- › Cycle times from 10 minutes
- › Integrated vacuum pump
- › Small and easily stackable: W 44cm x H 20cm x D 67cm
- › Built in water tank and automatic water filling - no external storage
- › For up to 22 instruments, offers complete traceability
- › 5 l. model includes internal vacuum pump
- › Includes integrated water tank

**3 liters version available**  
Requires compressed air



£ 5,550

17 l. RRP £8,467

£ 6,290

22 l. RRP £9,580

£ 3,400

5 l. RRP £4,900





Coming Soon



seethrough

Meet the Sales Team

### Richard Mansbridge

Area: London, South East, Anglia,  
East Midlands and Humber-side

**07775 792628**

richard.mansbridge@wh.com

### Caroline Caine

Area: West Midlands, Manchester, Merseyside  
South West England, Isle of Man, Channel  
Islands, Wales, Northern Ireland and Eire

**07786 515331**

caroline.caine@wh.com

### David Barrass

Area: Scotland, Humber-side,  
North of England and North West

**07795 315089**

david.barrass@wh.com

To place an order  
contact W&H UK on  
**+44 (0) 1727 874990**  
or contact your preferred  
dental dealer

### W&H (UK) Ltd

6 Stroud Wood Business Centre,  
Park Street, St Albans, AL2 2NJ

t +44 (0) 1727 874990

f +44 (0) 1727 874628

office.uk@wh.com

**www.wh.com**

Offers valid from 1<sup>st</sup> July to 23<sup>rd</sup> December 2025.  
Information is correct at time of publication and  
relates to W&H products and services available in  
the UK and Ireland. Prices are quoted in £ sterling  
and exclude VAT. Savings are based on RRP and  
rounded to the nearest £. Free goods are not  
transferrable and no cash equivalent is available.  
Terms and Conditions apply. E. & O. E.

Disinfection



## BePro

Professional Cleaning  
and Disinfection



### Cleaner TE

Tri-enzymatic detergent  
with neutral pH. No  
neutraliser required!

£ 45

Canister 5l RRP £56



### Rinse Aid

For improved  
drying phase.

£ 14

Bottle 1l RRP £18



### Soaking bath

Soaking bath for the  
immersion of dental  
instruments in the  
decontamination  
phase.

£ 38

RRP £46



### Disinfectant I

Powerful disinfectant  
concentrate for medical  
and dental instruments.

£ 38

Bottle 2l RRP £48



### Disinfectant R

Ready-to-use solution for the  
disinfection and cleaning of  
rotary dental instruments.

£ 13

Bottle 1l RRP £16

£ 44

Refill 5l RRP £53



### Disinfectant Wipes

120 wipes for the  
disinfection of non-  
invasive medical devices.

£ 9

Tube RRP £11

£ 7

Refill RRP £9



### Disinfectant Spray

Surface disinfectant for the  
rapid disinfection and cleaning  
of non-invasive medical devices.

£ 13

Spray 1l RRP £16

£ 46

Refill 5l RRP £56