

#whdentalwerk



wh.com



Piezomed
Module
Coming
Soon!

Show
casing
W&H | May-August
2022

KNOW&HOW

Turbines

All W&H turbines have hygienic head¹

anti-retraction systems to prevent cross infection.



Synea Vision TK-98 L

- › 5x Ring LED+
- › 5x Spray
- › MIDI Head Ø 11.5 mm
- › 24 W/360,000 rpm
- › Roto Quick® Connection

£709

£1,417 RRP



Synea Vision TK-97 L

- › 5x Ring LED+
- › 5x Spray
- › MINI Head Ø 10 mm
- › 21 W/400,000 rpm
- › Roto Quick® Connection

£709

£1,417 RRP



Synea Vision TK-94 L

- › LED+
- › 5x Spray
- › MICRO Head Ø 9 mm
- › 14 W/410,000 rpm
- › Roto Quick® Connection

£629

£1,258 RRP



Synea Fusion TG-98 L

- › LED+
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

£485

£969 RRP



Synea Fusion TG-97 L

- › LED+
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

£485

£969 RRP



Alegra TE-98

- › Without light
- › 3x Spray
- › STANDARD Head Ø 12.2 mm
- › 18 W/330,000 rpm
- › Roto Quick® Connection

£219

£475 RRP



Alegra TE-97

- › Without light
- › 3x Spray
- › MINI Head Ø 10.4 mm
- › 16 W/390,000 rpm
- › Roto Quick® Connection

£219

£475 RRP



Synea Fusion TG-98

- › Without light
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

£420

£841 RRP



Synea Fusion TG-97

- › Without light
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

£420

£841 RRP



Alegra TE-97

- › Without light
- › 3x Spray
- › MINI Head Ø 10.4 mm
- › 16 W/390,000 rpm
- › Roto Quick® Connection

£219

£475 RRP




RC Series RC-95 BC/RM

- › Without light
- › Single Spray
- › STANDARD Head Ø 12.2 mm
- › 11 W/330,000 rpm
- › Borden 2/3-hole/Ritter Midwest

£152

£240 RRP

50%
Off RRP



RC series are non washer disinfectable

¹Exclusive W&H system patented in Europe and the USA

Available for other fittings

Contra-angles

1:5 Speed increasing contra-angles, with FG turbine burs, are recommended for cavity and crown preparation. Due to the maximum speed of 200,000 rpm, aerosol creation is reduced.

1:1 contra-angles are recommended for finishing and polishing.



Synea Fusion WG-99 LT

- › Contra-angle 1:5
- › Optic
- › 4x Spray
- › Head Ø 9.5 mm

£613

£1,225 RRP



Synea Fusion WG-99 A

- › Contra-angle 1:5
- › Without light
- › 4x Spray
- › Head Ø 9.5 mm

£548

£1,095 RRP



Synea Fusion WG-56 LT

- › Contra-angle 1:1
- › Optic
- › Single Spray
- › Head Ø 9.5 mm

£485

£969 RRP



Synea Fusion WG-56 A

- › Contra-angle 1:1
- › Without light
- › Single Spray
- › Head Ø 9.5 mm

£387

£773 RRP



Alegra WE-56

- › Contra-angle 1:1
- › Without light
- › External Spray (Optional)
- › Head Ø 10.0 mm

£114

£237 RRP



Alegra HE-43

- › Straight Handpiece 1:1
- › Without light
- › External Spray (Optional)

£109

£238 RRP



RC Series RC-58

- › Contra-angle 1:1
- › Without light
- › External Spray (Optional)
- › Latch Grip

£89

£142 RRP



RC Series RC-43

- › Straight Handpiece 1:1
- › Without light
- › External Spray (Optional)

£81

£126 RRP

Electric Motor



Electric Motor EM-12 L

- LED light
- Internal spray

Add-on MF-100

2,000 – 40,000 rpm


£1,699

£3,012 RRP



2x
WG-99 LT
FREE

EM-12 L with 2xWG-99 LT
1:5 Speed increasing
contra-angle handpiece
Table-top, Plug and Play



RC series are non washer disinfectable

Roto Quick® Couplings

All W&H couplings have a **non-retraction valve** to prevent cross infection. Enjoy free 360° instrument rotation.



Roto Quick® RQ-24

- › Optic / Integrated spray
- › 6-hole connection

£145
£287 RRP



Roto Quick® RQ-04

- › Non-optic / Integrated spray
- › 4-hole Standard

£91
£174 RRP



Roto Quick® RQ-03

- › Non-optic / Integrated spray
- › (3)-hole Borden

£91
£174 RRP

Air Motors



Air Motor AM-25 LT

- › Optic with internal spray
- › Midwest for use with RQ-24 coupling

£490
£1,308 RRP



Air Motor AM-25 A BC/RM

- › Non-optic with internal spray
- › Borden / Midwest
- › Fixed Connection

£332
£870 RRP



Air Motor AM-25 BC/RM

- › Non-optic without internal spray
- › Borden / Midwest
- › Fixed Connection


£207
£374 RRP



Air Motor RC-20 BC/RM

- › Non-optic without spray
- › Borden / Midwest
- › Fixed Connection

£124
£168 RRP



Other models available

Why Choose W&H

In addition to our Accredited Service Partners, we have an in-house Technical Service Department and dedicated team of 10 factory-trained and qualified W&H Field Service Engineers covering the UK, providing state-of-the art service and support.

We at W&H, take great pride in our after sales service and understand the importance of taking care of our customers.



Joel has been at W&H for 12 years, of which 7 years were spent in the workshop servicing and repairing handpieces, Assistina's and surgical equipment. He decided to take the opportunity to train as a Field Service Engineer and has spent the past 5 years as an engineer on the road. Joel said 'he enjoys all aspects of the job but mostly travelling to parts of the country that he wouldn't otherwise get to see'. Outside of work Joel is a keen gardener and enjoys spending time in his garden, keeping it maintained the best he can. He has recently taken up golf as a new hobby and enjoys watching football and motorsport.

Joel McCann Field Service Engineer



Leah joined W&H in 2007 and has been an integral part of the Service Team ever since. Leah's role includes co-ordinating the Service Team, taking customer calls and making sure that an engineer is available when required amongst many other responsibilities. With her great wealth of knowledge, over the years Leah has become well known by many of W&H's Service Customers as the go to person for initial technical enquiries.

Out of work Leah and her husband enjoy walking the dog, who she says is spoilt rotten, shopping and loves to travel to warmer climes.

Leah Johnson Technical Service Advisor



Technical advice

Our technical team are happy to advise on the repair and maintenance of your W&H products.



Speedy Service

W&H guarantees that your equipment will be repaired normally within 24 hours, from acceptance of the estimate.



Free Transportation

Use our free handipak to send your handpieces in for repair or maintenance.



Quality Repairs

All of our technicians are factory trained to ensure top quality repairs and maintenance are carried out using only original spare parts.



Service Agreements

Gain added peace of mind with a W&H Service plan, ask for details.



ProService



Air Polishing System Proxeo Aura



Perio spray head
Subgingival pockets
> 3 mm

120° spray head
Supragingival,
sulcus, subgingival
pockets < 3 mm



Proxeo Aura Perio-Set

- › Air polishing handpiece
120° and standard perio
spray head with 40 perio
tips (subgingival)
- › Powder SOFT & SENSITIVE

£987

£1,518 RRP



Proxeo Aura Propy-Set

- › Air polishing handpiece
120° propy spray
head with Powder SOFT

£792

£1,218 RRP

Two Systems One Handpiece

The Proxeo Aura Air Polishing System is a real all-rounder. **Soft deposits and stains can be removed thoroughly and efficiently without changing the handpiece.** The powders, handpiece and spray heads work together perfectly and allow thorough and atraumatic supragingival and subgingival cleaning.



Polishing by rotation



Proxeo Twist LatchShort WP-66 W

- › Optimum access thanks to lower working height
- › Better visibility thanks to the smaller head
- › Change cup/brush at the push of a button
- › Sealing system prevents ingress of debris and fluids, ensuring longer working life and improved infection control
- › Ideal speed for gentle and efficient polishing/cleaning
- › Ease of maintenance and hygiene-friendly
- › Thermodisinfected and sterilisable

£199

£381 RRP



Proxeo Twist Cordless PL-40 HW

- › High quality polishing for professionals
- › Portable, simple to use solution
- › Cordless foot control
- › Premium quality rotary polishing system
- › Flexible cordless polishing
- › Ideal for limited mobility or domiciliary patients
- › Available for disposable heads or polishing cups

£795

£1,780 RRP



Proxeo WP-64 MU

- › Propy contra-angle 4:1
- › Without light
- › For universal screw-in and snap on (adaptor) cups and brushes.

£210

£483 RRP

Piezo Scaler Proxeo Ultra



**Compatible with
pacemakers**

Tested by
Graz University
of Technology

2
Years
Warranty



W&H Q-Link



W&H / EMS / Mectron



W&H / NSK / Satelec

The Proxeo Ultra is available with the unique **Q-Link quick connection system**. In addition, **handpieces are available to fit other manufacturers systems**. The tip changers with integrated torque control prevent overtightening and make changing the tip simple and safe.



**GERMAN
DESIGN
AWARD
WINNER
2020**

**Specialist implant
tip available**
for implant cleaning



Full range of tips

Tailored to meet the needs of all treatment indications including implantology

Proxeo Ultra PB-530

- Wireless foot control
- Cleaning program for hose reprocessing
- Compatibility with pacemakers
- LED ring for ultra-bright illumination
- Infinitely adjustable power setting
- Large tank (750 ml)

£2,115

£2,820 RRP

Proxeo Ultra PB-520

- Compatibility with pacemakers
- LED ring for ultra-bright illumination
- Infinitely adjustable power setting
- Large tank (750 ml)

£1,582

£2,434 RRP

Proxeo Ultra PB-510

- Compatibility with pacemakers
- LED ring for ultra-bright illumination
- Infinitely adjustable power setting

£1,080

£1,662 RRP


All Proxeo Ultra models are supplied with 3 tips (2x universal & 1x Perio)

Piezo Scaler Tips


The W&H piezo scaler tip range includes tips for prophylaxis and periodontology, endodontics, implantology and restoration & prosthetics. The W&H range of tips for piezo scalers including Tigon+, Proxeo Ultra and Pyon 2. Also compatible with other scaler brands.

Universal


Simple, gentle and efficient supragingival calculus removal in all quadrants.




Tip 1U **£56** £111 RRP



Tip 2U **£56** £111 RRP




Tip 3U **£64** £127 RRP



Tip 4U **£56** £111 RRP

Endodontics


W&H offers the perfect range of endo tips whether for flushing a root canal or complex revision in an endodontic patient.




Tip 1E **£83** £166 RRP



Tip 2E **£65** £129 RRP



Tip 3E **£83** £166 RRP



Tip 4E **£89** £177 RRP




Tip 5E **£81** £162 RRP




Tip 6E **£81** £162 RRP

Periodontology


Diamond-coated tips for thorough cleaning of the root surface during the initial treatment of periodontitis affected patients.



Tip 1P **£72** £143 RRP



Tip 3PI/Pr **£92** £184 RRP



Tip 5P **£111** £221 RRP



Tip 4P **£105** £173 RRP



Tip 5PI/Pr **£111** £221 RRP

Implant Set

Tips for the efficient and careful cleaning of implants and the attendant restorations.




Set 3x1I **£206** £411 RRP



Implantmed Plus SI-1023

- › Ultra short, lightweight motor with 6.2 Ncm torque
- › Speed range between 200 – 40,000 rpm
- › Precise torque between 5 – 80 Ncm
- › Automatic thread cutting function
- › Modular system for ultimate upgradeability




2
Years*
Warranty



Kit 01 SI-1023
Non-optic/Cable
› Implantmed SI-1023
› Non-optic motor
› Cable Foot Control

£2,662
£4,840 RRP



Kit 02 SI-1023
Optic/Cable
› Implantmed SI-1023
› Optic motor
› Cable Foot Control

£3,121
£5,674 RRP



Kit 03 SI-1023
Non-optic/Wireless
› Implantmed SI-1023
› Non-optic motor
› Wireless Foot Control

£2,983
£5,424 RRP



Kit 04 SI-1023
Optic/Wireless
› Implantmed SI-1023
› Optic motor
› Wireless Foot Control

£3,442
£6,268 RRP



*Warranty only for Implantmed Plus SI-1023 kits.

Implantmed Classic SI-923

- › Ultra short motor with 5.5 Ncm torque
- › Speed range between 300 – 40,000 rpm
- › Precise torque between 5 – 70 Ncm
- › Automatic thread cutting function



1
Year
Warranty




SI-923
Non-optic/Cable
› Implantmed SI-923
› Non-optic motor
› Cable Foot Control

£2,352
£4,277 RRP



Surgical Handpieces



- › Contra-angle 20:1
- › External single spray
- › Dismantleable

WS-75

- › Without light

£673
£1,345 RRP

WS-75 L

- › Mini LED+

£736
£1,471 RRP



- › Contra-angle 1:2,7 45°
- › 3x spray
- › Dismantleable

WS-91

- › Without light

£634
£1,267 RRP

WS-91 L

- › Mini LED+

£701
£1,402 RRP



- › Contra-angle 1:1
- › External single spray
- › Dismantleable

WS-56

- › Without light

£583
£1,166 RRP

WS-56 L

- › Mini LED+

£687
£1,373 RRP



- › Straight handpiece 1:1
- › External single spray
- › Dismantleable

S-11

- › Without light

£543
£1,086 RRP

S-11 L

- › Mini LED+

£765
£1,529 RRP

Surgical Handpieces with self-generating LED+ light are available for Non-optic units.
Call your local Representative for details.



S-16

- › Angled handpiece 1:2
- › External single spray
- › Without light

£750
£1,558 RRP

SZ-75

- › Zygoma handpiece 20:1
- › External single spray
- › Without light

£995
£2,181 RRP



Osstell Beacon

Implant stability measurement ISQ

Supported by over 1,200 scientific articles

1
Year
Warranty
WGH



Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the colour light indicator, depending on the ISQ value. Reduce treatment times and avoid rejection in high-risk patients.

£1,895

£4,424 RRP

Osstell SmartPegs

The SmartPeg used together with BEACON measures the resonance frequency when it is attached to an implant. FOR RELIABLE DIAGNOSIS AND PEACE OF MIND ALWAYS REQUEST THE OSSTELL SMARTPEG.

£75

£96 RRP

Pack of 5




The ISQ Scale

Osstell uses proprietary and evidenced-based RFA technology, based on more than 1,200 scientific studies. Visit osstell.com/clinical-guidelines to learn more.



Evolution of stability over time

As a result of osseointegration, the initial mechanical stability is supplemented and / or replaced by biological stability, the final value of implant stability being the sum of the two. Stability generally does not remain stable after implant insertion. There is an initial decrease in (mechanical) stability, followed by an increase, as the implant becomes biologically stable.



Analysis of the Frequency of Resonance (RFA) and ISQ

Method of measurement

The probe sends magnetic pulses at different frequencies so that the SmartPeg, connected to the implant, vibrates.

The probe measures the vibration of the SmartPeg to determine the frequency that gives the strongest vibration, this being the resonance frequency of the implant / SmartPeg assembly.

The more stable the implant, the higher the resonant frequency.



OsstellConnect

Osstell also gives you **free access to OsstellConnect**, where you can analyze your data, keep track of healing times and success rates over time, and compare with other Osstell users all over the world. **Simply a tool to help you excel.**

*Register online to claim your free Smartpegs when you buy an Osstell IDx or Beacon

Piezomed Classic



Automatic
instrument
detection

Piezomed, the ultimate force in ultrasonic bone surgery, offering extremely high cutting performance, whilst being gentle on soft tissue ensuring reduced patient discomfort and quicker healing.



Piezomed Optic/Wireless

- › Piezomed Classic
- › Handpiece with LED ring
- › Wireless Foot Control

£4,625

£8,409 RRP

Piezomed Instruments

The W&H Piezomed instrument range includes different instruments for bone surgery, periodontology and extraction, for retrograde endodontics and for lateral sinus lift.

Tips compatible with all Piezomed units.

Bone surgery

The fine saw tooth design – more teeth per instrument – sets new standards in the field of bone preparation.



Periodontology and extraction

Instruments for thorough cleaning in periodontitis treatment and for atraumatic tooth extraction.



Retrograde endodontics

Diamond-coated, delicate instruments for easy access during retrograde root canal preparation.



Lateral sinus lift

The instruments permit safe lateral preparation of the maxillary sinus wall and gentle mobilization of the Schneiderian membrane.



Preparation of the implant site and crestal sinus lift

The preparation instruments are specially matched to the bone qualities predominant in the maxilla.



Piezomed Module

This NEW modular, plug and play device is compatible with the new Implantmed SI-1023 versions. This modular system offers the possibility of adapting the device to the requirements of the individual user.

- › The ergonomically designed handpiece offers absolute comfort and safety.
- › Utilising an innovative design to direct the spray to the treatment area the Piezomed ensures ultra-efficient cooling of the operating site and avoids the need for a case fan.
- › Including the two technologies in one modular device, utilising only one foot control and one disposable coolant tube per treatment makes it the ultimate space and cost saving solution.



Piezomed Classic Module

- › Power 18 W
- › Handpiece without light SA-40
- › Instrument changer
- › 6 Irrigation tubes

£2,255

£4,100 RRP

Piezomed Plus Module

- › Power 24W
- › Handpiece with light SA-40 L
- › Enhanced automatic instrument detection
- › Instrument changer
- › Instrument library
- › 6 Irrigation tubes


£2,860

£5,200 RRP

Modular Surgery Solution

A single device for each treatment

The only system that offers oral surgery, implantology and piezosurgery in one device.



Documentation
of osseointegration



Beacon
for Stability
Measurement



Surgical
Instruments



Powerful
Electric Motor



Wireless
Foot Control

Chain of Hygiene & Sterilization

aims Advanced Infection Prevention Management Solutions

130 years of experience in research, development and production of dental solutions to prevent infections.



Dr Weigert high quality washer disinfector solutions available from W&H

Neodisher MediClean	Neodisher Z	Neodisher MediKlar Dental
› Cleaner › 5 litres	› Neutralizer › 5 litres	› Rinse Aid › 1 litre
£40 £46 RRP	£45 £56 RRP	£14 £19 RRP

Can be used in some other washer disinfectors – check with manufacturer



Step 1 Cleaning & Desinfection

The new Teon washer disinfector is easy to use and offers an exceptional standard of cleaning, disinfection and reprocessing, maximising efficiency and increasing the effectiveness of the decontamination protocol.



Step 2 Internal cleaning & lubrication

Increase lifetime and reduce costs for repairs!



Step 3 Packaging

Pack it!



Step 4 Sterilization

Always with high quality demineralised water!



Step 5 Traceability

Store and follow it!



Teon thermal washer disinfector

- › Validated cleaning of all instruments
- › Traceability and connectivity through Ethernet module
- › 1 Year Warranty

Standard Starter Set

£6,300

£7,816 RRP

Handpiece Injector Starter Set

£7,530

£8,840 RRP


Teon Chemical Starter Set

£98

£132 RRP

Assistina Twin

- › Two process chambers allowing loading of a 2nd handpiece while the 1st is processed.
- › Superior Maintenance and Lubrication Guaranteed
- › Intuitive and easy to use
- › A flexible range of handpiece adaptors available



Maintenance with Oil spray can



Maintenance with Assistina

£1,400

£2,684 RRP



Seal 2

- › Universal Pouch Sealing System
- › 12 mm weight sealing
- › 1 Pack of Seal Check FREE

£599

£1,042 RRP

Seal Check

- › Pack of 100 sealing check

£26

£51 RRP



Multidem C27

- › Produces good quality demineralised water guaranteed to:
- › Improve efficiency
- › Ensure consistent performance
- › Increase steriliser service life
- › Reduce service costs
- › Offer environmentally friendly solution, reducing plastic waste

Multidem C27


£390

£576 RRP

Refill Cartridge

£140

£295 RRP



Lisa Safe Label Printer & QR Barcode Reader

Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Lisa Safe Printer

£485

£668 RRP

Barcode Reader

£142

£195 RRP

Sterilization at its Best!



Unique upgrade
flexibility with
Activation Codes

- › B Cycle
- › Remote data storage
- › Traceability
- › All in One



Lyla. New concept sterilizer

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Mobile App available via upgrade for remote data storage
- › Ergonomic design, with modular feet for flexible positioning. Easy to clean surfaces and accessible water
- › Pre-programmable vacuum and helix test function to run automatically out of surgery hours
- › Automatic Water Fill function
- › Future-proof upgradability to a type B sterilizer using the new Activation Code concept

Lyla 17l.

£2,995

£3,996 RRP

Lyla 22l.

£3,338

£4,172 RRP



*Based on average 2kg load for immediate use



Unique upgrade
flexibility with
Activation Codes

- › Performance
- › Traceability
- › Fast Cycle
- › All in One



Lara B type Sterilizer

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Mobile App for Real Time Remote Monitoring
- › Ergonomic design, with modular feet for flexible positioning. Easy to clean surfaces and accessible water
- › Pre-programmable vacuum and helix test function to run automatically out of surgery hours
- › Automatic water fill function
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradeability to increase Performance and Traceability or to add a Fast Cycle

Lara 17l.

£4,121

£5,723 RRP

Lara 22l.

£4,630

£6,431 RRP



*with fast cycle upgrade



Lisa B type sterilizer

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Full extended traceability thanks to EliTrace
- › Mobile App for Real Time Remote Monitoring
- › Ergonomic design, with modular feet for flexible positioning. Easy to clean surfaces and accessible water
- › Pre-programmable vacuum and helix test function to run automatically out of surgery hours
- › Automatic water fill function
- › Eco Dry+ technology automatically adapts the drying time to the mass of the load
- › Elisense user interface generates status, temperature and other cycle information
- › Cycle times from 13 minutes for specific loads
- › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use

Lisa 17l.

£4,945

£7,074 RRP

Lisa 22l.

£5,590

£8,007 RRP



W&H UK Staying Connected

With an experienced, specialist sales team and highly skilled, qualified field service engineers, we are proud to offer you unparalleled support and advice on our full product range.

We at W&H, take great pride in our levels of service offering a high quality customer experience. We understand the importance of taking care of our customers.

To place an order
contact W&H UK on
+44 (0) 1727 874990
or contact your preferred
dental dealer

W&H (UK) Ltd
6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ
t +44 (0) 1727 874990
f +44 (0) 1727 874628
office.uk@wh.com
www.wh.com

Meet the Sales and Marketing Team



Kate Scheer

Marketing Executive

01727 876091

kate.scheer@wh.com

Having worked in the dental industry for over 30 years, Kate has a wealth of dental experience as well as a proven track record in marketing. She joined W&H as Marketing Executive in 2011 and works alongside the Sales Team and Management to increase the profile of the W&H brand and range. In addition to Kate's experience she is a qualified Chartered Marketer with a Masters in Strategic Marketing and is currently the Chairperson of the BDIA Statistics Committee. Kate is responsible for all W&H's marketing and communications.



In her spare time, Kate enjoys walking her 2 dogs, spending time with her teenage daughter, reading and travelling especially to France her home from home.



Richard Mansbridge

North, North West & East London,
Anglia, East Midlands, Yorkshire/
Humberside and North East England

07775 792628

richard.mansbridge@wh.com



Jon Bryant

South & West London,
South England, Scotland &
Channel Islands

07917 886326

jon.bryant@wh.com



Caroline Caine

West Midlands, North West and
South West England, Isle of Man,
Wales, Northern Ireland and Eire

07786 515331

caroline.caine@wh.com