

#whdentalwerk



video.wh.com



Show
casing
W&H

Special offers
July – December 2026



Click here to
discover the
Seethrough Flex

Turbines

High speed, low weight

All W&H turbines have a hygienic head¹ anti-retraction system to prevent cross infection.



Alegria
TE-97/98

- › Without light
- › 3x Spray
- › Head Ø 10.4 / Ø 12.2 mm
- › 16 / 18 W
- › 390,000 / 330,000 rpm
- › Roto Quick® Connection

£ 265

RRP £586



Synea Fusion
TG-97/98 - TG-97/98 L

- › Without light / LED+
- › 4x Spray
- › MINI Head Ø 10 / 11.5 mm
- › 18 W-390,000 rpm / 21 W-360,000 rpm
- › Roto Quick® Connection

TG-97/98 Non-optic

£ 519

RRP £1,037

TG-97/98 L Optic

£ 596

RRP £1,192



Special Offer Price
When you buy
3x TG 97/98

x3



3x TG-97/98 Non-optic

£ 1,399

RRP £3,111



3x TG-97/98 L Optic

£ 1,609

RRP £3,576

The TK-94 L
microturbine is
perfect for hard
to access areas
and the mouths
of children.



Synea Vision
TK-94 L

- › LED+
- › 5x Spray
- › MICRO Head Ø 9 mm
- › 14 W/410,000 rpm
- › Roto Quick® Connection

£ 774


RRP £1,548



Click here



¹Exclusive W&H system patented in Europe and the USA



Contra-angles

High torque and precision



Synea Fusion
WG-99 LT / WG-99 A

- › Contra-angle 1:5
- › Optic / Without light
- › 4x Spray
- › Head Ø 9.5 mm

£ 754

RRP £1,507

£ 674

RRP £1,348



Synea Fusion
WG-56 LT / WG-56 A

- › Contra-angle 1:1
- › Optic / Without light
- › Single Spray
- › Head Ø 9.5 mm

£ 596

RRP £1,192

£ 476

RRP £951



Alegria WE-56

- › Contra-angle 1:1
- › Without light
- › Optional External Spray
- › Head Ø 10.0 mm

£ 140

RRP £367



Synea Fusion Power Edition
WG-900 LT / WG-900 A

- › 1:4 high-speed contra-angle, designed to remove high-strength materials, like zirconia
- › Optic / Without light
- › 3x Spray
- › Head Ø 10.7 mm

£ 767

RRP £1,534

£ 688

RRP £1,375



Special Offer Price
When you buy
3x WG 56 LT / A

x3

3x WG 56 LT

£ 1,609

RRP £3,576

3x WG 56 A

£ 1,284

£2,853



Alegria HE-43

- › Straight Handpiece 1:1
- › Without light
- › Optional External Spray

£ 130

RRP £319

Air motors
& Couplings

Flexible – for use with all standardised turbine tubing. Full compatibility with all turbines, air motors and scalers with Roto Quick and Multiflex connections.



360° instrument
rotation
for easy handling
and to prevent the
tube from twisting.



Air motors
Available from

£ 256

RRP £462



Couplings
Available from

£ 112

RRP £217

Our turbines are available to fit other manufacturers couplings.

Implantmed Classic SI-923



Implantmed Classic SI-923 Non-optic/Cable

- › Non-optic motor
- › Cable Foot Control

FREE WI-75 E/KM

£ 2,967

RRP £6,594

Implantmed Plus SI-1023



Implantmed Plus SI-1023 Kit 01 Non-optic/Cable

- › Non-optic motor
- › Cable Foot Control

FREE WS-75

£ 3,421

RRP £7,602

Implantmeds available from £2,967

Oral Surgery & Implantology

Implantmed Plus II

Smart assistance for advanced implantology



**Save up to 10 minutes
per implant case
with the ScanEdition**



Implantmed Plus II Optic / Cable or Wireless

- › Optic motor
- › Cable or Wireless Foot Control

FREE WS-75 L

£ 4,235

RRP £7,700



Implantmed Plus II Optic / Wireless / Scan

- › Optic motor and Wireless Foot Control
- › Scan function for easy, efficient documentation

FREE WS-75 L

£ 4,648

RRP £8,450

Bundle Offer Implantmed Xpert Package

- › Implantmed Plus II SI-2102 Optic
- › With travel case
- › Piezomed Plus Module II SP-210

£ 8,038

RRP £15,833

Upgrade your Implantmed Plus/Plus II



Piezomed Plus Module

- › Power 24W. Handpiece with light SA-40 L
- › Automatic instrument detection
- › Instrument changer and instrument library
- › **FREE EXTRACTION SET INCLUDED**

£ 3,513

RRP £6,829



W&H are disrupting the industry with the Seethrough portfolio – it's genuinely one of the coolest imaging solutions I've seen. The clarity, workflow integration and clinical confidence it provides make a real difference in practice, and that's why it stands out."

Dr Aly Virani Clinical Director and Implant Dentist at North Cardiff Dental and Implants

Enjoy direct connectivity, linking Implantmed Plus II via Sicat to ioDent, making the Implantmed / Seethrough an intuitive synergistic solution.

Documentation ioDent®

ioDent® allows for uploading of digital treatment documentation including ISQ values, torque values and materials used, completing and saving the documentation flow.



Click here



Monitoring
Osstell Beacon
Osstell Connect

Treatment
Implantmed Plus II with
Piezomed Plus II

Acquisition
Seethrough Flex



Treatment Planning
Seethrough studio and Sicat

seethrough | studio
SICAT

Implant planning software which creates a direct link from Seethrough Studio Software to import and export planning information. Usable by over 13,000 implants in the library.

Enjoy a 3 month free trial of Sicat when you purchase an Implantmed Plus II

Piezomed Classic



Piezomed Classic Optic/Cable

- › Piezomed Classic
- › Handpiece with LED ring
- › Cable Foot Control
- › **FREE BONE SET**

£ 5,479

RRP £9,962



Click here to view the full
range of tips all with 50% off



Light



No light



Cable foot control



Wireless foot control



Scan function

Surgical Handpieces

High-quality stainless steel with a **scratch-resistant coating** makes W&H surgical handpieces robust and gives them a particularly long service life.



Self-generating LED Surgical Handpieces available.



Contra-angle 1:2.7 WS-91 / WS-91 L

- › Without light / Mini LED+
- › 45° head with 3x spray
- › Dismantleable

WS-91 Non-optic

£ 780

RRP £1,599

WS-91 L Optic

£ 863

RRP £1,725



Straight handpiece 1:1 S-11 / S-11 L

- › Without light / Mini LED+
- › External single spray
- › Dismantleable

S-11 Non-optic

£ 720

RRP £1,655

S-11 L Optic

£ 941

RRP £1,881



Contra-angle 20:1 WS-75 / WS-75

- › Without light / Mini LED+
- › External single spray
- › Dismantleable

WS-75 Non-optic

£ 828

RRP £1,655

WS-75 L Optic

£ 905

RRP £1,810



Contra-angle 20:1 WI-75 E/KM

- › Without light
- › External single spray
- › Not dismantable

£ 670

RRP £1,339



Handpiece 20:1 Zygora SZ-75

- › Without light
- › External single spray
- › Dismantleable

£ 897

RRP £2,682

Osstell

Measurement of implant stability



Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and avoid rejection in high-risk patients.

£ 1,700

RRP £5,436



Osstell Classic

Osstell Classic helps you to objectively determine implant stability and measure the degree of osseointegration.

£ 1,390

RRP £2,884



Click here

When you buy an Osstell Beacon register online at Osstell Connect to claim your 5 free single use SmartPeg*

*Free SmartPegs only available with the purchase of an Osstell Beacon and on registering for Osstell Connect



SINGLE-USE

Osstell Single Use SmartPegs

The SmartPeg measures the resonance frequency when it is attached to an implant.

Pack of 5

£ 88

RRP £121



REUSABLE

Osstell Reusable SmartPegs

Cost efficient ensuring accurate readings up to 20 uses. Soft aluminium avoids damage to the implant. Reusable up to 20 uses – needs sterilizing before initial use.

Sold singularly

£ 39

RRP £42

Two option: One Proven Method

NEW reusable SmartPegs provide a cost-efficient and ecofriendly solution for measuring implant stability. Designed to the same high quality as the single use version offering a cost-efficient, ecofriendly solution for measuring implant stability. **Surface coated to protect from wear and tear during a maximum of 20 sterilization cycles.**



see through

Discover the benefits of digital dentistry with proven W&H quality and see for yourself why this new must-see is a true must-have.



Click here



Seethrough Flex

High performance, small footprint: this is what Seethrough Flex stands for. With its compact design, intuitive operation and face-to-face positioning, this flexible CBCT device can be easily integrated into existing practice workflows.

£ POA

see through | studio

In addition to sophisticated hardware technologies, the revolutionary Seethrough Studio software is a decisive factor for the high standard of the entire W&H imaging portfolio. With its beneficial features and advanced algorithms, Seethrough Studio is an all-in-one solution that takes imaging quality and digital dentistry workflows to a new level.



Seethrough Studio

- › Free FOV selection on the scout image
- › Patient motion and position correction
- › Metal artefact removal (MAR)
- › AI-supported functions
- › Seamless integration
- › Implant planning tool integrated



Cleaning & Disinfection



U-son
› Available with 3 or 6 litre tank
› Stainless steel chamber
› Programmable temperature and timer
› Automatic 'de-gas' function

£ 535
3 liters RRP £824

£ 837
6 liters RRP £1,288



Soaking bath
Soaking bath for the immersion of dental instruments in the decontamination phase.

£ 40
RRP £48

BePro
Professional
Cleaning and
Disinfection



Click here

Disinfectant I
Powerful disinfectant concentrate for medical and dental instruments.

£ 41
Bottle 2l RRP £50

Disinfectant R
Ready-to-use solution for the disinfection and cleaning of rotary dental instruments.

£ 14
Bottle 1l RRP £17

£ 45
Refill 5l RRP £55

Disinfectant Wipes
120 wipes for the disinfection of non-invasive medical devices.

£ 10
Tube RRP £12

£ 8
Refill RRP £10

Disinfectant Spray
Surface disinfectant for the rapid disinfection and cleaning of non-invasive medical devices.

£ 14
Spray 1l RRP £17

£ 48
Refill 5l RRP £58

Sterilization, Hygiene & Maintenance



Internal cleaning & lubrication

Increase the life of your handpieces and reduce the costs of repairs!



Assistina Twin
› Two process chambers allowing loading of a 2nd handpiece while the 1st is processed
› Superior Maintenance & Lubrication
› Intuitive and easy to use
› A flexible range of adaptors available

£ 1,484
RRP £3,299

Assistina Twin and Assistina One includes Care set and Quick RM and ISO adaptors



Assistina One
Internal cleaning & lubrication of one handpiece in only 18 seconds, with Assistina One the aerosols stay in the HEPA filter thereby releasing clean air.

£ 1,051
RRP £2,235

Enjoy long term savings of 40% on handpiece lubrication with an Assistina

Accessories for sterilization



Lisa Safe Label Printer & QR Barcode Reader
Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Barcode reader £ 150 RRP £242

Lisa Safe printer £ 538 RRP £827



Osmo
Supplies multiple sterilizers simultaneously with top quality demineralized water.
› Automatic water filling
› Most cost-effective solution
› Increase sterilizer service life
› Up to 6 devices simultaneously

£ 954
RRP £1,590

ProService



PremiumCare

Gives you the flexibility to select a plan that best suits the individual needs of your practice.

The standard PremiumCare package offers a flexible plan allowing you to add optional extras as you please or the **PremiumCare Plus or Flex**** plan, which include fully comprehensive service packages and the required pressure vessel inspection and validation testing at no extra cost.

PremiumCare
Maintenance plan for all W&H sterilizer models, with a maximum limit of 4,000 cycles over the 3 year term*.

PremiumCare Plus & Flex**
A fully comprehensive maintenance plan for all W&H sterilizer models, including PVI certification, unlimited cycles and validation in line with manufacturer's or current HTM 01-05 guidance (model dependent).

Our Technical Service Team
are here to support you nationwide!



W&H
ProService

Services	Premium Care	Premium Care Plus & Flex**
All routine service visits for three years	✓	✓
All service parts and scheduled sundries	✓	✓
All breakdown calls	✓	✓
All breakdown parts including labour	✓	✓
All engineer travel in the UK Excluding Highlands, Islands and ferries etc,	✓	✓
Free on-site operator training	✓	✓
Free technical software upgrades	✓	✓
Unlimited Cycles - for the duration of the plan	⊗	✓
PVI Certification - in line with Pressure Systems Safety Regulations 2000:128	⊗	✓
Periodic Test Logbook	⊗	✓
Validation in line with manufacturer's or current HTM 01-05 guidance (model dependent)	⊗	✓
Installation and Commissioning in line with current guidelines	⊗	✓
Helix Test Kit***	⊗	✓
Extra Discount	⊗	✓

✓ Included ⊗ Not Included ⊕ Optional

*An additional cycle tariff will be levied if the 4,000 cycle limit is exceeded within the 3 year period ** PremiumCare Flex Plans are only available for certain models *** Only applies to B Type Sterilizers



Thermal washer
disinfector

The only system with two washing levels, for powerful cleaning, full flexibility, and faster turnaround in even the busiest clinics.



New!

Thed+

Includes Thed features, plus:

- › EliSense: device status indicator
- › Automatic rack and manifold injector recognition
- › Automatic program recommendation
- › Automatic adjustment of cycle time, water, process chemicals, and power consumption
- › Washing chamber with internal light

Offers include upper and lower rack

Thed+ Tabletop

£ 6,887

RRP £9,850

Thed+ Underbench

£ 7,505

RRP £10,468

Cleaner TE

Tri-enzymatic detergent with neutral pH. No neutraliser required!

£ 47

Canister 5l RRP £58



Rinse Aid

For improved drying phase.

£ 19

Bottle 1l RRP 23



Sterilization at its Best!

With a wide range of sterilizers, W&H is able to offer the perfect option to meet individual needs.



Click here



Lisa MINI 5 Litre B type sterilizer

- › Cycle times from 10 minutes
- › Integrated vacuum pump
- › Small and easily stackable: W44xH20xD67cm
- › Built in water tank and automatic water filling
- › For up to 22 instruments, complete traceability

£ 3,502

5 l. RRP £5,214

Special Offer Lisa Mini 5 liters & Assistina One/Assistina Twin

£ 4,420

RRP £7,449



ioDent® Cloud

Smart Digital Connection



Lisa B type sterilizer

- › Full extended traceability thanks to EliTrace
- › Mobile App for Real Time Remote Monitoring
- › Eco Dry+ technology automatically adapts the drying time to the mass of the load
- › Elisense user interface generates status, temperature and other cycle information
- › Cycle times from 13 minutes for specific loads
- › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use

£ 5,717

17 l. RRP £8,722

£ 6,479

22 l. RRP £9,868



Lara B type Sterilizer

- › Mobile App for Real Time Remote Monitoring
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradeability to increase: Performance and Traceability or to add a Fast Cycle
- › Available in 4 chamber capacities of 17 – 38 Litres

- Unique upgrade flexibility with **Activation Codes**
 - › Performance**
 - › Traceability
 - › Fast Cycle
 - › All in One** Only for Lara 17 and 22 litres

£ 4,790

17 l. RRP 7,063

£ 5,408

22 l. RRP £7,931

*with fast cycle upgrade

£ 5,923

Lara XL 28 l. RRP £8,747

£ 7,365

Lara XXL 38 l. RRP £10,682



Lyla S type Sterilizer
New concept sterilizer

- › Mobile App available via upgrade for remote data storage
- › Future-proof upgradeability to a type B sterilizer using the new Activation Code concept

- Unique upgrade flexibility with **Activation Codes**
 - › B Cycle
 - › Remote data storage
 - › Traceability
 - › All in One

£ 3,480

17 l. RRP £4,905

£ 3,860

22 l. RRP £5,159

*Based on average 2kg load for immediate use



Click here

- All W&H Sterilizers offer:
- › Intuitive user interface and touchscreen display for a smoother user experience
 - › Cycle traceability available thanks to Wi-fi connection and built in data logger
 - › Ergonomic design, with modular feet for flexible positioning
 - › Easy to clean surfaces and accessible water
 - › Automatic Water Fill function
 - › Pre-programmable vacuum and helix (for B Type Sterilizers) test function to run automatically out of surgery hours



Sales Team

Get in touch with us!
We will treat your request immediately.



Tom James

National Sales and
Marketing Manager

07876 527720

tom.james@wh.com



Richard Mansbridge

Area: London, South East, Anglia,
East Midlands and Humberside

07775 792628

richard.mansbridge@wh.com



David Barrass

Area: Scotland, Humberside,
North of England and North West

07795 315089

david.barrass@wh.com



Caroline Caine

Area: West Midlands, Manchester, Merseyside
South West England, Isle of Man, Channel
Islands, Wales, Northern Ireland and Eire

07786 515331

caroline.caine@wh.com

To place an order contact W&H UK on

+44 (0) 1727 874990

or contact your preferred dental dealer

W&H (UK) Ltd

6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ

office.uk@wh.com

www.wh.com



Click here

Seethrough Flex

Seethrough Flex stands for high
performance and small footprint.



Compact design
Intuitive operation
Face-to-face patient positioning
Easy to use
Flexibility
Impressive connectivity
Improved outcomes

seethrough

Offers valid from: 1st July to 31st December 2026. Information is correct at time of publication and relates to W&H products and services available in the UK and Ireland. Prices are quoted in £ sterling and exclude VAT. Savings are based on RRP and rounded to the nearest £. Free goods are not transferrable and no cash equivalent is available. Terms and Conditions apply. E. & O. E.