

#whdentalwerk



wh.com/na



Assistina Twin

Get the green light for more time, money and a healthier environment.



assistina^{TWIN}



- ✓ Thanks to smart process monitoring, the Assistina Twin always achieves thorough maintenance with minimum consumption.
- ✓ The Assistina Twin Care Set only needs to be changed after approximately 2,800 instruments (= 20 oil spray cans at 400 ml each).
- ✓ Your instrument maintenance is quick, reliable, cost-efficient, resource efficient and waste minimizing.
- ✓ Maintenance is also crucial to expand the lifespan of your handpieces and turbines, decreasing its carbon footprint over the total product lifecycle.



Advantage: 12 instruments in 2 min
thorough and precise oil service of up to 360 instruments per hour.



Advantage: HEPA filter
protects staff from harmful aerosols and prevents unpleasant odors in the practice air.



Advantage: less waste
less environmental harm and waste, more space in your practice.

Comparison of key advantages:

Assistina Twin versus oil spray cans

- 1 Faster:** 10 seconds per instrument
- 2 More efficient:** dual chambers for an uninterrupted workflow
- 3 More precise:** innovative oil nebulization for perfect cleaning and maintenance of even the farthest reaches of the instrument
- 4 Cost-effective:** considerably less expensive compared to oil spray cans
- 5 More flexible:** thanks to a variety of corresponding adaptors
- 6 Space saving:** no storage of bulky cans required
- 7 Safer:** no exposure to dangerous aerosol spray and no residue is left on hands
- 8 More environmentally friendly:** absolutely no emissions thanks to the HEPA filter and less harmful waste

Assistina Twin versus multiport systems

- 1 Faster:** 10 seconds per instrument, no unnecessary waiting times
- 2 More efficient:** dual chambers for an uninterrupted workflow
- 3 More precise:** innovative oil nebulization for perfect cleaning and maintenance of even the farthest reaches of the instrument
- 4 Cost-effective:** less expensive than other multiport systems
- 5 More user-friendly:** easy to operate, simple, tool-free cartridge replacement
- 6 More economical:** one Assistina Twin Care Set for up to 2,800 instruments
- 7 More practical:** instruments are available immediately
- 8 More space-saving:** Compact design with sliding window vs. drop down door


Innovative and reliable. The Assistina Twin is a device that has been developed based on W&H's **130+ years of experience** and produced according to the highest quality guidelines. The Assistina Twin comes with a 24-month warranty.

Tailored to your needs


Thanks to the integrated Quick Connect by W&H, the Assistina Twin fits into your processes more smoothly and makes them even more efficient. Why? Because you can equip the maintenance device with exactly the adaptors you need for your transmission instruments.



Quick RM
for turbine handpieces.



Quick ISO
for straight and contra-angle handpieces.



quick
connect
by W&H

Quick and simple adaptor changes optimize your processes and increase efficiency in your practice.

You decide which adaptors you actually need. Equip your Assistina Twin according to your needs and only use the right adaptors suitable for you in the two chambers.

W&H Impex Inc.
6490 Hawthorne Drive
Windsor, Ontario, N8T 1J9

t 1 800 265 6277 Within USA/Canada
t 1.519.944.6739 Outside USA/Canada
f 1.519.974.6121 **wh.com/na**