

#whdentalwerk



video.wh.com



Show casing W&H

Special offers January - April 2024



**New Synea
Power Edition**

For the hardest
materials.

Restoration & Prosthetics

Contra-angle Synea Power Edition

For the hardest materials

The Synea Power Edition is especially designed for heavy-duty applications, such as the removal of decayed materials, cavities, fillings, crowns, and bridges. They are the optimal solution for advanced ceramic material like **zirconia** that is increasingly being used in restorative dentistry.

synea



**Synea Fusion
WG-900 LT / WG-900 A**

- › Contra-angle 1:4
- › Optic / Without light
- › 3x Spray
- › Head Ø 10.7 mm

£ 720

RRP £1,445

£ 525

RRP £1,295



**Synea Vision
WK-900 LT**

- › Contra-angle 1:4
- › Optic
- › 4x Spray
- › Head Ø 10.7 mm

£ 849

RRP £1,699




Turbines

High speed, low weight

All W&H turbines have a hygienic head¹ anti-retraction system to prevent cross infection.

¹Exclusive W&H system patented in Europe and the USA



Our turbines are available to fit other manufacturers couplings.



**Synea Fusion
TG-97 / TG-97 L**

- › LED+ / Without light
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

£ 439

RRP £976

£ 505

RRP £1,123



**Synea Fusion
TG-98 / TG-98 L**


- › LED+ / Without light
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

£ 439

RRP £976

£ 505

RRP £1,123



**Synea Vision
TK-94 L**

- › LED+
- › 5x Spray
- › MICRO Head Ø 9 mm
- › 14 W/410,000 rpm
- › Roto Quick® Connection

£ 729

RRP £1,458



**Alegra
TE-97/98**

- › Without light
- › 3x Spray
- › Head Ø 10.4 / Ø 12.2 mm
- › 16 / 18 W
- › 390,000 / 330,000 rpm
- › Roto Quick® Connection

£ 239

RRP £551



**Synea Fusion
WG-99 LT / WG-99 A**

- › Contra-angle 1:5
- › Optic / Without light
- › 4x Spray
- › Head Ø 9.5 mm

£ 710

RRP £1,420

£ 508

RRP £1,269



**Synea Fusion
WG-56 LT / WG-56 A**

- › Contra-angle 1:1
- › Optic / Without light
- › Single Spray
- › Head Ø 9.5 mm

£ 449

RRP £1,123

£ 358

RRP £896



Alegra WE-56

- › Contra-angle 1:1
- › Without light
- › Optional External Spray
- › Head Ø 10.0 mm

£ 125

RRP £345




Alegra HE-43

- › Straight Handpiece 1:1
- › Without light
- › Optional External Spray

£ 120

RRP £300



Air motors & Couplings

The perfect connection



Air motors
Available from

£ 229

RRP £434



Couplings
Available from

£ 101

RRP £203

Prophylaxis & Periodontology

Proxeo Ultra

Scaling can be so easy!

The Proxeo Ultra is available with the unique **Q-Link quick** connection system. In addition, **handpieces are available to fit other manufacturers systems.** The tip changers with integrated torque control prevent overtightening and make changing the tip simple and safe.



Compatible with pacemakers
Tested by Graz University of Technology



GERMAN DESIGN AWARD WINNER 2020



Proxeo Ultra PB-530

- › Wireless foot control
- › Cleaning program for hose reprocessing
- › Compatibility with pacemakers
- › LED ring for ultra-bright illumination
- › Infinitely adjustable power setting
- › Large tank (750 ml)

£ **1,960**

RRP £3,266

Proxeo Ultra PB-520

- › Compatibility with pacemakers
- › LED ring for ultra-bright illumination
- › Infinitely adjustable power setting
- › Large tank (750 ml)

£ **1,255**

RRP £2,205

Proxeo Ultra PB-510


- › Compatibility with pacemakers
- › LED ring for ultra-bright illumination
- › Infinitely adjustable power setting

£ **903**

RRP £1,470

Full range of tips

The enhanced W&H piezo scaler tip range includes different tips for prophylaxis and periodontology, endodontics, implantology and restoration & prosthetics. The W&H range of tips were developed in close co-operation with end users. The range therefore fulfils the exact requirements of dental practice and implantology.



Tips from

£ **54**

RRP £98

Contact W&H for details



All Proxeo Ultra models are supplied with 3 tips (2x universal & 1x Perio)

Proxeo Air Scalers

Regular removal of calculus and plaque leads to a much greater sense of well-being and self-confidence amongst your patients.



Proxeo ZA-55

- › Without power regulation
- › Without light
- › Roto Quick® Connection
- › Tips 1AU, 2AU, 3AU, 1AP

£ **499**

RRP £1,249



Proxeo ZE-55 RM/BC

- › Without power regulation
- › Without light
- › Borden/Midwest Connection
- › 2 Tips 1AU

£ **341**

RRP £716

Polishing by rotation

Professional tooth cleaning is an effective method, giving your patients a beaming smile.



Proxeo WP-64 MU

- › Prophy contra-angle 4:1
- › Without light
- › For universal screw-in and snap on (adaptor) cups and brushes.

£ **230**

RRP £561



Proxeo HP-44 M

- › Handpiece 4:1
- › System Doriot
- › For disposable contra-angles with plastic shank

£ **199**

RRP £318



Dental Laboratory

Laboratory Devices

Good, better, perfecta

Based on decades of experience W&H has developed with the Perfecta 300 and Perfecta 600 a dental laboratory motor control unit which is available as a table control, knee control or foot control unit. In the Perfecta 600 laboratory handpiece the air-jet function is integrated, this is a world first for your dental laboratory.



Perfecta 600

- › Motor handpiece 7.5 Ncm
- › 1,000 - 50,000 rpm
- › Brushless motor
- › Tool Ø 2.35 mm
- › Air-jet function
- › Automatic: bistable/speed control operation
- › Forward/reverse drive

£ **1,392**

RRP £2,531
Foot

£ **1,552**

RRP £2,821
Knbee / Table

Perfecta 300

- › Motor handpiece 5.5 Ncm
- › 1,000 - 40,000 rpm
- › Brushless motor
- › Tool Ø 2.35 mm
- › Automatic: bistable mode
- › Forward/reverse drive

£ **1,025**

RRP £1,865
Foot

£ **1,185**

RRP £2,156
Knbee / Table

Oral Surgery & Implantology

Implantmed Classic SI-923

Easy to use, powerful and reliable

1
Year
Warranty



Implantmed SI-923

Non-optic/Cable

- › Ultra short motor with 5.5 Ncm torque
- › Speed range between 300 – 40,000 rpm
- › Precise torque between 5 – 70 Ncm
- › Automatic thread cutting function
- › Non-optic motor
- › Cable Foot Control

£ **2,724**


RRP £4,952

Implantmed Plus SI-1023

Workflow driven solution

W&H offer a synergistic solution for oral surgery guiding the clinician through the cycle from implementation to evaluation.

You can enjoy the benefits of utilising the modular Implantmed surgical motor and Piezomed piezosurgery unit alongside the Osstell Connect consultation tool and the Osstell Beacon, which ensures predictable assessment of ideal loading.




Implantmed SI-1023 Kit 01

Non-optic/Cable

- › Non-optic motor
- › Cable Foot Control

£ **3,082**

RRP £5,604




Implantmed SI-1023 Kit 03

Non-optic/Wireless

- › Non-optic motor
- › Wireless Foot Control

£ **3,455**

RRP £6,281



Implantmed SI-1023 Kit 02

Optic/Cable

- › Optic motor
- › Cable Foot Control

£ **3,613**

RRP £6,569



Implantmed SI-1023 Kit 04

Optic/Wireless

- › Optic motor
- › Wireless Foot Control

£ **3,985**


RRP £7,245

- Light
- No light
- Cable foot control
- Wireless foot control

New!

Piezomed Module
Plug & Play piezo surgery device for Implantmed Plus

Display menu and documentation in your local language



Motors EM-19 EM-19 LC
With and without light

Wireless foot control
for controlling multiple W&H devices

Straight & Contra-angle Handpieces



Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and avoid rejection in high-risk patients.

2
Years*
Warranty

Surgical Handpieces

High-quality stainless steel with a **scratch-resistant coating** makes the contra-angle handpieces robust and gives them a particularly long service life.



Contra-angle 20:1 WS-75 L / WS-75
› Mini LED+ / Without light
› External single spray
› Dismantleable

£ **852**

RRP £1,705

£ **779**

RRP £1,559



Contra-angle 1:2.7 WS-91 L / WS-91

- › Mini LED+ / Without light
- › 45° head with 3x spray
- › Dismantleable

£ **812**

RRP £1,625

£ **734**

RRP £1,468



Straight handpiece 1:1 S-11 L / S-11

- › Mini LED+ / Without light
- › External single spray
- › Dismantleable

£ **886**

RRP £1,772

£ **629**

RRP £1,259



Handpiece 20:1 Zygoma SZ-75

- › Without light
- › External single spray
- › Dismantleable

£ **1,260**

RRP £2,527



Self-generating LED
Surgical Handpieces available.

*Warranty only for Implantmed Plus SI-1023 kits.

Oral Surgery & Implantology

Piezomed Classic & Modules

"A real practice builder"

"I first purchased the Piezomed from W&H when I left general practice to commit to being a peripatetic implant surgeon in 2019. Since I began using the unit, I've noticed that surgical procedures that involve bone removal, block bone grafting and sinus surgery have become safer and more predictable. Patients also report less postoperative pain versus traditional methods. It has been a real practice builder and led to fewer complications which is so important when I am only occasionally travelling to some practices. More recently, I have become dependent on the kit for 'simple' extractions as the attachments have led to more atraumatic and faster tooth removal. It's a feature that I wish I had explored sooner!

I would say the biggest benefit is that the unit automatically detects the tips that are attached and changes the settings accordingly. During complicated implant surgeries I prefer not to have to worry about adjusting the settings. So, for extractions, sinus surgery, surgical extractions involving bone removal and block bone grafting, I just have to plug in and play. Ultimately, modern oral surgery shouldn't be performed without piezo bone removal. The ease of use of this unit in particular makes it so simple to work with – I can't recommend it enough."



Dr Aly Virani
North Cardiff Dental & Implants

**FREE
Extraction
Set
Included**

**Piezomed
Classic Module**

- Power 18 W
- Handpiece without light SA-40
- Instrument changer
- 6 Irrigation tubes

£ 3,086

RRP £4,748

**Piezomed
Plus Module**

- Power 24W
- Handpiece with light SA-40 L
- Automatic instrument detection
- Instrument changer
- Instrument library
- 6 Irrigation tubes

£ 3,902

RRP £6,004

**Piezomed Classic
Optic/Wireless**

- Piezomed Classic
- Handpiece with LED ring
- Wireless Foot Control

£ 5,842

RRP £9,736

Automatic
instrument
detection

**With FREE
Bone Set**
Worth £890

Piezomed Instruments

The W&H Piezomed instrument range includes different instruments for bone surgery, periodontology and extraction, for retrograde endodontics and for lateral sinus lift.

Tips compatible with all Piezomed units.

Bone surgery

The fine saw tooth design – more teeth per instrument – sets new standards in the field of bone preparation.



Periodontology and extraction

Instruments for thorough cleaning in periodontitis treatment and for atraumatic tooth extraction.



Retrograde endodontics

Diamond-coated, delicate instruments for easy access during retrograde root canal preparation.



Lateral sinus lift

The instruments permit safe lateral preparation of the maxillary sinus wall and gentle mobilization of the Schneiderian membrane.



Preparation of the implant site and crestal sinus lift

The preparation instruments are specially matched to the bone qualities predominant in the maxilla.



Osstell Beacon

Measurement of
implant stability

Osstell Beacon

Measure and evaluate the stability of the implant whilst monitoring its osseointegration.

£ 1,989

RRP £5,123

Osstell SmartPegs

The SmartPeg measures the resonance frequency when it is attached to an implant.

£ 83

RRP £113

What are the positive outcomes when using the Osstell Implant stability solutions?

- A positive experience for patients, encouraging trust through understanding
- Enhanced patient education builds confidence, resulting in increased treatment acceptance
- Reduction in treatment visits and failures with improved average healing times
- Improved outcomes for patients with risk factors
- Enhanced communication between specialists ensuring a more synergistic approach to treatment
- Improved treatment documenting


Osstell ISQ devices with the Osstell ISQ Scale help implantologists to make informed implant treatment decisions.

**Osstell Connect
FREE App for
Osstell Beacon**

Osstell Connect is designed to document quantified data which ensures traceability, whilst complimenting data in your Patient Management System.

Systems with graphical presentation of data and statistics which is great for patient communication. You can easily present treatments in a graphic and totally personalized way, even generating pdf reports, while having access to the world's largest online surgical portal for peace of mind. The patient consultation tool allows you to easily communicate your proposed implant.


*Register online at Osstell Connect to claim your free Smartpegs when you buy an Osstell Beacon



OSSTELL
A GUDIP COMPANY

Sterilization, Hygiene & Maintenance

Chain of Hygiene & Sterilization



Step 1 Cleaning & Disinfection

The new Teon washer disinfector is easy to use and offers an exceptional standard of cleaning, disinfection and reprocessing, maximising efficiency and increasing the effectiveness of the decontamination protocol.

- Teon thermal washer disinfector**
- › Validated cleaning of all instruments
 - › Traceability and connectivity through Ethernet module
 - › 1 Year Warranty

Teon with Standard Starter Set

£ 7,315
RRP £9,071

Teon with Handpiece Injector Starter Set

£ 8,738
RRP £10,257

BePro WDD Starter Pack

- £ 69
- › BePro WDD Cleaner TE – 5000 ml
 - › BePro WDD Rinse – 1000 ml
 - › Canister opener

New!

BePro Automated Reprocessing

Maximum performance in automated reprocessing. Discover the W&H BePro products for automated reprocessing of instruments, a holistic solution for thermal washer disinfectors like the W&H Teon+ or Teon.



BePro Cleaner TE
Tri-enzymatic non-foaming detergent for automatic cleaning in washer disinfector.

£ 44
RRP £54



BePro Rinse
For rapid, stain-free drying during automated reprocessing of dental instruments.

£ 14
RRP £17



Step 2 Internal cleaning & lubrication

Increase your handpieces lifetime and reduce the costs of repairs!



quick connect
by W&H

Assistina Twin

- › Two process chambers allowing loading of a 2nd handpiece while the 1st is processed
- › Superior Maintenance & Lubrication
- › Intuitive and easy to use
- › A flexible range of adaptors available

£ 1,554
RRP £3,108



quick connect
by W&H

Assistina One

Internal cleaning & lubrication of one handpiece in only 18 seconds, with Assistina One the aerosols stay in the HEPA filter thereby releasing clean air.

£ 1,100
RRP £2,105

Assistina One includes
Care set and Quick RM and ISO adaptors



Service Oil F1
Ask your dealer about W&H Handpiece Oil's



Step 3 Packaging

Seal2: Pack it!



Step 4 Sterilization

Osmo RO System: Always with high quality demineralized water!



Step 5 Traceability

Store and follow it!



Sterilization, Hygiene & Maintenance

Sterilization at its Best!

Saturated steam under pressure, is by far the most reliable medium, known for the destruction of all forms of microbial life. With the wide range of sterilizers W&H is able to offer the perfect sterilizer to meet individual needs.






All W&H Sterilizers offer:

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Ergonomic design, with modular feet for flexible positioning
- › Easy to clean surfaces and accessible water
- › Automatic Water Fill function
- › Pre-programmable vacuum and helix (for B Type Sterilizers) test function to run automatically out of surgery hours



Lyla S type Sterilizer New concept sterilizer

- › Mobile App available via upgrade for remote data storage
- › Future-proof upgradability to a type B sterilizer using the new Activation Code concept



Unique upgrade flexibility with **Activation Codes**

- › B Cycle
- › Remote data storage
- › Traceability
- › All in One

£ **3,318**

17 l. RRP £4,608

£ **3,696**

22 l. RRP £4,847



Lara B type Sterilizer

- › Mobile App for Real Time Remote Monitoring
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradeability to increase: Performance and Traceability or to add a Fast Cycle



Unique upgrade flexibility with **Activation Codes**

- › Performance
- › Traceability
- › Fast Cycle
- › All in One

£ **4,559**

17 l. RRP £6,642





£ **5,121**

22 l. RRP £7,461



Lara XL B type Sterilizer

- › With its 28 litre chamber, Lara XL is able to reprocess up to 7 kg of wrapped instruments or 9 kg containers
- › Mobile App for Real Time Remote Monitoring
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradeability to increase: Traceability or to add a Fast Cycle



Unique upgrade flexibility with **Activation Codes**

- › EliTrace
- › Fast Cycle
- › All in One

£ **5,648**

28 l. RRP £8,230



Lisa B type sterilizer

- › Full extended traceability thanks to EliTrace
- › Mobile App for Real Time Remote Monitoring
- › Eco Dry+ technology automatically adapts the drying time to the mass of the load
- › Elisense user interface generates status, temperature and other cycle information
- › Cycle times from 13 minutes for specific loads
- › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use

£ **5,467**

17 l. RRP £8,205

£ **6,178**

22 l. RRP £9,286

*Based on average 2kg load for immediate use

*with fast cycle upgrade

*with fast cycle upgrade for a 2 kg load

Sterilization, Hygiene & Maintenance

The best accessories for your sterilizer



Lisa Safe Label Printer & QR Barcode Reader
Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Barcode reader Lisa Safe printer
£ 142 £ 499
RRP £227 RRP £777



Osmo
Supplies multiple sterilizers simultaneously with top quality demineralized water.
• Automatic water filling
• Most cost-effective solution
• Increase sterilizer service life
• Up to 6 devices simultaneously
• Water free from minerals and contaminants

£ 899
RRP £1,498

Seal2 Starter Pack
• Universal Pouch Sealing System
• 12 mm weight sealing
• 3 x EliReel Rolls FREE

£ 656
RRP £1,208

Seal Check
• Pack of 100 sealing check
£ 29
RRP £60

W&H EliReel sealing pouches
5.5cm RRP £21 £13
7.5cm RRP £29 £17
10cm RRP £37 £22
£ 13 Per Roll
From RRP £21

BePro Professional Cleaning and Disinfection for All Surfaces



Soaking bath
Soaking basin for the immersion of dental instruments in the decontamination phase.

£ 37
RRP £44



Disinfectant R 1l.
Ready-to-use solution for the disinfection and cleaning of rotary dental instruments.

£ 13
RRP £15

Disinfectant I 2l.
Powerful disinfectant concentrate for medical and dental instruments.

£ 38
RRP £46

Disinfectant Spray 1l.
Surface disinfectant for the rapid disinfection and cleaning of non-invasive medical devices.

£ 13
RRP £15



Disinfectant Wipes
120 wipes for the disinfection of non-invasive medical devices.

£ 9 £ 7
Tube RRP £10 Refill RRP £8

PremiumCare

Gives you the flexibility to select a plan that best suits the individual needs of your practice.

The standard PremiumCare package offers a flexible plan allowing you to add optional extras as you please or the **PremiumCare Plus or Flex**** plan, which include fully comprehensive service packages and the required pressure vessel inspection and validation testing at no extra cost.

PremiumCare
Maintenance plan for all W&H sterilizer models, with a maximum limit of 4,000 cycles over the 3 year term*.

PremiumCare Plus & Flex**
A fully comprehensive maintenance plan for all W&H sterilizer models, including PVI certification, unlimited cycles and validation in line with manufacturer's or current HTM 01-05 guidance (model dependent).

Our Technical Service Team are here to support you nationwide!



W&H ProService

Services	Premium Care	Premium Care Plus & Flex**
All routine service visits for three years	✓	✓
All service parts and scheduled sundries	✓	✓
All breakdown calls	✓	✓
All breakdown parts including labour	✓	✓
All engineer travel in the UK Excluding Highlands, Islands and ferries etc,	✓	✓
Free on-site operator training	✓	✓
Free technical software upgrades	✓	✓
Unlimited Cycles - for the duration of the plan	⊖	✓
PVI Certification - in line with Pressure Systems Safety Regulations 2000:128	⊖	✓
Periodic Test Logbook	⊖	✓
Validation in line with manufacturer's or current HTM 01-05 guidance (model dependent)	⊗	✓
Installation and Commissioning in line with current guidelines	⊗	✓
Helix Test Kit***	⊗	✓
Extra Discount	⊗	✓

Optional
⊖
Not Included
⊗
Included
⊙

*An additional cycle tariff will be levied if the 4,000 cycle limit is exceeded within the 3 year period ** PremiumCare Flex Plans are only available for certain models *** Only applies to B Type Sterilizers

Meet the Sales Team



Jon Bryant

National Sales Manager

07917 886326

jon.bryant@wh.com



Richard Mansbridge

Area: London, South East, Anglia,
East Midlands and Humberside

07775 792628

richard.mansbridge@wh.com



Caroline Caine

Area: West Midlands, Manchester, Merseyside
South West England, Isle of Man, Channel Islands,
Wales, Northern Ireland and Eire

07786 515331

caroline.caine@wh.com



David Barrass

Area: Scotland, Humberside,
North of England and North West

07995 315089

david.barrass@wh.com

W&H would like to introduce David Barrass who has joined the Sales Team as Business Development Manager looking after Scotland and the North East of England. David has worked in the dental industry for 8 years, most recently as Sales Manager at Trigiene Dental, specializing in handpieces and small equipment. He has a wealth of experience in this area and is looking forward to working closely with his customers and local dealer representatives to ensure that they continue to receive a consistently high level of service, advice, and support. David is keen to maximize W&H's profile across Scotland to provide customers with the information and support they need to make informed purchasing decisions based on product awareness.

Outside of work David enjoys long hikes with his labradoodle Bernie, skiing holidays in the winter and beach holidays in the summer, with ticking off visits to all the Greek Islands being a high priority on his bucket list. To relax, David is happy to sit in front of the TV especially if he can watch Top Gear or Friends.



"Suits our needs"

Dr Stephen Young

Dentist

"We invested in a Lisa sterilizer from W&H based on our need for a vacuum autoclave. From the Lisa's touchscreen display panel, we can select a setting that suits our needs from a number of different sterilizing cycles. We would recommend the Lisa sterilizer as it's easy to use. The best feature of this system is that it automatically data logs all cycles, so there is no need to print out information on every cycle and file it for future reference. Moreover, we have always received very good customer service from W&H and if we have needed an engineer, we haven't had to wait long for one to visit our practice."



"Exceptional service from W&H"

Dr Rajen Vadher

Quadrant Dental Care

"I have been using handpieces from W&H since the turn of the century and have always found them reliable and easy-to-use. By chance, one day I was sent a pre-paid envelope by W&H that allowed me to receive a free quote/inspection before repair on the handpieces. W&H very quickly emailed me back regarding the handpieces, providing a detailed fault analysis and their recommendations. They followed this up with very useful quotes for exchange and repair, allowing me to invest in some new handpieces which I look forward to lasting me until retirement! I would definitely recommend W&H to other professionals – the customer care is exceptional and the team is always happy to help."

W&H (UK) Ltd

6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ

t +44 (0) 1727 874990

f +44 (0) 1727 874628

office.uk@wh.com

www.wh.com

**To place an order
contact W&H UK on
+44 (0) 1727 874990
or contact your preferred
dental dealer**

Offers valid until 30th April 2024. Information is correct at time of publication and relates to W&H products and services available in the UK and Ireland.

Prices are quoted in £ sterling and exclude VAT.

Savings are based on RRP and rounded to the nearest £.

Free goods are not transferrable and no cash equivalent is available.

Terms and Conditions apply. E. & O. E.