

## PRODUCT DATA SHEET

# FINE TIP PASTETTE®

### SPECIALIST TRANSFER PIPETTES FOR SMALL DROP WORK AND DIAGNOSTIC PROCEDURES.

With a soft responsive bulb, Fine Tip Pastettes provide a high level of control and dexterity for quick and easy use with exacting procedures.

Frequently referred to in National Standard Methods and NHS SOP, Fine Tip Pastettes are used in many laboratory and diagnostic procedures as an effective solution for specific liquid handling tasks.

#### Ideal for:

- Careful separation of supernatant in centrifuge tubes
- Plasma separation during blood chemistry analysis
- Complete aspiration of media from microwell plates
- Immuno-staining on slides and *in-situ* hybridisation
- Microbiology for cell staining, reagent handling and aseptic transfers and colony picking
- Point of care blood glucose monitoring for drop-wise application to test strips

#### Sterile Packs

Sterile packs are securely sealed and e-beam irradiated to give a sterility assurance level of 10<sup>-4</sup> to provide Pastettes suitable for aseptic transfers. Certificates of sterility are available by lot number. A choice of sterile pack sizes provides cost effective quantities to suit your sterile procedures.

#### Varied Designs

Fine Tip Pastettes are available in a choice of forms with alternative drop volumes, draw volumes and lengths to suit various tasks; from the Extended Fine Tip at 157mm in length to the Small Bulb Fine Tip at just 36mm long.

*See reverse for full range of Fine Tip Pastettes.*



	Extended Fine Tip Standard	Extended Fine Tip	Fine Tip Maxi	Fine Tip	Extended Fine Tip Mini	Extended Fine Tip Ultra Micro	Fine Tip Ultra Micro	Small Bulb Fine Tip
Length (mm)	157	153	147	147	104	63	51	36
Drops per ml	50	50	45	55	50	45	45	40

**Bulb Draw (ml)**

3.5   3.3   4.8   3.3   1.0   1.0   1.0   0.5

**Pastette® Pluses**

Fine Tips for Great Control

No Sharps, No Hazards

Known drop volume

**samples available**

**Buy in BULK!**  
Save even more!

All drawings depict actual size

**Drop Volume (µl)**

## Are you using genuine Pastettes®?

Often imitated, but never bettered, Alpha Laboratories offer the most comprehensive range of disposable transfer pipettes available to suit all clinical, experimental and routine liquid handling procedures.

Pastettes are manufactured from virgin low-density polyethylene (LDPE), safely providing precise transfer and delivery of liquids, reagents and cell suspensions. Graduated, Non-graduated, Fine Tip and Micro ranges feature designs specific for your applications.

Pastettes are ideally suited to sterile and non-sterile applications. Available in sterile packs of 10's, 20's or single wrapped as well as non-sterile bulk packaging.

Request the Alpha Laboratories' catalogue to see our complete range or request samples to try before you buy.

Non Sterile									Non-sterile
Cat No	LW4635	LW4233	LW4234	LW4060	LW4231	LW4244	LW4243	LW4697	
Pack Size	500	250	400	500	400	500	500	500	
Sterile 1's									1's
Cat No	LW4636	LW4238	LW4235	LW4061	LW4239				
Pack Size	400	400	400	400	400				
Sterile 20's									20's
Cat No	LW4638	LW4232	LW4236	LW4070	LW4237				
Pack Size	500	500	400	500	500				