


## DNA SEQUENCING UNIT

### VERTICAL GEL SYSTEMS



The Clarit-E DNA Sequencing Unit offers advanced features for enhancing gel resolution and ease of use, essential when handling gels of this size. Each unit contains ultra soft silicone seals for easy plate sealing and trouble-free runs, even over extended run times.

#### KEY FEATURES:

The Clarit-E DNA Sequencing Unit is ideal for a variety of large format vertical gel applications.

- α Resolution is enhanced by using an aluminium heat sink plate, which is essential for even sample migration.
- α Added convenience is provided by a removable lower buffer tank and upper buffer drainage tap.



Fan heater-sensor kit  
For the cooling down of gels

#### Tanks:

Product Code	Description
EL2560	20x50cm Sequencing Unit
EL2580	33x35cm Sequencing Unit

#### Glass Plates:

Product Code	Description
EL2562	Notched for 20cm Sequencing Unit
EL2564	Plain for 20cm Sequencing Unit
EL2582	Notched for 33cm Sequencing Unit
EL2584	Plain for 33cm Sequencing Unit

## Accessories:

Product Code	Description
EL2566	0.25mm Spacers for 20cm Sequencing Unit
EL2568	0.35mm Spacers for 20cm Sequencing Unit
EL2570	1mm Spacers for 20cm Sequencing Unit
EL2572	1.5mm Spacers for 20cm Sequencing Unit
EL2586	0.25mm Spacers for 33cm Sequencing Unit
EL2588	0.35mm Spacers for 33cm Sequencing Unit
EL2590	1mm Spacers for 33cm Sequencing Unit
EL2592	1.5mm Spacers for 33cm Sequencing Unit
EL2594	Fan Heater Sensor Kit


## Compatible Power Supplies and Adapters:

Product Code	Description
EL2707	Mini Power Supply 500V, 400mA, 120W
EL2705	Mini-300 Power Supply, 300V, 400mA, 60W
EL2702	Power Pack Adapters 4mm To 2mm
EL2752	Maxi Power Pack 300V, 3A, 300W
EL2732	Midi Power Pack 300V, 700mA, 150W
EL2714	Replacement Cables

## Compatible Combs:

Product Code	Description
EL2560-24-0.25	Comb 24 Sample, 0.25mm
EL2560-48-0.25	Comb 48 Sample, 0.25mm
EL2560-24-0.35	Comb 24 well, 0.35mm
EL2560-48-0.35	Comb 48 well, 0.35mm
EL2560-24-1	Comb 24 Sample, 1mm
EL2580-48-0.25	Comb 48 Sample, 0.25mm
EL2580-96-0.25	Comb 96 Sample, 0.25mm
EL2580-48-0.35	Comb 48 Sample, 0.35mm
EL2580-96-0.35	Comb 96 Sample, 0.35mm
EL2580-48-1	Comb 48 Sample, 1mm
EL2580-80-1	Comb 80 Sample, 1mm

# Clarit-E®



Visit [www.alphalabs.co.uk/electrophoresis](http://www.alphalabs.co.uk/electrophoresis) to find out more.