


MIDI

HORIZONTAL GEL SYSTEMS



The Clarit-E Midi is Designed for routine checks of a medium number of samples over short to medium gel lengths. They are ideal for quick checks of samples from PCR and cloning.

KEY FEATURES:

The compact Clarit-E-Midi system is perfectly sized for intermediate workloads, whilst still offering economical gel and buffer usage.

- α Available with 10 x 7cm, 10 x 10cm or with both gel trays
- α Run up to 100 samples
- α Low buffer volumes
- α Ideal for rapid electrophoresis



Casting Dams
Cast gels separately from the gel tank and save time while the tank is in use.



Buffer Saver Blocks
Save costs by physically reducing the volume of buffer required in the gel chamber.

Tanks:

Product Code	Description
EL1200	Gel Tank 10x7cm UV Tray
EL1250	Gel Tank 10x7 & 10x10cm UV Trays
EL1300	Gel Tank 10x10cm UV Tray

Trays:

Product Code	Description
EL1202	UV Tray 10x7cm
EL1302	UV Tray 10x10cm

Accessories:

Product Code	Description
EL1204	Casting Dams
EL1206	Electrode, positive
EL1208	Electrode, negative
EL1210	Loading Guides
EL1212	Viewing Platform
EL1214	Cool Pack & Platform
EL1216	Buffer Saver Blocks
EL1220	10cm UV Gel Scoop
EL1018	Flexi Caster
EL1026	Gel Levelling Table
EL1222	Lid
EL1224	Tank
EL1230	Horizontal Electrophoresis Package

Clarit-E®

Compatible Power Supplies and Adapters:

Product Code	Description
EL2707	Mini Power Supply 500V, 400mA, 120W
EL2705	Mini-300 Power Supply 300V, 400mA, 60W
EL2702	Power Pack Adapters 4mm to 2mm
EL2752	Maxi Power Pack 300V, 3A, 300W
EL2732	Midi Power Pack 300V, 700mA, 150W
EL2714	Replacement Cables

Compatible Combs:

Product Code	Description
EL1200-1-0.75	Comb Prep 1 Well, 1 Marker, 0.75mm
EL1200-1-1	Comb Prep 1 Well, 1 Marker, 1mm
EL1200-1-1.5	Comb Prep 1 Well, 1 Marker, 1.5mm
EL1200-1-2	Comb Prep 1 Well, 1 Marker, 2mm
EL1200-2-0.75	Comb Prep 2 Well, 2 Marker, 0.75mm
EL1200-2-1	Comb Prep 2 Well, 2 Marker, 1mm
EL1200-2-1.5	Comb Prep 2 Well, 2 Marker, 1.5mm
EL1200-2-2	multiSUB Midi Comb, Prep 2, Marker 2, 2mm
EL1200-4-0.75	Comb Prep 4 Well, 2 Marker, 0.75mm
EL1200-4-1	Comb Prep 4 Well, 2 Marker, 1mm
EL1200-4-1.5	Comb Prep 4 Well, 2 Marker, 1.5mm
EL1200-4-2	Comb Prep 4 well, 2mm
EL1200-8-0.75	Comb 8 Well, 0.75mm
EL1200-8-1	Comb 8 Well, 1mm
EL1200-8-1.5	Comb 8 Well, 1.5mm
EL1200-8-2	Comb 8 Well, 2mm
EL1200-10MC-0.75	Comb 10MC Well, 0.75mm
EL1200-10MC-1	Comb 10MC Well, 1mm

Product Code	Description
EL1200-10MC-1.5	Comb 10MC Well, 1.5mm
EL1200-10MC-2	Comb 10 Well, 2mm
EL1200-12-0.75	Comb 12 Well, 0.75mm
EL1200-12-1	Comb 12 Well, 1mm
EL1200-12-1.5	Comb 12 Well, 1.5mm
EL1200-12-2	Comb 12 Well, 2mm
EL1200-16-0-75	Comb 16 Well, 0.75mm
EL1200-16-1	Comb 16 Well, 1mm
EL1200-16-1.5	Comb 16 Well, 1.5mm
EL1200-16-2	Comb 16 Well, 2mm
EL1200-20MC-0.75	Comb 20 Well, 0.75mm
EL1200-20MC-1	Comb 20MC Well, 1mm
EL1200-20MC-1.5	Comb 20MC Well, 1.5mm
EL1200-20MC-2	Comb 20 Well, 2mm
EL1200-25-0.75	Comb 25 Well, 0.75mm
EL1200-25-1	Comb 25 Well, 1mm
EL1200-25-1.5	Comb 25 Well, 1.5mm
EL1200-25-2	Comb 25 Well, 2mm

Visit www.alphalabs.co.uk/electrophoresis to find out more.