


MINI

HORIZONTAL GEL SYSTEMS



The Clarit-E Mini is the smallest unit in the range, designed to maximise run economy without compromising versatility.

KEY FEATURES:

The Clarit-E Mini is designed for quick checks of low to medium numbers of samples.

- α Available with 7 x 7cm, 7 x 10cm or with both gel trays
- α Economic low gel and buffer volumes
- α Small lab bench footprint



Casting Dams

Cast gels separately from the gel tank and save time while the tank is in use.



Buffer Saver Blocks

Save costs by physically reducing the volume of buffer required in the gel chamber.

Tanks:

Product Code	Description
EL1000	Gel Tank 7x7cm UV Tray
EL1050	Gel Tank with 7x7 and 7x10cm UV Tray
EL1100	Gel Tank 7x10cm UV Tray


Trays:

Product Code	Description
EL1002	UV Tray 7x7cm
EL1102	UV Tray 7x10cm

Accessories:

Product Code	Description
EL1004	Casting Dams
EL1006	Electrode, positive
EL1008	Electrode, negative
EL1010	Loading Guides
EL1012	Viewing Platform
EL1014	Cool Pack & Platform
EL1016	Buffer Saver Blocks
EL1020	7cm UV Gel Scoop
EL1018	Flexi Caster
EL1026	Gel Levelling Table

Clarit-E[®]



Compatible Power Supplies and Adapters:

Product Code	Description
EL2707	Mini Power Supply 500V, 400mA, 120W
EL2705	Mini-300 Power Supply, 300V, 400mA, 60W
EL2702	Power Pack Adapters 4mm To 2mm
EL2752	Maxi Power Pack 300V, 3A, 300W
EL2732	Midi Power Pack 300V, 700mA, 150W
EL2714	Replacement Cables

Compatible Combs:

Product Code	Description
EL1000-1-0.75	Comb Prep 1 Well, 1 Marker, 0.75mm
EL1000-1-1	Comb Prep 1 Well, 1 Marker, 1mm
EL1000-1-1.5	Comb Prep 1 Well, 1 Marker, 1.5mm
EL1000-1-2	Comb Prep 1 Well, 1 Marker, 2mm
EL1000-2-0.75	Comb Prep 2 Well, 2 Markers, 0.75mm
EL1000-2-1	Comb Prep 2 Well, 2 Markers, 1mm
EL1000-2-1.5	Comb Prep 2 Well, 2 Markers, 1.5mm
EL1000-2-2	Comb Prep 2 Well, 2 Markers, 2mm
EL1000-4-0.75	Comb Prep 4 Well, 2 Markers, 0.75mm
EL1000-4-1	Comb Prep 4 Well, 2 Markers, 1mm
EL1000-4-1.5	Comb Prep 4 Well, 2 Markers, 1.5mm
EL1000-4-2	Comb Prep 4 Well, 2 Markers, 2mm
EL1000-8-0.75	Comb 8 Well, 0.75mm
EL1000-8-1	Comb 8 Well, 1mm

Product Code	Description
EL1000-8-1.5	Comb 8 Well, 1.5mm
EL1000-8-2	Comb 8 Well, 2mm
EL1000-10-0.75	Comb 10 Well, 0.75mm
EL1000-10-1	Comb 10 Well, 1mm
EL1000-10-1.5	Comb 10 Well, 1.5mm
EL1000-10-2	Comb 10 Well, 2mm
EL1000-12-0.75	Comb 12 Well, 0.75mm
EL1000-12-1	Comb 12 Well, 1mm
EL1000-12-1.5	Comb 12 Well, 1.5mm
EL1000-12-2	Comb 12 Well, 2mm
EL1000-16-0.75	Comb 16 Well, 0.75mm
EL1000-16-1	Comb 16 Well, 1mm
EL1000-16-1.5	Comb 16 Well, 1.5mm
EL1000-16-2	Comb 16 Well, 2mm

Visit www.alphalabs.co.uk/electrophoresis to find out more.